



SiriusT3

From Pion Inc.

For compound screening and
detailed PhysChem characterisation

SiriusT3 measures pK_a , log P, log D and solubility of ionisable drugs and small molecules, needing little amounts of sample. SiriusT3 instruments are used in labs that support drug discovery, crop science research, by universities, biochem and CROs.



Five miniature probes entering sample vial

Use it for

Early stage compound screening


- In synthetic chemistry, measuring pK_a helps in selecting the best structures for development
- Good compound selection reduces the risk of costly late stage failures
- Up to 80 pK_a samples a day; automation reduces operator time

PhysChem characterisation

- Understand the ionisation state of a drug, how it varies with pH and its effect on absorption
- Learn how the pK_a changes due to solvent interactions, partitioning, precipitation or binding
- Improve ADME predictions and evaluate how changes to the molecule can affect the log P for a desired application.

Pre-formulation

- Preparing detailed PhysChem profiles
- Learn about the molecule's capacity for supersaturation and LLPS (liquid-liquid phase separation) formation



Spectral data vs. pH for sample with two pK_a s

Importance of pK_a

Most small-molecule drugs ionize in aqueous solutions. Ionization occurs when the molecules acquire or lose hydrogen ions in response to changes in the pH. When pH goes up (e.g. in the GI Tract), bases convert from ionised to neutral and become more lipophilic and less soluble. These changes affect the rate of absorption. pK_a values affect many other properties, too. Chemists often improve drug candidates by making substitutions to change their pK_a .

Screening

SiriusT3 is ideal for high throughput pK_a measurement in support of drug discovery. In most labs, 80% of samples will succeed first time by a Fast UV/cosolvent method. Failures can be re-measured by pH-metric, or flagged for future attention if the compound gets through DMPK.

What customers say

"SiriusT3 integrates perfectly in our analytical workflow".

"It is a good fit, filled the gap we had in the company".

"...a pre-requisite for DMPK studies".

"It has become a part of our screening cascade to aid in selection of compounds to move forward into further PK testing".

"The SiriusT3 is an integral part of a pool of instruments for physicochemical compound characterisation. For pK_a and log P/D measurements, SiriusT3 is the Golden Standard".




Measuring pK_a

If a molecule has a pK_a between 2 and 12, we will find it. Many drugs absorb UV, and the intensity vs. wavelength changes when they ionise. Our UV-metric methods determine pK_a values from these changes. Around 80% of drug pK_a values can be successfully measured in about 15 minutes by our Fast UV method, which requires only 5 μ L of 10 mM stock solution (typically 10 μ g). Those pK_a s that are not associated with UV change can be measured by our pH-metric method, which requires less than 1 mg of solid sample.

Prediction and measurement

Scientists often use software to predict pK_a values. It's powerful and fast, but it doesn't always get things right. This graph shows measured vs. predicted pK_a for 115 pK_a s measured in our Applications Laboratory in 2013. We're confident in the measured values, but it would be risky to base development strategies on some of the predictions



log P and solubility – find out more

Pion pH-metric methods for measuring log P and solubility of ionisable compounds are supported by several validation studies. Methods are described on www.pion-inc.com, along with recordings of webinars about pK_a , log P and solubility.

TECHNICAL SPECIFICATION

SIRIUST3 TECHNICAL SPECIFICATION	
pH electrode	Ag/AgCl, double junction reference
pH range	1.8 - 12.2, standardised by Sirius Four-Plus™
Electrode storage	Home position for electrode storage pH7 buffer positions for calibration
Stirrer	Overhead, variable speed, computer controlled
Temperature	Measured by thermocouple at every datapoint. Controlled by Peltier, range 12°C to 70°C
Turbidity detection	Pass-through light detection – percentage light detected reported to 0.1%
Reagents	Water, Acid, Base, partition solvent, CoSolvent
Partition solvents	Supports 7 including Octanol, Dodecane, Toluene
CoSolvents	Supports 10 including Methanol, Acetonitrile, DMSO and MDM. Six can be installed, via automated valve.
Dispensers	Precision dispensers, 0.5 mL syringes, multi-tip capillary bundle, minimum volume of 0.024 L
Probe movement	Fully automated X-Y axis on titrator module
Washes	Wash station provides static and flowing washes
Autoloader	Holds up to 192 samples, automatic gripper arm for transfer to sample position
Purge gas	Two internal flow meters, nitrogen or argon supply required
Assay volume	1.0 to 3.5 mLs
Environment	For indoor use only Altitude ≤2000m Temp. 5°C to 40°C Max relative humidity 80% @ 31°C Mains voltage fluctuation +/- 10% Installation overvoltage Cat II Pollution degree 2

SIRIUST3 ASSAY SPECIFICATION	
Analytical methods	Templates supplied for all supported assays. Users can save customised templates if required. 1-2-3 template for setting up cosolvent pKa assays.
pK _a	Determine pK _a values to the highest commercially available accuracy <ul style="list-style-type: none"> • Results in 15 minutes by Fast UV for most samples • UV-metric method for pK_as range of 0.5 to 13.5 • pH-metric assays for samples without pH/UV activity • Cosolvent methods for poorly soluble samples • Yasuda-Shedlovsky plots for cosolvent pK_as • Generate distribution of species plots • Determine buffer capacity Standard assay conditions: 0.15 M KCl, 25 or 37°C. Sample requirement: UV, 10 µg; pH-metric, 1 mg
log P	<ul style="list-style-type: none"> • Determine log P and log D vs. pH • Study octanol/water & many other mixtures • Built-in methods for high, medium and low log P • No need for time-consuming shake-flask style experiments in multiple pH buffers. • Sample requirement: 1 - 2 mg
Solubility	<ul style="list-style-type: none"> • Determine intrinsic solubility and kinetic solubility • Derive pH vs. solubility profile • Determine supersaturation • Sample requirement: 1 -10 mg, depends on solubility
Assay Expert	Optional integrated software prediction module (powered by ACD) to analyse structure and identify ionizable groups, identify whether they are acidic or basic and estimate pK _a and log P values. Predicted properties are used to automatically optimize the design of experiments.
Reagents	<ul style="list-style-type: none"> • Minimal running costs - uses off the shelf reagents with no special "kits" required.
Getting samples to dissolve	<ul style="list-style-type: none"> • Built in ultra-sonic bath to aid measurements of insoluble compounds. Automated acid and base pre-dose can be programmed. Short automated procedures for cleaning electrode and probes.
Clean-up assays	Low to high and high to low pH with cosolvent present; surfactant wash at end.

PHYSICAL DIMENSIONS	WEIGHT	HEIGHT	WIDTH	DEPTH
DISPENSER MODULE	35kg/77.16lbs	700mm	260mm	460mm
TITRATOR	23kg/50.71lbs	490mm	215mm	460mm
AUTOLOADER	27kg/59.52lbs	490mm	350mm	460mm
TOTAL	85kg/187.39lbs	700mm	825mm	460mm

Please contact your appropriate region below to find out more about Pion products.

Pion stands behind the science. Let us prove it to you. Contact us today.

North America | Pion Inc. | 10 Cook St. | Billerica, MA 01821 | USA | +1-978-528-2020 | sales@pion-inc.com | www.pion-inc.com
Rest of World | Pion Inc (UK) Ltd | Forest Row Business Park | Station Road Forest Row | East Sussex RH18 5DW

