

## PREDICTABLY SAFE | DISTINCTLY INTUITIVE | CONSISTENTLY EFFECTIVE



Your Sound Calcium Strategy.

# Unique Mechanism of Action to Disrupt Calcium



The IVL catheter is delivered across a calcified lesion over an 0.014" wire and the integrated balloon is expanded to 4atm to facilitate efficient energy transfer



An electrical discharge from the emitters vaporizes the fluid within the balloon, creating a rapidly expanding & collapsing bubble that generates sonic pressure waves



The waves create a localized field effect that travels through soft arterial tissue, selectively cracking superficial and deep calcium within the vessel wall



After calcium modification, the integrated balloon may subsequently be used to dilate the lesion at low pressure in order to maximize luminal gain prior to stent placement

**0.96**mm 2.1% 0% 98% **IVL Crossing & Severe Dissections Abrupt Closure Calcium Length** Average Max **Calcium Thickness** Post-IVL Post-IVL Therapy Delivery **6.5**mm<sup>2</sup> 99% 02% 11.9% **Stent Delivery Stent Expansion** MSA **Residual Stenosis** Acute Gain

# Intravascular Lithotripsy System for Cardiovascular Calcium



### **IVL Generator**

- Portable and rechargeable
- No external connections
- Quick & easy setup with no settings

### **IVL Connector Cable**

- Simple magnetic connections
- Push-button activated

## SHOCKWAVE C<sup>2</sup> IVL Catheter

- Standard interventional technique
- Rapid exhange system
- 0.014" guidewire of choice





384 Patients at 47 Sites





0% No Reflow Post IVL

# **IVL MOA Under OCT**

Multi-plane & Longitudinal Calcium Fractures

#### **PRE-PROCEDURE**



SEVERE CALCIUM DETECTED<sup>1</sup>

#### **POST-IVL**



TRANSMURAL CALCIUM FRACTURES<sup>1</sup>

#### **POST-STENT**



STENT EXPANDS FRACTURES<sup>1</sup>

PREDICTABLY SAFE

### DISTINCTLY INTUITIVE

CONSISTENTLY EFFECTIVE Establish a new standard of care for safe calcium modification by reducing the risk of perforations and other complications to make procedures more predictable and efficient

Optimally modify calcium from the very first case via a unique mechanism of action in an intuitive platform

Achieve excellent lumen gains and stent expansion consistently despite the calcium morphology – whether superficial or deep, concentric or eccentric, and short or long lesions

## Learn more at ShockwavelVL.com

## Important Safety Information

Indications for Use— The Shockwave Intravascular Lithotripsy (IVL) System with the Shockwave C<sup>2</sup> Coronary IVL Catheter is indicated for lithotripsy-enabled, low-pressure balloon dilatation of severely calcified, stenotic *de novo* coronary arteries prior to stenting.

Contraindications— The Shockwave C<sup>2</sup> Coronary IVL System is contraindicated for the following: This device is not intended for stent delivery. This device is not intended for use in carotid or cerebrovascular arteries.

Warnings— Use the IVL Generator in accordance with recommended settings as stated in the Operator's Manual. The risk of a dissection or perforation is increased in severely calcified lesions undergoing percutaneous treatment, including IVL. Appropriate provisional interventions should be readily available. Balloon loss of pressure was associated with a numerical increase in dissection which was not statistically significant and was not associated with MACE. Analysis indicates calcium length is a predictor of dissection and balloon loss of pressure. IVL generates mechanical pulses which may cause atrial or ventricular capture in bradycardic patients. In patients with implantable pacemakers and defibrillators, the asynchronous capture may interact with the sensing capabilities. Monitoring of the electrocardiographic rhythm and continuous arterial pressure during IVL treatment is required. In the event of clinically significant hemodynamic effects, temporarily cease delivery of IVL therapy.

Precautions— Only to be used by physicians trained in angiography and intravascular coronary procedures. Use only the recommended balloon inflation medium. Hydrophilic coating to be wet only with normal saline or water and care must be taken with sharp objects to avoid damage to the hydrophilic coating. Appropriate anticoagulant therapy should be administered by the physician. Precaution should be taken when treating patients with previous stenting within 5mm of target lesion.

Potential adverse effects consistent with standard based cardiac interventions include – Abrupt vessel closure - Allergic reaction to contrast medium, anticoagulant and/or antithrombotic therapy-Aneurysm-Arrhythmia-Arteriovenous fistula-Bleeding complications-Cardiac tamponade or pericardial effusion-Cardiopulmonary arrest-Cerebrovascular accident (CVA)-Coronary artery/vessel occlusion, perforation, rupture or dissection-Coronary artery spasm-Death-Emboli (air, tissue, thrombus or atherosclerotic emboli)-Emergency or non-emergency coronary artery bypass surgery-Emergency or non-emergency percutaneous coronary intervention-Entry site complications-Fracture of the guide wire or failure/malfunction of any component of the device that may or may not lead to device embolism, dissection, serious injury or surgical intervention-Hematoma at the vascular access site(s)-Hemorrhage-Hypertension/Hypotension-Infection/sepsis /fever-Myocardial Infarction-Myocardial Ischemia or unstable angina-Pain-Peripheral Ischemia-Pseudoaneurysm-Renal failure/insufficiency-Restensis of the treated coronary artery leading to revascularization-Shock/pulmonary edema-Slow flow, no reflow, or abrupt closure of coronary artery-Stroke-Thrombus-Vessel closure, abrupt-Vessel injury requiring surgical repair-Vessel dissection, perforation, rupture, or spasm.

Risks identified as related to the device and its use: Allergic/immunologic reaction to the catheter material(s) or coating-Device malfunction, failure, or balloon loss of pressure leading to device embolism, dissection, serious injury or surgical intervention-Atrial or ventricular extrasystole-Atrial or ventricular capture.

Prior to use, please reference the Instructions for Use for more information on warnings, precautions and adverse events. www.shockwavemedical.com/IFU © 2021 Shockwave Medical, Inc. All Rights Reserved. IVL Overview | SPL 64337 Rev. B | 1. Disrupt CAD III Case