



PELVIS CONTROL HARNESS

PCH HARNESS DESIGN ALLOW FOR SAFE AND EASY TRANSFER OF PATIENT, IN OR OUT OF A WHEELCHAIR, WHEN OTHER BODY WEIGHT SYSTEMS OR LIFTS ARE NOT AVAILABLE.



PCH HARNESS PROVIDES NEEDED SECURITY OR ASSISTANCE DURING WALKING ACTIVITIES.



HORIZONTAL HANDLES

Provide support or assistance on the pelvis region during standing, sit to stand or stand to sit, tall kneeling, quadruped positions. Provide manual resistance by holding on to the plastic handles. The length can be adjusted when needed. Another option is to use the horizontal handles on the side or back of the harness.



PLASTIC HANDLES

PELVIS CONTROL HARNESS



Thigh straps male part of the buckle and tension adjustment ring provides the ability to secure thigh straps going from the gluteals area in the left and right groin and connect to the diagonally pointed female part of the buckle on both sides of the harness. This way the patient can be secured or spotted in the case of a fall and these elements prevent the harness from sliding up.



Sacrum loop on the harness can be used to add resistance on pelvis using cable machines, elastic bands with resistance.

QUICK REFERENCE GUIDE:

1. Harness can be placed on a patient in a sitting or supine position.
2. Make sure that leg loops are connected to the harness.
3. Place the harness behind the patient's back (lower back region).
4. Connect the 2" buckle and then, lifting patient's thigh slide each leg loop around the patient thigh into the groin and connect buckles on each side, adjust the tension to prevent harness from slipping up.
5. If you want max control on the patient's pelvis.
 - Top line of the harness in front should be aligned with ASIS.

SAFETY AND PRECAUTIONS:

1. This harness is not meant to completely unload or lift patient.
2. Before each use, inspect all components of the harness for any structural damage.
3. While adding resistance on pelvis through the sacrum loop, make sure to pad groin region and the stomach to avoid discomfort or skin damage.
4. Stop using harness if patient feel discomfort or pain.
5. Usage of harness differently then advocated may cause serious injury.

Fitting guide:

Best fits users with **32 to 42 inches** pelvis circumference.

CLEANING & MAINTENANCE INSTRUCTIONS:

To clean use the wipes or soap with warm water and let it dry.
Do not place in the sun. Do not use washing machine or dryer.

Warranty:

18 months. Do not use it if you see any damage of the material, buckles or webbing.

We recommend replacing this product every 5 years.

Material:

2mm black nylon foam/600D black cordura/2mm black nylon foam.

ORDER INFORMATION

EMAIL: dockust@ngulab.com | **PHONE NUMBER:** 708-203-4661

NEVERGIVEUPLABS.COM

MADE IN THE USA
U.S. PATENT PENDING