


VIAL2BAG

UNIVERSAL, NEEDLE-FREE ADMIXTURE DEVICE FOR IMMEDIATE USE IV ADMINISTRATION

The Vial2Bag Advanced™ 20mm admixture device enables reconstitution and transfer of a drug between a vial and an IV bag prior to administration to the patient. The innovative needle-free and easy-to-use design is ideal for immediate use and can optimize pharmacy clean room time, while assisting in standardizing admixture procedures to maximize efficiency of point of care nursing.

NEW DESIGN

- Dual channel design provides dedicated fluid pathways into and out of the IV bag
- Enhanced vial spike design improves the connection between vial and Vial2Bag Advanced™ 20mm adapter and the transfer of medication
- Needle-free device meets OSHA and NAPPSI sharps prevention guidance





Vial2Bag Advanced™20mm admixture device can be an effective solution when meeting the compounding provisions associated with USP <797>.*



Optimized Patient Care

- Provides a simple solution for just-in-time IV medication preparation, improving timeliness at point of care
- Easy to implement immediate use solution

Effective Waste Reduction

- Reduces medication and IV fluid waste
- Universal compatibility with all manufacturers' 50, 100, and 250mL IV bags
- Suitable for use with liquid and powder medications with standard 20mm closure vials

*Note: The use of the Vial2Bag Advanced™ 20mm admixture device in association with any provision of compounding should not be interpreted as modifying, extending or superseding a drug manufacturer's labeling recommendations for storage and expiration dating. Refer to drug manufacturer's recommendations and USP compounding guidelines for shelf life and sterility information

Distribution and Ordering



Vial2Bag Advanced™ **20mm Admixture Device** 50/Box Item # 36098150

Quantity 150/Case



DISTRIBUTED BY: Progressive Medical, Inc. 997 Horan Dr. • Fenton, MO 63026 **ORDERING INFORMATION:**

Phone: 314-961-5786 • Fax: 314-961-4535 cs@progressivemedinc.com



INNOVATION IN RECONSTITUTION AND TRANSFER SYSTEMS

Scan to learn more about our products and services

For more information visit www.westpharma.com

NORTH AMERICA 1-800-345-9800 OPTION 8 | SOUTH AMERICA +55 11 4055 6061 | EUROPE +49 (0) 2403-7960 | ASIA PACIFIC +65 6862 3400 WEST PHARMACEUTICAL SERVICES, INC. 530 HERMAN O. WEST DRIVE | EXTON, PA 19341