

ORTHOMERICA'S STAR® FAMILY OF CRANIAL REMOLDING ORTHOSES



STARband®
Side-Opening Band



STARband
Bi-Valve



STARlight®
Side-Opening Band



STARlight
Bi-Valve



STARlight PRO
Living Hinge



STARband Plus™
Living Hinge



STARband 3D™
Bi-Valve

General Information

Most widely prescribed headband design

Most used for scaphocephalic head shapes

Clear material allows visual inspection

Clear material allows visual inspection

For post-endoscopic craniectomy

Low profile design with foam liner

3D printed design

Used to Treat*

- Plagiocephaly
- Brachycephaly
- Scaphocephaly
- Post-operative craniosynostosis

- Plagiocephaly
- Brachycephaly
- Scaphocephaly
- Post-operative craniosynostosis

- Plagiocephaly
- Brachycephaly
- Scaphocephaly
- Post-operative craniosynostosis

- Plagiocephaly
- Brachycephaly
- Scaphocephaly
- Post-operative craniosynostosis

- Post-operative craniosynostosis

- Plagiocephaly
- Brachycephaly
- Scaphocephaly
- Post-operative craniosynostosis

- Plagiocephaly
- Brachycephaly
- Scaphocephaly

Measurements

Custom fabricated from scan (or plaster cast)

Custom fabricated from scan (or plaster cast)

Custom fabricated from scan (or plaster cast)

Custom fabricated from scan (or plaster cast)

Custom fabricated from scan (or plaster cast)

Custom fabricated from scan (or plaster cast)

Custom fabricated from scan

Designer

Orthomerica

Orthomerica

Orthomerica

Orthomerica

Drs. Jimenez & Barone

Orthomerica

Orthomerica

Materials

Plastic shell with foam liner

Plastic shells with foam liner

Clear plastic shell with foam pads

Clear plastic shells with foam pads

Clear plastic shell with foam pads

Plastic shell with foam liner

3D printed nylon shell with foam pads

Style

Side and top openings

Overlap design

Side and top openings

Overlap design

One piece living hinge

One piece living hinge

Top opening with side closures

*Selection of cranial remolding orthosis design and treatment protocol is subject to clinical evaluation and physician prescription.



starband
When nature needs a nudge™



Before Treatment



After STARband Treatment

ORTHOMERICA®
(877) 737-8444 | www.orthomerica.com

EC REP Emergo Europe
Prinsessegracht 20
2514 AP The Hague
The Netherlands