

A woman with long brown hair is captured in mid-air, performing a dynamic jump or dance move. She is wearing a black long-sleeved top with a large, gathered collar, black leggings, and black sneakers with white soles. Her hair is blowing in the wind. The background shows a modern building with a glass facade and a blue sky. A large, semi-transparent blue circle is overlaid on the left side of the image, containing the text.

PRODUCT CATALOG MEDICAL COMPRESSION STOCKINGS

2023

THE BAUERFEIND-SERVICE

RELIABLE SUPPORT AND QUALITY

**BAUERFEIND
ONLINESHOP**
experts.bauerfeind.com



SERVICE GUARANTEES

INTERNATIONAL CUSTOMER SERVICE



100 % AVAILABILITY

By telephone within our service times
(Mo. – Th. 7.00 am – 05.00 pm,
Fr. 7.00 am – 04.00 pm CET).
By fax, e-mail and electronic ordering
systems 24 hours per day.



ELECTRONIC ORDERING SYSTEM

Make use of our offers for a fast,
faultless and effort reduced data
transfer: online shop, EDI service partner.



DISPATCH HANDLING

At your request we take care of the
full processing of your shipment
(incl. shipment quotations, additional
documents) and we can bundle your
order flexibly.



PROCESSING OF RETURNS

Within 5 days after receipt of goods.



ONE FACE TO THE CUSTOMER

You have one dedicated contact person
for all of your desires and concerns.



PROACTIVE INFORMATION

We keep you posted – you will receive
every important information about your
order and delivery status.




RETURN OF GOODS

Redemption of faultless goods within
30 days after delivery with 100 % refund
or replacement (T&C's apply).

SYSTEMATIC APPROACH TO COMPRESSION FOR THE RIGHT CHOICE

The Bauerfeind systematic approach will help you find the right VenoTrain compression stocking for your customers. The right choice depends on the severity of the condition.



¹ Lurie F. et al. The 2020 update of the CEAP classification system and reporting standards. J Vasc Surg Venous Lymphat Disord, 2020;8(3):342-52.

² Rabe, E. et al.: Bonner Venenstudie der DGP, Phlebologie 1 / 2003

BAUERFEIND.

BY PEOPLE FOR PEOPLE.

FOR A BETTER QUALITY OF LIFE.

Bauerfeind AG is one of the world's leading manufacturers of medical aids. With our products, we enable people to be more active and thus make an important contribution toward maintaining and restoring health.

MEDICAL KNOW-HOW AND TECHNICAL EXPERTISE

Bauerfeind combines decades of expertise in technology and manufacturing with innovation: the tradition of a medium-sized, family-run business with a global mindset. Our specialists have worked very closely with physicians, scientists and developers right from the start. The result: patented, award-winning products and services.

EFFECTIVE VARIETY

Our medical VenoTrain compression stockings cover all uses/indications in modern compression therapy, offering medically proven effectiveness and unique material properties.

Thanks to the wide range of colors available, the compression products blend seamlessly into their users' everyday lives and can be combined with various fashion styles. For the prevention or treatment of lymphedema: our indication-specific solutions include a wide variety of circular-knit and flat-knit designs that are available as made-to-measure or standard products.

UNCOMPROMISING QUALITY

Bauerfeind products don't simply meet standards – they fulfill the most stringent requirements when it comes to quality, effectiveness and wearing comfort. This is confirmed by inspection procedures, for which our products are permitted to bear the seals and marks after successful testing.

All medical VenoTrain compression stockings are CE-certified, thus complying with current applicable statutory regulations. They also bear the RAL quality label as well as the OEKO-TEX® STANDARD 100 quality seal.

From the thread to the completed product, our medical compression stockings are developed and manufactured in Germany.





CONTENTS

COMPRESSION THERAPY THE RIGHT PRODUCT FOR EVERY INDICATION

OVERVIEWS

- 06 Fashion Colors 2023
- 08 Edging and top bands
- 10 BODYTRONIC®
- 12 Bauerfeind training app
- 13 Academy

PHLEBOLOGICAL TREATMENT RANGE

For mild to moderate
vein problems



- 15 Perfect-Fit size chart
- 16 VenoTrain® discrétion
- 18 VenoTrain® look
- 20 VenoTrain® micro
- 22 Batik and Design Edition
- 26 VenoTrain® business
- 28 VenoTrain® soft
- 30 VenoTrain® soft S
- 34 VenoTrain® pure

ACCESSORIES

- 56 Care accessories
- 58 VenoTrain® glider
- 59 VenoTrain® glider plus

**PHLEBOLOGICAL
TREATMENT RANGE**

For moderate to severe
vein problems



- 38** **VenoTrain® cocoon**
- 40** **VenoTrain® impuls**
V6+ size chart
- 44** **VenoTrain® ulcertec**
X5+ size chart

**LYMPHOLOGICAL
TREATMENT RANGE***

For lipedema and
lymphedema patients



- 47** **VenoTrain® soft arm**
- 47** **VenoTrain® impuls**
- 47** **VenoTrain® delight**
- 47** **VenoTrain® curaflow**

**LIFESTYLE
AND SPORT**

For beautiful and active
legs



- 50** **VenoTrain® act**
SHEER ELEGANCE
- 52** **VenoTrain® act**
COTTON
- 52** **VenoTrain® act**
MICROFIBER
- 54** **Compression Sock**
TRAINING
- 54** **Compression Sock**
PERFORMANCE

PRODUCT SELECTION AND ORDER

- 60** **Circular knit product range**
- 62** **Additional services | Alterations | Repairs**
- 64** **Additional Bauerfeind product categories**
- 65** **Order sheets**
- 67** **Terms and Conditions of Sale and Delivery**

* Detailed information about products for
lymphedema and lipedema can be found in
the separate Lymphology Product Catalog
(Item no. 271 4210 1001 001) as well as in
the current price list.

VenoTrain® FASHION COLORS 2023

Brighten up your life.

NEW!

DEEP PETROL

Classy and confident make for a
sporty and elegant look.





Two Fashion Colors that change every year turn your legs into an eye-catcher and help you pay attention to your health whilst making a fashion statement.

NEW!

Now also available for VenoTrain soft S and VenoTrain impuls!













































WILD BERRY

The trendy warm berry shade can be combined easily.

COLOR OVERVIEW

THE RIGHT COMPRESSION STOCKING FOR EVERY TYPE

STANDARD-COLORS	look	micro	soft	soft S	business	impuls
CARAMEL						
NAVY						
ANTHRACITE						
BLACK						
ADDITIONAL-COLORS						
CREAM						
AMARETTO						
ESPRESSO						
FASHION COLORS						
WILD BERRY						
 DEEP PETROL						



Visit our Bauerfeind online shop now for more information:
experts.bauerfeind.com/b2b/fashioncolors

EDGING AND TOP BANDS

PERFECT HOLD WITH A HIGH LEVEL OF COMFORT

Our edging and top bands make sure that compression stockings fit securely and don't slip. Their soft materials are breathable, gentle on the skin and very comfortable to wear. Learn more about the individual styles and different designs.



EDGING

SECURE AND COMFORTABLE



TOP BAND
For all knee-high compression stockings

The soft edging ensures wearing comfort and doesn't dig into the skin. Its neutral look works for both women and men.



COMFORT BANDS
For all pantyhose

This edging provides an excellent fit and a high level of comfort thanks to its higher elasticity and wide band. It sits comfortably without constricting.

Similar image



GARTER BELT

The garter belt without silicone coating improves the hold of the stockings and offers plenty of freedom of movement. It is suitable for sensitive skin and can be worn on one or both sides.

SILICONE TOP BANDS

SECURE AND INDIVIDUALIZED



SILICONE MICRO DOTS

The highly elastic silicone top band with its small, flat silicone dots is very inconspicuous, even when worn under tight clothing, and leaves no marks on the skin. The skin-friendly material with its soft surface perfectly adapts to the contours of the leg.



LACE SENSITIVE

The highly elastic silicone top band has an ultra-thin, continuous coating that provides a perfect hold. It adapts perfectly to the leg contour, even with soft connective tissue.



LACE

The curved, flat coating is breathable and provides a good hold. Its lace design is especially appealing to women.



SILICONE DOTS

Our classic silicone top band is mainly used for thigh-high stockings. It is gentle on the skin and ideal for men thanks to its neutral design.



SILICONE DOT PATTERN

This skin-friendly silicone top band with its elegant pattern is especially popular with women.

QUALITY THAT CAN BE MEASURED – AND MADE-TO-MEASURE.

BAUERFEIND MEASUREMENT TECHNOLOGY



BODYTRONIC® 200

HEALTHY LEGS CAN BE MEASURED

- A quick, non-invasive vein function measurement process using infrared sensors
- These measurements can give indications of impairment to leg vein hemodynamics, e.g. as the result of varicose veins, thrombosis and post-thrombotic syndrome



BODYTRONIC® 410

FLEXIBLE 3D SCANNER FOR THE LOWER EXTREMITIES

- Contact-free measurement for compression stockings, with no need for special measurement clothing
- 3D image of the measurement area up to the waist
- Quick and easy to set up and take down





BODYTRONIC® 610

FULL-BODY SCANNER INCLUDING PRODUCT
SUPPLY FOR ALL EXTREMITIES

- Contact-free measurements for medical compression stockings, knee and foot supports, knee orthoses as well as for back, arm and shoulder products
- Optional full-body scan, export of 3D data
- Integrated leg volume measurement
- Standing space of only 1.60 m²

BAUERFEIND TRAINING APP

TRAINING AS IDEAL THERAPY ADDITION.

The body can only stay mobile and recover properly with training and movement. The Bauerfeind training App supports patients to exercise with products from Bauerfeind. Our trainer targets different problems and motivates users to do something to improve their health.



Available free of charge for iOS and Android.
For more information, please see:
apps.bauerfeind.com/therapy




**THE TRAINING APP
IS A MEDICAL DEVICE**



DEVELOPED BY EXPERTS
The exercises have been compiled by physicians and physiotherapists to suit the relevant conditions.

HOW DOES THE TRAINING APP HELP?



PERSONALIZED TRAINING FOR EACH PATIENT

Two or three training sessions per week will increase the patient's strength and mobility. The exercises are chosen to suit the individual indication and are constantly adapted to the patient's fitness level and pain threshold. This brief workout of less than 30 minutes can easily be integrated into daily activities.



MONITORING IMPROVEMENT AND TRAINING PROGRESS

The trainer keeps asking about how the patient is feeling and optimizes the training plan accordingly. At the same time, he records the individual pain assessment, meaning users can always keep an eye on their healing process.

BAUERFEIND ACADEMY

COMPETENT. GOAL-ORIENTED. SUCCESSFUL.

Would you like customers to clearly favor your medical retail store and reward your team's expertise?

We will be happy to support you!

We offer you customized training with strong practical relevance. These training opportunities can also be conducted online.

Discuss the options with your sales manager.

For further information and registration forms, please contact your local Sales Representative.

BENEFITS AT A GLANCE

+ TRAINING WITH STRONG PRACTICAL RELEVANCE

The training sessions offered by our Academy will always include medical information to improve understanding and ensure increased technical expertise.

Specifically trained coaches, such as orthotists, physical therapists, and lymph therapists, will provide professional guidance through the different subjects, while always considering the relevant target group.

+ IN-HOUSE TRAINING SESSIONS ACCORDING TO YOUR SPECIFICATIONS

Are you interested in personalized training that is tailor-made to your specific requirements?

You can book your training session according to your specifications. You choose the subject, the location and the time. Just get in touch!

+ ONLINE TRAINING

Would you like to learn in a time-saving and cost-efficient way? You can book tailored online training and specify the contents and the subject according to your needs.





PHLEBOLOGICAL TREATMENT RANGE FOR MILD TO MODERATE VEIN PROBLEMS

The phlebological range is suitable **for mild to moderate venous disorders**. It also covers situations in which there is a particular risk for the development of vein problems, such as during pregnancy or when traveling.

VenoTrain compression stockings, which have been developed precisely for these cases, improve venous activity and combat venous leg problems.



A video with measuring instructions can be found at:
video.bauerfeind.com/

PERFECT-FIT STANDARD SIZE CHART FOR CIRCULAR KNIT PRODUCTS

VenoTrain® look, VenoTrain® micro, VenoTrain® business, VenoTrain® soft, VenoTrain® soft S, VenoTrain® cocoon, VenoTrain® pure



		CIRCUMFERENCE		XS for VenoTrain soft S	S		M		L		XL	
					NORMAL	PLUS	NORMAL	PLUS	NORMAL	PLUS	NORMAL	PLUS
		cT			up to 85		up to 90		up to 100		up to 120	
		cH			up to 95		up to 105		up to 120		up to 140	
		cG										
		Silicone top band, wide*			51–57		56–63		60–67		63–71	
		Silicone top band, regular		44–50	45–51		49–56		53–60		55–63	
		Top band, AT, silicone top band		44–50	45–57		49–63		53–67		55–71	
		cF		39–50	40–52		43–55		45–59		49–63	
		cE		32–38	33–40		35–43		38–48		42–53	
		cD		28–31	29–35		32–38		34–42		37–45	
		cC		29–33	31–37		34–40		37–45		40–48	
		cB1		23–26	25–29		27–32		29–35		32–38	
		cB		18–20	20–22		22–25		25–28		28–31	
		cY		27–31	29–33		31–35		34–40		37–44	
		cA		17–22	19–24		22–27		25–31		28–35	

*Wide silicone top band with silicone dots not available for VenoTrain micro
Sizing system information given in cm

For shoe sizes smaller than 36 or larger than 46 as well as in cases of particular anatomical features of foot, the use of custom-made products is recommended.



STANDARD SIZE CHART: VenoTrain discrétion

LEG CIRCUMFERENCE

Size	1	2	3	4
cG	42–56	49–63	56–70	63–77
cC	28–36	32–40	35–43	38–46
cB	18–21	21–24	24–27	27–31

LENGTH

Size	short	long
AG AT	< 70	> 70
AD	< 38	> 38

Sizing system information given in cm


VenoTrain® discrétion

THE TRANSPARENT STOCKING

The exceptionally fine material of these compression stockings gives them a stylish transparency and a smooth shine.

Available in various designs and versatile colors, these stockings will give you an extra sense of vitality and well-being in your work, leisure and everyday life.

COLORS



FEATURES

- highly elastic, making the stockings easy to put on
- shapes the legs without constricting them
- breathable fine knitting for a long-lasting wearing comfort

DETAILS

SIZING SYSTEM VenoTrain discrétion (I–IV)

COMPRESSION CLASSES

- Class 1 (10–15 mmHg)
- Class 2 (15–20 mmHg)

TOE closed, open

MATERIAL CONTENT





- Class 1: 75 % polyamide, 25 % elastane
- Class 2: 70 % polyamide, 30 % elastane

CARE INSTRUCTIONS



● STANDARD

○ MADE-TO-MEASURE SOLUTIONS

	 AD Knee stockings	 AG Thigh-high stockings with silicone lace band	 AG Thigh-high stockings with silicone dotted band	 AT Pantyhose
EDGING SILICONE TOP BAND				
Clair	●*	●*	●*	
Beige	●○	●○	●○	●○
Espresso	●○	●○	●○	●○
Black	●○	●○	●○	●○
Java	●	●	●	●

* only available in classe 2 (15–20 mmHg)
available in selected countries






VenoTrain® look

ATTRACTIVE AND GLOSSY

The VenoTrain look is a real highlight among medical compression stockings, thanks to its glossy look: for a stylish appearance on any occasion!

The VenoTrain look pantyhose enhances the figure with a push-up effect. Stretch-Comfort and its special trilobal fiber guarantee the perfect fit.

STANDARD COLORS



FEATURES

- Fashionable and elegant gloss effect
- Slightly translucent

DETAILS



SIZING SYSTEM PERFECT-FIT (S–XL)

COMPRESSION CLASSES

- Ccl 1 (18–21 mmHg)
- Ccl 2 (23–32 mmHg)

TOE closed

MATERIAL CONTENT




- Ccl 1: 60 % polyamide, 40 % elastane
- Ccl 2: 55 % polyamide, 45 % elastane

CARE INSTRUCTIONS



● STANDARD SIZES

○ MADE-TO-MEASURE SOLUTIONS

	 AD	 AG	 AT
TOP BAND EDGING			
Top band / Comfort band	● ○	○	● ○
Lace sensitive		● ○	
Silicone micro dots		● ○	
Silicone dot pattern		○	



C0

C1

C2

C3

C4

C5

C6



VenoTrain® micro

FASHIONABLE, MODERN AND VERSATILE

With a particularly high microfiber content that makes it comfortable to wear, the VenoTrain micro is one of the most popular compression stockings.

Available in many designs and a wide selection of colors, it highlights the individual style of your customers.

FEATURES

- More than 50 per cent microfiber content: breathable, skin-friendly, and easy-care
- Figure-shaping panty area (pantyhose) and comfort sole
- Semi-transparent fine stocking quality
- Also available as pantyhose for pregnant women with a special panty area that expands as the mother's bump grows (ATU)

DETAILS



SIZING SYSTEM PERFECT-FIT (S–XL)

COMPRESSION CLASSES

- Ccl 1 (18–21 mmHg)
- Ccl 2 (23–32 mmHg)

TOE closed, open

MATERIAL CONTENT

- Ccl 1: 65 % polyamide, 35 % elastane
- Ccl 2: 55 % polyamide, 45 % elastane


CARE INSTRUCTIONS



STANDARD COLORS







ADDITIONAL COLORS



FASHION COLORS



- **STANDARD SIZES**
- **MADE-TO-MEASURE SOLUTIONS**

	 AD	 AG	 AT	 ATU
TOP BAND EDGING				
Top band / Comfort band	● ○	● ○	● ○	● ○
Lace sensitive		● ○		
Silicone micro dots		● ○		
Silicone dots	● ○	● ○		
Silicone dot pattern		● ○		



VenoTrain® micro BATIK- AND DESIGN-EDITION

A true fashion statement:
the Design Edition and individually
dyed VenoTrain micro batik versions
create real highlights.

With a subtle Tango or Jive pattern,
or wildly textured in batik with navy or
black, your customers' legs will be real
eye-catchers.



COLORS

STANDARD COLORS

	TANGO*	JIVE*	BATIK NAVY	BATIK BLACK
CARAMEL				
ANTHRACITE				
NAVY				
BLACK				

ADDITIONAL COLORS

	TANGO*	JIVE*	BATIK NAVY	BATIK BLACK
CREAM				
AMARETTO				
ESPRESSO				

FASHION COLORS




	TANGO*	JIVE*	BATIK NAVY	BATIK BLACK
WILD BERRY				
DEEP PETROL				

NEW!

DEEP PETROL TANGO

DESIGN EDITION

● STANDARD SIZES ○ MADE-TO-MEASURE SOLUTIONS

	 AD	 AG	 AT
TOP BAND EDGING			
Top band / Comfort band	● ○		● ○
Lace sensitive		● ○	
Silicone micro dots		● ○	
Silicone dots	● ○	● ○	
Silicone dot pattern		● ○	

*The Design Edition is available only with closed toes

C0

C1

C2

C3

C4

C5

C6

DEEP PETROL
FASHION COLOR 2023



VenoTrain® micro

FASHION

COLORS 2023

Gentle support with a dash of color – with a microfiber content of more than 50%, VenoTrain micro, our most popular compression stocking, is particularly cozy and goes with any style.



FASHION COLORS

	STANDARD	TANGO	JIVE	BATIK NAVY	BATIK BLACK
WILD BERRY					
DEEP PETROL ^(NEW)					

● STANDARD SIZES ○ MADE-TO-MEASURE SOLUTIONS

TOP BAND EDGING	AD	AG	AT	ATU*
Top band / Comfort band	● ○		● ○	● ○
Lace sensitive		● ○		
Silicone micro dots		● ○		
Silicone dots*	● ○	● ○		
Silicone dot pattern		● ○		

* not available in Design Edition

C0	C1	C2	C3	C4	C5	C6



VenoTrain® business

TIMELESSLY ELEGANT

While traveling or at work: the VenoTrain business meets every requirement for a premium knee-high stocking. The combination of cotton and Timbrelle® high-tech fibers with thermal regulation makes it exceptionally comfortable to wear.


Highly elastic elastane and the extra-wide edging band ensure that the stocking is easy to put on and stays perfectly in place.

FEATURES

- Innovative mix of materials
- Reinforced sole
- Timbrelle® high-tech fiber providing long-lasting freshness and a healthy foot climate



STANDARD COLORS



DETAILS



SIZING SYSTEM PERFECT-FIT (S – XL)

COMPRESSION CLASSES

- Ccl 1 (18–21 mmHg)
- Ccl 2 (23–32 mmHg)

TOE closed

MATERIAL CONTENT

- Ccl 1: 80 % polyamide, 14 % elastane, 6 % cotton
- Ccl 2: 71 % polyamide, 24 % elastane, 5 % cotton

CARE INSTRUCTIONS



● STANDARD SIZES

○ MADE-TO-MEASURE SOLUTIONS

TOP BAND EDGING



AD

Top band / Comfort band



Silicone dots



Available in shoe sizes 36–41 (foot short) and 41–46 (foot long)



C0

C1

C2

C3

C4

C5

C6



VenoTrain® soft

GENTLE AND EFFECTIVE

The VenoTrain soft, always the right choice:
This versatile stocking can withstand plenty of
wear, yet it is very soft and gentle on the skin.

Available in many styles with a proven medical
effect and beneficial massage action.


FEATURES

- For men and women
- Particularly gentle on the skin
- Also available as a single stocking with a standard design in caramel
- Also available as pantyhose for pregnant women with a special panty area that expands as the mother's bump grows (ATU)


● STANDARD SIZES

○ MADE-TO-MEASURE SOLUTIONS

STANDARD COLORS



ADDITIONAL COLORS



FASHION COLORS



DETAILS



SIZING SYSTEM PERFECT-FIT (S–XL)

COMPRESSION CLASSES

- Ccl 1 (18–21 mmHg)
- Ccl 2 (23–32 mmHg)

TOE closed, open, no foot (made-to-measure),
stump (made-to-measure)

MATERIAL CONTENT

- Ccl 1: 70 % polyamide, 30 % elastane
- Ccl 2: 65 % polyamide, 35 % elastane

CARE INSTRUCTIONS



TOP BAND EDGING	AD	AG	AGHB	AT	ATU	ATH	ATX	ATE
Top band / Comfort band	● ○	● ○		● ○	○	○	○	○
Lace sensitive		● ○						
Silicone micro dots		● ○						
Silicone dots	● ○	● ○						
Silicone dot pattern		● ○						
Garter belt			● ○					



C0

C1

C2

C3

C4

C5

C6



VenoTrain® soft S

ALL-ROUNDERS THAT ARE EASY TO PUT ON

Thanks to its spacious heel and reduced knitted fabric in the instep area, our VenoTrain soft S offers particular foot comfort. The anatomically shaped 90° heel and the robust knitted fabric make putting on and taking off the product easier.

They are therefore perfectly suited for post-traumatic or post-operative use in edema therapy*, as well as for those who sit and stand a lot.

FEATURES

- Spacious Easy-heel facilitates donning the stocking
- Crease-free instep for a comfortable fit
- Breathable material with a revitalizing effect

● STANDARD SIZES

○ MADE-TO-MEASURE SOLUTIONS

TOP BAND EDGING	AD	AG
Top band / Comfort band	● ○	
Lace sensitive		● ○
Silicone micro dots		● ○
Silicone dots	● ○	● ○
Silicone dot pattern		● ○

STANDARD COLORS



ADDITIONAL COLORS



FASHION COLORS



DETAILS



SIZING SYSTEM PERFECT-FIT (XS – XL)

COMPRESSION CLASSES

- Ccl 1 (18–21 mmHg)
- Ccl 2 (23–32 mmHg)

TOE closed, open

MATERIAL CONTENT

- Ccl 1: 70 % polyamide, 30 % elastane
- Ccl 2: 65 % polyamide, 35 % elastane

CARE INSTRUCTIONS



* For a complete list of all indications, please see the VenoTrain instructions for use.



C0

C1

C2

C3

C4

C5

C6



**DEEP
PETROL**
FASHION
COLOR 2023



**WILD
BERRY**
FASHION
COLOR 2023

VenoTrain® soft / VenoTrain® soft S

FASHION COLORS 2023

Skin-friendly, beautiful, robust – the opaque VenoTrain soft compression stocking is available in our Fashion Colors and is the perfect stocking to withstand high demands whilst providing ideal wearing comfort.

The anatomically shaped heel of the VenoTrain soft S ensures even distribution of the knitted fabric in the instep and heel areas. The newly available Fashion Colors therefore unite style and comfort in this compression stocking as well.



**Crease-free
instep**

for a comfortable fit



Easy heel

easier to put on
than ever before

FASHION COLORS

WILD BERRY 



DEEP PETROL 











● STANDARD SIZES

○ MADE-TO-MEASURE SOLUTIONS

VenoTrain soft

VenoTrain soft S

TOP BAND EDGING								
Top band / Comfort band	● ○			● ○	○	○	○	○
Lace sensitive		● ○						
Silicone micro dots		● ○						
Silicone dots	● ○	● ○						
Silicone dot pattern		● ○						
Garter belt			● ○					

TOP BAND EDGING		
Top band / Comfort band	● ○	
Lace sensitive		● ○
Silicone micro dots		● ○
Silicone dots	● ○	● ○
Silicone dot pattern		● ○



VenoTrain® pure

BASIC RANGE

The VenoTrain pure compression stocking allows economical basic treatment.

It supports venous activity and alleviates tired, heavy legs. Thanks to its skin-friendly and breathable knitted fabric, it provides a high level of wearing comfort.

COLORS

BEIGE



BLACK



FEATURES

- For men and women
- Available in the tried-and-tested Perfect Fit sizing system
- Hard-wearing material

DETAILS

SIZING SYSTEM PERFECT-FIT (S–XL)

COMPRESSION CLASSES

- Ccl 1 (18–21 mmHg)
- Ccl 2 (23–32 mmHg)

TOE closed, open

MATERIAL CONTENT

- Ccl 1: 77 % polyamide, 23 % elastane
- Ccl 2: 74 % polyamide, 26 % elastane

CARE INSTRUCTIONS



● STANDARD SIZES

○ MADE-TO-MEASURE SOLUTIONS

TOP BAND EDGING	AD	AG	AT
Top band / Comfort band	● ○		● ○
Silicone dots		● ○	

C0

C1

C2

C3

C4

C5

C6





PHLEBOLOGICAL TREATMENT RANGE

FOR MODERATE TO SEVERE VEIN PROBLEMS

Some patients with vein problems have special requirements that call for compression stockings with individually tailored modes of action and properties. These requirements can be very diverse. That is why VenoTrain offers compression stockings particularly for patients with:

- **DRY SKIN**
VenoTrain cocoon
- **ADVANCED VENOUS PROBLEMS**
VenoTrain impuls
- **VENOUS LEG ULCERS**
VenoTrain ulcertec



SKIN REPAIR
EFFECT



CLINICALLY
TESTED SKIN-
FRIENDLINESS

VenoTrain[®] cocoon

FOR DRY AND SENSITIVE SKIN

The medical compression stocking with innovative dual function protects and soothes dry and sensitive skin.¹

High-quality skincare substances repair the natural skin barrier, making dry skin more supple.

COLORS



FEATURES

- Clinically tested skincare substances with cosmetic effect¹
- Clinically proven skin-friendliness¹
- With an innovative mix of natural fibers
- The lipid care complex improves skin hydration and counteracts itching and scaling
- For men and women

DETAILS



SIZING SYSTEM PERFECT-FIT (S–XL)

COMPRESSION CLASS

- Ccl 2 (23–32 mmHg)

TOE closed

MATERIAL CONTENT

- Ccl 2: 52 % polyamide, 26 % elastane, 11 % cotton, 11 % cellulose

CARE INSTRUCTIONS



STANDARD SIZES

**TOP BAND
EDGING**




AD

Top band / Comfort band



Available in shoe sizes 36–41 (foot short) and 41–46 (foot long)



¹ Westphal, Thea, et al. "Medical compression stockings on the skin moisture in patients with chronic venous disease." Vasa (2019). 48, pp. 502-508.





V6+ STANDARD SIZE CHART

		CIRCUM-FERENCE	I		II		III		IV		V		VI	
			NORMAL	PLUS	NORMAL	PLUS	NORMAL	PLUS	NORMAL	PLUS	NORMAL	PLUS	NORMAL	PLUS
Length ℓ G		cG	45-53	51-59	48-56	54-62	51-59	57-66	54-63	60-69	57-66	63-73	60-70	66-76
short	62-70	cF	39-45	43-50	42-49	46-53	45-52	49-57	48-56	52-60	51-59	55-64	54-63	58-67
long	70-80	cE	32-37	35-40	34-39	37-42	36-41	39-44	38-44	41-47	40-46	43-49	42-48	45-51
Type	AG													
Length ℓ D		cD	28-32	30-34	30-34	32-37	32-37	34-39	34-39	36-41	36-41	38-43	38-43	40-46
short	35-40	cC	29-33	31-35	31-35	33-38	33-38	35-40	35-40	37-42	37-42	39-44	39-44	41-47
long	40-44													
Type	AD	cB1	23-27	24-28	25-29	26-30	27-31	28-32	29-34	30-35	31-36	32-37	33-38	34-39
		cB	19-20,5		20,5-22		22-24		24-26		26-28,5		28,5-31	
		cY	25-29		26-31		28-33		30-35		32-38		35-41	
		cA	19-22		20-23		22-25		24-27		26-30		28-33	

Sizing system information given in cm

VenoTrain® impuls


POWERFUL SUPPORT FOR ADVANCED VENOUS PROBLEMS



The VenoTrain impuls produces a particularly deep action by exerting high working pressure with even the slightest muscle movements. Thanks to its microfiber content, it provides powerful support and a high level of wearing comfort at the same time. It is breathable and easy to put on.

Recommended in cases of advanced vein problems, post-thrombotic syndrome, follow-up care for healed ulcers and treatment of early-stage lymphedema (0 – 1).*

STANDARD COLORS



FASHION COLORS



FEATURES

- The short-stretch knitted fabric generates a high working pressure with even the slightest muscle movements
- Skin-friendly microfiber material makes the product breathable and ensures climate control
- Highly differentiated V6+ sizing system (I – VI)

DETAILS

SIZING SYSTEM V6+ (I – VI)

COMPRESSION CLASSES

- Ccl 2 (23–32 mmHg)
- Ccl 3 (34–46 mmHg)

TOE closed, open

MATERIAL CONTENT

- Ccl 2: 80 % polyamide, 20 % elastane
- Ccl 3: 75 % polyamide, 25 % elastane

CARE INSTRUCTIONS



● STANDARD SIZES

○ MADE-TO-MEASURE SOLUTIONS

TOP BAND EDGING	AD	AG	AGHB	AT	ATH	ATX
Top band / Comfort band	● ○			○	○	○
Silicone dots	● ○	● ○				
Garter belt			● ○			

* For a complete list of all indications, please see the VenoTrain instructions for use.





**DEEP
PETROL**
FASHION
COLOR 2023

VenoTrain® impuls

FASHION

COLORS 2023

Look good and feel great, even in cases of advanced vein problems. The powerful knitted fabric of the VenoTrain impuls is opaque and breathable – now also available in our Fashion Colors!










FASHION COLORS



● STANDARD SIZES

○ MADE-TO-MEASURE SOLUTIONS

						
TOP BAND EDGING	AD	AG	AGHB	AT	ATH	ATX
Top band / Comfort band	● ○			○	○	○
Silicone dots	● ○	● ○				
Garter belt			● ○			



X5+ STANDARD SIZE CHART

	CIRCUM-FERENCE	XS		S		M		L		XL	
		NORMAL	PLUS	NORMAL	PLUS	NORMAL	PLUS	NORMAL	PLUS	NORMAL	PLUS
Length lG	cG	45-57	52-66	49-62	55-68	53-67	58-74	58-73	65-82	63-80	69-90
short 61-70											
long 70-80											
Type AG											
Length lD	cC	24-29	28-34	27-32	31-38	29-36	35-42	33-40	38-47	36-44	43-52
short 35-40											
long 40-45											
Type AD											
	cB	20-21,5		21,5-23		23-25		25-28		28-32	

Sizing system information given in cm

VenoTrain® ulcertec

FOR VENOUS LEG ULCERS

The combination of an overstocking with a patented rhomboid knitted fabric specially developed by Bauerfeind and a skin-friendly liner for holding the wound dressing in place makes this compression stocking easy to put on without any need for assistance.

The liners are also ideal for providing a continuous low level of compression during the night. These are available to order separately in a pack of three and are suitable for washing at high temperatures.

FEATURES

- Coordinated compression of overstocking and liner
- Unique compression profile that accelerates the healing process while ensuring exceptional patient comfort¹
- Differentiated treatment based on two designs within one compression class

39 – moderate

- Compression class 3 (34–46 mmHg)
- Liner with gentle compression providing a compression pressure in the lower range of compression class 3 in combination with the overstocking

46 – strong


- Compression class 3 (34–46 mmHg)
- Liner with moderate compression providing a compression pressure in the upper range of compression class 3 in combination with the overstocking

COLORS

BEIGE
OVERSTOCKING



WHITE
LINER



DETAILS

SIZING SYSTEM X5+ (XS–XL)

COMPRESSION CLASSES

- 39–moderate: 34–46 mmHg
Liner with gentle compression
- 46–strong: 34–46 mmHg
Liner with moderate compression

TOE

- open (overstocking)
- closed (liner)

MATERIAL CONTENT

LINER AG*

- 39–moderate: 68 % polyamide, 32 % elastane

LINER AD

- 39–moderate: 87 % polyamide, 13 % elastane
- 46–strong: 69 % polyamide, 31 % elastane

OVERSTOCKING

- 39–moderate / 46–strong: 79 % polyamide, 21 % elastane

CARE INSTRUCTIONS

- Overstocking:
- Liner:

● STANDARD SIZES

○ MADE-TO-MEASURE SOLUTIONS

TOP BAND EDGING



AD



AG

Top band / Comfort band



Silicone sensitive



¹ Jünger, M., Wollina, U., Kohnen, R., & Rabe, E. (2004). Efficacy and tolerability of an ulcer compression stocking for therapy of chronic venous ulcer compared with a below-knee compression bandage: results from a prospective, randomized, multicentre trial. Current medical research and opinion, 20(10), 1613–1623.

*The AG liner additionally has a silicone top band on the upper edge of the stocking. It consists of polyamide, elastane, and silicone.

C0

C1

C2

C3

C4

C5

C6





Christina is wearing the VenoTrain curaflow as a glove and arm sleeve in the color royal.

LYMPHOLOGICAL TREATMENT RANGE

Bauerfeind has a wide range of compression products for lymphological conditions which take into account the patients' different requirements. Based on the high microfiber content of VenoTrain products, the compression garments provide maximum wearing comfort and breathability even during movement.



The circular knit **VenoTrain soft** is recommended for mild edemas* in the arm area. With the comfortable inner elbow zone, it provides greater freedom of movement, wearing comfort and reduced feeling of pressure.



The circular-knit **VenoTrain impuls** compression stocking combines a wide range of benefits: its breathable microfiber and firm compression thread offer unique wearing comfort with a high degree of effectiveness. Only very little movement is needed for the stocking to exert the required therapeutic pressure. Recommended in cases of early-stage lymphedema.



The flat-knit **VenoTrain delight** compression stocking is indicated for use in cases of advanced venous disorders and lymphedema. The breathable double-stretch knit is skin-friendly and easy to put on and take off. Additionally, the firm material offers a good shaping effect and is very durable.



Thanks to our made-to-measure **VenoTrain curaflow** flat knit products, we can offer patients suffering from all stages of lymphedema and lipedema perfectly tailored solutions with more than 1,000 different options for combining products and additional features. The beneficial therapeutic features of the range include the knitted fabric's high microfiber content, the flat sensitive seams and special functional zones, such as the extra-soft area behind the knee.

COLORS

STANDARD COLORS



ADDITIONAL COLORS



VERSIONS

LEG & FOOT

Ccl 2 – 3

Ccl 2 – 3

Ccl 1 – 4

ARM & HAND

Ccl 1 – 2**

Ccl 1 – 3

* For a complete list of all indications, please see the VenoTrain instructions for use.

** VenoTrain soft arm is available only as arm sleeves.

Comprehensive information and product details of our lymphology portfolio can be found in the **Lymphology Product Catalog**. You are welcome to order it using item number 271 4210 1001 001 or download it directly using the following link: bauerfeind.com/downloads

LIFESTYLE AND SPORT

GENTLE SUPPORT AT WORK, DURING LEISURE TIME AND EXERCISE

Support stockings are an effective weapon in the fight against tired and heavy legs. The high-quality knitted fabric creates a noticeably pleasant compression, making you feel revitalized when you wear them. New skin-friendly materials meet stylish designs and a wide choice of colors.







VenoTrain[®] act SHEER ELEGANCE

PROVIDING ELEGANT SUPPORT

In 70 or 140 denier, the VenoTrain act sheer elegance provides light to medium support in a modern translucent design that is available in a wide range of colors and styles.

DETAILS

TOE closed

MATERIAL CONTENT*

AG (70 DEN) 81 % polyamide, 19 % elastane

AG (140 DEN) 73 % polyamide, 27 % elastane



CARE INSTRUCTIONS



FEATURES


- For her
- Elegant and transparent
- Available in seven shades
- Available as a knee-high or thigh-high stocking with a silicone lace band and as pantyhose
- Also available as pantyhose for pregnant women with a special pantyhose top that expands as the mother's bump grows
- Calibrato pantyhose (AT CAL) for particularly feminine curves

STANDARD SIZES


							
			AD	AG	AT	ATU	AT CAL
DEN			70 / 140	70 / 140	70 / 140	140	70 / 140
COLORS		SAND	●	●	●		
		CAMEL	●	●	●	●	●
		AMARETTO	●	●	●		
		ESPRESSO	●	●	●		
		BLACK	●	●	●	●	●
		GRAPHITE	●	●	●		
		NAVY	●	●	●		

SIZE CHARTS

AG, AT, ATU



AT CAL



AD

SIZE	0	1	2	3
SHOE SIZE	37-39	39-41	41-43	43-45

VenoTrain act is not a medical device.

*We are happy to share example specifications, additional material contents with you on request.

No alterations, repairs, or additional services are offered for the VenoTrain act.
It cannot be returned in exchange for a credit note or another item.



MICROFIBER

COTTON

VenoTrain[®] act COTTON

VenoTrain[®] act MICROFIBER

PROVIDING GENTLE SUPPORT

VenoTrain act cotton, available as a men's stocking, and VenoTrain act microfiber, available for men and women, offer a high degree of support in 280 denier and are reliable companions when traveling, at work or during leisure activities.

DETAILS

TOE closed

MATERIAL CONTENT

- COTTON: 27 % polyamide, 10 % elastane, 63 % cotton
- MICROFIBER: 94 % polyamide, 6 % elastane

CARE INSTRUCTIONS



FEATURES






COTTON

- For him
- Particularly comfortable to wear, thanks to cotton
- With distinctive ribbed look
- Opaque
- Knee-high stocking with a high level of support (280 DEN)
- Available in four colors

MICROFIBER

- For men and women
- Particularly skin-friendly and soft, thanks to high microfiber content
- Fits snugly around the leg, while allowing the skin to breathe
- Knee-high stocking with a high level of support (280 DEN)
- Available in four colors

STANDARD SIZES

			
			AD
DEN			280
COLORS		BEIGE	●
		BLACK	●
		GRAPHITE	●
		NAVY	●

SIZE CHART

AD

COTTON (for him)

SIZE	1	2	3	4
SHOE SIZE	39–41	41–43	43–45	45–47

AD

MICROFIBER (for men and women)

SIZE	0	1	2	3	4
SHOE SIZE	37–39	39–41	41–43	43–45	45–47

VenoTrain act is not a medical device.

No alterations, repairs, or additional services are offered for the VenoTrain act. It cannot be returned in exchange for a credit note or another item.



STANDARD SIZE CHART

		CALF CIRCUMFERENCE	S	M	L	XL
		cC	31–40	34–44	37–47	40–51
Shoe size						
short	35–40					
long	41–46					
Type	AD					

COMPRESSION SOCK TRAINING

COMPRESSION SOCK PERFORMANCE

The **Compression Socks Training** and **Performance** were developed on the basis of extensive expertise in biomechanics and compression therapy. They offer noticeable support during continuous exertion, prevent tiring muscle vibrations and help the muscles to regenerate more quickly.

FEATURES

TRAINING

- Muscle Toning Zone supports the calf muscles
- Shock absorption zones protect the forefoot and heel areas
- Muscle stimulation, thanks to gentle compression
- The taping zone helps to prevent injuries and improves proprioception
- Especially suited for sports with stopping and jumping movements



PERFORMANCE

- Protection zone on the Achilles tendon
- Breathable microfiber
- Especially suitable for endurance sports not requiring many changes of direction

DETAILS

TOE closed
MATERIAL CONTENT
TRAINING 93 % polyamide, 7 % elastane
PERFORMANCE 75 % polyamide, 25 % elastane
CARE INSTRUCTIONS


TRAINING AVAILABLE IN



PERFORMANCE AVAILABLE IN



STANDARD SIZES

TOP BAND EDGING		 AD	TOE CLOSED
Top band	TRAINING	●	●
	PERFORMANCE	●	●

ACCESSORIES THE PERFECT ADDITION FOR EVERYDAY USE

These handy extras from the VenoTrain accessories range are useful for ensuring that your patients feel comfortable using and wearing their VenoTrain compression stockings. They provide effective treatment support, make it easier to handle the products, and care for the compression stockings as well as your customers' skin.





01

01 Lingerie bag

To protect the compression stockings during every machine wash

02 Bauerfeind care foam*

Contains panthenol to protect and soothe particularly stressed and sensitive skin; before putting on the stockings and in the evening after taking them off



02

03 Detergent*

Cleans gently, ensuring careful product treatment and providing lasting protection



03

04 Bauerfeind special gloves

Help with putting the compression stockings on correctly



04

05 VenoTrain® Lotion*

Developed in collaboration with Sebapharma; contains chamomile extract and is specially designed for people with sensitive skin; for use in the evening after taking off their stockings



05

* only available in selected countries

VenoTrain® glider

COMPACT DONNING AID

A video with donning instructions can be found at:
video.bauerfeind.com/



The VenoTrain glider is a practical aid for putting on medical compression stockings. It can be prescribed by a physician.

The lightweight, hard-wearing material makes the stocking glide effortlessly over the arm or leg. It is suitable for compression stockings with an open or closed toe.

Saving space and time, it is the ideal companion when you're out and about.

FEATURES

- Extremely lightweight, robust, and tear-resistant
- Easy to handle
- Ideal for travel
- Product is available on prescription



STEP 1



STEP 2



STEP 3



STEP 4



STEP 5

VenoTrain® glider plus

ROBUST DONNING AID

The VenoTrain glider plus is a robust donning aid with an ergonomic design for all types of medical compression stockings. It can be prescribed by a physician.

Folding handles make it easy to put on the stockings, avoiding the need to bend down.

FEATURES

- Compact, ergonomic design
- Light metal frame with folding handles
- Easy to handle without the need for any great effort
- Product is available on prescription



STEP 1



STEP 2



STEP 3



STEP 4



STEP 5



STEP 6


STANDARD PRODUCT RANGE

CIRCULAR-KNIT PRODUCTS


DESIGN		SIZE	CCL			EDGING / TOP BAND								TOE	
						Top band / Comfort band	Lace	Lace sensitive	Silicone micro dots	Silicone dots	Silicone dot pattern	Silicone sensitive	Garter belt	closed	open
VenoTrain			1	2	3										
discrétion	AD	I – IV	●*	●**		●								●	●
	AG	I – IV	●*	●**			●			●				●	●
	AT	I – IV	●*	●**		●								●	●
look	AD	S – XL	●	●		●								●	
	AG	S – XL	●	●				●	●					●	
	AT	S – XL	●	●		●								●	
micro	AD	S – XL	●	●		●				●				●	●
	AG	S – XL	●	●		●		●	●	●	●			●	●
	AT	S – XL	●	●		●								●	●
	ATU	S – XL	●	●		●								●	●
micro Design Edition (Jive, Tango)	AD	S – XL	●	●		●								●	
	AG	S – XL	●	●				●	●	●	●			●	
	AT	S – XL	●	●		●								●	
business	AD	S – XL	●	●		●								●	
soft	AD	S – XL	●	●		●				●				●	●
	AG	S – XL	●	●		●		●	●	●	●			●	●
	AGHB	S – XL	●	●									●	●	●
	AT	S – XL	●	●		●								●	●
pure	AD	S – XL	●	●		●								●	●
	AG	S – XL	●	●						●				●	●
	AT	S – XL	●	●		●								●	●
soft S	AD	XS – XL	●	●		●				●				●	●
	AG	XS – XL	●	●				●	●	●	●			●	●
cocoon	AD	S – XL		●		●								●	
impuls	AD	I – VI		●	●	●				●				●	●
	AG	I – VI		●	●	●				●				●	●
	AGHB	I – VI		●	●								●	●	●
ulcertec 39 moderate	AD	XS – XL	34 – 46 mmHg			●								Liner	Over- stocking
	AG									●					
	AD/AG					●									
ulcertec 46 strong	AD	XS – XL	34 – 46 mmHg			●							Liner	Over- stocking	

* 10 – 15 mmHg ** 15 – 20 mmHg


VERSIONS




AD
Knee-high
stocking




AF
Mid-
high
stocking




AG
Thigh-
high
stocking




AGHB
Thigh-high
compression
stocking with
garter belt




AT
Pantyhose




ATU
Maternity
pantyhose




ATH
Gent's
pantyhose
without fly




ATX
Pantyhose
with com-
pressive
pantyhose
top



ATE
Single-leg
pantyhose



Toe
closed



Toe
open

MADE-TO-MEASURE PRODUCT RANGE

CIRCULAR-KNIT PRODUCTS

DESIGN		CCL			EDGING / TOP BAND								TOE			
		1	2	3	Top band / Comfort band	Lace	Lace sensitive	Silicone micro dots	Silicone dots	Silicone dot pattern	Silicone sensitive	Garter belt	closed	open	No foot	Stump
VenoTrain																
	discrétion	AD	•*	•**	•								•	•		
		AG	•*	•**		•			•				•	•		
look		AT	•*	•**	•								•	•		
	AD		•	•	•								•			
	AG		•	•	•		•	•		•			•			
micro		AT	•	•	•								•			
	AD		•	•	•				•				•	•		
	AF		•	•	•		•	•	•	•			•	•		
	AG		•	•	•		•	•	•	•			•	•		
	ATU		•	•	•								•	•		
micro Design Edition (Jive, Tango)		AD	•	•	•								•			
	AG		•	•			•	•	•	•			•			
	AT		•	•	•								•			
business	AD		•	•	•				•				•			
soft	AD		•	•	•				•				•	•	•	•
	AF		•	•	•		•	•	•	•			•	•	•	•
	AG		•	•	•		•	•	•	•			•	•	•	•
	AGHB		•	•								•	•	•	•	•
	AT		•	•	•								•	•		•
	ATU		•	•	•								•	•		
	ATH		•	•	•								•	•		
	ATX		•	•	•								•	•		
	ATE		•	•	•								•	•	•	
	AD		•	•	•								•	•		
pure		AG	•	•					•				•	•		
		AT	•	•	•								•	•		
	AD		•	•	•				•				•	•		
soft S		AG	•	•			•	•	•	•			•	•		
	AD		•	•	•				•				•	•		
impuls		AD		•	•	•			•				•	•		
		AG		•	•	•			•				•	•		
		AGHB		•	•							•	•	•		
		AT		•	•	•							•	•		
		ATH		•	•	•							•	•		
		ATX		•	•	•							•	•		
	AD		•	•	•								•	•		
ulcertec 39 moderate	AD		34–46 mmHg		•								Liner	Over- stocking		
	AG										•					
ulcertec 46 strong	AD		34–46 mmHg		•								Liner	Over- stocking		

* 10–15 mmHg ** 15–20 mmHg

Please refer to the current price list for information on standard and made-to-measure designs. The **made-to-measure lymphology product range** can be found in the current Lymphology Product Catalog (Item no. 271 4210 1001 001) as well as in the current price list.



ADDITIONAL SERVICES ALTERATIONS | REPAIRS

INDIVIDUAL CUSTOMER SERVICE

ADDITIONAL SERVICES FOR NEW ORDERS

○ MADE-TO-MEASURE SOLUTIONS

	look	micro	micro DESIGN EDITION	business	soft	pure	soft S	impuls
Knitted toe area To be lengthened or shortened (3 – 10 cm)		○			○		○	○
Toe Seam on the outside	○	○	○	○	○	○	○	○
Toe Wider, narrower, anatomical		○	○	○	○	○	○	○
Panty area Front/back – two seams	○	○			○*			○
ATH with longitudinal slit					○			○
ATH with transverse opening					○			○
AT with extra-wide gusset for the fuller figure	○	○			○	○		○
Sew in a stoma opening (only for repeat order)					○			

* Please get in touch with the Bauerfeind Customer Service Team.

ALTERATIONS AND REPAIRS

● STANDARD SIZES

○ MADE-TO-MEASURE SOLUTIONS

		look	micro	micro DESIGN EDITION	business	soft	pure	soft S	cocoon	impuls	ulcertec
Top of stocking	Replace silicone top band	● ○	● ○	● ○	○	● ○	● ○	● ○		● ○	● ○
	Change top band to silicone top band	○	● ○		○	● ○		● ○		● ○	
	Silicone top band narrower or wider	● ○	● ○	● ○	○	● ○	● ○	● ○		● ○	● ○
	Remove hip section and sew on silicone top band					● ○				● ○	
AT	Shorten pantyhose top (except ATX)	● ○	● ○	● ○		● ○	● ○			○	
	Lengthen pantyhose top					○*				○	
	Make gusset narrower (excluding ATX)	○	○*			○*	○			○	
	Make gusset wider (excluding ATX)	○	○*			○*	○*			○	
	Extra-wide gusset for the fuller figure (excluding ATX)	○	○			○	○			○	
	Alter AT to AG with silicone top band	○	○	○		○	○			○	
AGHB	Replace hip strap					● ○				● ○	
	Shorten hip section					● ○				● ○	
	Lengthen / shorten hip strap					○				○	
Add a stoma opening or an open crotch afterwards (except ATX)						○					
Shorten stocking	AD by up to 3 cm	● ○	● ○	● ○		● ○	● ○	● ○		● ○	● ○
	AF by up to 4 cm		○			○					
	AG by up to 5 cm	● ○	● ○	● ○		● ○	● ○	● ○		● ○	● ○
Redye from light to dark**						○		○		○	
Toe alteration	Alter closed toe to open toe		○*			○*	○*	○*		○*	
	Shorten open toe (from 3 cm)		○			○	○	○		○	○
Repairs, if possible	Hole*	● ○	● ○		● ○	● ○		● ○	●	● ○	● ○
	Run*	● ○	● ○		● ○	● ○		● ○	●	● ○	● ○
	Panty area seam	● ○	● ○			● ○				○	
	Loose thread*	● ○	● ○		● ○	● ○		● ○	●	● ○	● ○

*Please get in touch with the Bauerfeind Customer Service Team.

**Please do not try on stockings – re dyeing is no longer possible after contact with the skin.

Please refer to the current price list for information about designs, item numbers and prices for standard and made-to-measure solutions.

The **made-to-measure lymphology product range** can be found in the current Lymphology Product Catalog (Item no. 271 4210 1001 001) as well as in the current price list. If you have any questions about our products, please contact the Bauerfeind Customer Service Team.

ADDITIONAL BAUERFEIND PRODUCT CATEGORIES

Medical aids
have been
proven to be
effective*



SUPPORTS AND ORTHOSES

by Bauerfeind help after injuries as well as in cases of acute and chronic orthopedic conditions. They improve the feeling of stability and provide confidence during everyday activities and exercise. The supports and orthoses are made from lightweight and breathable materials to perfectly underpin therapeutic exercise. They are easy to handle during everyday use, thanks to innovative solutions such as quick-fit fasteners and integrated donning aids, and will win you over with premium quality "Made in Germany".

COMPRESSION STOCKINGS FOR EDEMA TREATMENT

Every day is more pleasant when you're feeling good. That is why the flat knit, custom-made VenoTrain® curaflow compression stockings were developed. This versatile and colorful product perfectly adapts to the individual shape of the lymphedema or lipedema patient's legs, arms, hands, or thorax. VenoTrain® curaflow alleviates pain caused by tension and provides noticeable relief so patients can live their everyday lives feeling carefree, confident, and attractive.



A healthy body needs healthy feet. **ORTHOPEDIC FOOT ORTHOSES** simultaneously relieve pressure on the feet, knees, hips and back. They help the foot to function naturally and can therefore positively influence the body's biomechanics. As different feet need different types of support, Bauerfeind offers more than 300 different types of foot orthosis.

For an overview of all treatment areas, please visit:
experts.bauerfeind.com

*** Supports / orthoses, compression stockings and foot orthoses help patients get back on their feet quickly after an injury or in case of disorders.**

(Source: representative patient survey conducted by the Allensbach Institute for Public Opinion Research on behalf of eurocom e.V., 2019, n=1,295 users of supports/orthoses, compression stockings and foot orthoses)

BAUERFEIND AG TERMS AND CONDITIONS OF SALE AND DELIVERY (GTC)

Prices

The prices quoted are net prices in euros ex works and do not include statutory VAT.

Shipping

Shipping is at the buyer's cost and risk even if the buyer does not bear the transportation cost. Disruptions of operations affecting us or any of our suppliers, lack of suitable personnel, strike, lockout, and force majeure are causes for extending delivery and relieve us of any delivery obligation for the time of the obstruction.

Shipping costs

We calculate our shipping charges as a flat rate of EUR 7.00 for a net total order value of up to EUR 50.00, and EUR 4.00 for a net total order value of up to EUR 200.00. No shipping costs applied, when the order contains custom-made goods. Special shipping requirements will be charged at cost.

Returns

Returns outside the warranty terms are only accepted after prior advice and approval from us, provided the original delivery documents are submitted; it is supposed that the returned goods are in good order. We reserve the right of rejecting the acceptance of any returned consignment. Except where covered by warranty, we charge for processing returns. For returns within one month of the date of delivery, this charge is EUR 10.00 for each return consignment. For goods returned one month or more after the date of delivery, a handling fee of 20 percent of the net order value, or a maximum of EUR 25.00 per product and unit, is charged; for returns six to twelve months after delivery date, the fee is 40 percent of the net order value, or a maximum fee of EUR 50.00 for each product and unit.

Goods returned outside the warranty or 12 months or more after the date of delivery will not be accepted back. Custom made goods are non-returnable.

Payments

If not otherwise agreed individually in advance and in writing, Payments are due immediately following delivery. We offer a 4 percent discount for payments made by direct debit and a 3 percent discount for payments made within 10 days of the invoice date. Otherwise, the period allowed for payment is 30 days net. The buyer shall be considered in default once 30 days have passed since the invoice date. From this point on, any further deliveries will only be made against advance payment. From the time at which payment is delayed, we will charge default interest of eight percent above the applicable base interest rate. Payments by check, bank draft, etc. are deemed to have been made when the amount has finally been credited to our account.

Retention of title

We retain the title to all goods delivered until the buyer has paid all amounts outstanding under our business relations with him. The retention of title shall also apply proportionately to goods which have been processed or otherwise modified. If the goods are sold, all third-party claims are deemed to have been assigned to us without provision up to the total amount owing to us.

Notwithstanding the assignment, the buyer is entitled to the collection of these amounts. Our right to collect the amounts is not affected by the preceding provision; however, we undertake to refrain from collection unless the buyer defaults on payment. If the buyer defaults on payment we can demand that the buyer discloses the assigned claims and the names of the debtors, provides all information required for collection of the debt, delivers all required documents

and informs the debtors (third party) of the assignment. The buyer shall not be entitled to sell, pledge or assign the goods as security. Pledges, confiscation, or any other dispositions taken by a third party shall be promptly communicated to us by the buyer, and all information and documents necessary for us to defend our rights shall be provided. We undertake to release, on the buyer's demand, any security which exceeds a secured claim that has not been satisfied by more than 20 percent

Complaints

Warranty claims in respect of obvious defects shall be void if we do not receive immediate written notification of these claims in accordance with Section 377 of the German Commercial Code (HGB). Warranty claims shall also lapse one year after the delivery of the purchased goods. If the product is defective, we shall make a replacement delivery or repair it at our discretion. For reasons of hygiene, repairs will only be performed on products that have been washed or are completely hygienic. Upon receipt of unwashed goods, we reserve the right to claim a suitable sum for additional work. The customer may not demand a price reduction or withdraw from the contract provided that we fulfil our duties of subsequent performance, particularly with regard to the repair of defects or the delivery of goods without defect, unless our attempts at subsequent performance have failed. Defects which are caused by normal wear and tear, overuse, misuse, neglect of care, or noncompliance with the instructions for use are not covered by warranty.

Liability

Liability for damages is excluded unless there has been grossly negligent or willful conduct on our part or on the part of our representatives or agents. This limitation of liability shall not apply to damage due to injury to life, body or health, or the violation of a material contractual duty.

Place of jurisdiction

The place of performance for the delivery is our company headquarters. The place of jurisdiction for all disputes arising from the business relation is Gera, Germany.

Severability clause

The event that any individual provision of these terms and conditions of sale and delivery is or becomes invalid, will not affect the validity of the remaining provisions. The invalid provision shall be replaced with a valid provision which reflects the original economic intention of the parties as closely as possible.

Data protection clause / recording of conversations

For quality- and other purposes, we record calls, faxes, e-mails and letters between the buyer and our employees. All recordings serve legal provisions, which have to be met with medical products ordered. Recordings will be deleted automatically in line with legal provisions. Non of the recorded data will be disclosed to third parties, if not required by decision of court or on request of local tax or healthcare authorities. The buyer hereby consents to recordings such.

Please see www.bauerfeind.de/ds for further information on the latest version of Bauerfeind AG's General Data Protection Policy.

Last update: February 2022

DISCLAIMER:

All information given in this catalog is up to date at the day of first publishing. The catalog allows – depending on the information given – an initial overview and suggests some, but not all options of medical diagnosis and treatments to experts only. This catalog does not address patients. Although scientifically based, the catalog by no means claims to be complete or show all existing doctrinal. The entire content of this catalog – also in extracts, no matter how published, and especially concerning diagnostics and forms of treatment (in text, pictures, graphics, images or other material) has to be understood as non-binding support for the therapy to be selected finally by medical doctor, physician or other medical professionals. The catalog's intension is not to offer the full range of all medical aspects and latest scientific discussion that should be mandatorily considered before taking actions. We strongly recommend individual medical advice. If in doubt, please seek additional opinion of medical experts, who we would be delighted to provide.

MOTION IS LIFE WITH COMPRESSION BY BAUERFEIND

Do you have any questions, suggestions or requests?
Please contact us:

BAUERFEIND AG

Triebeser Straße 16
07937 Zeulenroda-Triebes
Germany
P +49 (0) 36628 66-1000
F +49 (0) 36628 66-1999
E info@bauerfeind.com

AUSTRIA

Bauerfeind Ges.m.b.H.
Hainburger Straße 33
1030 Wien
P +43 (0) 800 4430-130
F +43 (0) 800 4430-131
E info@bauerfeind.at

BENELUX

Bauerfeind Benelux B.V.
Waarderveldweg 1
2031 BK Haarlem

THE NETHERLANDS

P +31 (0) 23 531-9427
F +31 (0) 23 532-1970
E info@bauerfeind.nl

BELGIUM

P +32 (0) 2 527-4060
F +32 (0) 2 792-5345
E info@bauerfeind.be

BOSNIA AND HERZEGOVINA

Bauerfeind d.o.o.
Meše Selmovića 19
71000 Sarajevo
P +387 (0) 33 710-100
F +387 (0) 33 619-422
E info@bauerfeind.ba

CROATIA

Bauerfeind d.o.o.
Goleška 20
10020 Zagreb
P +385 (0) 1 6542-855
F +385 (0) 1 6542-860
E info@bauerfeind.hr

FRANCE

Bauerfeind France S.A.R.L.
B.P. 59258
95957 Roissy CDG Cedex
P +33 (0) 1 4863-2896
F +33 (0) 1 4863-2963
E info@bauerfeind.fr

ITALY

Bauerfeind Italia SRL
Via Cornaggia 58
20092 Cinisello Balsamo (MI)
P +39 02 8977 6310
F +39 02 8977 5900
E info@bauerfeind.it

NORDIC

Bauerfeind Nordic AB
Storgatan 14
114 55 Stockholm
P +46 (0) 774 100 020
E info@bauerfeind.se

REPUBLIC OF NORTH MACEDONIA

Bauerfeind Dooel Skopje
50 Divizija 24 a
1000 Skopje
P +389 (0) 2 3179-002
F +389 (0) 2 3179-004
E info@bauerfeind.mk

SERBIA

Bauerfeind d.o.o.
102 Omladinskih brigada
11070 Novi Beograd
P +381 (0) 11 2287-050
F +381 (0) 11 2287-052
E info@bauerfeind.rs

SINGAPORE

Bauerfeind Singapore Pte Ltd.
Blk 41 Cambridge Road #01-21
Singapore 210041
P +65 6396-3497
F +65 6295-5062
E info@bauerfeind.com.sg

SLOVENIA

Bauerfeind d.o.o.
Dolenjska cesta 242 b
1000 Ljubljana
P +386 (0) 1 4272-941
F +386 (0) 1 4272-951
E info@bauerfeind.si

SPAIN

Bauerfeind Ibérica, S.A.
C / San Vicente Mártir,
nº 71 – 4º - 7ª
46007 Valencia
P +34 96 385-6633
F +34 96 385-6699
E info@bauerfeind.es

SWITZERLAND

Bauerfeind AG
Vorderi Böde 5
5452 Oberrohrdorf
P +41 (0) 56 485-8242
F +41 (0) 56 485-8259
E info@bauerfeind.ch

UNITED ARAB EMIRATES

Bauerfeind Middle East FZ-LLC
Dubai Healthcare City
Building 40, Box 505116
Dubai
P +971 4 4335-684
F +971 4 4370-344
E info@bauerfeind.ae

UNITED KINGDOM

Bauerfeind UK
85 Tottenham Court Road
London
W1T 4TQ
P +44 (0) 121 446-5353
F +44 (0) 121 446-5454
E info@bauerfeind.co.uk

USA

Bauerfeind USA, Inc.
75 14th St NE
Suite 2350
Atlanta, GA 30309
P +1 800 423-3405
P +1 404 201-7800
F +1 404 201-7839
E info@bauerfeindusa.com