



People.Health.Care.

VENOSAN® Compression stockings

Complete catalogue



Contents


Medical compression stockings and pantyhose	04	How to take the correct measurements Information on the guidelines
	06	Product overview, colours and models
	08	Product overview: Compression stockings
	10	VENOSAN 3000 COTTON NATURAL LOOK
	12	VENOSAN 4000 CLIMATE EFFECT
	17	VENOSAN 5000 THE MAGIC OF SILVER
	22	VENOSAN 6000 CLASSIC
	24	VENOSAN 7000 NATURAL TOUCH
Compression armsleeves	30	How to take the correct measurements
	31	VENOSAN 4000 CLIMATE EFFECT
	32	VENOSAN 7000 NATURAL TOUCH
Medical support stockings	34	Product overview: Support stockings
	36	VENOSAN LEGLINE 20/30 ELEGANT LUXURY
	40	VENOSAN SILK SUPPORT TRAVEL & LEISURE
	42	VENOSAN COTTON SUPPORT TRAVEL & LEISURE
Travel socks	44	VENOSAN JET LEGS TRAVEL WITH COMFORT SOCKS

How to take the correct measurements

Medical compression stockings

When? Take the measurements in the morning if swelling tends to be severe.

How? Standing, with a tape measure.



Sizing charts

Length

Size	Version	Length (cm)
lD	AD / AD with top band, short	35–40
	AD / AD with top band, long	37–42
lG	AGH short	63–72
	AGH long	68–78
lK	AGG / AT / AT Materna short	67–74
	AGG / AT / AT Materna long	74–85

Circumference

Size	S (cm)	M (cm)	L (cm)	XL (cm)
cB	20–23	23–26	26–29	29–31
cC	29–39	33–42	36–45	39–47
cD	28–38	32–41	35–44	38–46
cG	48–63	54–67	58–72	63–74

Fit

Sizes

It is important that the stockings provide the pressure prescribed by the doctor with the correct pressure gradient. The pressure exerted on the leg is determined by two factors:

- Laplace's law of physics, which states that pressure decreases when the circumference increases at the same tension and vice versa.
- The greater the degree of stretch (lateral extension), the greater the tension and vice versa.

To some extent, one effect compensates for the other, i.e. larger circumference = loss of pressure due to the circumference, but an increase in pressure due to the material being stretched more.

Length

Our VENOSAN products have soft, high lengthwise elasticity and readily adapt to the individual leg length without slipping. Furthermore, our comparative measurements demonstrate that the pressure only changes 3–5% with ± 3 cm of length.





Where should I measure?



For AD below-knee stockings

Measure

- the leg length from the lower edge of the heel to below the knee ℓD
- the ankle circumference cB
- the calf circumference cC
- the circumference below the knee cD



For AGH thigh-length stockings with stay-up top band, AGG thigh-high stockings with belt, AT pantyhose and AT Materna maternity pantyhose

Measure

- the leg length ℓG for AG and AGH;
- the full leg length ℓK for AGG, AT and AT Materna
- the ankle circumference cB
- the calf circumference cC
- the circumference below the knee cD
- the thigh circumference cG


Then check the sizing chart to

- find the appropriate size for the measured ankle circumference cB (Small, Medium, Large or X-large)
- and for the other circumferences cC , cD and cG

Product overview

Medical compression stockings and pantyhose

Models:





Line	CCL	Ankle compression (mmHg)	AD	AD with top band	AG	AG with top band	AG with lace top band	AGG	AT	AT Materna	ATH
VENOSAN 3000 COTTON NATURAL LOOK	CCL 2	23–32	●	●	○	●	○	○	○		○
	CCL 3	34–46	●	○	○	○	○	○	○		○
VENOSAN 4000 CLIMATE EFFECT	CCL 1	18–21	●	●	○	●	●	○	●	●	
	CCL 2	23–32	●	●	○	●	●	●	●	●	
VENOSAN 5000 THE MAGIC OF SILVER	CCL 1	18–21	●	●	○	●	●	●	●	○	
	CCL 2	23–32	●	●	○	●	●	●	●	○	
	CCL 3	34–46	●	○	○	○	○	○	○	○	
VENOSAN 6000 CLASSIC	CCL 1	18–21	⊙			⊙			⊙		
	CCL 2	23–32	⊙			⊙		⊙	⊙		
	CCL 3	34–46	⊙			⊙			⊙		
VENOSAN 7000 NATURAL TOUCH	CCL 1	18–21	●	●	○	●	○	●	●		○
	CCL 2	23–32	●	●	○	●	○	●	●		○
	CCL 3	34–46	●	●	○	●	○	○	●		○

available in standard sizes
 available in standard and custom sizes
 available in custom sizes only

Colours and models

Medical compression stockings and pantyhose

Colours:

Line	Colours	Actual product colours may vary slightly from the colours in this catalogue due to the printing process.
VENOSAN 3000 COTTON NATURAL LOOK	Beige	
VENOSAN 4000 CLIMATE EFFECT	Mexico	
	Morocco	
	Sand	
	Marine	
	Anthracite	
	Black	
VENOSAN 5000 THE MAGIC OF SILVER	Mexico	
	Morocco	
	Sand	
	Black	
VENOSAN 6000 CLASSIC	Beige	
	Black	
VENOSAN 7000 NATURAL TOUCH	Beige	
	Black	


Model overview:



AD
Below-knee
stocking




**AD with
top band**
Below-knee
stocking with
dotted top band




AG
Thigh-length
stocking




**AG with
top band**
Thigh-length
stocking with
dotted top band




**AG with lace
top band**
Thigh-length
stocking with lace
top band




AGG
Thigh-length
stocking with belt



AT
Pantyhose



ATH
Pantyhose
for men



AT Materna
Maternity
pantyhose

Product overview

Medical compression stockings and pantyhose

VENOSAN® 3000 COTTON

NATURAL LOOK

- Durable
- 100% cotton on the skin
- Sturdy
- Easy to don

Indications

- CVI (chronic venous insufficiency) 2–3
- Lipoedema stage 2–3
- Following venous surgery
- PTS (post-thrombotic syndrome)
- Mild lymphoedema



CCL 2 CCL 3

Composition

CCL 2 49% cotton
27% nylon polyamide
24% elastane
CCL 3 40% cotton
27% nylon polyamide
33% elastane

VENOSAN® 4000

CLIMATE EFFECT

- Climate effect
- Semi-transparent
- Comfortable
- Extremely elastic

Indications

- CVI (chronic venous insufficiency) 1
- Pregnancy-induced varicose veins
- Milder symptoms
- Following venous surgery



CCL 1 CCL 2

Composition

71% TACTEL® polyamide
29% elastane

VENOSAN® 5000

THE MAGIC OF SILVER

- Odour-inhibiting
- Anti-static
- Temperature-equalising

Indications

- CVI (chronic venous insufficiency) 1–2
- Following venous surgery
- Mild skin conditions
- Lipoedema stage 1



CCL 1 CCL 2 CCL 3

Composition

40% TACTEL® polyamide
30% elastane
20% cotton
10% Silber X-Static®

VENOSAN® 6000

CLASSIC

- Slight transparency
- Made of highly elastic yarn for a perfect fit
- Economical stocking which nevertheless satisfies all requirements

Indications

- CVI (chronic venous insufficiency) 1–2
- Varicose veins
- Varicosis during pregnancy
- After vein surgery



Composition

70% Nylon Polyamide
30% Elastane

VENOSAN® 7000

NATURAL TOUCH

- High stiffness (resistance)
- Sturdy
- Skin-friendly
- Temperature-equalising

Indications

- CVI (chronic venous insufficiency) 2–3
- Lipoedema stage 2–3
- Following venous surgery
- Following acute thrombosis (following dressing)
- PTS (post-thrombotic syndrome)
- Mild lymphoedema
- Marked orthostatic stress



Composition

62% TACTEL® polyamide
26% elastane
8% cotton
4% SeaCell® (algae fibre)

NATURAL LOOK

VENOSAN® 3000 COTTON

Cotton – natural and soft

- **Natural**

The natural cotton on the skin increases wearing comfort and is ideal for wearers who are allergic to synthetic fibres. 100% cotton on the skin.

- **Sturdy**

VENOSAN 3000 COTTON impresses wearers with its durable, hard-wearing quality, making it ideal for long-term compression therapy.

- **Easy to don**

The elastic cotton knit fabric makes the compression stockings easy to don.

Composition

CCL 2: 49% cotton, 27% nylon polyamide, 24% elastane

CCL 3: 40% cotton, 27% nylon polyamide, 33% elastane



Standard product range

VENOSAN® 3000 COTTON

Compression class 2: Ankle compression 23–32 mmHg

CCL
2

Item	Colour	Length	Toe	Top band	REF/Size				
					XS	S	M	L	XL
AD	Beige	long	open	–		SB 32001	SB 32002	SB 32003	SB 32004
		long	closed	–		SB 32011	SB 32012	SB 32013	SB 32014
		short	open	–		SB 32021	SB 32022	SB 32023	SB 32024
		short	closed	–		SB 32031	SB 32032	SB 32033	SB 32034
AD with top band	Beige	long	open	dotted		SB 32041	SB 32042	SB 32043	SB 32044
		long	closed	dotted		SB 32051	SB 32052	SB 32053	SB 32054
		short	open	dotted		SB 32061	SB 32062	SB 32063	SB 32064
		short	closed	dotted		SB 32071	SB 32072	SB 32073	SB 32074
AGH	Beige	long	open	dotted		SB 32601	SB 32602	SB 32603	SB 32604
		short	open	dotted		SB 32621	SB 32622	SB 32623	SB 32624
AT	Beige	long	open	–					
		short	open	–					

VENOSAN® 3000 COTTON

Compression class 3: Ankle compression 34–46 mmHg

CCL
3

Item	Colour	Length	Toe	Top band	REF/Size				
					XS	S	M	L	XL
AD	Beige	long	open	–		SB 33001	SB 33002	SB 33003	SB 33004
		short	open	–		SB 33021	SB 33022	SB 33023	SB 33024
AT	Beige	long	open	–					
		short	open	–					

Custom production

VENOSAN® 3000 COTTON

Custom sizes

CCL
2

CCL
3

Item	Colour	Toe		
AD	Beige	open/closed	●	●
AD with top band		open/closed	●	●
AG		open/closed	●	●
AGH		open/closed	●	●
AGG		open/closed	●	●
AT		open/closed	●	●
ATH		open/closed	●	●

● = available

0004

CLIMATE EFFECT VENOSAN® 4000

Microfibre – fine and stylish

■ **Climate effect**

With its refined sandwich construction, VENOSAN 4000 is the first compression stocking to feature the TACTEL® climate effect label. This makes it more comfortable to wear.

■ **Extremely elastic**

New technology allows the force/stretch curve to be much flatter. It is very easy to don.

■ **Comfortable**

The premium quality materials ensure high wearing comfort.

Composition

71% TACTEL® polyamide, 29% elastane



Standard product range

VENOSAN® 4000

Compression class 1: Ankle compression 18–21 mmHg



Item	Colour	Length	Toe	Top band	REF / Size			
					S	M	L	XL
AD	Mexico	long	open	–	SG 41001	SG 41002	SG 41003	SG 41004
		long	closed	–	SG 41011	SG 41012	SG 41013	SG 41014
		short	open	–	SG 41021	SG 41022	SG 41023	SG 41024
		short	closed	–	SG 41031	SG 41032	SG 41033	SG 41034
	Morocco	long	open	–	SM 41001	SM 41002	SM 41003	SM 41004
		long	closed	–	SM 41011	SM 41012	SM 41013	SM 41014
		short	open	–	SM 41021	SM 41022	SM 41023	SM 41024
		short	closed	–	SM 41031	SM 41032	SM 41033	SM 41034
	Sand	long	open	–	SD 41001	SD 41002	SD 41003	SD 41004
		long	closed	–	SD 41011	SD 41012	SD 41013	SD 41014
		short	open	–	SD 41021	SD 41022	SD 41023	SD 41024
		short	closed	–	SD 41031	SD 41032	SD 41033	SD 41034
	Black	long	open	–	SN 41001	SN 41002	SN 41003	SN 41004
		long	closed	–	SN 41011	SN 41012	SN 41013	SN 41014
		short	open	–	SN 41021	SN 41022	SN 41023	SN 41024
		short	closed	–	SN 41031	SN 41032	SN 41033	SN 41034
AD with top band	Mexico	long	open	dotted	SG 41041	SG 41042	SG 41043	SG 41044
		long	closed	dotted	SG 41051	SG 41052	SG 41053	SG 41054
		short	open	dotted	SG 41061	SG 41062	SG 41063	SG 41064
		short	closed	dotted	SG 41071	SG 41072	SG 41073	SG 41074
	Morocco	long	open	dotted	SM 41041	SM 41042	SM 41043	SM 41044
		long	closed	dotted	SM 41051	SM 41052	SM 41053	SM 41054
		short	open	dotted	SM 41061	SM 41062	SM 41063	SM 41064
		short	closed	dotted	SM 41071	SM 41072	SM 41073	SM 41074
	Black	long	open	dotted	SN 41041	SN 41042	SN 41043	SN 41044
		long	closed	dotted	SN 41051	SN 41052	SN 41053	SN 41054
		short	open	dotted	SN 41061	SN 41062	SN 41063	SN 41064
		short	closed	dotted	SN 41071	SN 41072	SN 41073	SN 41074
AGH	Mexico	long	closed	lace	SG 41601	SG 41602	SG 41603	SG 41604
		long	open	dotted	SG 41611	SG 41612	SG 41613	SG 41614
		long	closed	dotted	SG 41621	SG 41622	SG 41623	SG 41624
		short	closed	lace	SG 41631	SG 41632	SG 41633	SG 41634
		short	open	dotted	SG 41641	SG 41642	SG 41643	SG 41644
		short	closed	dotted	SG 41651	SG 41652	SG 41653	SG 41654
	Morocco	long	closed	lace	SM 41601	SM 41602	SM 41603	SM 41604
		long	open	dotted	SM 41611	SM 41612	SM 41613	SM 41614
		long	closed	dotted	SM 41621	SM 41622	SM 41623	SM 41624
		short	closed	lace	SM 41631	SM 41632	SM 41633	SM 41634
		short	open	dotted	SM 41641	SM 41642	SM 41643	SM 41644
		short	closed	dotted	SM 41651	SM 41652	SM 41653	SM 41654
	Sand	long	closed	lace	SD 41601	SD 41602	SD 41603	SD 41604
		long	open	dotted	SD 41611	SD 41612	SD 41613	SD 41614
		long	closed	dotted	SD 41621	SD 41622	SD 41623	SD 41624
		short	closed	lace	SD 41631	SD 41632	SD 41633	SD 41634
		short	open	dotted	SD 41641	SD 41642	SD 41643	SD 41644
		short	closed	dotted	SD 41651	SD 41652	SD 41653	SD 41654

Standard product range

VENOSAN® 4000

Compression class 1: Ankle compression 18–21 mmHg

CCL
1

Item	Colour	Length	Toe	Top band	REF / Size			
					S	M	L	XL
AGH	Black	long	open	dotted	SN 41611	SN 41612	SN 41613	SN 41004
		long	closed	dotted	SN 41621	SN 41622	SN 41623	SN 41624
		long	closed	lace	SN 41601	SN 41602	SN 41603	SN 41604
		short	open	dotted	SN 41641	SN 41642	SN 41643	SN 41644
		short	closed	dotted	SN 41651	SN 41652	SN 41653	SN 41654
		short	closed	lace	SN 41631	SN 41632	SN 41633	SN 41634
AT	Mexico	long	closed	–	SG 41221	SG 41222	SG 41223	SG 41224
		short	closed	–	SG 41231	SG 41232	SG 41233	SG 41234
	Morocco	long	closed	–	SM 41221	SM 41222	SM 41223	SM 41224
		short	closed	–	SM 41231	SM 41232	SM 41233	SM 41234
	Sand	long	closed	–	SD 41221	SD 41222	SD 41223	SD 41224
		short	closed	–	SD 41231	SD 41232	SD 41233	SD 41234
	Black	long	closed	–	SN 41221	SN 41222	SN 41223	SN 41224
		short	closed	–	SN 41231	SN 41232	SN 41233	SN 41234
AT Materna	Mexico	long	closed	–	SG 41301	SG 41302	SG 41303	SG 41304
	Black	short	closed	–	SN 41301	SN 41302	SN 41303	SN 41304

VENOSAN® 4000

Compression class 2: Ankle compression 23–32 mmHg

CCL
2

Item	Colour	Length	Toe	Top band	REF / Size			
					S	M	L	XL
AD	Mexico	long	open	–	SG 42001	SG 42002	SG 42003	SG 42004
		long	closed	–	SG 42011	SG 42012	SG 42013	SG 42014
		short	open	–	SG 42021	SG 42022	SG 42023	SG 42024
		short	closed	–	SG 42031	SG 42032	SG 42033	SG 42034
	Morocco	long	open	–	SM 42001	SM 42002	SM 42003	SM 42004
		long	closed	–	SM 42011	SM 42012	SM 42013	SM 42014
		short	open	–	SM 42021	SM 42022	SM 42023	SM 42024
		short	closed	–	SM 42031	SM 42032	SM 42033	SM 42034
	Sand	long	open	–	SD 42001	SD 42002	SD 42003	SD 42004
		long	closed	–	SD 42011	SD 42012	SD 42013	SD 42014
		short	open	–	SD 42021	SD 42022	SD 42023	SD 42024
		short	closed	–	SD 42031	SD 42032	SD 42033	SD 42034
	Black	long	open	–	SN 42001	SN 42002	SN 42003	SN 42004
		long	closed	–	SN 42011	SN 42012	SN 42013	SN 42014
		short	open	–	SN 42021	SN 42022	SN 42023	SN 42024
		short	closed	–	SN 42031	SN 42032	SN 42033	SN 42034
AD with top band	Mexico	long	open	dotted	SG 42041	SG 42042	SG 42043	SG 42044
		long	closed	dotted	SG 42051	SG 42052	SG 42053	SG 42054
		short	open	dotted	SG 42061	SG 42062	SG 42063	SG 42064
		short	closed	dotted	SG 42071	SG 42072	SG 42073	SG 42074

VENOSAN® 4000

Compression class 2: Ankle compression 23–32 mmHg



Item	Colour	Length	Toe	Top band	REF/Size			
					S	M	L	XL
AD with top band	Morocco	long	open	dotted	SM 42041	SM 42042	SM 42043	SM 42044
		long	closed	dotted	SM 42051	SM 42052	SM 42053	SM 42054
		short	open	dotted	SM 42061	SM 42062	SM 42063	SM 42064
		short	closed	dotted	SM 42071	SM 42072	SM 42073	SM 42074
	Black	long	open	dotted	SN 42041	SN 42042	SN 42043	SN 42044
		long	closed	dotted	SN 42051	SN 42052	SN 42053	SN 42054
		short	open	dotted	SN 42061	SN 42062	SN 42063	SN 42064
		short	closed	dotted	SN 42071	SN 42072	SN 42073	SN 42074
AGH	Mexico	long	closed	lace	SG 42601	SG 42602	SG 42603	SG 42604
		long	open	dotted	SG 42611	SG 42612	SG 42613	SG 42614
		long	closed	dotted	SG 42621	SG 42622	SG 42623	SG 42624
		short	closed	lace	SG 42631	SG 42632	SG 42633	SG 42634
		short	open	dotted	SG 42641	SG 42642	SG 42643	SG 42644
		short	closed	dotted	SG 42651	SG 42652	SG 42653	SG 42654
	Morocco	long	closed	lace	SM 42601	SM 42602	SM 42603	SM 42604
		long	open	dotted	SM 42611	SM 42612	SM 42613	SM 42614
		long	closed	dotted	SM 42621	SM 42622	SM 42623	SM 42624
		short	closed	lace	SM 42631	SM 42632	SM 42633	SM 42634
		short	open	dotted	SM 42641	SM 42642	SM 42643	SM 42644
		short	closed	dotted	SM 42651	SM 42652	SM 42653	SM 42654
	Sand	long	closed	lace	SD 42601	SD 42602	SD 42603	SD 42604
		long	open	dotted	SD 42611	SD 42612	SD 42613	SD 42614
		long	closed	dotted	SD 42621	SD 42622	SD 42623	SD 42624
		short	closed	lace	SD 42631	SD 42632	SD 42633	SD 42634
		short	open	dotted	SD 42641	SD 42642	SD 42643	SD 42644
		short	closed	dotted	SD 42651	SD 42652	SD 42653	SD 42654
	Black	long	closed	lace	SN 42601	SN 42602	SN 42603	SN 42604
		long	open	dotted	SN 42611	SN 42612	SN 42613	SN 42614
		long	closed	dotted	SN 42621	SN 42622	SN 42623	SN 42624
		short	closed	lace	SN 42631	SN 42632	SN 42633	SN 42634
		short	open	dotted	SN 42641	SN 42642	SN 42643	SN 42644
		short	closed	dotted	SN 42651	SN 42652	SN 42653	SN 42654
AGG (left)	Mexico	long	open	–	SG 42151	SG 42152	SG 42153	SG 42154
AGG (right)		long	open	–	SG 42171	SG 42172	SG 42173	SG 42174
AT	Mexico	long	closed	–	SG 42221	SG 42222	SG 42223	SG 42224
		short	closed	–	SG 42231	SG 42232	SG 42233	SG 42234
	Morocco	long	closed	–	SM 42221	SM 42222	SM 42223	SM 42224
		short	closed	–	SM 42231	SM 42232	SM 42233	SM 42234
	Sand	long	closed	–	SD 42221	SD 42222	SD 42223	SD 42224
		short	closed	–	SD 42231	SD 42232	SD 42233	SD 42234
	Black	long	closed	–	SN 42221	SN 42222	SN 42223	SN 42224
		short	closed	–	SN 42231	SN 42232	SN 42233	SN 42234
AT Materna	Mexico	long	closed	–	SG 42301	SG 42302	SG 42303	SG 42304
		short	closed	–	SG 42311	SG 42312	SG 42313	SG 42314
	Black	long	closed	–	SN 42301	SN 42302	SN 42303	SN 42304
		short	closed	–	SN 42311	SN 42312	SN 42313	SN 42314

Custom production

VENOSAN® 4000

Custom sizes

CCL
1

CCL
2

Item	Colour	Toe		
AD	CCL 1 Mexico, Black	open/closed	●	●
AD with top band		open/closed	●	●
AG	CCL 2 Mexico, Morocco, Sand, Black	open/closed	●	●
AGH		open/closed	●	●
AGG+		open/closed	●	●
AT		open/closed	●	●
AT Materna		open/closed	●	●

+ available in Mexico and black only

● = available

THE MAGIC OF SILVER

VENOSAN® 5000

Silver & cotton fibre

■ Pure silver

Contains silver-coated X-Static® fibre. No reduction of the properties of silver even after 150 washing cycles.

■ Temperature-equalising

The fine mesh stocking has an equalising action and transports moisture from the inside to the outside of the stocking. This ensures an ideal temperature sensation and optimal wearing comfort.

■ Anti-static

The silver fibres conduct away the static electricity associated with everyday clothing and increase well-being.

Composition

40% TACTEL® polyamide, 30% elastane, 20% cotton, 10% Silver X-Static®



Standard product range

VENOSAN® 5000

Compression class 1: Ankle compression 18–21 mmHg

CCL
1

Item	Colour	Length	Toe	Top band	REF / Size			
					S	M	L	XL
AD	Mexico	long	open	–	SG 45001	SG 45002	SG 45003	SG 45004
		long	closed	–	SG 45011	SG 45012	SG 45013	SG 45014
		short	open	–	SG 45021	SG 45022	SG 45023	SG 45024
		short	closed	–	SG 45031	SG 45032	SG 45033	SG 45034
	Morocco	long	open	–	SM 45001	SM 45002	SM 45003	SM 45004
		long	closed	–	SM 45011	SM 45012	SM 45013	SM 45014
		short	open	–	SM 45021	SM 45022	SM 45023	SM 45024
		short	closed	–	SM 45031	SM 45032	SM 45033	SM 45034
	Sand	long	open	–	SD 45001	SD 45002	SD 45003	SD 45004
		long	closed	–	SD 45011	SD 45012	SD 45013	SD 45014
		short	open	–	SD 45021	SD 45022	SD 45023	SD 45024
		short	closed	–	SD 45031	SD 45032	SD 45033	SD 45034
	Black	long	open	–	SN 45001	SN 45002	SN 45003	SN 45004
		long	closed	–	SN 45011	SN 45012	SN 45013	SN 45014
		short	open	–	SN 45021	SN 45022	SN 45023	SN 45024
		short	closed	–	SN 45031	SN 45032	SN 45033	SN 45034
AD with top band	Mexico	long	open	dotted	SG 45041	SG 45042	SG 45043	SG 45044
		long	closed	dotted	SG 45051	SG 45052	SG 45053	SG 45054
		short	open	dotted	SG 45061	SG 45062	SG 45063	SG 45064
		short	closed	dotted	SG 45071	SG 45072	SG 45073	SG 45074
	Morocco	long	open	dotted	SM 45041	SM 45042	SM 45043	SM 45044
		long	closed	dotted	SM 45051	SM 45052	SM 45053	SM 45054
		short	open	dotted	SM 45061	SM 45062	SM 45063	SM 45064
		short	closed	dotted	SM 45071	SM 45072	SM 45073	SM 45074
	Black	long	open	dotted	SN 45041	SN 45042	SN 45043	SN 45044
		long	closed	dotted	SN 45051	SN 45052	SN 45053	SN 45054
		short	open	dotted	SN 45061	SN 45062	SN 45063	SN 45064
		short	closed	dotted	SN 45071	SN 45072	SN 45073	SN 45074
AGH	Mexico	long	closed	lace	SG 45601	SG 45602	SG 45603	SG 45604
		long	open	dotted	SG 45611	SG 45612	SG 45613	SG 45614
		short	closed	lace	SG 45631	SG 45632	SG 45633	SG 45634
		short	open	dotted	SG 45641	SG 45642	SG 45643	SG 45644
	Morocco	long	closed	lace	SM 45601	SM 45602	SM 45603	SM 45604
		long	open	dotted	SM 45611	SM 45612	SM 45613	SM 45614
		short	closed	lace	SM 45631	SM 45632	SM 45633	SM 45634
		short	open	dotted	SM 45641	SM 45642	SM 45643	SM 45644
	Sand	long	closed	lace	SD 45601	SD 45602	SD 45603	SD 45604
		long	open	dotted	SD 45611	SD 45612	SD 45613	SD 45614
		short	closed	lace	SD 45631	SD 45632	SD 45633	SD 45634
		short	open	dotted	SD 45641	SD 45642	SD 45643	SD 45644
	Black	long	closed	lace	SN 45601	SN 45602	SN 45603	SN 45604
		long	open	dotted	SN 45611	SN 45612	SN 45613	SN 45614
		short	closed	lace	SN 45631	SN 45632	SN 45633	SN 45634
		short	open	dotted	SN 45641	SN 45642	SN 45643	SN 45644

Standard product range

VENOSAN® 5000

Compression class 1: Ankle compression 18–21 mmHg

CCL
1

Item	Colour	Length	Toe	Top band	REF/Size			
					S	M	L	XL
AGG (left)	Mexico	long	open	–	SG 45151	SG 45152	SG 45153	SG 45154
AGG (right)		long	open	–	SG 45171	SG 45172	SG 45173	SG 45174
AT	Mexico	long	closed	–	SG 45221	SG 45222	SG 45223	SG 45224
		short	closed	–	SG 45231	SG 45232	SG 45233	SG 45234
	Morocco	long	closed	–	SM 45221	SM 45222	SM 45223	SM 45224
		short	closed	–	SM 45231	SM 45232	SM 45233	SM 45234
	Sand	long	closed	–	SD 45221	SD 45222	SD 45223	SD 45224
		short	closed	–	SD 45231	SD 45232	SD 45233	SD 45234
	Black	long	closed	–	SN 45221	SN 45222	SN 45223	SN 45224
		short	closed	–	SN 45231	SN 45232	SN 45233	SN 45234

VENOSAN® 5000

Compression class 2: Ankle compression 23–32 mmHg

CCL
2

Item	Colour	Length	Toe	Top band	REF/Size			
					S	M	L	XL
AD	Mexico	long	open	–	SG 46001	SG 46002	SG 46003	SG 46004
		long	closed	–	SG 46011	SG 46012	SG 46013	SG 46014
		short	open	–	SG 46021	SG 46022	SG 46023	SG 46024
		short	closed	–	SG 46031	SG 46032	SG 46033	SG 46034
	Morocco	long	open	–	SM 46001	SM 46002	SM 46003	SM 46004
		long	closed	–	SM 46011	SM 46012	SM 46013	SM 46014
		short	open	–	SM 46021	SM 46022	SM 46023	SM 46024
		short	closed	–	SM 46031	SM 46032	SM 46033	SM 46034
	Sand	long	open	–	SD 46001	SD 46002	SD 46003	SD 46004
		long	closed	–	SD 46011	SD 46012	SD 46013	SD 46014
		short	open	–	SD 46021	SD 46022	SD 46023	SD 46024
		short	closed	–	SD 46031	SD 46032	SD 46033	SD 46034
	Black	long	open	–	SN 42001	SN 42002	SN 42003	SN 42004
		long	closed	–	SN 46011	SN 46012	SN 46013	SN 46014
		short	open	–	SN 46021	SN 46022	SN 46023	SN 46024
		short	closed	–	SN 46031	SN 46032	SN 46033	SN 46034
AD with top band	Mexico	long	open	dotted	SG 46041	SG 46042	SG 46043	SG 46044
		long	closed	dotted	SG 46051	SG 46052	SG 46053	SG 46054
		short	open	dotted	SG 46061	SG 46062	SG 46063	SG 46064
		short	closed	dotted	SG 46071	SG 46072	SG 46073	SG 46074
	Morocco	long	open	dotted	SM 46041	SM 46042	SM 46043	SM 46044
		long	closed	dotted	SM 46051	SM 46052	SM 46053	SM 46054
		short	open	dotted	SM 46061	SM 46062	SM 46063	SM 46064
		short	closed	dotted	SM 46071	SM 46072	SM 46073	SM 46074
	Black	long	open	dotted	SN 46041	SN 46042	SN 46043	SN 46044
		long	closed	dotted	SN 46051	SN 46052	SN 46053	SN 46054
		short	open	dotted	SN 46061	SN 46062	SN 46063	SN 46064
		short	closed	dotted	SN 46071	SN 46072	SN 46073	SN 46074

Standard product range

VENOSAN® 5000

Compression class 2: Ankle compression 23–32 mmHg

CCL
2

Item	Colour	Length	Toe	Top band	REF / Size			
					S	M	L	XL
AGH	Mexico	long	closed	lace	SG 46601	SG 46602	SG 46603	SG 46604
		long	open	dotted	SG 46611	SG 46612	SG 46613	SG 46614
		long	closed	dotted	SG 46621	SG 46622	SG 46623	SG 46624
		short	closed	lace	SG 46631	SG 46632	SG 46633	SG 46634
		short	open	dotted	SG 46641	SG 46642	SG 46643	SG 46644
		short	closed	dotted	SG 46651	SG 46652	SG 46653	SG 46654
	Morocco	long	closed	lace	SM 46601	SM 46602	SM 46603	SM 46604
		long	open	dotted	SM 46611	SM 46612	SM 46613	SM 46614
		long	closed	dotted	SM 46621	SM 46622	SM 46623	SM 46624
		short	closed	lace	SM 46631	SM 46632	SM 46633	SM 46634
		short	open	dotted	SM 46641	SM 46642	SM 46643	SM 46644
		short	closed	dotted	SM 46651	SM 46652	SM 46653	SM 46654
	Sand	long	closed	lace	SD 46601	SD 46602	SD 46603	SD 46604
		long	open	dotted	SD 46611	SD 46612	SD 46613	SD 46614
		long	closed	dotted	SD 46621	SD 46622	SD 46623	SD 46624
		short	closed	lace	SD 46631	SD 46632	SD 46633	SD 46634
		short	open	dotted	SD 46641	SD 46642	SD 46643	SD 46644
		short	closed	dotted	SD 46651	SD 46652	SD 46653	SD 46654
	Black	long	closed	lace	SN 46601	SN 46602	SN 46603	SN 46604
		long	open	dotted	SN 46611	SN 46612	SN 46613	SN 46614
		long	closed	dotted	SN 46621	SN 46622	SN 46623	SN 46624
		short	closed	lace	SN 46631	SN 46632	SN 46633	SN 46634
		short	open	dotted	SN 46641	SN 46642	SN 46643	SN 46644
		short	closed	dotted	SN 46651	SN 46652	SN 46653	SN 46654
AGG (left)	Mexico	long	open	–	SG 46151	SG 46152	SG 46153	SG 46154
AGG (right)		long	open	–	SG 46171	SG 46172	SG 46173	SG 46174
AT	Mexico	long	closed	–	SG 46221	SG 46222	SG 46223	SG 46224
		short	closed	–	SG 46231	SG 46232	SG 46233	SG 46234
	Morocco	long	closed	–	SM 46221	SM 46222	SM 46223	SM 46224
		short	closed	–	SM 46231	SM 46232	SM 46233	SM 46234
	Sand	long	closed	–	SD 46221	SD 46222	SD 46223	SD 46224
		short	closed	–	SD 46231	SD 46232	SD 46233	SD 46234
	Black	long	closed	–	SN 46221	SN 46222	SN 46223	SN 46224
		short	closed	–	SN 46231	SN 46232	SN 46233	SN 46234

Standard product range

VENOSAN® 5000

Compression class 3: Ankle compression 34–46 mmHg

CCL
3

Item	Colour	Length	Toe	Top band	REF/Size			
					S	M	L	XL
AD	Mexico	long	open	–	SG 47001	SG 47002	SG 47003	SG 47004
		long	closed	–	SG 47011	SG 47012	SG 47013	SG 47014
	Morocco	long	open	–	SM 47001	SM 47002	SM 47003	SM 47004
		long	closed	–	SM 47011	SM 47012	SM 47013	SM 47014
	Black	long	open	–	SN 47001	SN 47002	SN 47003	SN 47004
		long	closed	–	SN 47011	SN 47012	SN 47013	SN 47014

Custom production

VENOSAN® 5000

Custom sizes

CCL
1

CCL
2

Item	Colour	Toe		
AD	CCL 1 Mexico, Black	open/closed	●	●
AD with top band		open/closed	●	●
AG		open/closed	●	●
AGH	CCL 2 Mexico, Morocco, Sand, Black	open/closed	●	●
AGG+		open/closed	●	●
AT		open/closed	●	●
AT Materna		open/closed	●	●

+ available in Mexico and black only

● = available

CCL
3

Item	Colour	Toe	
AD	Mexico	open/closed	●
AD with top band		open/closed	●
AG		open/closed	●
AGH		open/closed	●
AGG		open/closed	●
AT		open/closed	●
AT Materna		open/closed	●

● = available

CLASSIC VENOSAN® 6000

Economical & reliable

- Slight transparency
- Made of highly elastic yarn for a perfect fit
- Economical stocking which nevertheless satisfies all requirements

Material composition

70% Nylon Polyamide
30% Elastane



Standard product range

VENOSAN® 6000

Compression class 1: pressure at the ankle 18–21 mm Hg

CCL
1

Item	Colour	Length	Toe	Self supporting top	REF/Size				
					XS	S	M	L	XL
AD	Beige	long	open	–	CG61000	CG61001	CG61002	CG61003	CG61004
		short			CG61020	CG61021	CG61022	CG61023	CG61024
		long	closed	–	CG61010	CG61011	CG61012	CG61013	CG61014
		short			CG61030	CG61031	CG61032	CG61033	CG61034
	Black	long	open	–	CN61000	CN61001	CN61002	CN61003	CN61004
		short			CN61020	CN61021	CN61022	CN61023	CN61024
		long	closed	–	CN61010	CN61011	CN61012	CN61013	CN61014
		short			CN61030	CN61031	CN61032	CN61033	CN61034
AGH	Beige	long	open	plain	CG61600	CG61601	CG61602	CG61603	CG61604
		short			CG61620	CG61621	CG61622	CG61623	CG61624
		long	closed	plain	CG61610	CG61611	CG61612	CG61613	CG61614
		short			CG61630	CG61631	CG61632	CG61633	CG61634
	Black	long	open	plain	CN61600	CN61601	CN61602	CN61603	CN61604
		short			CN61620	CN61621	CN61622	CN61623	CN61624
		long	closed	plain	CN61610	CN61611	CN61612	CN61613	CN61614
		short			CN61630	CN61631	CN61632	CN61633	CN61634
AG	Beige	long	open	–					
		short							
		long	closed	–					
		short							
	Black	long	open	–					
		short							
		long	closed	–					
		short							
AT	Beige	long	open	–	CG61200	CG61201	CG61202	CG61203	CG61204
		short			CG61220	CG61221	CG61222	CG61223	CG61224
		long	closed	–	CG61210	CG61211	CG61212	CG61213	CG61214
		short			CG61230	CG61231	CG61232	CG61233	CG61234
	Black	long	open	–	CN61200	CN61201	CN61202	CN61203	CN61204
		short			CN61220	CN61221	CN61222	CN61223	CN61224
		long	closed	–	CN61210	CN61211	CN61212	CN61213	CN61214
		short			CN61230	CN61231	CN61232	CN61233	CN61234

Standard product range

VENOSAN® 6000

Compression class 2: pressure at the ankle 23–32 mm Hg

CCL
2

Item	Colour	Length	Toe	Self supporting top	REF/Size				
					XS	S	M	L	XL
AD	Beige	long	open	–	CG62000	CG62001	CG62002	CG62003	CG62004
		short			CG62020	CG62021	CG62022	CG62023	CG62024
		long	closed	–	CG62010	CG62011	CG62012	CG62013	CG62014
		short			CG62030	CG62031	CG62032	CG62033	CG62034
	Black	long	open	–	CN62000	CN62001	CN62002	CN62003	CN62004
		short			CN62020	CN62021	CN62022	CN62023	CN62024
		long	closed	–	CN62010	CN62011	CN62012	CN62013	CN62014
		short			CN62030	CN62031	CN62032	CN62033	CN62034
AGH	Beige	long	open	plain	CG62600	CG62601	CG62602	CG62603	CG62604
		short			CG62620	CG62621	CG62622	CG62623	CG62624
		long	closed	plain	CG62610	CG62611	CG62612	CG62613	CG62614
		short			CG62630	CG62631	CG62632	CG62633	CG62634
	Black	long	open	plain	CN62600	CN62601	CN62602	CN62603	CN62604
		short			CN62620	CN62621	CN62622	CN62623	CN62624
		long	closed	plain	CN62610	CN62611	CN62612	CN62613	CN62614
		short			CN62630	CN62631	CN62632	CN62633	CN62634
AG	Beige	long	open	–					
		short							
		long	closed	–					
		short							
	Black	long	open	–					
		short							
		long	closed	–					
		short							
AGGL	Beige	long	open	left	–				
		short			–				
AGGR		long	closed	right	–				
		short			–				
AT	Beige	long	open	–	CG62200	CG62201	CG62202	CG62203	CG62204
		short			CG62220	CG62221	CG62222	CG62223	CG62224
		long	closed	–	CG62210	CG62211	CG62212	CG62213	CG62214
		short			CG62230	CG62231	CG62232	CG62233	CG62234
	Black	long	open	–	CN62200	CN62201	CN62202	CN62203	CN62204
		short			CN62220	CN62221	CN62222	CN62223	CN62224
		long	closed	–	CN62210	CN62211	CN62212	CN62213	CN62214
		short			CN62230	CN62231	CN62232	CN62233	CN62234

Standard product range

VENOSAN® 6000

Compression class 3: pressure at the ankle 34–46 mmHg

CCL
3

Item	Colour	Length	Toe	Self supporting top	REF/Size				
					XS	S	M	L	XL
AD	Beige	long	open	–	CG63000	CG63001	CG63002	CG63003	CG63004
		short			CG63020	CG63021	CG63022	CG63023	CG63024
AGH	Beige	long	open	plain	CG63600	CG63601	CG63602	CG63603	CG63604
		short			CG63620	CG63621	CG63622	CG63623	CG63624
AT	Beige	long	open	–	CG63200	CG63201	CG63202	CG63203	CG63204
		short			CG63220	CG63221	CG63222	CG63223	CG63224

NATURAL TOUCH VENOSAN® 7000

Strong natural fibres – high stiffness

■ **High stiffness**

With its high stiffness – the increase in pressure related to the increase in circumference due to muscle movement – VENOSAN 7000 also exerts a pronounced massaging action on the leg.

■ **Skin-friendly**

The composition ensures both high resistance and optimal wearing comfort.

Composition

62% TACTEL® polyamide, 26% elastane, 8% cotton
4% SeaCell® (algae fibre)



Standard product range

VENOSAN® 7000

Compression class 1: Ankle compression 18–21 mmHg

CCL
1

Item	Colour	Length	Toe	Top band	REF/Size			
					S	M	L	XL
AD	Beige	long	open	–	SG71001	SG71002	SG71003	SG71004
		short			SG71021	SG71022	SG71023	SG71024
		long	closed	–	SG71011	SG71012	SG71013	SG71014
		short			SG71031	SG71032	SG71033	SG71034
	Black	long	open	–	SN71001	SN71002	SN71003	SN71004
		short			SN71021	SN71022	SN71023	SN71024
		long	closed	–	SN71011	SN71012	SN71013	SN71014
		short			SN71031	SN71032	SN71033	SN71034
ADH	Beige	long	open	plain	SG71041	SG71042	SG71043	SG71044
		short			SG71061	SG71062	SG71063	SG71064
		long	closed	plain	SG71051	SG71052	SG71053	SG71054
		short			SG71071	SG71072	SG71073	SG71074
	Black	long	open	plain	SN71041	SN71042	SN71043	SN71044
		short			SN71061	SN71062	SN71063	SN71064
		long	closed	plain	SN71051	SN71052	SN71053	SN71054
		short			SN71071	SN71072	SN71073	SN71074
AGH	Beige	long	open	plain	SG71601	SG71602	SG71603	SG71604
		short			SG71621	SG71622	SG71623	SG71624
		long	closed	plain	SG71611	SG71612	SG71613	SG71614
		short			SG71631	SG71632	SG71633	SG71634
	Black	long	open	plain	SN71601	SN71602	SN71603	SN71604
		short			SN71621	SN71622	SN71623	SN71624
		long	closed	plain	SN71611	SN71612	SN71613	SN71614
		short			SN71631	SN71632	SN71633	SN71634
AGGL	Beige	long	open	–	SG71151	SG71152	SG71153	SG71154
short		SG71161			SG71162	SG71163	SG71164	
AGGR		long	closed	–	SG71171	SG71172	SG71173	SG71174
		short			SG71181	SG71182	SG71183	SG71184
AT	Beige	long	open	–				
		short						
		long	closed	–				
		short						
	Black	long	open	–				
		short						
		long	closed	–				
		short						

Standard product range

VENOSAN® 7000

Compression class 2: Ankle compression 23–32 mmHg



Item	Colour	Length	Toe	Top band	REF/Size			
					S	M	L	XL
AD	Beige	long	open	–	SG 72001	SG 72002	SG 72003	SG 72004
		long	closed	–	SG 72011	SG 72012	SG 72013	SG 72014
		short	open	–	SG 72021	SG 72022	SG 72023	SG 72024
		short	closed	–	SG 72031	SG 72032	SG 72033	SG 72034
	Black	long	open	–	SN 72001	SN 72002	SN 72003	SN 72004
		long	closed	–	SN 72011	SN 72012	SN 72013	SN 72014
		short	open	–	SN 72021	SN 72022	SN 72023	SN 72024
		short	closed	–	SN 72031	SN 72032	SN 72033	SN 72034
AD with top band	Beige	long	open	dotted	SG 72041	SG 72042	SG 72043	SG 72044
		long	closed	dotted	SG 72051	SG 72052	SG 72053	SG 72054
		short	open	dotted	SG 72061	SG 72062	SG 72063	SG 72064
		short	closed	dotted	SG 72071	SG 72072	SG 72073	SG 72074
	Black	long	open	dotted	SN 72041	SN 72042	SN 72043	SN 72044
		long	closed	dotted	SN 72051	SN 72052	SN 72053	SN 72054
		short	open	dotted	SN 72061	SN 72062	SN 72063	SN 72064
		short	closed	dotted	SN 72071	SN 72072	SN 72073	SN 72074
AGH	Beige	long	open	dotted	SG 72601	SG 72602	SG 72603	SG 72604
		long	closed	dotted	SG 72611	SG 72612	SG 72613	SG 72614
		short	open	dotted	SG 72621	SG 72622	SG 72623	SG 72624
		short	closed	dotted	SG 72631	SG 72632	SG 72633	SG 72634
	Black	long	open	dotted	SN 72601	SN 72602	SN 72603	SN 72604
		long	closed	dotted	SN 72611	SN 72612	SN 72613	SN 72614
		short	open	dotted	SN 72621	SN 72622	SN 72623	SN 72624
		short	closed	dotted	SN 72631	SN 72632	SN 72633	SN 72634
AGG (left)	Beige	long	open	–	SG 72151	SG 72152	SG 72153	SG 72154
		short	open	–	SG 72161	SG 72162	SG 72163	SG 72164
AGG (right)		long	open	–	SG 72171	SG 72172	SG 72173	SG 72174
		short	open	–	SG 72181	SG 72182	SG 72183	SG 72184
AT	Beige	long	open	–				
		short	open	–				
		long	closed	–				
		short	closed	–				
	Black	long	open	–				
		short	open	–				
		long	closed	–				
		short	closed	–				

Standard product range

VENOSAN® 7000

Compression class 3: Ankle compression 34–46 mmHg

CCL
3

Item	Colour	Length	Toe	Top band	REF/Size			
					S	M	L	XL
AD	Beige	long	open	–	SG 73001	SG 73002	SG 73003	SG 73004
		long	closed	–	SG 73011	SG 73012	SG 73013	SG 73014
		short	open	–	SG 73021	SG 73022	SG 73023	SG 73024
		short	closed	–	SG 73031	SG 73032	SG 73033	SG 73034
	Black	long	open	–	SN 73001	SN 73002	SN 73003	SN 73004
		long	closed	–	SN 73011	SN 73012	SN 73013	SN 73014
		short	open	–	SN 73021	SN 73022	SN 73023	SN 73024
		short	closed	–	SN 73031	SN 73032	SN 73033	SN 73034
AD with top band	Beige	long	open	dotted	SG 73041	SG 73042	SG 73043	SG 73044
		long	closed	dotted	SG 73051	SG 73052	SG 73053	SG 73054
		short	open	dotted	SG 73061	SG 73062	SG 73063	SG 73064
		short	closed	dotted	SG 73071	SG 73072	SG 73073	SG 73074
	Black	long	open	dotted	SN 73041	SN 73042	SN 73043	SN 73044
		long	closed	dotted	SN 73051	SN 73052	SN 73053	SN 73054
		short	open	dotted	SN 73061	SN 73062	SN 73063	SN 73064
		short	closed	dotted	SN 73071	SN 73072	SN 73073	SN 73074
AGH	Beige	long	open	dotted	SG 73601	SG 73602	SG 73603	SG 73604
		long	closed	dotted	SG 73611	SG 73612	SG 73613	SG 73614
		short	open	dotted	SG 73621	SG 73622	SG 73623	SG 73624
		short	closed	dotted	SG 73631	SG 73632	SG 73633	SG 73634
	Black	long	open	dotted	SN 73601	SN 73602	SN 73603	SN 73604
		long	closed	dotted	SN 73611	SN 73612	SN 73613	SN 73614
		short	open	dotted	SN 73621	SN 73622	SN 73623	SN 73624
		short	closed	dotted	SN 73631	SN 73632	SN 73633	SN 73634
AT	Beige	long	open	–				
		short						
		long	closed	–				
		short						
	Black	long	open	–				
		short						
		long	closed	–				
		short						

Custom production

VENOSAN® 7000

Custom sizes

CCL
1

CCL
2

CCL
3


Item	Colour	Toe	CCL 1	CCL 2	CCL 3
AD	Beige, Black	open/closed	●	●	●
AD with top band		open/closed	●	●	●
AG		open/closed	●	●	●
AGH		open/closed	●	●	●
AGG		open/closed	●	●	●
AT		open/closed	● = available	●	●
ATH		open/closed	●	●	●

How to take the correct measurements

Compression armsleeves

When? Take the measurements in the morning if swelling tends to be severe.

How? With a tape measure.



Compression armsleeves sizing chart

Wrist circumference (cm)				
	S	M	L	XL
cC	14–17	17–20	20–23	23–26


Forearm circumference (cm)				
	S	M	L	XL
cD	22–26	26–30	30–34	34–38

Upper arm circumference (cm)				
	S	M	L	XL
cF	24–30	30–36	36–42	42–48

Length (cm)		short	long
Wrist (C) to upper arm (G)		35–45	40–50

Measure the length with the arm outstretched.


Models, series and size:




CG
without
hand attachment




CG-H
without
hand attachment,
with dotted top band




ZG
with
hand attachment



ZG-H
with
hand attachment,
with dotted top band



CH
without
hand attachment,
with shoulder piece



ZH
with
hand attachment,
with shoulder piece

CLIMATE EFFECT

VENOSAN® 4000

Microfibre – fine and stylish

■ Climate effect

With its refined sandwich construction, VENOSAN 4000 is the first compression stocking to feature the TACTEL® climate effect label. This makes it more comfortable to wear.

■ Extremely elastic

New technology allows the force/stretch curve to be much flatter. It is very easy to don.

■ Comfortable

The premium quality materials ensure high wearing comfort.

Composition

71% TACTEL® polyamide, 29% elastane



VENOSAN® 4000

Compression class 2: Compression at wrist 23–32 mmHg



Item	Colour	Hand attachment	REF/Size				Custom sizes
			S	M	L	XL	
CG	Mexico	open	SG 42501	SG 42502	SG 42503	SG 42504	•
CG	Black	open	SN 42501	SN 42502	SN 42503	SN 42504	•

• = available

NATURAL TOUCH VENOSAN® 7000

Strong natural fibres – high stiffness

■ High stiffness

With its high stiffness – the increase in pressure related to the increase in circumference due to muscle movement – VENOSAN 7000 also exerts a pronounced massaging action on the arm.

■ Skin-friendly

The composition ensures both high resistance and optimal wearing comfort.

Composition

62% TACTEL® polyamide, 26% elastane, 8% cotton, 4% SeaCell®



VENOSAN® 7000

Compression class 1: Compression at wrist 18–21 mmHg

CCL
1

Item	Colour	Length	Hand attachment	Top band	REF/Size			
					S	M	L	XL
ZG	Beige	long	with		SG71541	SG71542	SG71543	SG71544
		short			SG71551	SG71552	SG71553	SG71554
	Black	long			SN71541	SN71542	SN71543	SN71544
		short			SN71551	SN71552	SN71553	SN71554
ZG-H	Beige	long	with	dotted	SG71561	SG71562	SG71563	SG71564
		short			SG71571	SG71572	SG71573	SG71574
	Black	long		dotted	SN71561	SN71562	SN71563	SN71564
		short			SN71571	SN71572	SN71573	SN71574
ZH	Beige	long	with		SG71581	SG71582	SG71583	SG71584
		short			SG71591	SG71592	SG71593	SG71594
CG	Beige	long	without		SG71501	SG71502	SG71503	SG71504
		short			SG71511	SG71512	SG71513	SG71514
	Black	long			SN71501	SN71502	SN71503	SN71504
		short			SN71511	SN71512	SN71513	SN71514
CG-H	Beige	long	with	dotted	SG71521	SG71522	SG71523	SG71524
		short			SG71531	SG71532	SG71533	SG71534
	Black	long		dotted	SN71521	SN71522	SN71523	SN71524
		short			SN71531	SN71532	SN71533	SN71534

Standard product range

VENOSAN® 7000

Compression class 2: Compression at wrist 23–32 mmHg

CCL
2

Item	Colour	Length	Hand attachment	Top band	REF/Size			
					S	M	L	XL
CG	Beige	long	without		SG 72501	SG 72502	SG 72503	SG 72504
		short	without		SG 72511	SG 72512	SG 72513	SG 72514
CG	Black	long	without		SN 72501	SN 72502	SN 72503	SN 72504
		short	without		SN 72511	SN 72512	SN 72513	SN 72514
CG-H	Beige	long	without	dotted	SG 72521	SG 72522	SG 72523	SG 72524
		short	without	dotted	SG 72531	SG 72532	SG 72533	SG 72534
CG-H	Black	long	without	dotted	SN 72521	SN 72522	SN 72523	SN 72524
		short	without	dotted	SN 72531	SN 72532	SN 72533	SN 72534
ZG	Beige	long	with		SG 72541	SG 72542	SG 72543	SG 72544
		short	with		SG 72551	SG 72552	SG 72553	SG 72554
ZG	Black	long	with		SN 72541	SN 72542	SN 72543	SN 72544
		short	with		SN 72551	SN 72552	SN 72553	SN 72554
ZG-H	Beige	long	with	dotted	SG 72561	SG 72562	SG 72563	SG 72564
		short	with	dotted	SG 72571	SG 72572	SG 72573	SG 72574
ZG-H	Black	long	with	dotted	SN 72561	SN 72562	SN 72563	SN 72564
		short	with	dotted	SN 72571	SN 72572	SN 72573	SN 72574
ZH	Beige	long	with		SG 72581	SG 72582	SG 72583	SG 72584
		short	with		SG 72591	SG 72592	SG 72593	SG 72594

Custom production

VENOSAN® 7000

Custom sizes

Item	Colour	Hand attachment	7002 CCL 2	7003 CCL 3
CG	Beige & Black	without	●	●
CG-H	Beige & Black	without	●	●
CH	Beige & Black	without	●	●
ZG	Beige & Black	with	●	●
ZG-H	Beige & Black	with	●	●
ZH	Beige & Black	with	●	●

● = available

CCL
1

Compression class 1:

Compression at wrist
18–21mmHg

CCL
2

Compression class 2:

Compression at wrist
23–32 mmHg

CCL
3

Compression class 3:

Compression at wrist
34–46 mmHg

Product overview

Medical support stockings

VENOSAN® LEGLINE 20

ELEGANT LUXURY

- For non-heavy legs in everyday situations
- Retains shape, ensures maximum freedom of movement
- Sheer look
- Optimal wearing comfort

Use

- To prevent varicose veins and heavy legs
- For people who sit or stand a lot at work
- Long journeys by plane, train and car



20 mmHg
140 den

Composition

78% polyamide
22% LYCRA® elastane


VENOSAN® LEGLINE 30

ELEGANT LUXURY

- For non-heavy legs in everyday situations
- Retains shape, ensures maximum freedom of movement
- Sheer look
- Optimal wearing comfort

Use

- To prevent varicose veins and heavy legs
- For people who sit or stand a lot at work
- Long journeys by plane, train and car



30 mmHg
280 den

Composition

AD / AGH:
70% polyamide
30% LYCRA® elastane
AT:
78% polyamide
22% LYCRA® elastane

VENOSAN® SILK SUPPORT

TRAVEL & LEISURE

- For non-heavy legs for everyday situations and travel
- Promotes blood circulation
- Efficient thrombosis prophylaxis
- Soft, smooth silk fibres on the skin
- Temperature-equalising

Use

- For people who sit or stand a lot at work
- Long journeys by plane, train and car
- For sensitive and dry skin



15–20
mmHg

Composition

65% polyamide
20% elastane
15% silk

VENOSAN® COTTON SUPPORT

TRAVEL & LEISURE

- For non-heavy legs for everyday situations and travel
- Promotes blood circulation
- Efficient thrombosis prophylaxis
- 100% cotton on the skin
- Can be worn with jeans and suits

Use

- For people who sit or stand a lot at work
- Long journeys by plane, train and car



15–20
mmHg

Composition
65% polyamide
20% cotton
15% elastane

Travel socks

JET LEGS® by VENOSAN®

TRAVEL WITH COMFORT SOCKS

- Gradually decreasing compression (15–20 mmHg)
- Promotes blood circulation
- Efficient thrombosis prophylaxis

Use

- Wear for plane, train and car travel
- Long journeys



Composition
86% nylon polyamide
14% LYCRA® elastane


Legline 20

Support Stockings 140 den

20 mmHg
140 den

HEEL AND TOE
REINFORCEMENT

Supporting socks



AD
below-knee
stockings



AGH
thigh-length
stockings



AT*
pantyhose
*Comfort top



AT Materna*
Maternity pantyhose
*Comfort top




Gradually decreasing compression


AD	AGH, AT, AT Materna
60% 12 mmHg	cF
80% 16 mmHg	cD
	cC
100% 20 mmHg	cB

Circumference (cm)

XS	S	M	L	XL
42–44	45–47	48–50	51–53	54–60
29–31	32–34	35–37	38–40	41–43
31–33	34–36	37–39	40–42	43–45
20–21	22–23	24–25	26–27	28–29



Effective compression
Prevents swollen legs



Everyday elegance
Sheer look



Optimal comfort
For optimal wearing comfort,
does not slip down

Lace top band



Outside



Inside

Size

	100	105	110	115	120	125	130	135	140	145	150	155	160	165	170	175	180	185	190	195	
4'11"																					150
5'0"																					152
5'1"																					155
5'2"																					157
5'3"																					160
5'4"																					163
5'5"																					165
5'6"																					168
5'7"																					170
5'8"																					173
5'9"																					175
5'10"																					178
5'11"																					180
6'0"																					182
	45	48	50	52	54	57	59	61	64	66	68	70	73	75	77	79	81	84	86	88	
	Weight (kg)																				

Height (cm)

Item

Type	Colours*	XS	S	M	L	XL
AD	Nude	VL 303 NU	VL 304 NU	VL 305 NU	VL 306 NU	–
	Sahara	VL 303 SA	VL 304 SA	VL 305 SA	VL 306 SA	–
	Black	VL 303 BL	VL 304 BL	VL 305 BL	VL 306 BL	–
AGH	Nude	VL 323 NU	VL 324 NU	VL 325 NU	VL 326 NU	–
	Sahara	VL 323 SA	VL 324 SA	VL 325 SA	VL 326 SA	–
	Black	VL 323 BL	VL 324 BL	VL 325 BL	VL 326 BL	–
AT	Nude	VL 333 NU	VL 334 NU	VL 335 NU	VL 336 NU	VL 337 NU
	Sahara	VL 333 SA	VL 334 SA	VL 335 SA	VL 336 SA	VL 337 SA
	Black	VL 333 BL	VL 334 BL	VL 335 BL	VL 336 BL	VL 337 BL
AT Materna	Sahara	VL 343 SA	VL 344 SA	VL 345 SA	VL 346 SA	VL347 SA
	Black	VL 343 BL	VL 344 BL	VL 345 BL	VL 346 BL	VL 347 BL

*Colours



Nude



Sahara




Black

Actual product colours may vary slightly from the colours in this catalogue due to the printing process.

Composition

78% polyamide | 22% LYCRA® elastane




Legline 30

Support stockings 280 den

30 mmHg
280 den

Support stockings




Gradually decreasing compression


AD	AGH, AT
60% 18 mmHg	cF
80% 24 mmHg	cD
	cC
100% 30 mmHg	100% 30 mmHg
	cB

Circumference (cm)

XS	S	M	L	XL
42–44	45–47	48–50	51–53	54–60
29–31	32–34	35–37	38–40	41–43
31–33	34–36	37–39	40–42	43–45
20–21	22–23	24–25	26–27	28–29



Effective compression
Prevents swollen legs



Everyday elegance
Sheer look



Optimal comfort
For optimal wearing comfort,
does not slip down

Lace top band



Outside



Inside

Size

	100	105	110	115	120	125	130	135	140	145	150	155	160	165	170	175	180	185	190	195	
4'11"																					150
5'0"						8½															152
5'1"							XS														155
5'2"								9													157
5'3"									S												160
5'4"										9½											163
5'5"											M										165
5'6"												10									168
5'7"													L			10½					170
5'8"														L							173
5'9"															L						175
5'10"																					178
5'11"																					180
6'0"																					182
	45	48	50	52	54	57	59	61	64	66	68	70	73	75	77	79	81	84	86	88	
	Weight (kg)																				

Height (cm)

Item

Type	Colours*	XS	S	M	L	XL
AD	Nude	VL 403 NU	VL 404 NU	VL 405 NU	VL 406 NU	–
	Sahara	VL 403 SA	VL 404 SA	VL 405 SA	VL 406 SA	–
	Black	VL 403 BL	VL 404 BL	VL 405 BL	VL 406 BL	–
AGH	Nude	VL 423 NU	VL 424 NU	VL 425 NU	VL 426 NU	–
	Sahara	VL 423 SA	VL 424 SA	VL 425 SA	VL 426 SA	–
	Black	VL 423 BL	VL 424 BL	VL 425 BL	VL 426 BL	–
AT	Nude	VL 433 NU	VL 434 NU	VL 435 NU	VL 436 NU	VL 437 NU
	Sahara	VL 433 SA	VL 434 SA	VL 435 SA	VL 436 SA	VL 437 SA
	Black	VL 433 BL	VL 434 BL	VL 435 BL	VL 436 BL	VL 437 BL

*Colours



Nude



Sahara



Black

Actual product colours may vary slightly from the colours in this catalogue due to the printing process.

Composition

AD|AGH: 70% polyamide | 30% LYCRA® elastane

AT: 78% polyamide | 22% LYCRA® elastane




Silk Support

Support stockings · Travel & leisure


15–20
mmHg

Support stockings



For non-heavy legs on long journeys

Gradually decreasing compression




Sizes


	S	M	L
EUR	36–39	39–42	42–45




**SWISS
COMFORT**



Effective compression
Prevents swollen legs



Soft band
Does not constrict




Satin touch
For optimal
wearing comfort,
does not slip down

Item

Colours*	S	M	L
White	SG 58800	SG 58801	SG 58802
Black	SG 58820	SG 58821	SG 58822
Jeans	SG 58870	SG 58871	SG 58872
Marine	SG 58810	SG 58811	SG 58812
Silver	SG 58880	SG 58881	SG 58882
Beige	SG 58840	SG 58841	SG 58842
Apricot	SG 58900	SG 58901	SG 58902
Sisko (Violet)	SG 58910	SG 58911	SG 58912

*Colours



DURABLE COLOURS

Actual product colours may vary slightly from the colours in this catalogue due to the printing process.

Composition

65% polyamide | 20% elastane | 15% silk




Cotton Support

Support stockings · Travel & leisure


15–20
mmHg

Support stockings




For non-heavy legs on long journeys

Gradually decreasing compression



Sizes

	S	M	L	XL
EUR	36–39	39–42	42–45	45–48



Effective compression
Prevents swollen legs

Soft band
Does not constrict

Anatomically knitted
For optimal wearing
comfort, does not slip down


Item

Colours*	S	M	L	XL
White	SG 59800	SG 59801	SG 59802	SG 59803
Navy	SG 59810	SG 59811	SG 59812	SG 59813
Black	SG 59820	SG 59821	SG 59822	SG 59823
Beige	SG 59840	SG 59841	SG 59842	SG 59843
Anthracite	SG 59850	SG 59851	SG 59852	SG 59853
Olive	SG 59860	SG 59861	SG 59862	SG 59863
Jeans	SG 59870	SG 59871	SG 59872	SG 59873
Silver	SG 59880	SG 59881	SG 59882	SG 59883
Wood	SG 59890	SG 59891	SG 59892	SG 59893

*Colours




White




Navy



Black




Beige



Anthracite



Olive



Jeans



Silver



Wood

Actual product colours may vary slightly from the colours in this catalogue due to the printing process.

DURABLE COLOURS

Composition

65% polyamide | 20% cotton | 15% elastane



S
 G
 E
 L
 E
 G
 S

JET LEGS®

by VENOSAN®

The ideal travel companions

JET LEGS by VENOSAN are the ideal travel companions. Their gradually decreasing compression provides efficient thrombosis prevention.

■ Compact product range

1 colour – 2 sizes

■ Above average

Benefit-to-cost ratio

Composition

86% nylon polyamide, 14% LYCRA® elastane



JET LEGS® by VENOSAN®

Ankle compression 15–20 mmHg

Item	Colour	Toe	REF/Size	
			Standard (36–40)	Large (41–45)
AD	Black	closed	SG 59610	SG 59611

[illegible]



People.Health.Care.



Lohmann & Rauscher GmbH & Co. KG · Westerwaldstraße 4 · 56579 Rengsdorf · Germany
Lohmann & Rauscher GmbH · Johann-Schorsch-Gasse 4 · 1141 Vienna · Austria

VENOSAN® is a registered trademark of
SWISSLASTIC AG ST. GALLEN, Switzerland – an  Company

2002683 0420 e