

# SCOPE BUDDY™ PLUS

Endoscope Flushing Aid



# A REVOLUTIONARY SOLUTION TO AN EVOLUTIONARY PROCESS

SCOPE BUDDY™ PLUS Endoscope Flushing Aid helps create a consistent, reliable, and safe means to flush endoscope channels. It provides disposable tubing, increased data monitoring, and automated flushing modes to properly flush endoscopes.

Manual cleaning\* is the most important step in reprocessing endoscopes, ultimately helping to reduce the risk of Healthcare-Associated Infections (HAIs).<sup>1</sup>

Over 1.7 million contracted HAIs/year  
Over 99,000 HAI deaths/year<sup>2</sup>

\* SCOPE BUDDY™ PLUS Endoscope Flushing Aid is only intended to be used as a Flushing Aid during the Manual Cleaning step of endoscope reprocessing. It is not intended to manually clean endoscopes.



## THE COMPLETE CIRCLE OF PROTECTION

As the global vanguard in infection prevention, **only Cantel Medical delivers the Complete Circle of Protection**, a full-value, proactive partnership dedicated to helping you remove risk, streamline operational efficiencies and optimize your success.

**MANUAL CLEAN** Manual cleaning is essential in removing bioburden and preparing the endoscope for high-level disinfection. For optimal results, Cantel's manual cleaning products and detergents are designed to be a critical first step to efficiently and effectively eliminating bacteria on endoscopes prior to reprocessing.



## BEST PRACTICE IS OUR PROMISE

*Industry Standards of Infection Prevention in Reprocessing Flexible GI Endoscopes confirm manual cleaning of endoscopes is necessary prior to high-level disinfection.\**

- Standardizing the process of cleaning endoscopes by utilizing both manual and automated methods will help ensure optimal cleaning while complying with the AAMI, Regulatory and Society Standards and Guidelines
- Recent investigations conclude that manual cleaning is clearly not being performed consistently or correctly resulting in serious negative patient outcomes<sup>3,4</sup>
- Performing meticulous manual cleaning is of paramount importance and implies that it should have the same attention to quality control as high-level disinfection and sterilization<sup>5</sup>

# Progressively Automated Thoroughly Validated

## INCREASING SAFETY

- **Reduce cross-contamination risks** by using 24-hour disposable tubing instead of reusable tubing
- **Continuously monitor** the endoscope manual flushing process with the Track and Trace features
- **Validated 1:35 default flushing time** to meet worst case flushing times for all major endoscope brands

## IMPROVING EFFICIENCY

- **Automatic and Manual Modes** provide easy to use steps for a dependable flushing process
- **Compatible** with all major endoscope brands
- **Improves staff efficiency** by simultaneously flushing multiple endoscope channels

## IMPACTING LIVES

- **Best practice flushing standards** help reduce the risk of hospital acquired infection breakouts for patients
- **Consistent, hands-free flushing** replaces repetitive and labor-intensive syringing of channels



**Continuing  
our 40 year  
legacy in  
infection  
prevention**





## THE MOST EXTENSIVE MANUAL CLEANING PORTFOLIO

- INTERCEPT™ Bedside Kit
- INTERCEPT™ Wipes
- VERISCAN™ LT Automated Endoscope Leak Detection System
- PULL THRU™ Channel Cleaning Device
- SCOPE BUDDY™ PLUS Endoscope Flushing Aid
- INTERCEPT™ Detergent



# SCOPE BUDDY™ PLUS Endoscope Flushing Aid Ordering Information

## SCOPE BUDDY™ PLUS ENDOSCOPE FLUSHING AID

### MODEL NUMBER DESCRIPTION

SBP-1000	USA, Canada
SBP-1001	INT

## DETERGENT

### PART NUMBER DESCRIPTION

ML02-0106	INTERCEPT™ Detergent
-----------	----------------------

## ACCESSORIES

SB-110-HU1015	24-Hour Multi-Use Pump Tube
SB-110-HU1007	24-Hour Multi-Use Hookup for OLYMPUS® Endoscopes
SB-110-HU1008	24-Hour Multi-Use Hookup for PENTAX® Endoscopes
SB-110-HU1009	24-Hour Multi-Use Hookup for FUJIFILM™ Endoscopes
SB-110-HU1010	24-Hour Multi-Use DSD Hookup Adapter
SB-110-HU1016	24-Hour Multi-Use Universal Hookup
SB-110-HU1012	Auxiliary Detergent Out Tube (24-Hour Multi-Use)
SB-110-HU1011	Detergent In Tube
SB-110-HU1017	OLYMPUS® Air/Water Blank Cap

## SCOPE BUDDY™ PLUS Endoscope Flushing Aid

### Unit Dimensions

10.3 H x 7.5 W x 11.5 D (inches)  
26.2 H x 19 W x 29.2 D (cm)

### Unit Weight

10 lbs (4.5 kg)

### Warranty

15 Months from the date of shipment, or one year from date of installation, whichever comes first.

### Electrical Requirements

100 – 240 V, 50/60 Hz

- 1 Rutala, W.A., Weber, D. J., and the Healthcare Infection Control Practices Advisory Committee (2008). *Guideline for Disinfection and Sterilization in Healthcare Facilities* (Last update: February 15, 2017). Retrieved from CDC: <https://www.cdc.gov/infectioncontrol/pdf/guidelines/disinfection-guidelines.pdf>
- 2 Centers for Disease Control, *Preventing Healthcare-Associated Infections*. Retrieved from CDC: <https://www.cdc.gov/washington/~cdcatWork/pdf/infections.pdf>
- 3 FDA Safety Communication Design of Endoscopic Retrograde Cholangiopancreatography (ERCP) Duodenoscopes May Impede Effective Cleaning. March 2015. <http://www.fda.gov/medicaldevices/safety/alertsandnotices/ucm434871.htm>.
- 4 Zachary L. Smith, Young S. Oh, Kia Saeian, Charles E. Edmiston Jr, Abdul H. Khan, Benson T. Massey, Kulwinder S. Dua., 2015. Transmission of carbapenem-resistant Enterobacteriaceae during ERCP: time to revisit the current reprocessing guidelines. *Gastrointestinal Endoscopy* 81(4) pp.1041-1045.
- 5 Infection Control during GI Endoscopy. ASGE Standard of Practice Committee. *Gastrointestinal Endoscopy*. 2008. Vol. 67. No. 6 pp.781-790.

\*SCOPE BUDDY™ PLUS Flushing Aid is only intended to be used as a Flushing Aid during the Manual Cleaning step of endoscope reprocessing. It is not intended to manually clean endoscopes.

SCOPE BUDDY™, INTERCEPT™, VERISCAN™, and MEDICLEAN EZ™ are trademarks of Medivators Inc. All other company and product names are trademarks of their respective owner.

© 2017-2020 Cantel. All Rights Reserved. The Complete Circle of Protection diagram is a trademark owned by Cantel Medical, and is subject to registered copyright protection in the U.S. and elsewhere under international treaties. U.S. Copyright Reg. Nos. VA0002197834 & VA0002201399.

## MEDIVATORS IS NOW CANTEL

[www.medivators.com](http://www.medivators.com)

### TO PLACE AN ORDER

p: 1.800.328.3340 (customer service) f: 1.800.686.8493 e: [custserv@medivators.com](mailto:custserv@medivators.com)

© 2021 Medivators Inc. LIT-00383/A

