



Laboratory identification and tracking solutions

Product catalog

Laboratory identification and tracking solutions



When it comes to reducing errors in your lab to improve productivity and efficiency, you need durable labels and an inventory management system you can trust. That's where we can help. With Brady, you get a reliable partner with full solutions in lab identification and specimen tracking, so you can maintain the integrity of your work while improving efficiencies.

Experience and support you can trust

Brady has been offering identification solutions to the global laboratory market for more than 15 years. Our experience is enhanced by the dedication of our highly qualified support team. This group is on hand to help you find a stock or customized labeling solution that aligns with your needs.

Proven lab products

Our identification solutions are backed by proven laboratory label materials, printers, software and scanners designed specifically for the needs and challenges of your laboratory environment, including harsh chemicals, autoclaves, liquid nitrogen and more.

▶ Let us show you what we have to offer



Table of contents

Laboratory overview.....	4
Printers, software and scanners.....	8
Printer comparison	8
Printers	10
Software	18
Barcode scanners	20
Label dispenser.....	21
Label material guide	22
Microbiology labels	26
Cell culture flasks.....	26
Petri dishes	28
Histology labels.....	30
Microscope slides	30
Tissue cassettes	31
Cryogenic labels	32
Cryo boxes and freezer racks	32
Cryo canes.....	33
Straws and rods.....	33
Cryo tubes.....	34
Sample collection and processing labels.....	40
Blood tubes.....	40
Specimen cups	40
Tamper-resistant seals	41
Sample bags	41
Wells and microplates.....	42
Centrifuge and conical tubes	43
Microcentrifuge and PCR tubes.....	46
General labels.....	52
Bottles, flasks and dishes	52
Water dissolvable	55
Autoclave indicating.....	55
Lean lab / visual workflow	56
GHS HazCom	58
Biohazard identification.....	60
Lab signs	61
Absorbents and spill kits.....	62
Index.....	64
Labels by size index.....	64
Part number index	70
Product index.....	74

Laboratory overview

Laboratory identification challenges

Loss means different things across various types of labs, from a minor inconvenience to major devastation:

- In clinical settings, sample loss can lead to misdiagnoses, delays in diagnoses or even fatalities and lawsuits
- In pharmaceutical environments, it can postpone or halt drug development and production
- In academic laboratories, it can impair or even void your findings
- In biotechnology settings, it can compromise your quality control, data integrity or drug safety. It also opens the door for the possibility of intellectual property loss or a delay in technology and product development

In any of these situations, a single sample loss can mean critical data or information will have to be reproduced — putting millions of dollars at risk.



Handwritten labels vs. printed labels



Handwritten labels

- Difficult to read
- Smears when handled
- Limited writing space
- Fades over time when stored
- Erased by chemicals
- Results in lost or misinterpreted information



Printed Brady labels

- Clear, crisp and readable
- Will not smear or fade
- Barcodes hold large amounts of information
- Durable and chemical resistant

Stunning stats

We surveyed 350 scientists about their sample management processes. Nearly **60% reported** having lost samples due to label failure and almost **50% reported** the loss impacted more than **2% of their samples**.

Sample management regulators

Numerous regulators, funders and agencies have expanded their policies on sample management. The new regulations are centered on one primary goal — to help your lab mitigate or eliminate sample loss.

Regulatory Agency	Policy
The Joint Commission	National Patient Safety Goals (NPSG 01.01.01)
College of American Pathologists	Quality Management Plan
Clinical and Laboratory Standards Institute (CLSI)	AUTO12-A
Clinical Laboratory Improvement Amendment (CLIA)	42 USC 263A
Centers for Disease Control (CDC)	Laboratory Best Practices ISBN 1-56238-748-0
National Cancer Institute (NCI)	Best Practices for Biospecimen (ID B.6.2)
Food and Drug Administration (FDA)	GLP 21. CFR Part 56

Brady solutions for today's leading labs

Brady is trusted throughout the industry for reliable, first-rate identification and chain-of-custody tracking systems. No matter what your laboratory studies, Brady labeling solutions can help you maintain effective laboratory practices.

- Label materials that endure freezer storage, liquid nitrogen, autoclave, hot water baths and solvents including xylene, DMSO and ethanol
- Continuous and die-cut labels that match existing sizes of vials, slides, well plates and more
- Printed labels that remain legible throughout sample processing and storage
- Label software and printers that allow barcoding for extensive information storage on even the smallest samples
- Easy-to-use printers and software that save valuable time during training and printing processes
- Ability to seamlessly integrate with your Laboratory Information Management System (LIMS) through Brady's Application Engineering Team



Label printers, software and scanners



Materials that withstand extreme lab conditions without smearing, fading or falling off

Improve productivity with color

Color helps to create a visual language that employees are able to process before they are even aware of it. When color is used as part of a standardized process, it improves efficiency and reduces errors. Brady provides a variety of color solutions that:

- Differentiate samples — find what you're looking for and move through your work faster
- Reduce errors — get context that can help reduce errors and help enhance memory performance
- Labels for every environment — protect your research with colors and labels that survive the harshest lab environments



Color makes it easier to visually process multiple objects at once

Try Brady labels with a free sample pack

Brady's labels are designed for harsh lab environments — test them for yourself.

Request your sample pack at BradyID.com/lab



Laboratory overview

Laboratory overview

The best lab is one that's custom built

You love that your lab is growing. But you also need to stay on top of your growing volume by designing a better workspace. With Brady, it's possible. From label printers and scanners to signs and floor marking tape, Brady offers many ways to design a space that is effective, safe and 100% your own.

General identification

Absorbent spill kits

Microbiology labels

Portable printers and software

Benchtop printers and scanners



Laboratory overview



Printers, software and scanners

Printer comparison

Your labeling solution starts with clear and readable print, and Brady has the easy-to-use printers that give you just that. Create custom labels on-site and on demand for lab applications such as cryo tubes, tissue cassettes, slides and more.

				
	BMP®21-LAB Label Printer Page 10	BMP®51 Label Printer Page 11	BradyPrinter M611 Mobile Label Printer Page 12	BBP®12 Label Printer Page 13
Description	Handheld printer for low volume lab labeling where die-cut labels are not needed.	Entry-level printer with the added versatility of die-cut labels.	Peripheral printer, powered by Brady Workstation mobile or desktop software, allows you to stay connected to your data.	A compact entry-level desktop label printer designed to combine efficiency, versatility and performance.
Applications	Tubes and vials, slides, bottles, conicals, lab storage and general purpose	Cryo samples, tubes and vials, bottles, conicals, slides, tissue cassettes and general purpose	Cryo samples, tubes and vials, bottles, conicals, petri dishes, slides and general purpose	Cryo samples, tubes and vials, bottles, conicals, petri dishes and slides
Max labels per day	100 labels per day	500 labels per day	750 labels per day	1,000 labels per day
Print resolution	203 dpi	300 dpi	300 dpi	300 dpi
Print technology	Thermal transfer	Thermal transfer	Thermal transfer	Thermal transfer and direct thermal
Memory	16 MB, 10 labels (2 MB available for storage)	48 MB / up to 25,000 labels	512 MB RAM / 256 MB flash	2 MB flash memory
Standalone printing	Yes	Yes	No	No
PC connectivity	No	Yes	Yes	PC-connect only
Connectivity options	Not applicable	USB, Wi-Fi, LAN	Bluetooth, Wi-Fi, USB 2.0	USB, ethernet, serial, parallel
Software compatibility	Not applicable	Brady Workstation	Brady Workstation Express Labels Mobile App (IOS / Android) or Brady Workstation Desktop	Brady Workstation
File transfer from PC to printer	Not applicable	Yes, with software	Yes	Yes
Wireless connectivity	No	Yes	Yes	No
Color capability	Single print color on a label	Single print color on a label	Single print color on a label	Single print color on a label
Keyboard	ABC or Cyrillic	QWERTY	Not applicable	Not applicable
Display	2 in. x 1 in., user adjustable backlit	2.7 in. x 1.7 in., backlit	3.5 in. (88.9 mm) color, user adjustable backlit	Not applicable
Text sizes	6-40 pt with autosizing	4 to 103 pt + autosizing	Yes, via label creation software	Yes, via label creation software
Multi-line print	Yes	Yes	Yes, via label creation software	Yes, via label creation software
Print rotation	No	Yes	Yes, via label creation software	Yes, via label creation software
Serialization	Yes	Yes	Yes, via label creation software	Yes, via label creation software
Barcode	39 and 128	39 and 128	Yes, via label creation software	Yes, via label creation software
Symbols	104	450	Yes, via label creation software	Yes, via label creation software
Portability	Yes	Yes	Yes	No
Power supply	6 AA batteries, AC adapters, or Li-Ion battery	8 AA Alkaline batteries, AC power, or Li-Ion battery	Rechargeable Li-Ion, AC power / battery charger / car charger	AC power
Weight	1.66 lbs. (0.75 kg)	2.6 lbs. (1.2 kg)	3.52 lbs. (1.6 kg)	5.51 lbs. (2.5 kg)
Warranty	2 years	2 years	2 years	5 years
Material and ribbons				
Continuous or die-cut labels	Continuous labels only	Continuous and die cut	Continuous and die cut	Continuous and die cut
Tape widths	0.250 to 0.750 in. (6.4 to 19.1 mm)	0.235 to 1.500 in. (6.4 to 38.1 mm)	0.812 to 2.000 in. (20.6 to 50.8 mm)	1.000 to 4.250 in. (25.4 to 108 mm)
Max label width	0.750 in. (19.1 mm)	1.500 in. (38.1 mm)	2.000 in. (50.8 mm)	4.250 in. (108 mm)
Materials supported	M21-series cartridges (Label and ribbon in one cartridge)	M-, MC-, MC1-series cartridges (Label and ribbon in one cartridge)	BPTL, PTL, PSHT, PTLEP, PTS, PTLFP, PSPT, PTPSL, PSFR, PTPS, PTLSL, BM71, M61, M61C, BM61C labels and ribbon cartridges	1-in. small core label rolls with 0.5-in. core ribbon 3-in. core label rolls with use of an external media holder

Printers, software and scanners



**BradyPrinter i3300
Industrial Label Printer**
Page 14

**BradyJet J2000
Color Label Printer**
Page 15

**BradyPrinter i5100
Industrial Label Printer**
Page 16

**BradyPrinter i7100
Industrial Label Printer**
Page 17

Description	Small footprint desktop printer featuring hassle-free printing where smart-chip supplies automatically set up the label, eliminating calibration and label waste.	Full-color, high-resolution inkjet printer for creating photo-quality color labels while reducing inventory of labeling supplies.	High-performance printing with smart-cell supplies allowing for printer, materials and ribbon to communicate. Makes set up and changeovers significantly faster, easier and less wasteful.	Reliable and heavy-duty label printer designed to print thousands of labels per day, every day. High-accuracy printing on the widest range of materials, including labels as small as 0.125-in. wide.
Applications	Cryo samples, tubes and vials, bottles, conicals, petri dishes, slides, GHS and general purpose	Cryo samples, tubes and vials, bottles, conicals, petri dishes, slides, GHS and general purpose	Cryo samples, tubes and vials, bottles, conicals, slides, tissue cassettes and general purpose	Cryo samples, tubes and vials, bottles, conicals, slides, tissue cassettes and general purpose
Max labels per day	3,000 labels per day	3,000 labels per day	5,000 labels per day	7,000 labels per day
Print resolution	300 dpi	Up to 4,800 dpi	300 dpi or 600 dpi	300 dpi or 600 dpi
Print technology	Thermal transfer	Thermal inkjet CMY	Thermal transfer and direct thermal	Thermal transfer and direct thermal
Memory	512 MB	Not applicable	256 MB	256 MB
Standalone printing	No	No	No	No
PC connectivity	PC connect only	PC connect only	PC connect only	PC connect only
Connectivity options	USB, ethernet, Wi-Fi	USB, ethernet	USB, ethernet, serial	USB, ethernet, serial
Software compatibility	Brady Workstation	Brady Workstation	Brady Workstation	Brady Workstation
File transfer from PC to printer	Yes	Yes	Yes	Yes
Wireless connectivity	Yes	No	No	No
Color capability	Single print color on a label	Full color CMY — 16 million color combinations	Single print color on a label	Single print color on a label
Keyboard	Not applicable	Not applicable	Not applicable	Not applicable
Display	Touchscreen color LCD	Not applicable	Touchscreen color LCD	Touchscreen color LCD
Text sizes	Yes, via label creation software	Yes, via label creation software	Yes, via label creation software	Yes, via label creation software
Multi-line print	Yes, via label creation software	Yes, via label creation software	Yes, via label creation software	Yes, via label creation software
Print rotation	Yes, via label creation software	Yes, via label creation software	Yes, via label creation software	Yes, via label creation software
Serialization	Yes, via label creation software	Yes, via label creation software	Yes, via label creation software	Yes, via label creation software
Barcode	Yes, via label creation software	Yes, via label creation software	Yes, via label creation software	Yes, via label creation software
Symbols	Yes, via label creation software	Yes, via label creation software	Yes, via label creation software	Yes, via label creation software
Portability	No	No	No	No
Power supply	AC power	AC power	AC power	AC power
Weight	11.5 lbs. (5.2 kg)	7 lbs. (3.2 kg)	14.6 lbs. (6.6 kg)	22.05 lbs. (10.0 kg)
Warranty	5 years	1 year	5 years	5 years
Material and ribbons				
Continuous or die-cut labels	Continuous and die cut	Continuous and die cut	Continuous and die cut	Continuous and die cut
Tape widths	0.710 to 4.250 in. (18.03 to 108.0 mm)	1.000 to 4.250 in. (25.4 to 107.95 mm)	0.360 to 4.400 in. (9.14 to 111.8 mm)	0.360 to 4.490 in. (9.14 to 114.0 mm)
Max label width	4.250 in. (108.0 mm)	4.250 in. (108 mm)	4.330 in. (109.9 mm)	4.160 in. (105.7 mm)
Materials supported	B30-, B33-series labels and ribbons	J20-series labels and ink cartridges	3-in. core label rolls with IP-enabled ribbon	3-in. core label rolls with ribbon

Printers, software and scanners

Printers, software and scanners



BMP[®]21-LAB LABEL PRINTER

Where pens fail, printers succeed

Once you experience the power of the BMP21-LAB printer, you'll never use a pen again.

- Vial label shortcut key automatically fits labels to vial and tube sizes
- Greek symbols and time and date stamps for enhanced sample tracking
- Built-in serialization and barcode symbologies provide more information in a smaller footprint
- Drop-lock-and-print cartridges for quick printing with no set up required
- Easy to use and train others to use
- Variety of power options (rechargeable battery, AA batteries or AC power supply — all sold separately)

Learn more at BradyID.com/bmp21lab



Printer and accessories

Catalog #	Description
BMP21-LAB	BMP21-LAB Label Printer*
BMP21-LAB-KIT1	BMP21-LAB Label Printer Kit**
M-AC-110937	AC Power Adapter / Charger
BMP21-PLUS-BATT	Li-Ion Battery Pack
BMP21-LAB-BATTCOV	Spare Battery Cover
BMP21-MAGNET	Magnet Accessory
BMP21-PLUS-HC	Hard Carrying Case

*Printer includes M21-750-7425 label cartridge and Quick Start Guide (batteries not included).

**Printer kit includes printer, M21-750-7425 label cartridge, Li-Ion battery, AC power / charger and Quick Start Guide.

An affordable laboratory printer designed for all your lab stations




Vial labels



Cryo tube labels



Slide labels



Petri dishes



General lab ID



Lean lab / visual workflow

Printers, software and scanners



BRADY 360
2 YEAR
WARRANTY

BMP[®]51 LABEL PRINTER

Made to get the job done

A powerful tool at your fingertips! Produce legible, professional-looking labels for your lab samples while drastically reducing the time it takes to label.

- Continuous and die-cut label styles, including round tops, self-laminating labels and two-label combos
- Ribbon and labels in one easy-to-load cartridge
- Fast print speed — prints 1 in. per second
- 1.5-in. wide continuous output
- Material recognition and automatic label formatting
- Large, backlit graphical display
- Durable built-in cutter with label retention
- Advanced label creation capabilities with Brady Workstation software
- Tough, durable label materials that won't fade, smear or fall off

Learn more at
BradyID.com/bmp51



Printer and accessories

Catalog #	Description
BMP51	BMP51 Printer with AC Adapter / Battery Charger and Li-Ion Battery*
BMP51-KIT-AC	BMP51 Printer with AC Adapter / Battery Charger*
BMP51-KIT-LS	BMP51 Laboratory Starter Kit**
BMP-UBP-AC	Li-Ion Battery Pack with AC Adapter / Battery Charger
BMP-UBP	Li-Ion Battery Pack
M-AC-143110	AC Adapter / Battery Charger
M50-BATT-TRAY	Spare AA Battery Tray
142267	Brady Network Card — Ethernet only
142268	Brady Network Card — Wi-Fi + Ethernet
M50-CABLE	USB Cable
M50-HC	Hard Carrying Case
BMP-SC-1	Soft Carrying Case

*Printer ships with hard case, Quick Start Guide, USB cable, Brady Workstation Basic Design Suite and one MC-1500-595-WT-BK cartridge. (Network card sold separately.)

**Printer kit ships with Li-Ion battery, AC adapter / battery charger, hard case, Quick Start Guide, USB cable, Brady Workstation Basic Design Suite, and one each of M-103-498, M-120-492, M-156-492, M-81-488 and MC-1500-595-WT-BK cartridges. (Network card sold separately.)

U.S. Patent Nos. 6,732,619; 7,070,347; 8,542,155

Get the job done faster, easier and with confidence



Tube and vial labels



Cryo can labels



Slide labels



Well and microplate labels



Lean lab / visual workflow



And more

Printers, software and scanners

Printers, software and scanners



Learn more at BradyID.com/m611



Your durable print-from-everywhere solution



Tube and vial labels



Cryo tube labels



Slide labels



Petri dishes



General lab ID



Lean lab / visual workflow

BradyPrinter M611

MOBILE LABEL PRINTER

Keep your data traceable

Your sample is only as valuable as the data behind it. Don't lose data integrity with poor identification — stay connected to your systems with the BradyPrinter M611 Mobile Label Printer and Brady Workstation Mobile or Desktop Software.

- Seamlessly get your data from your LIMS using data automation and the Brady Workstation Express Labels Mobile App
- High-performance printing in a small footprint when bench space is a premium
- Long-lasting identification with materials designed to survive your harsh lab conditions
- Print from your mobile devices via Bluetooth or Wi-Fi connectivity
- Edit your label data on the go with the Express Labels Mobile App
- Enter data through a form entry window or by scanning barcodes using the device's camera
- Use with Brady Workstation Lab Suite for streamlined template creation

Printer and accessories

Catalog #	Description
M611-AM-BT-W	BradyPrinter M611 Bluetooth / Wi-Fi AC*
M611-BC-KIT1	BradyPrinter M611 with (1) M611-R6010 ribbon and (2) rolls of B-432 1-in. x 2-in. labels*
M-AC-143110	AC adapter and battery charger
BMP-MAGNET-2	Magnet for BradyPrinter M611
BMP-HC-2	Hard case for BradyPrinter M611
BMP-SC-2	Soft case for BradyPrinter M611
BMP-STRAP-1	Carrying strap†
PCK-5	Cleaning kit 5-pack

*Printer ships with AC adapter / battery charger, 2.0 cable, utility strap, PTL-97-488 label material, M611-R4310 ribbon and USB flash drive with product information.

Print ribbons

For best print performance, use the ribbon that corresponds to your material. All ribbons are 2.0 in. x 75 ft.

Catalog #	Color	For use with label materials
M611-R4310	Black	B-403, B-461, B-488, B-490, B-7425
M611-R6010	Black	B-439
M611-R6210	Black	B-499
M611-R6410	Black	B-481, B-492, B-494,



Learn more at
BradyID.com/bbp12



BBP[®]12 LABEL PRINTER

Efficiency, versatility and performance packed into one multi-use printer

With high-speed printing and a simple, user-friendly design, the BBP[®]12 Label Printer packs a lot of punch into one compact tool.

From slide labels to cryo tube labels, the BBP[®]12 Label Printer gives you the versatility to print many of your labels from one machine — saving costs and reducing downtime.

- Equipped with an Ethernet connection and three communication ports (USB, LPT and serial)
- Label materials up to 4-in. wide
- Prints up to 1,000 labels per day

Printer and accessories

Catalog #	Description
BBP12-NA-PWID	BBP12 with Brady Workstation Product and Wire ID Software Suite*
BBP12-US	BBP12 Printer*
BBP12-U	BBP12 External Media Holder for 3-in. core label rolls**

*Printer ships with drivers CD, power cord, printer, Quick Start Guide, ribbon paper core, USB cable.

**While the BBP12 Printer primarily uses 1-in. core materials, the external media holder (#BBP12-U) is available for printing 3-in. core rolls in limited quantities.

Print ribbons


For best print performance, use the ribbon that corresponds to your material. Ribbons below are also compatible with the BBP[®]11 Label Printer. All ribbons are 4.0 in. x 242 ft.

Catalog #	Color	For use with materials
R4313	Black	B-424, B-449, B-461, B-488, B-490
R4913	Black	B-499
R6013	Black	B-533
R6213	Black	B-498
R6413	Black	B-481, B-482, B-492, B-7425

Versatility to print all your labels



Tube and vial labels



Slide labels



Cryogenic labels



Sample bag labels



Freezer rack labels



And much more

Printers, software and scanners

Printers, software and scanners



BRADY 360
5 YEAR
WARRANTY



Learn more at
BradyID.com/i3300



BradyPrinter i3300

INDUSTRIAL LABEL PRINTER

No waste. More labels. Less effort.

Get your ideal label printing scenario with the BradyPrinter i3300 Industrial Label Printer — labels and ribbons that snap into place with no calibration and no wasted blank labels.

- Print smart technology allows the printer to print on the first label — even small labels
- Changeover in a flash — foolproof drop-in label rolls and ribbons snap into place in 20 seconds
- Talk about options, this printer prints on a variety of lab-specific materials
- Expand the possibilities — add Brady Workstation Software to your solution and make label creation even easier
- Get auto calibration in the printer and auto setup in the software with print smart technology that allows the printer, label materials, ribbons and software to work together

Printer and accessories

Catalog #	Description
150642	BradyPrinter i3300 with Brady Workstation Laboratory ID Suite*
149552	BradyPrinter i3300 with Brady Workstation Basic Suite*
150639	BradyPrinter i3300 with Brady Workstation Scan and Print and Scanner*
150645	BradyPrinter i3300 with Brady Workstation GHS App*
120994	Padded rolling case with handle and supply pockets
B31-WIPER	Extra media wipers (5 per pack)
B31-CCT	Cutting tool
PCK-6	IPA printer cleaning swabs 50/pk

*Printer ships with cleaning kit, cutter cleaning tool, power cord, product CD, Quick Start Guide and USB cable.

More than 300 label parts designed specifically for the laboratory



Tube and vial labels



Cryogenic labels



Well and microplate labels



Lean lab / visual workflow



GHS chemical labels



And much more

Print ribbons

For best print performance, use the ribbon that corresponds to your material. Ribbons below are also compatible with the BBP®33 Printer. All ribbons are 4.330 in. x 200 ft.

Catalog #	Color	For use with materials
B30-R4300	Black	B-403, B-424, B-427, B-461, B-490, B-499, B-7425
B30-R6000	Black	B-430, B-481, B-499
B30-R6200	Black	B-498
B30-R6400	Black	B-481, B-482, B-492, B-494
B30-R10000	Black	B-595, B-7569



Learn more at
BradyID.com/j2000



BradyJet J2000

COLOR LABEL PRINTER

Improve productivity with lab-durable color

Harsh environments and limited space can significantly narrow your identification options. With the BradyJet J2000, you get colorful identification that improves productivity and reduces errors while surviving even the harshest lab environments.

- Quickly and easily differentiate your samples with multiple color options
- Protect your research with colors and labels that survive the harshest lab environments
- Get clear identification without taking up too much bench space
- Take the guesswork out of label setup with Intelligent Printing (IP) enabled materials that communicate with the printer and software so you can get back to work faster

Printer and accessories

Catalog #	Description
J2000-BWLAB	BradyJet J2000 with Brady Workstation Laboratory ID Suite ^{††}
J2000	BradyJet J2000 Inkjet Label Printer ^{††}
J20-CUTTER	J20 Replacement Cutter Blade
J20-PAD	J20 Ink Pads

**Printer ships with power cable, USB cable, inkjet cartridge and roll of print head calibration material.*

Ink cartridge

Catalog #	Description
J20-CMY	J20 Ink Cartridge

Colorful identification in harsh environments



Tube and vial labels



Cryogenic labels



Slide labels



Boxes, plates and glassware



GHS chemical labels



Lab equipment labels and more

[†] **WARNING:** Cancer and Reproductive Harm www.P65Warnings.ca.gov



Learn more at BradyID.com/i5100



BradyPrinter i5100

INDUSTRIAL LABEL PRINTER

The perfect mix of performance and smarts

When you're looking for a printer with the fundamental printing capabilities for the work you need to get done every day, look to the BradyPrinter i5100 Industrial Label Printer. You'll get the perfect combination of high-volume performance and smart printing technology to make high-mix printing easier throughout your lab. The BradyPrinter i5100 makes printing easier with:

- An intuitive touchscreen that allows you to quickly and easily change settings on your printer
- Intelligent Printing (IP) takes the guesswork out of changeovers since materials and software communicate with the printer
- Changing out ribbons and label materials is easier than ever with the specially designed ribbon spindles, guides and roll holders that lock into place without any clumsy extra steps

Printer and accessories

Catalog #	Description
149452	BradyPrinter i5100, 300 dpi*
149455	BradyPrinter i5100, 600 dpi*
149453	BradyPrinter i5100, 300 dpi with auto cutter*
149456	BradyPrinter i5100, 600 dpi with auto cutter*
PCK-6	IPA printer cleaning swabs 50/pkg

*Printer ships with Brady Workstation Product and Wire ID Software Suite, one RFID-enabled media roll holder, one ribbon take-up core, power cord, USB cable and CD that contains drivers, user and configuration manuals.

Print ribbons

For best print performance, use the ribbon that corresponds to your material and roll width below. Ribbons below are also compatible with the Brady IP Printer. All ribbons are 984 ft. long.

Catalog #	Color	For use with materials	For use with label width
IP-R4306	Black	B-403, B-424, B-427, B-449, B-461, B-488,	Up to 1.57-in. width label
IP-R4302	Black	B-490, B-494, B-7425, B-7425-AC, B-813	Up to 2.36-in. width label
IP-R4300	Black		Up to 3.27-in. width label
IP-R4307	Black		Up to 4.33-in. width label
IP-R6002	Black	B-430, B-499, B-533	Up to 2.36-in. width label
IP-R6000	Black		Up to 3.27-in. width label
IP-R6200	Black	B-351, B-498	Up to 3.27-in. width label
IP-R6207	Black		Up to 4.33-in. width label
IP-R6402	Black	B-481, B-482, B-492	Up to 2.36-in. width label
IP-R6400	Black		Up to 3.27-in. width label

High-performance labeling solutions



Tube and vial labels



Tamper-resistant seals



Tissue cassette and slide labels



Specimen cup labels



Cell culture flasks and dishes



And much more

Printers, software and scanners



Learn more at
BradyID.com/i7100

BradyPrinter i7100

INDUSTRIAL LABEL PRINTER

Heavy duty meets high tech

No matter the demands of your labeling jobs — from higher volumes to precise printing — this reliable, heavy-duty and high-accuracy label printer gets it done.

- Print on a wide range of materials with confidence
- Ultra-sharp, high-accuracy printing and center alignment allow you to print reliably on labels as small as 0.125 in.
- Adjust printer settings quickly with an intuitive color touchscreen
- With print speeds up to 11.8 inches per second and the “peel” configured printer, you save time by automatically peeling labels from the liner



Printer and accessories

Catalog #	Description
149050	BradyPrinter i7100, 300 dpi*
149056	BradyPrinter i7100, 600 dpi*
149045	BradyPrinter i7100, peel version 300 dpi*
149053	BradyPrinter i7100, peel version 600 dpi*
149048	Rubber-coated print roller for media widths up to 25 mm
149055	Rubber-coated print roller for media widths up to 50 mm
149054	Rubber-coated print roller for media widths up to 80 mm
149043	Rubber-coated print roller for media widths up to 120 mm
149044	Silicone-coated print roller for media widths up to 114 mm
149051	Rewind assist roller
PCK-6	IPA printer cleaning swabs 50/pk
149208	Spare hex key tool to change rollers and install accessories
149069	Auto cutter with tray

Additional cutter, label presentation and I/O interface accessories available at BradyID.com.

**Printer ships with Brady Workstation Product and Wire ID Software Suite, two ribbon take-up cores, 120 mm rubber-coated print roller, hex key tool, USB cable, power cord, owners manual (English) and product CD with driver, configuration / settings manual and owners manual in multiple languages.*

Identification even on difficult-to-print materials



Tube and vial labels



Cryo labels



Slide labels



Dish, tray and well plate labels



Sample bag labels



General lab equipment labels

Print ribbons

For best print performance, use the ribbon that corresponds to the material and label width below. All ribbons are 984 ft. long.

Catalog #	Color	For use with materials	For use with label width
R4306	Black	B-403, B-424, B-427,	Up to 1.57-in. width label
R4302	Black	B-449, B-461, B-488,	Up to 2.36-in. width label
R4300	Black	B-490, B-494, B-7425,	Up to 3.27-in. width label
R4307	Black	B-7425-AC, B-813	Up to 4.33-in. width label
R6002	Black	B-430, B-499, B-533	Up to 2.36-in. width label
R6000	Black		Up to 3.27-in. width label
R6200	Black		Up to 3.27-in. width label
R6207	Black	B-351, B-498	Up to 4.33-in. width label
R6402	Black		Up to 2.36-in. width label
R6400	Black	B-481, B-482, B-492	Up to 3.27-in. width label

Printers, software and scanners

Printers, software and scanners



Brady Workstation Laboratory ID Suite

This comprehensive software suite is specifically designed to help your lab run more efficiently. Advanced design features allow you to create custom labels that fit the needs of your laboratory. With this collection of apps, you can implement a new identification process that will help eliminate costs and errors and keep things running smoothly.

Brady Workstation Laboratory ID Suite includes the following apps:



Custom Designer
Complete, customized label creation and editing.



Excel Import
Easily import Excel data right onto your labels.



Advanced Import
Import data from sources beyond Excel onto your labels.



Sequence
Labels based on numeric and alphanumeric sequences.



Advanced Sequence
Create dependent and complex sequences for label creation.



Template Editor
Set elements on labels to be editable or locked.



Print Partner
Quickly create labels by filling in a template and printing.



Batch Print
Group files to print together in a batch



Label Settings
Set default horizontal and vertical printing adjustments for specific supplies and printers.



Manage Files
Manage saved files from one screen, with options to preview, categorize and delete.



Network Printers
Connect to your printers via IP address or DNS address.

▶ **Free 30-day trial at**
Workstation.BradyID.com/lab

Catalog #	Description
BWS-LABS-EM	Email activation code

Label design software for your entire lab



Tube and vial labels



Slide labels



Vial top labels



Tissue cassette labels



Equipment labels



And more

Additional software for your identification needs

The following suites are designed to improve and automate your printing and labeling process:



Scan and Print Suite

This suite can solve your workflow needs by consuming data from a barcode scanner, keyboard or custom script to populate your label templates as part of a labeling workflow.

- Custom Designer
- Text Labels
- Template Editor
- Scan and Print
- Label Settings
- Manage Files
- Network Printers

Catalog #	Description
BWS-SPS-EM	Email activation
BWS-SPS-CD	CD*



Print Partner Suite

This suite includes everything you need to fill in your templates and print, including importing data from a variety of supported data sources or automatically generating sequential data for your labels.

- Excel Import
- Advanced Import
- Sequence
- Advanced Sequence
- Print Partner
- Network Printers

Catalog #	Description
BWS-PPS-EM	Email activation
BWS-PPS-EM-VOL	Email with volume license activation
BWS-PPS-CD	CD*
BWS-PPS-CD-VOL	Volume license CD*



Basic Label Design Suite

Use this free software option to get started designing your labels. You get basic formatting capabilities, graphics and simple barcode libraries.

- Custom Designer LITE
- Excel Import
- Sequence
- Batch Print
- Label Settings
- Manage Files
- Network Printers

Free download at Workstation.BradyID.com/FREE



GHS Chemical Identification App

We've streamlined the label creation process so you can easily enter the data once from your SDS sheets, then create and print GHS labels. You can also auto-populate your GHS labels using the built-in chemical database.

Catalog #	Description
BWRK-GHS-DWN	Email activation



Printers, software and scanners

Barcode scanners

Brady has partnered with Code® Corp., a technology leader in image-based barcode systems, in order to offer several different barcode scanner solutions to meet your data collection and tracking needs. All of the scanners are capable of decoding both 1D (linear) and 2D (data matrix) barcodes, allowing flexibility in label design.



Code Reader™ 950

- High-speed, omnidirectional imaging of 1D and 2D barcodes
- Manual or automatic triggering
- Wide area image sensor
- 2 programmable indicators: LED and tone
- IP54 housing


Catalog #	Description
149561	CR950 with 6-ft. USB cable and stand



Code Reader™ 1500

- High-speed, omnidirectional reading of 1D and 2D barcodes
- Dual-field optics, high density and wide field
- Manual or automatic triggering
- Glare reduction technology for reading codes on shiny surfaces
- Data editing and parsing with JavaScript
- Disinfectant-ready housing

Catalog #	Description
150508	Light gray CR1500 with 6-ft. USB cable and stand
150507	Dark gray CR1500 with 6-ft. USB cable



Code Reader™ 2600

- High-speed, omnidirectional reading of 1D and 2D barcodes
- Patented dual-field optics, high density and wide field
- Unequaled performance on rounded, curved and shiny surfaces
- Multiple programmable buttons for customized workflow processes
- Data management capability with JavaScript
- Disinfectant-ready housing

Catalog #	Description
Scanner with handle	
133080	CR2600 with handle, Bluetooth, battery, charging station with US power supply and M3 modem with USB cable
133081	CR2600 with handle, Bluetooth and battery
Scanner in palm configuration	
133082	CR2600 palm with Bluetooth, battery, charging station with US power supply and M3 modem with USB cable
133083	CR2600 palm with Bluetooth and battery
Accessories for CR2600	
133084	Charging station with US power supply
133085	Quad bay B5 battery charging station with US power supply
133086	B5 battery
133087	Weighted base / adapter plate
133088	6-ft. USB cable (USB to 10PIN RJ45)

Simple tools for simple tasks

Brady Material Dispenser with Liner Remover

Your one-off identification tasks are quick and uncomplicated — keep it that way with the Brady Material Dispenser.

Perfectly designed to fit your 3-in. core Brady materials, this material dispenser is easy to use with one hand. The narrow design means space is saved for your more important work projects. Get bonus green points by reusing your old material rolls as a liner take up roll. Other helpful features include:

- Sturdy weighted base and rubber feet keep the dispenser in place during your on-the-go labeling
- Minimal adhesive exposure to the elements so your labels last even in the harshest environments
- Brady's complete 3-in. core continuous material offering designed to stand up to your toughest conditions
- Options for custom-printed materials to fit your unique needs



▶ **Learn more at
BradyID.com**

Catalog #	Description
DISPENSER	Material Dispenser
DSP-403KIT	Material Dispenser with five rolls of B-403 water-dissolvable paper (1 in. x 100 ft.)

Use with pre-printed labels or load with blank material for one-off identification



Brady identification materials

Products that perform. Quality that lasts.

At Brady, we have a steadfast commitment to providing our customers with unrivaled products and an unrivaled buying experience. We use customer feedback to continually design, improve and innovate our products to serve all of your needs, from general identification solutions to custom engineered labels.

With more than 20 label material patents, Brady works to provide customers with innovative products of the highest technological quality. Our Tobey Research Center material science team develops and rigorously tests all of our products to ensure high performance, industry compliance, durability and reliability.

Our world-class label series are engineered to be long-lasting and to withstand the harshest environments and work for a variety of applications. When it comes to identifying and tracking specimen samples, Brady is the company you can trust when performance matters most.



Label materials for laboratories

Brady offers a variety of thermal transfer label materials engineered to survive and perform in the harsh chemicals and range of temperatures of the lab environment.

These durable labels remain legible and in place when put to the test making them ideal for the biotech, pharmaceutical, academic, industrial and health care industries.

Custom lab labels

We understand that each lab is unique and you may need custom label solutions that are specifically designed and compatible with your existing equipment or print systems.

Custom label solutions include:

- Blank labels in a variety of sizes and colors to easily organize and visualize batches of samples
- Custom text, graphics, colors and serialized barcodes
- Custom print and apply labels for automation equipment
- Pre-printed labels in the sizes and materials needed
- Fully printed custom labels or partially pre-printed labels so you can add additional text as needed



Test Brady materials yourself with a free sample pack

Request a free sample pack of materials and put them to the test. Get your samples at BradyID.com/lab

View Technical Data Sheets for all materials at BradyID.com/techdata



Label material guide



Microbiology labels

Page 26

In microbiology, cells are grown under controlled conditions to support the formulation of culture media or sera. Depending on the cell type, culture conditions and procedures may vary widely, from incubators and autoclaves to freezers, liquid nitrogen and more.

Labels used on flasks, petri dishes, tubes, centrifuge tubes and slides in microbiology should be designed with the durability to withstand these conditions in order to reduce lab errors and improve efficiency.



Histology labels

Page 30

Histology labs require clear, accurate and reliable specimen identification and tracking to help reduce errors and improve productivity. Brady slide labels can be applied before staining, and they resist chemicals (including xylene and H&E), withstand most soaks and stains, and remain durable when faced with paraffin, processing, long-term storage and more.



Cryogenic labels

Page 32

Brady offers cryo-specific label materials that survive in liquid nitrogen or ultra-low temperature freezers. Whether you're applying a label to a tube at room temperature or already frozen, your label material needs to withstand exposure to LN2, liquid or vapor phases, nitrogen freezers, thaw cycles, autoclaves, chemicals and solvents, such as DMSO, methanol or isopropyl alcohol.



Sample collection and processing labels

Page 40

Proper sample management is critical to the reliability of testing, accuracy of data and patient safety. When it comes to effective sample handling, Brady labels remain intact and legible through sample collection, processing and storage.



General labels

Page 52

Label options range from cost effective paper to more durable choices, such as nylon cloth that will last in harsh environments and cleanup. Removable and repositionable adhesives are available so you can easily reuse your glassware or container. Permanent adhesives are also available for rough surface applications.



Biohazard identification

Page 60

The most essential step toward safety is the clear and instantly readable identification of all biohazard areas and waste products. Brady has the identification products you need, including biohazard tape, disposal bags and signs, in order to clearly communicate the presence of infectious substances.

Label material guide

Label materials	Descriptions	Common applications					Common environments					
		Microbiology	Histology	Cryogenic	Sample collection	General lab ID	Liquid nitrogen (196° C)	Autoclave (121° C)	Hot water bath (100° C)	Freezer (-80° C)	Ambient temp only	Dissolves in water
B-351 vinyl	Tamper-resistant vinyl designed to fracture easily to prevent one-piece removal. Use caution when removing from the liner.				●						✓	
B-403 paper	Paper material completely dissolves in warm water within 30 seconds, leaving no adhesive residue behind.					●	✓					✓
B-424 paper	Economical paper material with good contrast and smear resistance.					●				✓		
B-427 vinyl	Durable vinyl with clear overlaminate is good for wrapping and conformability. Must be wrapped upon itself when used in liquid nitrogen.			●	●		✓	✓	✓			
B-430 polyester	Clear high-performance polyester material.	●		●	●	●	✓	✓	✓	✓		
B-432 polyester	Clear polyester withstands exposure to numerous solvents and exhibits good adhesion to metals and a variety of other surfaces.	●			●	●						
B-449 polypropylene	Removeable polypropylene for temporary labeling applications requiring clean removability. Also has write-on capabilities.			●		●		✓		✓		
B-461 polyester	Ultra-thin polyester with a clear overlaminate tail that protects the print and allows liquid levels to be visible. Also available in flat-white format without the clear overlaminate.			●	●		✓	✓	✓	✓		
B-481 polyester	StainerBondz™ polyester with extreme chemical resistance that withstands the rigors of typical deparaffinization, including ethanol, xylene and H&E staining process.	●	●					✓		✓		
B-488 polyester	Polyester with moderate chemical resistance and short-term exposure to xylene, DMSO and ethanol.	●	●	●	●			✓		✓		
B-490 polyester	FreezerBondz™ polyester ideal for flat, frozen surfaces. Adheres to frozen tubes and vials when wrapped around itself by at least 1/8 inch.			●	●		✓	✓	✓	✓		
B-492 polyester	FreezerBondz™ ultra-thin polyester for labeling frosted, frozen or room-temperature samples. Superior chemical and abrasion resistance.			●	●		✓			✓		
B-494 color polyester	Color polyester for quick differentiation between samples. Withstands harsh environments, including liquid nitrogen, chemicals and more.	●		●	●	●	✓		✓	✓		
B-498 vinyl cloth	Repositionable vinyl cloth for re-identification of post-process tissue cassettes or general lab applications. Can be applied directly onto the cassettes or on top of the previous label.		●		●	●				✓		
B-499 nylon cloth	Nylon cloth with permanent adhesive for easy handling and conformability.	●		●	●	●	✓	✓		✓		
B-533 polyester	Removable polyester for applications requiring clean removal.					●					✓	
B-595 vinyl	Vinyl with aggressive adhesive for drawers, shelves, cabinets and general lab labeling. Available in a variety of colors.					●					✓	
B-813 direct thermal paper	Paper material offering good adhesion to corrugated surfaces, glass and many plastics. For direct thermal printing.					●						
B-7425 polypropylene	Economical polypropylene with moderate chemical resistance, ideal for labeling dry, room-temperature specimen containers. Must be applied to containers prior to being frozen.	●		●	●	●	✓	✓	✓	✓		
B-7425J polypropylene	Inkjet printable polypropylene ideal for full-color identification and barcodes where resistance to liquid nitrogen, autoclave conditions and harsh chemicals are required. Must be applied to containers prior to being frozen.	●		●	●	●	✓	✓	✓	✓		
B-7425-AC polypropylene	Steam sterilization indicating polypropylene with a pre-printed header that displays STERILIZED when exposed to autoclave conditions.	●		●	●	●	✓	✓	✓	✓		

Label material guide

Brady products are sold with the understanding that the buyers will test them in actual use and determine for themselves their adaptability to their intended uses. Brady warrants to the buyers that its products are free from defects in material and workmanship, but limits its obligation under this warranty to replacement of the product shown to Brady's satisfaction to have been defective at the time Brady sold it. This warranty does not extend to any persons obtaining the product from the buyers. This warranty is in lieu of any other warranty, express or implied, including, but not limited to, any implied warranty of merchantability or fitness for a particular purpose and of any other obligations or liability on Brady's part. Under no circumstances will Brady be liable for any loss, damage, expense, or consequential damages of any kind arising in connection with the use, or inability to use, Brady's products.

Microbiology labels

Cell culture flask labels

Cell culture flask labels need to withstand varying storage conditions and thaw cycles while allowing an easy view of any growth in media. These labels are designed to withstand both temperature changes and chemicals so they stay intact and clearly legible throughout lab processes.

- B-430 – Clear high-performance polyester label
- B-432 – Clear polyester withstands exposure to solvents and exhibits good adhesion to metal and other surfaces
- B-488 – White polyester with moderate chemical resistance and short-term exposure to xylene, DMSO and ethanol
- B-494 – Glossy polyester with a printed color bar that allows for quick visual differentiation between samples. Withstands liquid nitrogen, a variety of chemicals and more. Label dimensions listed below include the color bar area. Vial tops are solid color
- B-7425 / B-7425J – Economical white polypropylene label with moderate chemical resistance that is ideal for labeling dry, room-temperature specimen containers. Must be applied to containers prior to being frozen









Figure	Catalog #	Material	Color	Label width A inch (mm)	Label length B inch (mm)	Secondary label C inch (mm)	Labels per row	Labels per roll	Recommended ribbon
BMP®21-LAB Label Printer cartridges									
Fig. 33b	M21-500-430	B-430	Black on clear	0.500 (12.7)	21 ft. (6.4 m)	-	1	1 roll (cont.)	In cartridge
Fig. 33b	M21-500-7425	B-7425	Black on white	0.500 (12.7)	21 ft. (6.4 m)	-	1	1 roll (cont.)	In cartridge
Fig. 33b	M21-750-430	B-430	Black on clear	0.750 (19.1)	21 ft. (6.4 m)	-	1	1 roll (cont.)	In cartridge
Fig. 33b	M21-750-7425	B-7425	Black on white	0.750 (19.1)	21 ft. (6.4 m)	-	1	1 roll (cont.)	In cartridge
BMP®51 Label Printer cartridges									
Fig. 2	M-17-432-CL-BK	B-432	Black on clear	1.000 (25.4)	0.500 (12.7)	-	1	280	In cartridge
Fig. 2	M-49-494-BL	B-494	Black on blue / white	1.000 (25.4)	1.000 (25.4)	-	1	260	In cartridge
Fig. 2	M-49-494-GN	B-494	Black on green / white	1.000 (25.4)	1.000 (25.4)	-	1	260	In cartridge
Fig. 2	M-49-494-OR	B-494	Black on orange / white	1.000 (25.4)	1.000 (25.4)	-	1	260	In cartridge
Fig. 2	M-49-494-PK	B-494	Black on pink / white	1.000 (25.4)	1.000 (25.4)	-	1	260	In cartridge
Fig. 2	M-49-494-RD	B-494	Black on red / white	1.000 (25.4)	1.000 (25.4)	-	1	260	In cartridge
Fig. 2	M-49-494-YL	B-494	Black on yellow / white	1.000 (25.4)	1.000 (25.4)	-	1	260	In cartridge
Fig. 2	M-81-488	B-488	Black on white	1.900 (48.3)	0.250 (6.4)	-	1	160	In cartridge
BradyPrinter M611 label rolls									
Fig. 2	M61-18-7425	B-7425	White	1.000 (25.4)	0.750 (19.1)	-	1	250	M61-R4310
Fig. 2	M61-19-494-BL	B-494	Blue / white	1.000 (25.4)	1.000 (25.4)	-	1	250	M61-R6410
Fig. 2	M61-19-494-GN	B-494	Green / white	1.000 (25.4)	1.000 (25.4)	-	1	250	M61-R6410
Fig. 2	M61-19-494-OR	B-494	Orange / white	1.000 (25.4)	1.000 (25.4)	-	1	250	M61-R6410
Fig. 2	M61-19-494-PK	B-494	Pink / white	1.000 (25.4)	1.000 (25.4)	-	1	250	M61-R6410
Fig. 2	M61-19-494-RD	B-494	Red / white	1.000 (25.4)	1.000 (25.4)	-	1	250	M61-R6410
Fig. 2	M61-19-494-YL	B-494	Yellow / white	1.000 (25.4)	1.000 (25.4)	-	1	250	M61-R6410

Cell culture flask labels continue on next page.

Cell culture flask labels (continued)

Figure	Catalog #	Material	Color	Label width A inch (mm)	Label length B inch (mm)	Secondary label C inch (mm)	Labels per row	Labels per roll	Recommended ribbon
BradyPrinter i3300 and BBP®33 Printer label rolls									
Fig. 24	B33-251-7425	B-7425	White	1.000 (25.4)	0.500 (12.7)	0.440 (11.2) dia.	1	1,500	B30-R4300
Fig. 25	B33-88-7425	B-7425	White	1.000 (25.4)	0.750 (19.1)	-	2	2,000	B30-R4300
Fig. 2	B33-179-494-BL	B-494	Blue / white	1.000 (25.4)	1.000 (25.4)	-	1	1,500	B30-R6400
Fig. 2	B33-179-494-GN	B-494	Green / white	1.000 (25.4)	1.000 (25.4)	-	1	1,500	B30-R6400
Fig. 2	B33-179-494-OR	B-494	Orange / white	1.000 (25.4)	1.000 (25.4)	-	1	1,500	B30-R6400
Fig. 2	B33-179-494-PK	B-494	Pink / white	1.000 (25.4)	1.000 (25.4)	-	1	1,500	B30-R6400
Fig. 2	B33-179-494-RD	B-494	Red / white	1.000 (25.4)	1.000 (25.4)	-	1	1,500	B30-R6400
Fig. 2	B33-179-494-YL	B-494	Yellow / white	1.000 (25.4)	1.000 (25.4)	-	1	1,500	B30-R6400
Fig. 2	B33-179-7425	B-7425	White	1.000 (25.4)	1.000 (25.4)	-	1	1,500	B30-R4300
Fig. 2	B33-256-7425	B-7425	White	1.300 (33.0)	0.600 (15.2)	-	1	1,500	B30-R4300
Fig. 2	B33-235-7425	B-7425	White	1.500 (38.1)	1.000 (25.4)	-	1	1,250	B30-R4300
Fig. 2	B33-17-430	B-430	Clear	2.000 (50.8)	1.000 (25.4)	-	1	1,500	B30-R6000
Fig. 2	B33-137-7425	B-7425	White	2.000 (50.8)	1.000 (25.4)	-	1	1,500	B30-R4300
BradyJet J2000 Printer label rolls									
Fig. 24	J20-251-7425J	B-7425J	White	1.000 (25.4)	0.500 (12.7)	0.440 (11.2) dia.	2	1,500	J20-CMY
Fig. 25	J20-88-7425J	B-7425J	White	1.000 (25.4)	0.750 (19.1)	-	2	2,000	J20-CMY
Fig. 2	J20-179-7425J	B-7425J	White	1.000 (25.4)	1.000 (25.4)	-	1	1,500	J20-CMY
Fig. 2	J20-256-7425J	B-7425J	White	1.300 (33.0)	0.600 (15.2)	-	1	1,500	J20-CMY
Fig. 2	J20-235-7425J	B-7425J	White	1.500 (38.1)	1.000 (25.4)	-	1	1,250	J20-CMY
Fig. 2	J20-137-7425J	B-7425J	White	2.000 (50.8)	1.000 (25.4)	-	1	1,500	J20-CMY
BBP®12 Printer and other industry printers using 1-inch small core label rolls									
Fig. 25	THT-140-488-1.5-SC	B-488	White	0.900 (22.9)	0.750 (19.1)	-	2	1,500	R4313
Fig. 25	THT-88-7425-2-SC	B-7425	White	1.000 (25.4)	0.750 (19.1)	-	2	2,000	R6413
Fig. 24	THT-251-7425-2-SC	B-7425	White	1.000 (25.4)	0.500 (12.7)	0.440 (11.2) dia.	1	2,000	R6413
Fig. 2	THT-179-7425-2-SC	B-7425	White	1.000 (25.4)	1.000 (25.4)	-	1	2,000	R6413
Fig. 2	THT-256-7425-2-SC	B-7425	White	1.300 (33.0)	0.600 (15.2)	-	1	2,000	R6413
Fig. 2	THT-235-7425-2-SC	B-7425	White	1.500 (38.1)	1.000 (25.4)	-	1	2,000	R6413
Fig. 2	THT-137-7425-2-SC	B-7425	White	2.000 (50.8)	1.000 (25.4)	-	1	2,000	R6413
BradyPrinter i5100, i7100, Brady IP Printers and other industry printers using 3-inch core label rolls									
Fig. 25	THT-140-488-3	B-488	White	0.900 (22.9)	0.750 (19.1)	-	2	3,000	R4302 / IP-R4302
Fig. 1	THT-5-488-10	B-488	White	1.000 (25.4)	0.500 (12.7)	-	3	10,000	R4307 / IP-R4307
Fig. 24	THT-251-7425-2	B-7425	White	1.000 (25.4)	0.500 (12.7)	0.440 (11.2) dia.	1	2,000	R4302 / IP-R4302
Fig. 25	THT-88-7425-2	B-7425	White	1.000 (25.4)	0.750 (19.1)	-	2	2,000	R4302 / IP-R4302
Fig. 2	THT-179-494-BL	B-494	Blue / white	1.000 (25.4)	1.000 (25.4)	-	1	3,000	R6402 / IP-R6402
Fig. 2	THT-179-494-GN	B-494	Green / white	1.000 (25.4)	1.000 (25.4)	-	1	3,000	R6402 / IP-R6402
Fig. 2	THT-179-494-OR	B-494	Orange / white	1.000 (25.4)	1.000 (25.4)	-	1	3,000	R6402 / IP-R6402
Fig. 2	THT-179-494-PK	B-494	Pink / white	1.000 (25.4)	1.000 (25.4)	-	1	3,000	R6402 / IP-R6402
Fig. 2	THT-179-494-RD	B-494	Red / white	1.000 (25.4)	1.000 (25.4)	-	1	3,000	R6402 / IP-R6402
Fig. 2	THT-179-494-YL	B-494	Yellow / white	1.000 (25.4)	1.000 (25.4)	-	1	3,000	R6402 / IP-R6402
Fig. 2	THT-179-7425-2	B-7425	White	1.000 (25.4)	1.000 (25.4)	-	1	2,000	R4302 / IP-R4302
Fig. 2	THT-256-7425-2	B-7425	White	1.300 (33.0)	0.600 (15.2)	-	1	2,000	R4302 / IP-R4302
Fig. 2	THT-235-7425-2	B-7425	White	1.500 (38.1)	1.000 (25.4)	-	1	2,000	R4302 / IP-R4302
Fig. 2	THT-15-488-2.5	B-488	White	2.000 (50.8)	0.250 (6.4)	-	1	2,500	R4302 / IP-R4302
Fig. 2	THT-137-7425-2	B-7425	White	2.000 (50.8)	1.000 (25.4)	-	1	2,000	R4302 / IP-R4302
Fig. 2	THT-17-430-3	B-430	Clear	2.000 (50.8)	1.000 (25.4)	-	1	3,000	R6000 / IP-R6000

Petri dish labels

When labeling petri dishes, label materials need to withstand long culture periods, have a low profile for lab automation processes and resist chemicals and varying temperatures.

- B-430 – Clear high-performance polyester label
- B-432 – Clear polyester withstands exposure to solvents and exhibits good adhesion to metal and other surfaces
- B-488 – White polyester with moderate chemical resistance and short-term exposure to xylene, DMSO and ethanol
- B-499 – White nylon cloth with permanent adhesive for easy handling and conformability
- B-7425 / B-7425J – Economical white polypropylene label with moderate chemical resistance that is ideal for labeling dry, room-temperature specimen containers. Must be applied to containers prior to being frozen









Figure	Catalog #	Material	Color	Label width A inch (mm)	Label length B inch (mm)	Labels per row	Labels per roll	Recommended ribbon
BMP®21-LAB Label Printer cartridges								
Fig. 33b	M21-250-430	B-430	Black on clear	0.250 (6.4)	21 ft. (6.4 m)	1	1 roll (cont.)	In cartridge
Fig. 33b	M21-375-430	B-430	Black on clear	0.375 (9.5)	21 ft. (6.4 m)	1	1 roll (cont.)	In cartridge
Fig. 33b	M21-375-7425	B-7425	Black on white	0.375 (9.5)	21 ft. (6.4 m)	1	1 roll (cont.)	In cartridge
Fig. 33b	M21-500-430	B-430	Black on clear	0.500 (12.7)	21 ft. (6.4 m)	1	1 roll (cont.)	In cartridge
Fig. 33b	M21-500-7425	B-7425	Black on white	0.500 (12.7)	21 ft. (6.4 m)	1	1 roll (cont.)	In cartridge
BMP®51 Label Printer cartridges								
Fig. 2	M-17-432-CL-BK	B-432	Black on clear	1.000 (25.4)	0.500 (12.7)	1	280	In cartridge
Fig. 2	M-81-488	B-488	Black on white	1.900 (48.3)	0.250 (6.4)	1	160	In cartridge
BradyPrinter M611 label rolls								
Fig. 82	M61-2-7425	B-7425	White	0.250 (6.4)	2.000 (50.8)	1	100	M61-R4310
BradyPrinter i3300 and BBP®33 Printer label rolls								
Fig. 25	B33-152-499	B-499	White	1.000 (25.4)	0.375 (9.5)	2	1,500	B30-R6000
Fig. 2	B33-59-7425	B-7425	White	1.000 (25.4)	0.500 (12.7)	1	2,000	B30-R4300
Fig. 25	B33-37-7425	B-7425	White	1.500 (38.1)	0.500 (12.7)	2	1,500	B30-R4300
Fig. 2	B33-151-490	B-490	White	2.000 (50.8)	0.250 (6.4)	1	2,500	B30-R4300
Fig. 2	B33-151-499	B-499	White	2.000 (50.8)	0.250 (6.4)	1	1,500	B30-R6000
Fig. 2	B33-151-7425	B-7425	White	2.000 (50.8)	0.250 (6.4)	1	2,000	B30-R4300
Fig. 2	B33-53-7425	B-7425	White	2.000 (50.8)	0.500 (12.7)	1	2,000	B30-R4300
Fig. 2	B33-17-430	B-430	Clear	2.000 (50.8)	1.000 (25.4)	1	1,500	B30-R6000

Petri dish labels continue on next page.

Petri dish labels (continued)

Figure	Catalog #	Material	Color	Label width A inch (mm)	Label length B inch (mm)	Labels per row	Labels per roll	Recommended ribbon
BradyJet J2000 Printer label rolls								
Fig. 2	J20-59-7425J	B-7425J	White	1.000 (25.4)	0.500 (12.7)	1	2,000	J20-CMY
Fig. 25	J20-37-7425J	B-7425J	White	1.500 (38.1)	0.500 (12.7)	2	1,500	J20-CMY
Fig. 2	J20-151-7425J	B-7425J	White	2.000 (50.8)	0.250 (6.4)	1	2,000	J20-CMY
Fig. 2	J20-53-7425J	B-7425J	White	2.000 (50.8)	0.500 (12.7)	1	2,000	J20-CMY
BBP®12 Printer and other industry printers using 1-inch small core label rolls								
Fig. 2	THT-59-7425-2-SC	B-7425	White	1.000 (25.4)	0.500 (12.7)	1	2,000	R6413
Fig. 25	THT-37-7425-2-SC	B-7425	White	1.500 (38.1)	0.500 (12.7)	2	2,000	R6413
Fig. 2	THT-151-7425-2-SC	B-7425	White	2.000 (50.8)	0.250 (6.4)	1	2,000	R6413
Fig. 2	THT-53-7425-2-SC	B-7425	White	2.000 (50.8)	0.500 (12.7)	1	2,000	R6413
BradyPrinter i5100, i7100, Brady IP Printers and other industry printers using 3-inch core label rolls								
Fig. 1	THT-5-488-10	B-488	White	1.000 (25.4)	0.500 (12.7)	3	10,000	R4307 / IP-R4307
Fig. 2	THT-59-7425-2	B-7425	White	1.000 (25.4)	0.500 (12.7)	1	2,000	R4302 / IP-R4302
Fig. 25	THT-37-7425-2	B-7425	White	1.500 (38.1)	0.500 (12.7)	2	2,000	R4307 / IP-R4307
Fig. 2	THT-15-488-2.5	B-488	White	2.000 (50.8)	0.250 (6.4)	1	2,500	R4302 / IP-R4302
Fig. 2	THT-151-7425-2	B-7425	White	2.000 (50.8)	0.250 (6.4)	1	2,000	R4302 / IP-R4302
Fig. 2	THT-53-7425-2	B-7425	White	2.000 (50.8)	0.500 (12.7)	1	2,000	R4302 / IP-R4302
Fig. 2	THT-17-430-3	B-430	Clear	2.000 (50.8)	1.000 (25.4)	1	3,000	R6000 / IP-R6000

Custom lab labels

Brady understands that each lab is unique and may need custom labels that are specifically designed and compatible with your existing equipment or print systems.

Custom label solutions include:

- Blank labels in any size or color to easily organize and visualize batches of samples
- Custom text, graphics, colors and serialized barcodes
- Custom print and apply labels for automation equipment
- Pre-printed labels in the sizes and materials needed
- Fully printed custom labels or partially pre-printed labels so you can add additional text as needed

Order your custom labels today by calling 1-888-272-3946.

Histology labels

Microscope slide labels

Brady StainerBondz™ Labels are designed specifically to identify microscope slides and to resist exposure to harsh chemicals during the slide staining process – including hematoxylin and eosin.

- B-481 – White polyester with extreme chemical resistance that withstands the rigors of typical deparaffinization, including ethanol, xylene and H&E staining process.
- B-488 – White polyester with moderate chemical resistance for short-term exposure to xylene, DMSO and ethanol




Figure 1




Figure 2




Figure 25




Figure 33b

Figure	Catalog #	Material	Color	Label width A inch (mm)	Label length B inch (mm)	Labels per row	Labels per roll	Recommended ribbon
BMP®21-LAB Label Printer cartridges								
Fig. 33b	M21-375-488	B-488	Black on white	0.375 (9.5)	21 ft. (6.4 m)	1	1 roll (cont.)	In cartridge
Fig. 33b	M21-500-488	B-488	Black on white	0.500 (12.7)	21 ft. (6.4 m)	1	1 roll (cont.)	In cartridge
Fig. 33b	M21-750-488	B-488	Black on white	0.750 (19.1)	21 ft. (6.4 m)	1	1 roll (cont.)	In cartridge
BMP®51 Label Printer cartridges								
Fig. 2	M-195-481	B-481	Black on white	0.900 (22.9)	0.375 (9.5)	1	320	In cartridge
Fig. 2	M-194-481	B-481	Black on white	0.900 (22.9)	0.750 (19.1)	1	320	In cartridge
Fig. 2	M-97-481	B-481	Black on white	0.900 (22.9)	0.900 (22.9)	1	280	In cartridge
Fig. 2	M-97-488	B-488	Black on white	0.900 (22.9)	0.900 (22.9)	1	280	In cartridge
BradyPrinter M611 label rolls								
Fig. 2	PTL-80-488	B-488	White	0.750 (19.1)	0.900 (22.9)	1	250	M61-R4310
Fig. 2	M61-195-481	B-481	White	0.900 (22.9)	0.375 (9.5)	1	500	M61-R6410
Fig. 2	M61-97-481	B-481	White	0.900 (22.9)	0.900 (22.9)	1	250	M61-R6410
Fig. 2	M61-194-481	B-481	White	0.750 (19.1)	0.750 (19.1)	1	250	M61-R4310
Fig. 2	PTL-97-488	B-488	White	0.900 (22.9)	0.900 (22.9)	1	250	M61-R4310
BradyPrinter i3300 and BBP®33 Printer label rolls								
Fig. 2	B33-195-481	B-481	White	0.900 (22.9)	0.375 (9.5)	1	3,000	B30-R6000
Fig. 2	B33-194-481	B-481	White	0.900 (22.9)	0.750 (19.1)	1	2,000	B30-R6400
Fig. 25	B33-141-481	B-481	White	0.900 (22.9)	0.900 (22.9)	2	1,500	B30-R6400
BBP®12 Printer and other industry printers using 1-inch small core label rolls								
Fig. 2	THT-194-481-2.5-SC	B-481	White	0.900 (22.9)	0.750 (19.1)	1	2,500	R6413
Fig. 25	THT-141-481-1.5-SC	B-481	White	0.900 (22.9)	0.900 (22.9)	2	1,500	R6413
Fig. 2	THT-164-481-2.5-SC	B-481	White	0.900 (22.9)	0.900 (22.9)	1	2,500	R6413
Fig. 1	THT-207-481-2.5-SC	B-481	White	0.900 (22.9)	0.900 (22.9)	4	2,500	R6413
Fig. 2	THT-164-488-1-SC	B-488	White	0.900 (22.9)	0.900 (22.9)	1	1,000	R4313
Fig. 2	THT-142-481-2.5-SC	B-481	White	2.500 (63.5)	0.200 (5.1)	1	2,500	R6413
BradyPrinter i5100, i7100, Brady IP Printers and other industry printers using 3-inch core label rolls								
Fig. 2	THT-195-481-3	B-481	White	0.900 (22.9)	0.375 (9.5)	1	3,000	R6400 / IP-R6400
Fig. 2	THT-194-481-3	B-481	White	0.900 (22.9)	0.750 (19.1)	1	3,000	R6400 / IP-R6400
Fig. 2	THT-164-481-3	B-481	White	0.900 (22.9)	0.900 (22.9)	1	3,000	R6400 / IP-R6400
Fig. 25	THT-141-481-3	B-481	White	0.900 (22.9)	0.900 (22.9)	2	3,000	R6400 / IP-R6400
Fig. 2	THT-164-488-1	B-488	White	0.900 (22.9)	0.900 (22.9)	1	1,000	R4302 / IP-R4302
Fig. 25	THT-141-488-3	B-488	White	0.900 (22.9)	0.900 (22.9)	2	3,000	R4302 / IP-R4302
Fig. 2	THT-142-481-2.5	B-481	White	2.500 (63.5)	0.200 (5.1)	1	2,500	R6402 / IP-R6402

Tissue cassette labels

Tissue cassette labels are designed for post-process identification. They have high durability and maintain print quality for accurate, identifiable specimen labeling.

- B-498 – Removable white vinyl cloth label for re-identification of post-process tissue cassettes or general lab applications




Figure 2



Figure	Catalog #	Material	Color	Label width A inch (mm)	Label length B inch (mm)	Labels per row	Labels per roll	Recommended ribbon
BMP®51 Label Printer cartridges								
Fig. 2	M-11-498	B-498	Black on white	0.750 (19.1)	0.500 (12.7)	1	220	In cartridge
BradyPrinter i5100, i7100, Brady IP Printers and other industry printers using 3-inch core label rolls								
Fig. 2	THT-252-498-3	B-498	White	0.750 (19.1)	0.250 (6.4)	1	3,000	R6200 / IP-R6200
Fig. 2	THT-254-498-3	B-498	White	0.750 (19.1)	0.750 (19.1)	1	3,000	R6200 / IP-R6200



Test the materials yourself with a free sample pack

Request a free sample pack of materials and put them to the test. Get your samples at BradyID.com/lab

View the Technical Data Sheets

Brady offers online access to Technical Data Sheets for all materials at BradyID.com/techdata

Cryo box and freezer rack labels

Brady labels adhere to either frozen or room-temperature boxes and racks and are designed to remain legible and intact when faced with harsh lab conditions.

- B-490 – FreezerBondz™ white polyester is ideal for flat, frozen surfaces. Must be applied to containers prior to being frozen
- B-492 – FreezerBondz ultra-thin white polyester with superior chemical and abrasion resistance for labeling frosted, frozen or room-temperature boxes and racks
- B-499 – White nylon cloth with permanent adhesive for easy handling and conformability
- B-7425 / B-7425J – Economical white polypropylene label with moderate chemical resistance that is ideal for labeling dry, room-temperature specimen containers. Must be applied to containers prior to being frozen




Figure 2




Figure 3




Figure 24




Figure 25




Figure 33b




Figure 82




Figure 83

Figure	Catalog #	Material	Color	Label width A inch (mm)	Label length B inch (mm)	Secondary label C inch (mm)	Labels per row	Labels per roll	Recommended ribbon
BMP®21-LAB Label Printer cartridges									
Fig. 33b	M21-500-499	B-499	Black on white	0.500 (12.7)	16 ft. (4.8 m)	-	1	1 roll (cont.)	In cartridge
Fig. 33b	M21-500-7425	B-7425	Black on white	0.500 (12.7)	21 ft. (6.4 m)	-	1	1 roll (cont.)	In cartridge
Fig. 33b	M21-750-499	B-499	Black on white	0.750 (19.1)	16 ft. (4.8 m)	-	1	1 roll (cont.)	In cartridge
Fig. 33b	M21-750-7425	B-7425	Black on white	0.750 (19.1)	21 ft. (6.4 m)	-	1	1 roll (cont.)	In cartridge
BMP®51 Label Printer cartridges									
Fig. 82	M-124-490	B-490	Black on white	0.500 (12.7)	1.650 (41.9)	-	1	180	In cartridge
Fig. 83	M-95-492	B-492	Black on white	0.600 (15.2)	1.625 (41.3)	0.375 (9.5) dia.	1	200	In cartridge
Fig. 2	M-132-499	B-499	Black on white	1.500 (38.1)	1.250 (31.7)	-	1	160	In cartridge
BradyPrinter M611 label rolls									
Fig. 82	M61-124-490	B-490	White	0.500 (12.7)	1.650 (41.9)	-	1	500	M61-R4310
Fig. 2	M61-132-499	B-499	White	1.500 (38.1)	1.250 (31.8)	-	1	250	M61-R6210
Fig. 2	PTL-76-490	B-490	White	1.625 (41.3)	0.600 (15.2)	-	1	250	M61-R4310
Fig. 2	M61-121-7425	B-7425	White	1.900 (48.3)	1.500 (38.1)	-	1	100	M61-R4310
BradyPrinter i3300 and BBP®33 Printer label rolls									
Fig. 2	B33-59-7425	B-7425	White	1.000 (25.4)	0.500 (12.7)	-	1	2,000	B30-R4300
Fig. 25	B33-37-7425	B-7425	White	1.500 (38.1)	0.500 (12.7)	-	2	1,500	B30-R4300
Fig. 2	B33-136-499	B-499	White	1.500 (38.1)	0.750 (19.1)	-	1	1,500	B30-R6000
Fig. 2	B33-235-7425	B-7425	White	1.500 (38.1)	1.000 (25.4)	-	1	1,250	B30-R4300
Fig. 2	B33-155-490	B-490	White	1.625 (41.3)	0.600 (15.2)	-	1	1,500	B30-R4300
Fig. 24	B33-190-492	B-492	White	1.625 (41.3)	0.600 (15.2)	0.375 (9.5) dia.	1	1,500	B30-R6400
Fig. 2	B33-137-499	B-499	White	2.000 (50.8)	1.000 (25.4)	-	1	1,100	B30-R6000
BradyJet J2000 Printer label rolls									
Fig. 2	J20-59-7425J	B-7425J	White	1.000 (25.4)	0.500 (12.7)	-	1	2,000	J20-CMY
Fig. 3	J20C-1125-7425J	B-7425J	White	1.125 (28.6)	100 ft. (30.5 m)	-	1	1 roll (cont.)	J20-CMY
Fig. 25	J20-37-7425J	B-7425J	White	1.500 (38.1)	0.500 (12.7)	-	2	1,500	J20-CMY
Fig. 2	J20-235-7425J	B-7425J	White	1.500 (38.1)	1.000 (25.4)	-	1	1,250	J20-CMY
Fig. 3	J20C-2250-7425J	B-7425J	White	2.250 (57.2)	100 ft. (30.5 m)	-	1	1 roll (cont.)	J20-CMY
BBP®12 Printer and other industry printers using 1-inch small core label rolls									
Fig. 2	THT-59-7425-2-SC	B-7425	White	1.000 (25.4)	0.500 (12.7)	-	1	2,000	R6413
Fig. 25	THT-37-7425-2-SC	B-7425	White	1.500 (38.1)	0.500 (12.7)	-	2	2,000	R6413
Fig. 2	THT-235-7425-2-SC	B-7425	White	1.500 (38.1)	1.000 (25.4)	-	1	2,000	R6413

Cryo box and freezer rack labels continue on next page.

Cryo box and freezer rack labels (continued)

Figure	Catalog #	Material	Color	Label width A inch (mm)	Label length B inch (mm)	Secondary label C inch (mm)	Labels per row	Labels per roll	Recommended ribbon
Fig. 2	THT-155-490-1.5-SC	B-490	White	1.625 (41.3)	0.600 (15.2)	-	1	1,500	R4313
Fig. 2	THT-156-490-1.5-SC	B-490	White	1.900 (48.3)	1.000 (25.4)	-	1	1,500	R4313
Fig. 2	THT-156-492-1.5-SC	B-492	White	1.900 (48.3)	1.000 (25.4)	-	1	1,500	R6413
Fig. 2	THT-137-499-1.5-SC	B-499	White	2.000 (50.8)	1.000 (25.4)	-	1	1,500	R4913
BradyPrinter i5100, i7100, Brady IP Printers and other industry printers using 3-inch core label rolls									
Fig. 2	THT-125-499-3	B-499	White	0.900 (22.9)	0.500 (12.7)	-	1	3,000	R6000 / IP-R6000
Fig. 2	THT-59-7425-2	B-7425	White	1.000 (25.4)	0.500 (12.7)	-	1	2,000	R4302 / IP-R4302
Fig. 25	THT-37-7425-2	B-7425	White	1.500 (38.1)	0.500 (12.7)	-	2	2,000	R4307 / IP-R4307
Fig. 2	THT-235-7425-2	B-7425	White	1.500 (38.1)	1.000 (25.4)	-	1	2,000	R4302 / IP-R4302
Fig. 2	THT-155-490-3	B-490	White	1.625 (41.3)	0.600 (15.2)	-	1	3,000	R4302 / IP-R4302
Fig. 2	THT-190-492-3	B-492	White	1.625 (41.3)	0.600 (15.2)	0.375 (9.5) dia.	1	3,000	R6402 / IP-R6402
Fig. 2	THT-185-499-3	B-499	White	1.750 (44.5)	1.000 (25.4)	-	1	3,000	R6000 / IP-R6000
Fig. 2	THT-156-490-3	B-490	White	1.900 (48.3)	1.000 (25.4)	-	1	3,000	R4302 / IP-R4302
Fig. 2	THT-156-492-3	B-492	White	1.900 (48.3)	1.000 (25.4)	-	1	3,000	R6400 / IP-R6400
Fig. 2	THT-137-499-3	B-499	White	2.000 (50.8)	1.000 (25.4)	-	1	3,000	R6000 / IP-R6000
Fig. 2	THT-19-499-1	B-499	White	3.000 (76.2)	2.000 (50.8)	-	1	1,000	R6000 / IP-R6000

Cryo cane labels

Cryo canes typically have a tab on the top to hold a label that can be read while the vials are in cryo drawers or freezers. Brady labels are designed to adhere to the canes and withstand the ultra-low temperatures.

- B-492 – FreezerBondz™ ultra-thin white polyester label with superior chemical and abrasion resistance for labeling of frosted, frozen or room-temperature canes




Figure 66



Figure	Catalog #	Material	Color	Label width A inch (mm)	Label length B inch (mm)	Labels per row	Labels per roll	Recommended ribbon
BMP®51 Cartridge								
Fig. 66	M-155-492	B-492	Black on white	1.250 (31.8)	0.250 (6.4)	1	180	In cartridge

Straw and rod labels

Brady labels are engineered to effectively wrap around 2-3 mm diameter straws and rods with an overlaminates so that they withstand freezing temperatures.

- B-427 – Durable vinyl label with clear overlaminates is good for wrapping and conformability. Must be wrapped upon itself when used in liquid nitrogen




Figure 2

Figure 14

Figure 14c




Figure	Catalog #	Material	Color	Label width A inch (mm)	Label length B inch (mm)	Print-on area C inch (mm)	Labels per row	Labels per roll	Recommended ribbon
BMP®51 Label Printer cartridges									
Fig. 2	M-135-427-AW	B-427	Black on white	0.250 (6.4)	1.100 (27.9)	-	1	240	In cartridge
Fig. 14	M-133-427	B-427	Black on white	1.750 (44.5)	1.000 (25.4)	0.375 (9.5)	1	140	In cartridge
Fig. 14	M-134-427	B-427	Black on white	1.750 (44.5)	1.500 (38.1)	0.500 (12.7)	1	140	In cartridge
BradyPrinter i3300 and BBP®33 Printer label rolls									
Fig. 14	B33-120-427	B-427	White	1.750 (44.5)	1.000 (25.4)	0.375 (9.5)	1	1,250	B30-R4300
Fig. 14	B33-121-427	B-427	White	1.750 (44.5)	1.500 (38.1)	0.500 (12.7)	1	1,000	B30-R4300
BradyPrinter i5100, i7100, Brady IP Printers and other industry printers using 3-inch core label rolls									
Fig. 14c	THT-116-427-5	B-427	White	1.000 (25.4)	1.000 (25.4)	0.375 (9.5)	2	5,000	R4300 / IP-R4300

Cryo tube labels

FreezerBondz™ labels are designed to remain legible and stick to the tubes during processing, even when immersed in liquid or vapor-phase nitrogen and ultra-low temperature freezer storage.

- B-461 – All white, ultra-thin polyester label. Also available in a self-laminating format
- B-490 – White polyester ideal for flat frozen surfaces. Adheres to frozen tube when wrapped around itself by at least 0.125 inch
- B-492 – Ultra-thin white polyester with superior chemical and abrasion resistance for labeling frosted, frozen or room-temperature tubes and vials. Must have small overlap if applied to already frozen containers
- B-494 – Glossy polyester with a printed color bar that allows for quick visual differentiation between samples. Withstands liquid nitrogen, a variety of chemicals and more. Label dimensions listed below include the color bar area. Vial tops are solid color
- B-499 – White nylon cloth with permanent adhesive for easy handling and conformability
- B-7425 / B-7425J – Economical white polypropylene label with moderate chemical resistance that is ideal for labeling dry, room-temperature specimen containers. Must be applied to containers prior to being frozen



Fits tube size*	Figure	Catalog #	Material	Color	Label width A inch (mm)	Label length B inch (mm)	Secondary label C inch (mm)	Labels per row	Labels per roll	Recommended ribbon
BMP®21-LAB Label Printer cartridges										
1.2 mL	Fig. 33b	M21-375-461	B-461	Black on white	0.375 (9.5)	21 ft. (6.4 m)	-	1	1 roll (cont.)	In cartridge
	Fig. 33b	M21-375-7425	B-7425	Black on white	0.375 (9.5)	21 ft. (6.4 m)	-	1	1 roll (cont.)	In cartridge
	Fig. 33b	M21-500-461	B-461	Black on white	0.500 (12.7)	21 ft. (6.4 m)	-	1	1 roll (cont.)	In cartridge
	Fig. 33b	M21-500-7425	B-7425	Black on white	0.500 (12.7)	21 ft. (6.4 m)	-	1	1 roll (cont.)	In cartridge
2 mL	Fig. 33b	M21-750-461	B-461	Black on white	0.750 (19.1)	21 ft. (6.4 m)	-	1	1 roll (cont.)	In cartridge
	Fig. 33b	M21-750-7425	B-7425	Black on white	0.750 (19.1)	21 ft. (6.4 m)	-	1	1 roll (cont.)	In cartridge
BMP®51 Label Printer cartridges										
1.2 mL	Fig. 82	M-130-492	B-492	Black on white	0.375 (9.5)	0.825 (21.0)	-	1	220	In cartridge
	Fig. 82	M-130-499	B-499	Black on white	0.375 (9.5)	0.825 (21.0)	-	1	220	In cartridge
	Fig. 33	MC-375-499	B-499	Black on white	0.375 (9.5)	16 ft. (4.8 m)	-	1	1 roll (cont.)	In cartridge
	Fig. 33	MC-375-461	B-461	Black on white	0.375 (9.5)	25 ft. (7.6 m)	-	1	1 roll (cont.)	In cartridge
	Fig. 82	M-131-492	B-492	Black on white	0.500 (12.7)	1.000 (25.4)	-	1	180	In cartridge
	Fig. 82	M-131-494-BL	B-494	Black on blue / white	0.500 (12.7)	1.000 (25.4)	-	1	180	In cartridge
	Fig. 82	M-131-494-GN	B-494	Black on green / white	0.500 (12.7)	1.000 (25.4)	-	1	180	In cartridge
	Fig. 82	M-131-494-OR	B-494	Black on orange / white	0.500 (12.7)	1.000 (25.4)	-	1	180	In cartridge
	Fig. 82	M-131-494-PK	B-494	Black on pink / white	0.500 (12.7)	1.000 (25.4)	-	1	180	In cartridge
	Fig. 82	M-131-494-RD	B-494	Black on red / white	0.500 (12.7)	1.000 (25.4)	-	1	180	In cartridge
	Fig. 82	M-131-494-YL	B-494	Black on yellow / white	0.500 (12.7)	1.000 (25.4)	-	1	180	In cartridge
	Fig. 82	M-131-499	B-499	Black on white	0.500 (12.7)	1.000 (25.4)	-	1	180	In cartridge
	Fig. 82	M-124-492	B-492	Black on white	0.500 (12.7)	1.650 (41.9)	-	1	180	In cartridge
	Fig. 33	MC-500-499	B-499	Black on white	0.500 (12.7)	16 ft. (4.8 m)	-	1	1 roll (cont.)	In cartridge
	Fig. 33	MC-500-461	B-461	Black on white	0.500 (12.7)	25 ft. (7.6 m)	-	1	1 roll (cont.)	In cartridge
	Fig. 33	MC-750-499	B-499	Black on white	0.750 (19.1)	16 ft. (4.8 m)	-	1	1 roll (cont.)	In cartridge

Cryo tube labels continue on next page.

*Tube sizes may vary by brand.

Cryo tube labels (continued)


Fits tube size*	Figure	Catalog #	Material	Color	Label width A inch (mm)	Label length B inch (mm)	Secondary label C inch (mm)	Labels per row	Labels per roll	Recommended ribbon
1.2 mL	Fig. 24	M-118-492	B-492	Black on white	1.000 (25.4)	0.375 (9.5)	0.375 (9.5) dia.	2	240	In cartridge
	Fig. 24	M-118-494-BL	B-494	Black on blue / white	1.000 (25.4)	0.375 (9.5)	0.375 (9.5) dia.	1	240	In cartridge
	Fig. 24	M-118-494-GN	B-494	Black on green / white	1.000 (25.4)	0.375 (9.5)	0.375 (9.5) dia.	1	240	In cartridge
	Fig. 24	M-118-494-OR	B-494	Black on orange / white	1.000 (25.4)	0.375 (9.5)	0.375 (9.5) dia.	1	240	In cartridge
	Fig. 24	M-118-494-PK	B-494	Black on pink / white	1.000 (25.4)	0.375 (9.5)	0.375 (9.5) dia.	1	240	In cartridge
	Fig. 24	M-118-494-RD	B-494	Black on red / white	1.000 (25.4)	0.375 (9.5)	0.375 (9.5) dia.	1	250	In cartridge
	Fig. 24	M-118-494-YL	B-494	Black on yellow / white	1.000 (25.4)	0.375 (9.5)	0.375 (9.5) dia.	1	240	In cartridge
	Fig. 24	M-118-499	B-499	Black on white	1.000 (25.4)	0.375 (9.5)	0.375 (9.5) dia.	2	240	In cartridge
	Fig. 24	M-118-7425	B-7425	Black on white	1.000 (25.4)	0.375 (9.5)	0.375 (9.5) dia.	1	240	In cartridge
	Fig. 24	M-117-492	B-492	Black on white	1.000 (25.4)	0.500 (12.7)	0.438 (11.1) dia.	2	240	In cartridge
	Fig. 24	M-117-499	B-499	Black on white	1.000 (25.4)	0.500 (12.7)	0.438 (11.1) dia.	2	240	In cartridge
	Fig. 24	M-120-494-BL	B-494	Black on blue / white	1.000 (25.4)	0.500 (12.7)	0.375 (9.5) dia.	1	240	In cartridge
	Fig. 24	M-120-494-GN	B-494	Black on green / white	1.000 (25.4)	0.500 (12.7)	0.375 (9.5) dia.	1	240	In cartridge
	Fig. 24	M-120-494-OR	B-494	Black on orange / white	1.000 (25.4)	0.500 (12.7)	0.375 (9.5) dia.	1	240	In cartridge
	Fig. 24	M-120-494-PK	B-494	Black on pink / white	1.000 (25.4)	0.500 (12.7)	0.375 (9.5) dia.	1	240	In cartridge
	Fig. 24	M-120-494-RD	B-494	Black on red / white	1.000 (25.4)	0.500 (12.7)	0.375 (9.5) dia.	1	240	In cartridge
	Fig. 24	M-120-494-YL	B-494	Black on yellow / white	1.000 (25.4)	0.500 (12.7)	0.375 (9.5) dia.	1	240	In cartridge
	Fig. 24	M-120-492	B-492	Black on white	1.000 (25.4)	0.500 (12.7)	0.375 (9.5) dia.	1	240	In cartridge
Fig. 24	M-120-499	B-499	Black on white	1.000 (25.4)	0.500 (12.7)	0.375 (9.5) dia.	1	240	In cartridge	
Fig. 2	M-112-492	B-492	Black on white	1.250 (31.8)	0.375 (9.5)	-	1	360	In cartridge	
1.2 mL top only	Fig. 8	M-82-492	B-492	Black on white	0.375 (9.5) dia.	-	-	1	240	In cartridge
	Fig. 8	M-82-499	B-499	Black on white	0.375 (9.5) dia.	-	-	1	240	In cartridge
	Fig. 8	M-82-499-OR-BK	B-499	Black on orange	0.375 (9.5) dia.	-	-	1	240	In cartridge
	Fig. 8	M-82-499-RD-BK	B-499	Black on red	0.375 (9.5) dia.	-	-	1	240	In cartridge
	Fig. 8	M-82-499-YL-BK	B-499	Black on yellow	0.375 (9.5) dia.	-	-	1	240	In cartridge
Fig. 8	M-83-492	B-492	Black on white	0.500 (12.7) dia.	-	-	1	240	In cartridge	
2 mL	Fig. 2	M-156-492	B-492	Black on white	1.000 (25.4)	0.900 (22.9)	-	1	180	In cartridge
	Fig. 2	M-126-490	B-490	Black on white	1.800 (45.7)	0.600 (15.2)	-	1	180	In cartridge
5 mL	Fig. 2	M-49-494-BL	B-494	Black on blue / white	1.000 (25.4)	1.000 (25.4)	-	1	260	In cartridge
	Fig. 2	M-49-494-GN	B-494	Black on green / white	1.000 (25.4)	1.000 (25.4)	-	1	260	In cartridge
	Fig. 2	M-49-494-OR	B-494	Black on orange / white	1.000 (25.4)	1.000 (25.4)	-	1	260	In cartridge
	Fig. 2	M-49-494-PK	B-494	Black on pink / white	1.000 (25.4)	1.000 (25.4)	-	1	260	In cartridge
	Fig. 2	M-49-494-RD	B-494	Black on red / white	1.000 (25.4)	1.000 (25.4)	-	1	260	In cartridge
	Fig. 2	M-49-494-YL	B-494	Black on yellow / white	1.000 (25.4)	1.000 (25.4)	-	1	260	In cartridge
	Fig. 33	MC-1000-461	B-461	Black on white	1.000 (25.4)	25 ft. (7.6 m)	-	1	1 roll (cont.)	In cartridge
	Fig. 33	MC-1250-461	B-461	Black on white	1.250 (31.8)	25 ft. (7.6 m)	-	1	1 roll (cont.)	In cartridge
	Fig. 2	M-31-7425	B-7425	Black on white	1.500 (38.1)	1.000 (25.4)	-	1	250	In cartridge
	Fig. 2	M-125-490	B-490	Black on white	1.750 (44.5)	1.000 (25.4)	-	1	180	In cartridge
Fig. 33	M-125-492	B-492	Black on white	1.750 (44.5)	1.000 (25.4)	-	1	180	In cartridge	
15 mL	Fig. 82	M-113-490	B-490	Black on white	1.500 (38.1)	2.200 (55.9)	-	1	140	In cartridge
BradyPrinter M611 label rolls										
1.2 mL	Fig. 33	M61C-375-499	B-499	White	0.375 (9.5)	30 ft. (9.1 m)	-	1	1 roll (cont.)	M61-R6210
	Fig. 33	M61C-375-461	B-461	White	0.375 (9.5)	50 ft. (15.2 m)	-	1	1 roll (cont.)	M61-R4310
	Fig. 82	M61-47-499	B-499	White	0.500 (12.7)	1.000 (25.4)	-	1	500	M61-R6210
	Fig. 82	M61-124-492	B-492	White	0.500 (12.7)	1.650 (41.9)	-	1	500	M61-R6410
	Fig. 33	M61C-500-499	B-499	White	0.500 (12.7)	30 ft. (9.1 m)	-	1	1 roll (cont.)	M61-R6210
	Fig. 33	M61C-500-461	B-461	White	0.500 (12.7)	50 ft. (15.2 m)	-	1	1 roll (cont.)	M61-R4310
	Fig. 33	M61C-750-499	B-499	White	0.750 (19.1)	30 ft. (9.1 m)	-	1	1 roll (cont.)	M61-R6210
	Fig. 2	M61-130-492	B-492	White	0.825 (21.0)	0.375 (9.5)	-	1	500	M61-R6410
	Fig. 2	M61-130-499	B-499	White	0.825 (21.0)	0.375 (9.5)	-	1	500	M61-R6210
	Fig. 24	M61-117-492	B-492	White	0.900 (22.9)	0.500 (12.7)	0.500 (12.7) dia.	1	500	M61-R6410
	Fig. 24	PTL-117-499	B-499	White	0.900 (22.9)	0.500 (12.7)	0.500 (12.7) dia.	1	500	M61-R6210
	Fig. 24	M61-98-492	B-492	White	1.000 (25.4)	0.375 (9.5)	0.375 (9.5) dia.	1	500	M61-R6410
	Fig. 24	M61-98-494-BL	B-494	Blue / white	1.000 (25.4)	0.375 (9.5)	0.375 (9.5) dia.	1	500	M61-R6410
	Fig. 24	M61-98-494-GN	B-494	Green / white	1.000 (25.4)	0.375 (9.5)	0.375 (9.5) dia.	1	500	M61-R6410
	Fig. 24	M61-98-494-OR	B-494	Orange / white	1.000 (25.4)	0.375 (9.5)	0.375 (9.5) dia.	1	500	M61-R6410
	Fig. 24	M61-98-494-PK	B-494	Pink / white	1.000 (25.4)	0.375 (9.5)	0.375 (9.5) dia.	1	500	M61-R6410
	Fig. 24	M61-98-494-RD	B-494	Red / white	1.000 (25.4)	0.375 (9.5)	0.375 (9.5) dia.	1	500	M61-R6410
	Fig. 24	M61-98-494-YL	B-494	Yellow / white	1.000 (25.4)	0.375 (9.5)	0.375 (9.5) dia.	1	500	M61-R6410
	Fig. 24	PTL-98-499	B-499	White	1.000 (25.4)	0.375 (9.5)	0.375 (9.5) dia.	1	500	M61-R6210
	Fig. 24	M61-98-7425	B-7425	White	1.000 (25.4)	0.375 (9.5)	0.375 (9.5) dia.	1	500	M61-R4310

Cryo tube labels continue on next page.

*Tube sizes may vary by brand.

Cryogenic labels

Cryo tube labels (continued)



Fits tube size*	Figure	Catalog #	Material	Color	Label width A inch (mm)	Label length B inch (mm)	Secondary label C inch (mm)	Labels per row	Labels per roll	Recommended ribbon
1.2 mL	Fig. 2	M61-17-492	B-492	White	1.000 (25.4)	0.500 (12.7)	-	1	250	M61-R6410
	Fig. 2	M61-17-494-BL	B-494	Blue / white	1.000 (25.4)	0.500 (12.7)	-	1	500	M61-R6410
	Fig. 2	M61-17-494-GN	B-494	Green / white	1.000 (25.4)	0.500 (12.7)	-	1	500	M61-R6410
	Fig. 2	M61-17-494-OR	B-494	Orange / white	1.000 (25.4)	0.500 (12.7)	-	1	500	M61-R6410
	Fig. 2	M61-17-494-PK	B-494	Pink / white	1.000 (25.4)	0.500 (12.7)	-	1	500	M61-R6410
	Fig. 2	M61-17-494-RD	B-494	Red / white	1.000 (25.4)	0.500 (12.7)	-	1	500	M61-R6410
	Fig. 2	M61-17-494-YL	B-494	Yellow / white	1.000 (25.4)	0.500 (12.7)	-	1	500	M61-R6410
	Fig. 2	M61-17-7425	B-7425	White	1.000 (25.4)	0.500 (12.7)	-	1	500	M61-R4310
	Fig. 24	M61-120-492	B-492	White	1.000 (25.4)	0.500 (12.7)	0.375 (9.5) dia.	1	250	M61-R6410
	Fig. 24	M61-120-499	B-499	White	1.000 (25.4)	0.500 (12.7)	0.375 (9.5) dia.	1	250	M61-R6210
Fig. 2	M61-199-7425	B-7425	White	1.300 (33.0)	0.600 (15.2)	-	1	250	M61-R4310	
1.2 mL top only	Fig. 8	PTL-82-499	B-499	White	0.375 (9.5) dia.	-	-	1	500	M61-R6210
	Fig. 8	M61-82-492	B-492	White	0.375 (9.5) dia.	-	-	1	500	M61-R6410
	Fig. 8	PTL-82-499-YL	B-499	Yellow	0.375 (9.5) dia.	-	-	1	500	M61-R6210
	Fig. 8	M61-82-499-OR	B-499	Orange	0.375 (9.5) dia.	-	-	1	500	M61-R6210
	Fig. 8	M61-82-499-RD	B-499	Red	0.375 (9.5) dia.	-	-	1	500	M61-R6210
	Fig. 8	PTL-83-499	B-499	White	0.500 (12.7) dia.	-	-	1	500	M61-R6210
	Fig. 8	M61-83-492	B-492	White	0.500 (12.7) dia.	-	-	1	500	M61-R6410
	Fig. 8	PTL-83-499-YL	B-499	Yellow	0.500 (12.7) dia.	-	-	1	500	M61-R6210
	Fig. 8	M61-83-499-OR	B-499	Orange	0.500 (12.7) dia.	-	-	1	500	M61-R6210
	Fig. 8	M61-83-499-RD	B-499	Red	0.500 (12.7) dia.	-	-	1	500	M61-R6210
2 mL	Fig. 82	M61-156-492	B-492	White	0.900 (22.9)	1.000 (25.4)	-	1	250	M61-R6410
	Fig. 2	M61-18-7425	B-7425	White	1.000 (25.4)	0.750 (19.1)	-	1	250	M61-R4310
	Fig. 2	M61-126-490	B-490	White	1.800 (45.7)	0.600 (15.2)	-	1	250	M61-R4310
5 mL	Fig. 2	M61-19-492	B-492	White	1.000 (25.4)	1.000 (25.4)	-	1	250	M61-R6410
	Fig. 2	M61-19-494-BL	B-494	Blue / white	1.000 (25.4)	1.000 (25.4)	-	1	250	M61-R6410
	Fig. 2	M61-19-494-GN	B-494	Green / white	1.000 (25.4)	1.000 (25.4)	-	1	250	M61-R6410
	Fig. 2	M61-19-494-OR	B-494	Orange / white	1.000 (25.4)	1.000 (25.4)	-	1	250	M61-R6410
	Fig. 2	M61-19-494-PK	B-494	Pink / white	1.000 (25.4)	1.000 (25.4)	-	1	250	M61-R6410
	Fig. 2	M61-19-494-RD	B-494	Red / white	1.000 (25.4)	1.000 (25.4)	-	1	250	M61-R6410
	Fig. 2	M61-19-494-YL	B-494	Yellow / white	1.000 (25.4)	1.000 (25.4)	-	1	250	M61-R6410
	Fig. 82	M61-203-490	B-490	White	1.000 (25.4)	1.750 (44.5)	-	1	250	M61-R4310
Fig. 33	M61C-1000-461	B-461	White	1.000 (25.4)	50 ft. (15.2 m)	-	1	1 roll (cont.)	M61-R4310	
15 mL	Fig. 2	M61-31-7425	B-7425	White	1.500 (38.1)	1.000 (25.4)	-	1	250	M61-R4310
	Fig. 83	M61-205-492	B-492	White	0.600 (15.2)	1.625 (41.3)	0.375 (9.5) dia.	1	250	M61-R6410
	Fig. 82	M61-203-492	B-492	White	1.000 (25.4)	1.750 (44.5)	-	1	250	M61-R6410
	Fig. 82	PTL-113-490	B-490	White	1.500 (38.1)	2.200 (55.9)	-	1	100	M61-R4310
	BradyPrinter i3300 and BBP®33 Printer label rolls									
1.2 mL	Fig. 1	B33-68-499	B-499	White	0.500 (12.7)	0.750 (19.1)	-	4	1,500	B30-R6000
	Fig. 2	B33-257-7425	B-7425	White	0.825 (21.0)	0.250 (6.4)	-	1	1,500	B30-R4300
	Fig. 2	B33-125-499	B-499	White	0.900 (22.9)	0.500 (12.7)	-	1	1,500	B30-R6000
	Fig. 24	B33-126-499	B-499	White	0.900 (22.9)	0.500 (12.7)	0.500 (12.7) dia.	2	1,500	B30-R4300
	Fig. 25	B33-152-492	B-492	White	1.000 (25.4)	0.375 (9.5)	-	2	2,500	B30-R6400
	Fig. 25	B33-152-494-BL	B-494	Blue / white	1.000 (25.4)	0.375 (9.5)	-	2	1,500	B30-R6400
	Fig. 25	B33-152-494-GN	B-494	Green / white	1.000 (25.4)	0.375 (9.5)	-	2	1,500	B30-R6400
	Fig. 25	B33-152-494-OR	B-494	Orange / white	1.000 (25.4)	0.375 (9.5)	-	2	1,500	B30-R6400
	Fig. 25	B33-152-494-PK	B-494	Pink / white	1.000 (25.4)	0.375 (9.5)	-	2	1,500	B30-R6400
	Fig. 25	B33-152-494-RD	B-494	Red / white	1.000 (25.4)	0.375 (9.5)	-	2	1,500	B30-R6400
	Fig. 25	B33-152-494-YL	B-494	Yellow / white	1.000 (25.4)	0.375 (9.5)	-	2	1,500	B30-R6400
	Fig. 25	B33-152-7425	B-7425	White	1.000 (25.4)	0.375 (9.5)	-	2	2,000	B30-R4300
	Fig. 24	B33-249-7425	B-7425	White	1.000 (25.4)	0.375 (9.5)	0.375 (9.5) dia.	1	1,500	B30-R4300
	Fig. 2	B33-59-492	B-492	White	1.000 (25.4)	0.500 (12.7)	-	1	2,500	B30-R6400
	Fig. 2	B33-59-7425	B-7425	White	1.000 (25.4)	0.500 (12.7)	-	1	2,000	B30-R4300

Cryo tube labels continue on next page.

*Tube sizes may vary by brand.

Cryo tube labels (continued)

Fits tube size*	Figure	Catalog #	Material	Color	Label width A inch (mm)	Label length B inch (mm)	Secondary label C inch (mm)	Labels per row	Labels per roll	Recommended ribbon
1.2 mL	Fig. 24	B33-181-492	B-492	White	1.000 (25.4)	0.500 (12.7)	0.375 (9.5) dia.	1	1,500	B30-R6400
	Fig. 24	B33-250-494-BL	B-494	Blue / white	1.000 (25.4)	0.500 (12.7)	0.375 (9.5) dia.	1	1,500	B30-R6400
	Fig. 24	B33-250-494-GN	B-494	Green / white	1.000 (25.4)	0.500 (12.7)	0.375 (9.5) dia.	1	1,500	B30-R6400
	Fig. 24	B33-250-494-OR	B-494	Orange / white	1.000 (25.4)	0.500 (12.7)	0.375 (9.5) dia.	1	1,500	B30-R6400
	Fig. 24	B33-250-494-PK	B-494	Pink / white	1.000 (25.4)	0.500 (12.7)	0.375 (9.5) dia.	1	1,500	B30-R6400
	Fig. 24	B33-250-494-RD	B-494	Red / white	1.000 (25.4)	0.500 (12.7)	0.375 (9.5) dia.	1	1,500	B30-R6400
	Fig. 24	B33-250-494-YL	B-494	Yellow / white	1.000 (25.4)	0.500 (12.7)	0.375 (9.5) dia.	1	1,500	B30-R6400
	Fig. 24	B33-182-492	B-492	White	1.000 (25.4)	0.500 (12.7)	0.440 (11.2) dia.	1	1,500	B30-R6400
1.2 mL top only	Fig. 24	B33-251-7425	B-7425	White	1.000 (25.4)	0.500 (12.7)	0.440 (11.2) dia.	1	1,500	B30-R4300
	Fig. 2	B33-256-7425	B-7425	White	1.300 (33.0)	0.600 (15.2)	-	1	1,500	B30-R4300
	Fig. 5a	B33-144-490	B-490	White	0.500 (12.7) dia.	-	-	3	2,500	B30-R4300
2 mL	Fig. 5a	B33-144-492	B-492	White	0.500 (12.7) dia.	-	-	3	1,500	B30-R6400
	Fig. 5a	B33-144-499	B-499	White	0.500 (12.7) dia.	-	-	3	1,500	B30-R6000
	Fig. 2	B33-180-492	B-492	White	0.900 (22.9)	1.000 (25.4)	-	1	1,500	B30-R6400
5 mL	Fig. 25	B33-88-7425	B-7425	White	1.000 (25.4)	0.750 (19.1)	-	2	2,000	B30-R4300
	Fig. 2	B33-179-492	B-492	White	1.000 (25.4)	1.000 (25.4)	-	1	1,900	B30-R6400
	Fig. 2	B33-179-494-BL	B-494	Blue / white	1.000 (25.4)	1.000 (25.4)	-	1	1,500	B30-R6400
	Fig. 2	B33-179-494-GN	B-494	Green / white	1.000 (25.4)	1.000 (25.4)	-	1	1,500	B30-R6400
	Fig. 2	B33-179-494-OR	B-494	Orange / white	1.000 (25.4)	1.000 (25.4)	-	1	1,500	B30-R6400
	Fig. 2	B33-179-494-PK	B-494	Pink / white	1.000 (25.4)	1.000 (25.4)	-	1	1,500	B30-R6400
	Fig. 2	B33-179-494-RD	B-494	Red / white	1.000 (25.4)	1.000 (25.4)	-	1	1,500	B30-R6400
	Fig. 2	B33-179-494-YL	B-494	Yellow / white	1.000 (25.4)	1.000 (25.4)	-	1	1,500	B30-R6400
	Fig. 2	B33-179-7425	B-7425	White	1.000 (25.4)	1.000 (25.4)	-	1	1,500	B30-R4300
	Fig. 2	B33-156-490	B-490	White	1.900 (48.3)	1.000 (25.4)	-	1	1,500	B30-R4300
15 mL	Fig. 25	B33-37-7425	B-7425	White	1.500 (38.1)	0.500 (12.7)	-	2	1,500	B30-R4300
	Fig. 2	B33-155-492	B-492	White	1.625 (41.3)	0.600 (15.2)	-	1	1,500	B30-R6400
	Fig. 2	B33-157-490	B-490	White	2.200 (55.9)	1.500 (38.1)	-	1	1,000	B30-R4300
BradyJet J2000 Printer label rolls										
1.2 mL	Fig. 2	J20-257-7425J	B-7425J	White	0.825 (21.0)	0.250 (6.4)	-	1	1,500	J20-CMY
	Fig. 25	J20-152-7425J	B-7425J	White	1.000 (25.4)	0.375 (9.5)	-	2	2,000	J20-CMY
	Fig. 24	J20-249-7425J	B-7425J	White	1.000 (25.4)	0.375 (9.5)	0.375 (9.5) dia.	2	1,500	J20-CMY
	Fig. 2	J20-59-7425J	B-7425J	White	1.000 (25.4)	0.500 (12.7)	-	1	2,000	J20-CMY
	Fig. 24	J20-251-7425J	B-7425J	White	1.000 (25.4)	0.500 (12.7)	0.440 (11.2) dia.	2	1,500	J20-CMY
	Fig. 2	J20-256-7425J	B-7425J	White	1.300 (33.0)	0.600 (15.2)	-	1	1,500	J20-CMY
2 mL	Fig. 25	J20-88-7425J	B-7425J	White	1.000 (25.4)	0.750 (19.1)	-	2	2,000	J20-CMY
5 mL	Fig. 2	J20-179-7425J	B-7425J	White	1.000 (25.4)	1.000 (25.4)	-	1	1,500	J20-CMY
15 mL	Fig. 25	J20-37-7425J	B-7425J	White	1.500 (38.1)	0.500 (12.7)	-	2	1,500	J20-CMY
BBP®12 Printer and other industry printers using 1-inch small core label rolls										
1.2 mL	Fig. 2	THT-257-7425-2-SC	B-7425	White	0.825 (21.0)	0.250 (6.4)	-	1	2,000	R6413
	Fig. 2	THT-125-499-1.5-SC	B-499	White	0.900 (22.9)	0.500 (12.7)	-	1	1,500	R4913
	Fig. 25	THT-152-492-2.5-SC	B-492	White	1.000 (25.4)	0.375 (9.5)	-	2	2,500	R6413
	Fig. 25	THT-152-499-1.5-SC	B-499	White	1.000 (25.4)	0.375 (9.5)	-	2	1,500	R4913
	Fig. 25	THT-152-7425-2-SC	B-7425	White	1.000 (25.4)	0.375 (9.5)	-	2	2,000	R6413
	Fig. 24	THT-249-492-1.5-SC	B-492	White	1.000 (25.4)	0.375 (9.5)	0.375 (9.5) dia.	1	1,500	R6413
	Fig. 24	THT-249-7425-2-SC	B-7425	White	1.000 (25.4)	0.375 (9.5)	0.375 (9.5) dia.	1	2,000	R6413
	Fig. 24	THT-249-499-1.5-SC	B-499	White	1.000 (25.4)	0.375 (9.5)	0.375 (9.5) dia.	1	1,500	R4913
	Fig. 2	THT-59-492-2.5-SC	B-492	White	1.000 (25.4)	0.500 (12.7)	-	1	2,500	R6413
	Fig. 2	THT-59-7425-2-SC	B-7425	White	1.000 (25.4)	0.500 (12.7)	-	1	2,000	R6413
	Fig. 24	THT-250-492-1.5-SC	B-492	White	1.000 (25.4)	0.500 (12.7)	0.375 (9.5) dia.	1	1,500	R6413
	Fig. 24	THT-251-492-1.5-SC	B-492	White	1.000 (25.4)	0.500 (12.7)	0.440 (11.2) dia.	1	1,500	R6413
	Fig. 24	THT-251-7425-2-SC	B-7425	White	1.000 (25.4)	0.500 (12.7)	0.440 (11.2) dia.	1	2,000	R6413
	Fig. 2	THT-154-490-1.5-SC	B-490	White	1.250 (31.8)	0.375 (9.5)	-	1	1,500	R4313
	Fig. 2	THT-256-7425-2-SC	B-7425	White	1.300 (33.0)	0.600 (15.2)	-	1	2,000	R6413
1.2 mL top only	Fig. 5a	THT-144-492-1.5-SC	B-492	White	0.500 (12.7) dia.	-	-	3	1,500	R6413
	Fig. 5a	THT-144-499-1.5-SC	B-499	White	0.500 (12.7) dia.	-	-	3	1,500	R4913
2 mL	Fig. 25	THT-88-492-2.5-SC	B-492	White	1.000 (25.4)	0.750 (19.1)	-	2	2,500	R6413
	Fig. 25	THT-88-7425-2-SC	B-7425	White	1.000 (25.4)	0.750 (19.1)	-	2	2,000	R6413
5 mL	Fig. 2	THT-179-492-2.5-SC	B-492	White	1.000 (25.4)	1.000 (25.4)	-	1	2,500	R6413
	Fig. 2	THT-179-7425-2-SC	B-7425	White	1.000 (25.4)	1.000 (25.4)	-	1	2,000	R6413
	Fig. 2	THT-156-490-1.5-SC	B-490	White	1.900 (48.3)	1.000 (25.4)	-	1	1,500	R4313
15 mL	Fig. 25	THT-37-7425-2-SC	B-7425	White	1.500 (38.1)	0.500 (12.7)	-	2	2,000	R6413
	Fig. 2	THT-157-490-1.5-SC	B-490	White	2.200 (55.9)	1.500 (38.1)	-	1	1,500	R4313

Cryo tube labels continue on next page.

*Tube sizes may vary by brand.

Cryogenic labels

Cryo tube labels (continued)




Figure 2




Figure 5a




Figure 24




Figure 25

Fits tube size*	Figure	Catalog #	Material	Color	Label width A inch (mm)	Label length B inch (mm)	Secondary label C inch (mm)	Labels per row	Labels per roll	Recommended ribbon
BradyPrinter i5100, i7100, Brady IP Printers and other industry printers using 3-inch core label rolls										
1.2 mL	Fig. 2	THT-257-7425-2	B-7425	White	0.825 (21.0)	0.250 (6.4)	-	1	2,000	R4302 / IP-R4302
	Fig. 2	THT-125-499-3	B-499	White	0.900 (22.9)	0.500 (12.7)	-	1	3,000	R6000 / IP-R6000
	Fig. 25	THT-152-492-3	B-492	White	1.000 (25.4)	0.375 (9.5)	-	2	3,000	R6402 / IP-R6402
	Fig. 25	THT-152-499-3	B-499	White	1.000 (25.4)	0.375 (9.5)	-	2	3,000	R6000 / IP-R6000
	Fig. 25	THT-152-7425-2	B-7425	White	1.000 (25.4)	0.375 (9.5)	-	2	2,000	R4302 / IP-R4302
	Fig. 24	THT-249-7425-2	B-7425	White	1.000 (25.4)	0.375 (9.5)	0.375 (9.5) dia.	1	2,000	R4302 / IP-R4302
1.2 mL	Fig. 25	THT-152-494-BL	B-494	Blue / white	1.000 (25.4)	0.375 (9.5)	-	2	3,000	R6402 / IP-R6402
	Fig. 25	THT-152-494-GN	B-494	Green / white	1.000 (25.4)	0.375 (9.5)	-	2	3,000	R6402 / IP-R6402
	Fig. 25	THT-152-494-OR	B-494	Orange / white	1.000 (25.4)	0.375 (9.5)	-	2	3,000	R6402 / IP-R6402
	Fig. 25	THT-152-494-PK	B-494	Pink / white	1.000 (25.4)	0.375 (9.5)	-	2	3,000	R6402 / IP-R6402
	Fig. 25	THT-152-494-RD	B-494	Red / white	1.000 (25.4)	0.375 (9.5)	-	2	3,000	R6402 / IP-R6402
	Fig. 25	THT-152-494-YL	B-494	Yellow / white	1.000 (25.4)	0.375 (9.5)	-	2	3,000	R6402 / IP-R6402
	Fig. 2	THT-59-492-10	B-492	White	1.000 (25.4)	0.500 (12.7)	-	1	2 rolls of 5,000	R6400 / IP-R6400
	Fig. 2	THT-59-7425-2	B-7425	White	1.000 (25.4)	0.500 (12.7)	-	1	2,000	R4302 / IP-R4302
	Fig. 24	THT-250-494-BL	B-494	Blue / white	1.000 (25.4)	0.500 (12.7)	0.375 (9.5) dia.	1	3,000	R6402 / IP-R6402
	Fig. 24	THT-250-494-GN	B-494	Green / white	1.000 (25.4)	0.500 (12.7)	0.375 (9.5) dia.	1	3,000	R6402 / IP-R6402
	Fig. 24	THT-250-494-OR	B-494	Orange / white	1.000 (25.4)	0.500 (12.7)	0.375 (9.5) dia.	1	3,000	R6402 / IP-R6402
	Fig. 24	THT-250-494-PK	B-494	Pink / white	1.000 (25.4)	0.500 (12.7)	0.375 (9.5) dia.	1	3,000	R6402 / IP-R6402
	Fig. 24	THT-250-494-RD	B-494	Red / white	1.000 (25.4)	0.500 (12.7)	0.375 (9.5) dia.	1	3,000	R6402 / IP-R6402
	Fig. 24	THT-250-494-YL	B-494	Yellow / white	1.000 (25.4)	0.500 (12.7)	0.375 (9.5) dia.	1	3,000	R6402 / IP-R6402
	Fig. 24	THT-251-7425-2	B-7425	White	1.000 (25.4)	0.500 (12.7)	0.440 (11.2) dia.	1	2,000	R4302 / IP-R4302
Fig. 2	THT-154-490-3	B-490	White	1.250 (31.8)	0.375 (9.5)	-	1	3,000	R4302 / IP-R4302	
Fig. 2	THT-256-7425-2	B-7425	White	1.300 (33.0)	0.600 (15.2)	-	1	2,000	R4302 / IP-R4302	
1.2 mL top only	Fig. 5a	THT-144-492-3	B-492	White	0.500 (12.7) dia.	-	-	3	3,000	R6402 / IP-R6402
	Fig. 5a	THT-144-499-3	B-499	White	0.500 (12.7) dia.	-	-	3	3,000	R6000 / IP-R6000
2 mL	Fig. 2	THT-180-492-3	B-492	White	0.900 (22.9)	1.000 (25.4)	-	1	3,000	R6400 / IP-R6400
	Fig. 25	THT-88-492-5	B-492	White	1.000 (25.4)	0.750 (19.1)	-	2	5,000	R6402 / IP-R6402
	Fig. 25	THT-88-7425-2	B-7425	White	1.000 (25.4)	0.750 (19.1)	-	2	2,000	R4302 / IP-R4302
5 mL	Fig. 2	THT-179-492-3	B-492	White	1.000 (25.4)	1.000 (25.4)	-	1	3,000	R6400 / IP-R6400
	Fig. 2	THT-179-494-BL	B-494	Blue / white	1.000 (25.4)	1.000 (25.4)	-	1	3,000	R6402 / IP-R6402
	Fig. 2	THT-179-494-GN	B-494	Green / white	1.000 (25.4)	1.000 (25.4)	-	1	3,000	R6402 / IP-R6402
	Fig. 2	THT-179-494-OR	B-494	Orange / white	1.000 (25.4)	1.000 (25.4)	-	1	3,000	R6402 / IP-R6402
	Fig. 2	THT-179-494-PK	B-494	Pink / white	1.000 (25.4)	1.000 (25.4)	-	1	3,000	R6402 / IP-R6402
	Fig. 2	THT-179-494-RD	B-494	Red / white	1.000 (25.4)	1.000 (25.4)	-	1	3,000	R6402 / IP-R6402
	Fig. 2	THT-179-494-YL	B-494	Yellow / white	1.000 (25.4)	1.000 (25.4)	-	1	3,000	R6402 / IP-R6402
	Fig. 2	THT-179-7425-2	B-7425	White	1.000 (25.4)	1.000 (25.4)	-	1	2,000	R4302 / IP-R4302
	Fig. 2	THT-156-490-3	B-490	White	1.900 (48.3)	1.000 (25.4)	-	1	3,000	R4302 / IP-R4302
15 mL	Fig. 25	THT-37-7425-2	B-7425	White	1.500 (38.1)	0.500 (12.7)	-	2	2,000	R4307 / IP-R4307
	Fig. 2	THT-155-492-3	B-492	White	1.625 (41.3)	0.600 (15.2)	-	1	3,000	R6402 / IP-R6402
	Fig. 2	THT-157-490-3	B-490	White	2.200 (55.9)	1.500 (38.1)	-	1	3,000	R4300 / IP-R4300

Self-laminating cryo tube labels

Cryo tube labels are designed to remain legible and stick to the tubes during processing, even when immersed in liquid or vapor-phase nitrogen and ultra-low temperature freezer storage.

- B-461 – Ultra-thin polyester with clear overlaminates tail that protects print and allows liquid levels to be visible




Figure 14




Figure 14b




Figure 14c




Figure 15



Fits tube size*	Figure	Catalog #	Material	Color	Label width A inch (mm)	Label length B inch (mm)	Print-on area C inch (mm)	Labels per row	Labels per roll	Recommended ribbon
BMP®51 Label Printer cartridges										
1.2 mL	Fig. 15	M-72-461	B-461	Black on white	1.750 (44.5)	0.500 (12.7)	0.750 (19.0)	1	180	In cartridge
	Fig. 15	M-126-461	B-461	Black on white	1.800 (45.7)	0.600 (15.2)	1.000 (25.4)	1	180	In cartridge
	Fig. 15	M-74-461	B-461	Black on white	2.200 (55.9)	0.500 (12.7)	0.750 (19.0)	1	140	In cartridge
	Fig. 15	M-86-461	B-461	Black on white	2.625 (66.7)	0.600 (15.2)	1.000 (25.4)	1	120	In cartridge
5 mL	Fig. 15	M-125-461	B-461	Black on white	1.750 (44.5)	1.000 (25.4)	1.000 (25.4)	1	180	In cartridge
	Fig. 15	M-75-461	B-461	Black on white	2.625 (66.7)	1.000 (25.4)	1.000 (25.4)	1	120	In cartridge
BradyPrinter M611 label rolls										
1.2 mL	Fig. 14	M61-201-461	B-461	White	0.600 (15.2)	1.800 (45.7)	1.000 (25.4)	1	100	M61-R4310
	Fig. 14	PTL-86-461	B-461	White	0.600 (15.2)	2.625 (66.7)	1.000 (25.4)	1	100	M61-R4310
	Fig. 15	PTL-108-461	B-461	White	1.500 (38.1)	0.187 (4.75)	0.700 (17.8)	1	750	M61-R4310
	Fig. 15	PTL-73-461	B-461	White	1.800 (45.7)	0.375 (9.5)	0.750 (19.1)	1	500	M61-R4310
	Fig. 15	PTL-72-461	B-461	White	1.800 (45.7)	0.500 (12.7)	0.750 (19.1)	1	500	M61-R4310
	Fig. 14	PTL-74-461	B-461	White	0.500 (12.7)	2.200 (55.9)	0.500 (12.7)	1	100	M61-R4310
5 mL	Fig. 14	PTL-79-461	B-461	White	0.750 (19.1)	3.750 (95.3)	0.750 (19.1)	1	100	M61-R4310
	Fig. 14	M61-203-461	B-461	White	1.000 (25.4)	1.750 (44.5)	1.000 (25.4)	1	100	M61-R4310
	Fig. 14	PTL-75-461-YL	B-461	Yellow	1.000 (25.4)	2.625 (66.7)	1.000 (25.4)	1	100	M61-R4310
	Fig. 14	M61-157-461	B-461	White	1.500 (38.1)	2.500 (63.5)	1.500 (38.1)	1	100	M61-R4310
	Fig. 14	PTL-75-461	B-461	White	1.000 (25.4)	2.625 (66.7)	1.000 (25.4)	1	100	M61-R4310
15 mL	Fig. 14	PTL-77-461	B-461	White	1.500 (38.1)	3.500 (88.9)	1.250 (31.8)	1	250	M61-R4310
BradyPrinter i3300 and BBP®33 Printer label rolls										
1.2 mL	Fig. 15	B33-123-461	B-461	White	1.800 (45.7)	0.375 (9.5)	0.750 (19.1)	1	1,500	B30-R4300
	Fig. 15	B33-122-461	B-461	White	1.800 (45.7)	0.500 (12.7)	0.750 (19.1)	1	1,500	B30-R4300
	Fig. 14b	B33-127-461	B-461	White	2.200 (55.9)	0.500 (12.7)	0.750 (19.1)	3	1,500	B30-R4300
5 mL	Fig. 14c	B33-276-461	B-461	White	2.625 (66.7)	1.000 (25.4)	1.000 (25.4)	2	1,500	B30-R4300
15 mL	Fig. 14c	B33-139-461	B-461	White	1.000 (25.4)	5.250 (133.4)	1.500 (38.1)	2	650	B30-R4300
BBP®12 Printer and other industry printers using 1-inch small core label rolls										
1.2 mL	Fig. 15	THT-122-461-1.5-SC	B-461	White	0.500 (12.7)	1.800 (45.7)	0.750 (19.1)	1	1,500	R4313
	Fig. 15	THT-127-461-1.5-SC	B-461	White	0.500 (12.7)	2.200 (55.9)	0.750 (19.1)	1	1,500	R4313
	Fig. 15	THT-153-461-1.5-SC	B-461	White	0.600 (15.2)	2.625 (66.7)	1.000 (25.4)	1	1,500	R4313
	Fig. 15	THT-123-461-1.5-SC	B-461	White	0.750 (19.1)	1.800 (45.7)	0.750 (19.1)	1	1,500	R4313
BradyPrinter i5100, i7100, Brady IP Printers and other industry printers using 3-inch core label rolls										
1.2 mL	Fig. 15	THT-123-461-3	B-461	White	0.375 (9.5)	1.800 (45.7)	0.750 (19.1)	1	3,000	R4302 / IP-R4302
	Fig. 15	THT-122-461-3	B-461	White	0.500 (12.7)	1.800 (45.7)	0.750 (19.1)	1	3,000	R4302 / IP-R4302
	Fig. 15	THT-127-461-3	B-461	White	0.500 (12.7)	2.200 (55.9)	0.750 (19.1)	1	3,000	R4300 / IP-R4300
	Fig. 15	THT-153-461-3	B-461	White	0.600 (15.2)	2.625 (66.7)	1.000 (25.4)	1	3,000	R4300 / IP-R4300
5 mL	Fig. 15	THT-183-461-1	B-461	White	1.000 (25.4)	1.750 (44.5)	1.000 (25.4)	1	1,000	R4302 / IP-R4302
	Fig. 15	THT-131-461-3	B-461	White	1.000 (25.4)	2.625 (66.7)	1.000 (25.4)	1	3,000	R4300 / IP-R4300

*Tube sizes may vary by brand.

Sample collection and processing labels

Blood tube labels

Brady's labels keep text and barcodes legible without falling off in the lab while also not sticking to gloves when you're applying the label.

- B-424 – Paper label with good contrast and smear resistance
- B-813 – Direct thermal paper labels with good adhesion to corrugated surfaces, glass and many plastics
- B-7425 / B-7425J – Economical white polypropylene label with moderate chemical resistance that is ideal for labeling dry, room-temperature specimen containers




Figure	Catalog #	Material	Color	Label width A inch (mm)	Label length B inch (mm)	Labels per row	Labels per roll	Recommended ribbon
BradyPrinter i3300 and BBP®33 Printer label rolls								
Fig. 2	B33-17-424	B-424	White	2.000 (50.8)	1.000 (25.4)	1	1,500	B30-R4300
Fig. 2	B33-137-7425	B-7425	White	2.000 (50.8)	1.000 (25.4)	1	1,500	B30-R4300
BradyJet J2000 Printer label rolls								
Fig. 2	J20-137-7425J	B-7425J	White	2.000 (50.8)	1.000 (25.4)	1	1,500	J20-CMY
BBP®12 Printer and other industry printers using 1-inch small core label rolls								
Fig. 2	THT-17-424-1.5-SC	B-424	White	2.000 (50.8)	1.000 (25.4)	1	1,500	R4313
Fig. 2	THT-137-7425-2-SC	B-7425	White	2.000 (50.8)	1.000 (25.4)	1	2,000	R6413
BradyPrinter i5100, i7100, Brady IP Printers and other industry printers using 3-inch core label rolls								
Fig. 2	THT-17-424-3	B-424	White	2.000 (50.8)	1.000 (25.4)	1	3,000	R4300 / IP-R4300
Fig. 2	DT-17-813-3	B-813	White	2.000 (50.8)	1.000 (25.4)	1	3,000	(direct thermal)
Fig. 2	THT-137-7425-2	B-7425	White	2.000 (50.8)	1.000 (25.4)	1	2,000	R4302 / IP-R4302

Specimen cup labels

Brady labels easily adhere to collection cups providing clear identification of samples.

- B-424 – Paper label with good contrast and smear resistance
- B-813 – Direct thermal paper labels with good adhesion to corrugated surfaces, glass and many plastics
- B-7425 / B-7425J – Economical white polypropylene label with moderate chemical resistance that is ideal for labeling dry, room-temperature specimen containers




Figure	Catalog #	Material	Color	Label width A inch (mm)	Label length B inch (mm)	Labels per row	Labels per roll	Recommended ribbon
BradyPrinter i3300 and BBP®33 Printer label rolls								
Fig. 2	B33-17-424	B-424	White	2.000 (50.8)	1.000 (25.4)	1	1,500	B30-R4300
Fig. 2	B33-137-7425	B-7425	White	2.000 (50.8)	1.000 (25.4)	1	1,500	B30-R4300
Fig. 2	B33-18-424	B-424	White	3.000 (76.2)	1.000 (25.4)	1	1,500	B30-R4300
BradyJet J2000 Printer label rolls								
Fig. 2	J20-137-7425J	B-7425J	White	2.000 (50.8)	1.000 (25.4)	1	1,500	J20-CMY
BBP®12 Printer and other industry printers using 1-inch small core label rolls								
Fig. 2	THT-137-7425-2-SC	B-7425	White	2.000 (50.8)	1.000 (25.4)	1	2,000	R6413
Fig. 2	THT-18-424-1.5-SC	B-424	White	3.000 (76.2)	1.000 (25.4)	1	1,500	R4313
BradyPrinter i5100, i7100, Brady IP Printers and other industry printers using 3-inch core label rolls								
Fig. 2	DT-17-813-3	B-813	White	2.000 (50.8)	1.000 (25.4)	1	3,000	(direct thermal)
Fig. 2	THT-137-7425-2	B-7425	White	2.000 (50.8)	1.000 (25.4)	1	2,000	R4302 / IP-R4302
Fig. 2	THT-18-424-3	B-424	White	3.000 (76.2)	1.000 (25.4)	1	3,000	R4307 / IP-R4307
Fig. 2	DT-18-813-3	B-813	White	3.000 (76.2)	1.000 (25.4)	1	3,000	(direct thermal)

Sample collection and processing labels

Tamper-resistant seals

Protect samples with irreversible, visible evidence of tampering. These labels tear if the container has been opened or tampered with in any way.

- B-351 – White vinyl designed to fracture easily to prevent one-piece removal. Use caution when removing from liner as material is fragile




Figure 1




Figure 2




Figure 25



Figure	Catalog #	Material	Color	Label width A inch (mm)	Label length B inch (mm)	Labels per row	Labels per roll	Recommended ribbon
BradyPrinter i5100, i7100, Brady IP Printers and other industry printers using 3-inch core label rolls								
Fig. 1	THT-5-351-10	B-351	White	1.000 (25.4)	0.500 (12.7)	3	10,000	R6207 / IP-R6207
Fig. 25	THT-37-351-10	B-351	White	1.500 (38.1)	0.500 (12.7)	2	10,000	R6200 / IP-R6200
Fig. 2	THT-53-351-3	B-351	White	2.000 (50.8)	0.500 (12.7)	1	3,000	R6200 / IP-R6200
Fig. 2	THT-17-351-3	B-351	White	2.000 (50.8)	1.000 (25.4)	1	3,000	R6200 / IP-R6200
Fig. 2	THT-18-351-3	B-351	White	3.000 (76.2)	1.000 (25.4)	1	3,000	R6200 / IP-R6200
Fig. 2	THT-19-351-1	B-351	White	3.000 (76.2)	2.000 (50.8)	1	1,000	R6200 / IP-R6200

Sample bag labels

Brady's general identification labels easily adhere to plastic bags used in the lab. These labels offer high performance and durability to withstand freezers, coolers and ice packs.

- B-499 – White nylon cloth with permanent adhesive for easy handling and conformability
- B-7425 / B-7425J – Economical white polypropylene label with moderate chemical resistance that is ideal for labeling dry, room-temperature specimen containers




Figure 2




Figure	Catalog #	Material	Color	Label width A inch (mm)	Label length B inch (mm)	Labels per row	Labels per roll	Recommended ribbon
BradyPrinter i3300 and BBP®33 Printer label rolls								
Fig. 2	B33-179-7425	B-7425	White	1.000 (25.4)	1.000 (25.4)	1	1,500	B30-R4300
Fig. 2	B33-258-7425	B-7425	White	1.500 (38.1)	1.750 (44.5)	1	900	B30-R4300
Fig. 2	B33-137-7425	B-7425	White	2.000 (50.8)	1.000 (25.4)	1	1,500	B30-R4300
BradyJet J2000 Printer label rolls								
Fig. 2	J20-179-7425J	B-7425J	White	1.000 (25.4)	1.000 (25.4)	1	1,500	J20-CMY
Fig. 2	J20-258-7425J	B-7425J	White	1.500 (38.1)	1.750 (44.5)	1	900	J20-CMY
Fig. 2	J20-137-7425J	B-7425J	White	2.000 (50.8)	1.000 (25.4)	1	1,500	J20-CMY
BBP®12 Printer and other industry printers using 1-inch small core label rolls								
Fig. 2	THT-179-7425-2-SC	B-7425	White	1.000 (25.4)	1.000 (25.4)	1	2,000	R6413
Fig. 2	THT-136-499-1.5-SC	B-499	White	1.500 (38.1)	0.750 (19.1)	1	1,500	R4913
Fig. 2	THT-258-7425-2-SC	B-7425	White	1.750 (44.5)	1.500 (38.1)	1	2,000	R6413
Fig. 2	THT-137-499-1.5-SC	B-499	White	2.000 (50.8)	1.000 (25.4)	1	1,500	R4913
Fig. 2	THT-137-7425-2-SC	B-7425	White	2.000 (50.8)	1.000 (25.4)	1	2,000	R6413
BradyPrinter i5100, i7100, Brady IP Printers and other industry printers using 3-inch core label rolls								
Fig. 2	THT-179-7425-2	B-7425	White	1.000 (25.4)	1.000 (25.4)	1	2,000	R4302 / IP-R4302
Fig. 2	THT-136-499-3	B-499	White	1.500 (38.1)	0.750 (19.1)	1	3,000	R6000 / IP-R6000
Fig. 2	THT-258-7425-2	B-7425	White	1.750 (44.5)	1.500 (38.1)	1	2,000	R4302 / IP-R4302
Fig. 2	THT-137-499-3	B-499	White	2.000 (50.8)	1.000 (25.4)	1	3,000	R6000 / IP-R6000
Fig. 2	THT-137-7425-2	B-7425	White	2.000 (50.8)	1.000 (25.4)	1	2,000	R4302 / IP-R4302

Sample collection and processing labels

Well and microplate labels

Brady's well and microplate labels can withstand varying temperature processes and have a thin profile, short height and long width to fit the necessary information on the plate.

- B-488 – White polyester with moderate chemical resistance and short-term exposure to xylene, DMSO and ethanol
- B-499 – White nylon cloth with permanent adhesive for easy handling and conformability
- B-7425 / B-7425J – Economical white polypropylene label with moderate chemical resistance that is ideal for labeling dry, room-temperature specimen containers




Figure 2




Figure 25




Figure 82


Figure	Catalog #	Material	Color	Label width A inch (mm)	Label length B inch (mm)	Labels per row	Labels per roll	Recommended ribbon
BMP®51 Label Printer cartridges								
Fig. 2	M-81-488	B-488	Black on white	1.900 (48.3)	0.250 (6.4)	1	160	In cartridge
BradyPrinter M611 label rolls								
Fig. 82	M61-2-7425	B-7425	White	0.250 (6.4)	2.000 (50.8)	1	100	M61-R4310
BradyPrinter i3300 and BBP®33 Printer label rolls								
Fig. 25	B33-152-499	B-499	White	1.000 (25.4)	0.375 (9.5)	2	1,500	B30-R6000
Fig. 2	B33-151-490	B-490	White	2.000 (50.8)	0.250 (6.4)	1	2,500	B30-R4300
Fig. 2	B33-151-499	B-499	White	2.000 (50.8)	0.250 (6.4)	1	1,500	B30-R6000
Fig. 2	B33-151-7425	B-7425	White	2.000 (50.8)	0.250 (6.4)	1	2,000	B30-R4300
Fig. 2	B33-53-7425	B-7425	White	2.000 (50.8)	0.500 (12.7)	1	2,000	B30-R4300
BradyJet J2000 Printer label rolls								
Fig. 2	J20-151-7425J	B-7425J	White	2.000 (50.8)	0.250 (6.4)	1	2,000	J20-CMY
Fig. 2	J20-53-7425J	B-7425J	White	2.000 (50.8)	0.500 (12.7)	1	2,000	J20-CMY
BBP®12 Printer and other industry printers using 1-inch small core label rolls								
Fig. 2	THT-151-499-1.5-SC	B-499	White	2.000 (50.8)	0.250 (6.4)	1	1,500	R4913
Fig. 2	THT-151-7425-2-SC	B-7425	White	2.000 (50.8)	0.250 (6.4)	1	2,000	R6413
Fig. 2	THT-53-7425-2-SC	B-7425	White	2.000 (50.8)	0.500 (12.7)	1	2,000	R6413
BradyPrinter i5100, i7100, Brady IP Printers and other industry printers using 3-inch core label rolls								
Fig. 2	THT-15-488-2.5	B-488	White	2.000 (50.8)	0.250 (6.4)	1	2,500	R4302 / IP-R4302
Fig. 2	THT-151-499-3	B-499	White	2.000 (50.8)	0.250 (6.4)	1	3,000	R6000 / IP-R6000
Fig. 2	THT-151-7425-2	B-7425	White	2.000 (50.8)	0.250 (6.4)	1	2,000	R4302 / IP-R4302
Fig. 2	THT-53-7425-2	B-7425	White	2.000 (50.8)	0.500 (12.7)	1	2,000	R4302 / IP-R4302

Sample collection and processing labels

Centrifuge and conical tube labels

Designed with ultra-thin material that prevents vials from jamming in centrifuges or racks. These labels withstand harsh temperatures and chemicals without falling off, fading or smearing.

- B-461 – All-white, ultra-thin polyester label. Also available in a self-laminating format
- B-427 – Durable all-white vinyl label is good for wrapping and conformability. Must be wrapped upon itself when used in liquid nitrogen
- B-432 – Clear polyester withstands exposure to solvents and exhibits good adhesion to metal and other surfaces
- B-490 – White polyester ideal for flat frozen surfaces. Adheres to frozen tubes and vials when wrapped around itself by at least 0.125 inch
- B-492 – Ultra-thin white polyester with superior chemical and abrasion resistance for labeling frosted, frozen or room-temperature tubes and vials. Must have small overlap if applied to already frozen containers
- B-494 – Glossy polyester with a printed color bar that allows for quick visual differentiation between samples. Withstands liquid nitrogen, a variety of chemicals and more. Label dimensions listed below include the color bar area. Vial tops are solid color
- B-499 – White nylon cloth with permanent adhesive for easy handling and conformability
- B-7425 / B-7425J – Economical white polypropylene label with moderate chemical resistance that is ideal for labeling dry, room-temperature specimen containers. Must be applied to containers prior to being frozen



Fits tube size*	Figure	Catalog #	Material	Color	Label width A inch (mm)	Label length B inch (mm)	Secondary label C inch (mm)	Labels per row	Labels per roll	Recommended ribbon
BMP®21-LAB Label Printer cartridges										
15 mL	Fig. 33b	M21-750-499	B-499	Black on white	0.750 (19.1)	16 ft. (4.8 m)	-	1	1 roll (cont.)	In cartridge
BMP®51 Label Printer cartridges										
15 mL	Fig. 83	M-95-492	B-492	Black on white	0.600 (15.2)	1.625 (41.3)	0.375 (9.5) dia.	1	200	In cartridge
	Fig. 2	M-156-492	B-492	Black on white	1.000 (25.4)	0.900 (22.9)	-	1	180	In cartridge
	Fig. 2	M-49-494-BL	B-494	Black on blue / white	1.000 (25.4)	1.000 (25.4)	-	1	260	In cartridge
	Fig. 2	M-49-494-GN	B-494	Black on green / white	1.000 (25.4)	1.000 (25.4)	-	1	260	In cartridge
	Fig. 2	M-49-494-OR	B-494	Black on orange / white	1.000 (25.4)	1.000 (25.4)	-	1	260	In cartridge
	Fig. 2	M-49-494-PK	B-494	Black on pink / white	1.000 (25.4)	1.000 (25.4)	-	1	260	In cartridge
	Fig. 2	M-49-494-RD	B-494	Black on red / white	1.000 (25.4)	1.000 (25.4)	-	1	260	In cartridge
	Fig. 2	M-49-494-YL	B-494	Black on yellow / white	1.000 (25.4)	1.000 (25.4)	-	1	260	In cartridge
	Fig. 82	M-78-432-CL-BK	B-432	Black on clear	1.000 (25.4)	1.900 (48.3)	-	1	100	In cartridge
	Fig. 82	M-78-492	B-492	Black on white	1.000 (25.4)	1.900 (48.3)	-	1	100	In cartridge
	Fig. 82	M-128-499	B-499	Black on white	1.000 (25.4)	1.900 (48.3)	-	1	100	In cartridge
	Fig. 33	MC-1000-499	B-499	Black on white	1.000 (25.4)	16 ft. (4.8 m)	-	1	1 roll (cont.)	In cartridge
	Fig. 33	MC-1000-461	B-461	Black on white	1.000 (25.4)	25 ft. (7.6 m)	-	1	1 roll (cont.)	In cartridge
	Fig. 33	MC-1250-461	B-461	Black on white	1.250 (31.8)	25 ft. (7.6 m)	-	1	1 roll (cont.)	In cartridge
50 mL	Fig. 82	M-113-490	B-490	Black on white	1.500 (38.1)	2.200 (55.9)	-	1	140	In cartridge
	Fig. 82	M-114-490	B-490	Black on white	1.500 (38.1)	3.750 (95.3)	-	1	80	In cartridge

Centrifuge and conical tube labels continue on next page.

*Tube sizes may vary by brand.

Sample collection and processing labels

Centrifuge and conical tube labels (continued)




Figure 2




Figure 3




Figure 5a




Figure 33




Figure 82

Fits tube size*	Figure	Catalog #	Material	Color	Label width A inch (mm)	Label length B inch (mm)	Secondary label C inch (mm)	Labels per row	Labels per roll	Recommended ribbon
BradyPrinter M611 label rolls										
15 mL	Fig. 82	M61-156-492	B-492	White	0.900 (22.9)	1.000 (25.4)	-	1	250	M61-R6410
	Fig. 2	M61-19-492	B-492	White	1.000 (25.4)	1.000 (25.4)	-	1	250	M61-R6410
	Fig. 2	M61-19-494-BL	B-494	Blue / white	1.000 (25.4)	1.000 (25.4)	-	1	250	M61-R6410
	Fig. 2	M61-19-494-GN	B-494	Green / white	1.000 (25.4)	1.000 (25.4)	-	1	250	M61-R6410
	Fig. 2	M61-19-494-OR	B-494	Orange / white	1.000 (25.4)	1.000 (25.4)	-	1	250	M61-R6410
	Fig. 2	M61-19-494-PK	B-494	Pink / white	1.000 (25.4)	1.000 (25.4)	-	1	250	M61-R6410
	Fig. 2	M61-19-494-RD	B-494	Red / white	1.000 (25.4)	1.000 (25.4)	-	1	250	M61-R6410
	Fig. 2	M61-19-494-YL	B-494	Yellow / white	1.000 (25.4)	1.000 (25.4)	-	1	250	M61-R6410
	Fig. 82	M61-203-490	B-490	White	1.000 (25.4)	1.750 (44.5)	-	1	250	M61-R4310
	Fig. 82	M61-203-492	B-492	White	1.000 (25.4)	1.750 (44.5)	-	1	250	M61-R6410
	Fig. 82	M61-204-492	B-492	White	1.000 (25.4)	1.900 (48.3)	-	1	250	M61-R6410
	Fig. 33	M61C-1000-499	B-499	White	1.000 (25.4)	30 ft. (9.1 m)	-	1	1 roll (cont.)	M61-R6210
	Fig. 33	M61C-1000-461	B-461	All White	1.000 (25.4)	50 ft. (15.2 m)	-	1	1 roll (cont.)	M61-R4310
	Fig. 33	M61C-1250-461	B-461	All White	1.250 (31.7)	50 ft. (15.2 m)	-	1	1 roll (cont.)	M61-R4310
	Fig. 2	PTL-78-490	B-490	White	1.900 (48.3)	1.000 (25.4)	-	1	250	M61-R4310
50 mL	Fig. 82	PTL-113-490	B-490	White	1.500 (38.1)	2.200 (55.9)	-	1	100	M61-R4310
	Fig. 82	PTL-114-490	B-490	White	1.500 (38.1)	3.875 (98.4)	-	1	100	M61-R4310
BradyPrinter i3300 and BBP®33 Printer label rolls										
15 mL	Fig. 2	B33-179-492	B-492	White	1.000 (25.4)	1.000 (25.4)	-	1	1,900	B30-R6400
	Fig. 2	B33-179-494-BL	B-494	Blue / white	1.000 (25.4)	1.000 (25.4)	-	1	1,500	B30-R6400
	Fig. 2	B33-179-494-GN	B-494	Green / white	1.000 (25.4)	1.000 (25.4)	-	1	1,500	B30-R6400
	Fig. 2	B33-179-494-OR	B-494	Orange / white	1.000 (25.4)	1.000 (25.4)	-	1	1,500	B30-R6400
	Fig. 2	B33-179-494-PK	B-494	Pink / white	1.000 (25.4)	1.000 (25.4)	-	1	1,500	B30-R6400
	Fig. 2	B33-179-494-RD	B-494	Red / white	1.000 (25.4)	1.000 (25.4)	-	1	1,500	B30-R6400
	Fig. 2	B33-179-494-YL	B-494	Yellow / white	1.000 (25.4)	1.000 (25.4)	-	1	1,500	B30-R6400
	Fig. 2	B33-179-7425	B-7425	White	1.000 (25.4)	1.000 (25.4)	-	1	1,500	B30-R4300
	Fig. 2	B33-136-499	B-499	White	1.500 (38.1)	0.750 (19.1)	-	1	1,500	B30-R6000
	Fig. 2	B33-136-7425	B-7425	White	1.500 (38.1)	0.750 (19.1)	-	1	1,500	B30-R4300
	Fig. 2	B33-258-7425	B-7425	White	1.500 (38.1)	1.750 (44.5)	-	1	900	B30-R4300
	Fig. 2	B33-155-490	B-490	White	1.625 (41.3)	0.600 (15.2)	-	1	1,500	B30-R4300
	Fig. 2	B33-156-490	B-490	White	1.900 (48.3)	1.000 (25.4)	-	1	1,500	B30-R4300
	Fig. 2	B33-156-492	B-492	White	1.900 (48.3)	1.000 (25.4)	-	1	1,500	B30-R6400
	50 mL	Fig. 2	B33-137-499	B-499	White	2.000 (50.8)	1.000 (25.4)	-	1	1,100
Fig. 2		B33-137-7425	B-7425	White	2.000 (50.8)	1.000 (25.4)	-	1	1,500	B30-R4300
Fig. 2		B33-157-490	B-490	White	2.200 (55.9)	1.500 (38.1)	-	1	1,000	B30-R4300
Fig. 2		B33-158-490	B-490	White	3.750 (95.3)	1.500 (38.1)	-	1	1,000	B30-R4300
BradyJet J2000 Printer label rolls										
15 mL	Fig. 2	J20-179-7425J	B-7425J	White	1.000 (25.4)	1.000 (25.4)	-	1	1,500	J20-CMY
	Fig. 3	J20C-1125-7425J	B-7425J	White	1.125 (28.6)	100 ft. (30.5 m)	-	1	1 roll (cont.)	J20-CMY
	Fig. 2	J20-136-7425J	B-7425J	White	1.500 (38.1)	0.750 (19.1)	-	1	1,500	J20-CMY
	Fig. 2	J20-258-7425J	B-7425J	White	1.500 (38.1)	1.750 (44.5)	-	1	900	J20-CMY
	Fig. 3	J20C-2250-7425J	B-7425J	White	2.250 (57.2)	100 ft. (30.5 m)	-	1	1 roll (cont.)	J20-CMY
50 mL	Fig. 2	J20-137-7425J	B-7425J	White	2.000 (50.8)	1.000 (25.4)	-	1	1,500	J20-CMY
	Fig. 2	J20-18-7425J	B-7425J	White	3.000 (76.2)	1.000 (25.4)	-	1	1,000	J20-CMY

Centrifuge and conical tube labels continue on next page.




Sample collection and processing labels

Centrifuge and conical tube labels (continued)

Fits tube size*	Figure	Catalog #	Material	Color	Label width A inch (mm)	Label length B inch (mm)	Secondary label C inch (mm)	Labels per row	Labels per roll	Recommended ribbon
BBP®12 Printer and other industry printers using 1-inch small core label rolls										
15 mL	Fig. 2	THT-179-492-2.5-SC	B-492	White	1.000 (25.4)	1.000 (25.4)	-	1	2,500	R6413
	Fig. 2	THT-179-7425-2-SC	B-7425	White	1.000 (25.4)	1.000 (25.4)	-	1	2,000	R6413
	Fig. 2	THT-136-7425-2-SC	B-7425	White	1.500 (38.1)	0.750 (19.1)	-	1	2,000	R6413
	Fig. 2	THT-155-490-1.5-SC	B-490	White	1.625 (41.3)	0.600 (15.2)	-	1	1,500	R4313
	Fig. 2	THT-258-7425-2-SC	B-7425	White	1.750 (44.5)	1.500 (38.1)	-	1	2,000	R6413
	Fig. 2	THT-156-490-1.5-SC	B-490	White	1.900 (48.3)	1.000 (25.4)	-	1	1,500	R4313
	Fig. 2	THT-156-492-1.5-SC	B-492	White	1.900 (48.3)	1.000 (25.4)	-	1	1,500	R6413
15 mL top	Fig. 5a	THT-143-492-1.5-SC	B-492	White	0.375 (9.5) dia.	-	-	4	1,500	R6413
50 mL	Fig. 2	THT-137-7425-2-SC	B-7425	White	2.000 (50.8)	1.000 (25.4)	-	1	2,000	R6413
	Fig. 2	THT-157-490-1.5-SC	B-490	White	2.200 (55.9)	1.500 (38.1)	-	1	1,500	R4313
	Fig. 2	THT-158-490-1.5-SC	B-490	White	3.750 (95.3)	1.500 (38.1)	-	1	1,500	R4313
BradyPrinter i5100, i7100, Brady IP Printers and other industry printers using 3-inch core label rolls										
15 mL	Fig. 2	THT-135-427-3-AW	B-427	All White	0.750 (19.1)	1.500 (38.1)	-	1	3,000	R4302 / IP-R4302
	Fig. 2	THT-179-492-3	B-492	White	1.000 (25.4)	1.000 (25.4)	-	1	3,000	R6400 / IP-R6400
	Fig. 2	THT-179-494-BL	B-494	Blue / white	1.000 (25.4)	1.000 (25.4)	-	1	3,000	R6402 / IP-R6402
	Fig. 2	THT-179-494-GN	B-494	Green / white	1.000 (25.4)	1.000 (25.4)	-	1	3,000	R6402 / IP-R6402
	Fig. 2	THT-179-494-OR	B-494	Orange / white	1.000 (25.4)	1.000 (25.4)	-	1	3,000	R6402 / IP-R6402
	Fig. 2	THT-179-494-PK	B-494	Pink / white	1.000 (25.4)	1.000 (25.4)	-	1	3,000	R6402 / IP-R6402
	Fig. 2	THT-179-494-RD	B-494	Red / white	1.000 (25.4)	1.000 (25.4)	-	1	3,000	R6402 / IP-R6402
	Fig. 2	THT-179-494-YL	B-494	Yellow / white	1.000 (25.4)	1.000 (25.4)	-	1	3,000	R6402 / IP-R6402
	Fig. 2	THT-179-7425-2	B-7425	White	1.000 (25.4)	1.000 (25.4)	-	1	2,000	R4302 / IP-R4302
	Fig. 2	THT-136-7425-2	B-7425	White	1.500 (38.1)	0.750 (19.1)	-	1	2,000	R4302 / IP-R4302
	Fig. 2	THT-155-490-3	B-490	White	1.625 (41.3)	0.600 (15.2)	-	1	3,000	R4302 / IP-R4302
	Fig. 2	THT-258-7425-2	B-7425	White	1.750 (44.5)	1.500 (38.1)	-	1	2,000	R4302 / IP-R4302
	Fig. 2	THT-156-490-3	B-490	White	1.900 (48.3)	1.000 (25.4)	-	1	3,000	R4302 / IP-R4302
	Fig. 2	THT-156-492-3	B-492	White	1.900 (48.3)	1.000 (25.4)	-	1	3,000	R6400 / IP-R6400
50 mL	Fig. 2	THT-137-7425-2	B-7425	White	2.000 (50.8)	1.000 (25.4)	-	1	2,000	R4302 / IP-R4302
	Fig. 2	THT-137-427-3-AW	B-427	All White	2.000 (50.8)	1.000 (25.4)	-	1	3,000	R4302 / IP-R4302
	Fig. 2	THT-157-490-3	B-490	White	2.200 (55.9)	1.500 (38.1)	-	1	3,000	R4300 / IP-R4300
	Fig. 2	THT-158-490-3	B-490	White	3.750 (95.3)	1.500 (38.1)	-	1	3,000	R4307 / IP-R4307

Self-laminating centrifuge and conical tube labels

- B-461 – Ultra-thin polyester with clear overlaminat tail that protects print and allows liquid levels to be visible



Fits tube size*	Figure	Catalog #	Material	Color	Label width A inch (mm)	Label length B inch (mm)	Print-on area C inch (mm)	Labels per row	Labels per roll	Recommended ribbon
BMP®51 Label Printer cartridges										
15 mL	Fig. 15	M-129-461	B-461	Black on white	2.250 (57.2)	1.500 (38.1)	1.250 (31.7)	1	140	In cartridge
	Fig. 15	M-119-461	B-461	Black on white	3.750 (95.3)	1.500 (38.1)	1.250 (31.7)	1	80	In cartridge
BradyPrinter M611 label rolls										
15 mL	Fig. 14	PTL-77-461	B-461	White	3.500 (88.9)	1.500 (38.1)	1.250 (31.8)	1	250	M61-R4310
BradyPrinter i3300 and BBP®33 Printer label rolls										
15 mL	Fig. 14c	B33-276-461	B-461	White	2.625 (66.7)	1.000 (25.4)	1.000 (25.4)	2	1,500	B30-R4300
	Fig. 15	B33-134-461	B-461	White	3.500 (88.9)	1.500 (38.1)	1.250 (31.8)	1	1,000	B30-R4300
BradyPrinter i5100, i7100, Brady IP Printers and other industry printers using 3-inch core label rolls										
15 mL	Fig. 14c	THT-139-461-2	B-461	White	1.000 (25.4)	5.250 (133.4)	1.500 (38.1)	2	2,000	R4300 / IP-R4300
	Fig. 15	THT-134-461-3	B-461	White	1.500 (38.1)	3.500 (88.9)	1.250 (31.8)	1	3,000	R4307 / IP-R4307


*Tube sizes may vary by brand.

Sample collection and processing labels

Microcentrifuge and PCR tube labels

Designed to fit the small sizes of microcentrifuge and PCR tubes, these labels withstand harsh temperatures and chemicals without falling off, fading or smearing. Once in the microcentrifuge, vial top labels allow you to quickly extract specific vials.

- B-430 – Clear high-performance polyester label
- B-461 – All-white, ultra-thin polyester label. Also available in a self-laminating format
- B-490 – White polyester is ideal for flat frozen surfaces. Adheres to frozen tubes and vials when wrapped around itself by at least 0.125 inch
- B-492 – Ultra-thin white polyester with superior chemical and abrasion resistance for labeling frosted, frozen or room-temperature tubes and vials. Must have small overlap if applied to already frozen containers
- B-494 – Glossy polyester with a printed color bar that allows for quick visual differentiation between samples. Withstands liquid nitrogen, a variety of chemicals and more. Label dimensions listed below include the color bar area. Vial tops are solid color
- B-499 – White nylon cloth with permanent adhesive for easy handling and conformability
- B-7425 / B-7425J – Economical white polypropylene label with moderate chemical resistance that is ideal for labeling dry, room-temperature specimen containers. Must be applied to containers prior to being frozen



Sample collection labels

Fits tube size*	Figure	Catalog #	Material	Color	Label width A inch (mm)	Label length B inch (mm)	Secondary label C inch (mm)	Labels per row	Labels per roll	Recommended ribbon
BMP®21-LAB Label Printer cartridges										
2 mL	Fig. 33b	M21-750-499	B-499	Black on white	0.750 (19.1)	16 ft. (4.8 m)	-	1	1 roll (cont.)	In cartridge
	Fig. 33b	M21-750-430	B-430	Black on clear	0.750 (19.1)	21 ft. (6.4 m)	-	1	1 roll (cont.)	In cartridge
	Fig. 33b	M21-750-7425	B-7425	Black on white	0.750 (19.1)	21 ft. (6.4 m)	-	1	1 roll (cont.)	In cartridge
	Fig. 33b	M21-750-461	B-461	Black on white	0.750 (19.1)	21 ft. (6.4 m)	-	1	1 roll (cont.)	In cartridge
0.5 mL	Fig. 33b	M21-250-430	B-430	Black on clear	0.250 (6.4)	21 ft. (6.4 m)	-	1	1 roll (cont.)	In cartridge
	Fig. 33b	M21-375-499	B-499	Black on white	0.375 (9.5)	16 ft. (4.8 m)	-	1	1 roll (cont.)	In cartridge
	Fig. 33b	M21-375-430	B-430	Black on clear	0.375 (9.5)	21 ft. (6.4 m)	-	1	1 roll (cont.)	In cartridge
	Fig. 33b	M21-375-7425	B-7425	Black on white	0.375 (9.5)	21 ft. (6.4 m)	-	1	1 roll (cont.)	In cartridge
	Fig. 33b	M21-375-461	B-461	Black on white	0.375 (9.5)	21 ft. (6.4 m)	-	1	1 roll (cont.)	In cartridge
	Fig. 33b	M21-500-499	B-499	Black on white	0.500 (12.7)	16 ft. (4.8 m)	-	1	1 roll (cont.)	In cartridge
	Fig. 33b	M21-500-430	B-430	Black on clear	0.500 (12.7)	21 ft. (6.4 m)	-	1	1 roll (cont.)	In cartridge
	Fig. 33b	M21-500-7425	B-7425	Black on white	0.500 (12.7)	21 ft. (6.4 m)	-	1	1 roll (cont.)	In cartridge
Fig. 33b	M21-500-461	B-461	Black on white	0.500 (12.7)	21 ft. (6.4 m)	-	1	1 roll (cont.)	In cartridge	
BMP®51 Label Printer cartridges										
0.5 mL	Fig. 82	M-130-499	B-499	Black on white	0.375 (9.5)	0.825 (21.0)	-	1	220	In cartridge
	Fig. 33	MC-375-499	B-499	Black on white	0.375 (9.5)	16 ft. (4.8 m)	-	1	1 roll (cont.)	In cartridge
	Fig. 33	MC-375-461	B-461	Black on white	0.375 (9.5)	25 ft. (7.6 m)	-	1	1 roll (cont.)	In cartridge
	Fig. 33	MC-500-499	B-499	Black on white	0.500 (12.7)	16 ft. (4.8 m)	-	1	1 roll (cont.)	In cartridge
	Fig. 33	MC-500-492	B-492	Black on white	0.500 (12.7)	25 ft. (7.6 m)	-	1	1 roll (cont.)	In cartridge
	Fig. 33	MC-500-461	B-461	Black on white	0.500 (12.7)	25 ft. (7.6 m)	-	1	1 roll (cont.)	In cartridge

Microcentrifuge and PCR tube labels continue on next page.

*Tube sizes may vary by brand.

Sample collection and processing labels

Microcentrifuge and PCR tube labels (continued)

Fits tube size*	Figure	Catalog #	Material	Color	Label width A inch (mm)	Label length B inch (mm)	Secondary label C inch (mm)	Labels per row	Labels per roll	Recommended ribbon
1.5 mL	Fig. 82	M-131-492	B-492	Black on white	0.500 (12.7)	1.000 (25.4)	-	1	180	In cartridge
	Fig. 82	M-131-494-BL	B-494	Black on blue / white	0.500 (12.7)	1.000 (25.4)	-	1	180	In cartridge
	Fig. 82	M-131-494-GN	B-494	Black on green / white	0.500 (12.7)	1.000 (25.4)	-	1	180	In cartridge
	Fig. 82	M-131-494-OR	B-494	Black on orange / white	0.500 (12.7)	1.000 (25.4)	-	1	180	In cartridge
	Fig. 82	M-131-494-PK	B-494	Black on pink / white	0.500 (12.7)	1.000 (25.4)	-	1	180	In cartridge
	Fig. 82	M-131-494-RD	B-494	Black on red / white	0.500 (12.7)	1.000 (25.4)	-	1	180	In cartridge
	Fig. 82	M-131-494-YL	B-494	Black on yellow / white	0.500 (12.7)	1.000 (25.4)	-	1	180	In cartridge
	Fig. 82	M-131-499	B-499	Black on white	0.500 (12.7)	1.000 (25.4)	-	1	180	In cartridge
	Fig. 82	M-124-492	B-492	Black on white	0.500 (12.7)	1.650 (41.9)	-	1	180	In cartridge
	Fig. 24	M-118-492	B-492	Black on white	1.000 (25.4)	0.375 (9.5)	0.375 (9.5) dia.	2	240	In cartridge
	Fig. 24	M-118-494-BL	B-494	Black on blue / white	1.000 (25.4)	0.375 (9.5)	0.375 (9.5) dia.	1	240	In cartridge
	Fig. 24	M-118-494-GN	B-494	Black on green / white	1.000 (25.4)	0.375 (9.5)	0.375 (9.5) dia.	1	240	In cartridge
	Fig. 24	M-118-494-OR	B-494	Black on orange / white	1.000 (25.4)	0.375 (9.5)	0.375 (9.5) dia.	1	240	In cartridge
	Fig. 24	M-118-494-PK	B-494	Black on pink / white	1.000 (25.4)	0.375 (9.5)	0.375 (9.5) dia.	1	240	In cartridge
	Fig. 24	M-118-494-YL	B-494	Black on red / white	1.000 (25.4)	0.375 (9.5)	0.375 (9.5) dia.	1	240	In cartridge
	Fig. 24	M-118-494-RD	B-494	Black on yellow / white	1.000 (25.4)	0.375 (9.5)	0.375 (9.5) dia.	1	250	In cartridge
	Fig. 24	M-118-499	B-499	Black on white	1.000 (25.4)	0.375 (9.5)	0.375 (9.5) dia.	2	240	In cartridge
	Fig. 24	M-118-7425	B-7425	Black on white	1.000 (25.4)	0.375 (9.5)	0.375 (9.5) dia.	1	240	In cartridge
	Fig. 24	M-120-494-BL	B-494	Black on blue / white	1.000 (25.4)	0.500 (12.7)	0.375 (9.5) dia.	1	240	In cartridge
	Fig. 24	M-120-494-GN	B-494	Black on green / white	1.000 (25.4)	0.500 (12.7)	0.375 (9.5) dia.	1	240	In cartridge
	Fig. 24	M-120-494-OR	B-494	Black on orange / white	1.000 (25.4)	0.500 (12.7)	0.375 (9.5) dia.	1	240	In cartridge
	Fig. 24	M-120-494-PK	B-494	Black on pink / white	1.000 (25.4)	0.500 (12.7)	0.375 (9.5) dia.	1	240	In cartridge
	Fig. 24	M-120-494-RD	B-494	Black on red / white	1.000 (25.4)	0.500 (12.7)	0.375 (9.5) dia.	1	240	In cartridge
	Fig. 24	M-120-494-YL	B-494	Black on yellow / white	1.000 (25.4)	0.500 (12.7)	0.375 (9.5) dia.	1	240	In cartridge
	Fig. 24	M-120-492	B-492	Black on white	1.000 (25.4)	0.500 (12.7)	0.375 (9.5) dia.	1	240	In cartridge
	Fig. 24	M-120-499	B-499	Black on white	1.000 (25.4)	0.500 (12.7)	0.375 (9.5) dia.	1	240	In cartridge
	Fig. 24	M-117-492	B-492	Black on white	1.000 (25.4)	0.500 (12.7)	0.438 (11.1) dia.	2	240	In cartridge
	Fig. 24	M-117-499	B-499	Black on white	1.000 (25.4)	0.500 (12.7)	0.438 (11.1) dia.	2	240	In cartridge
	Fig. 2	M-112-490	B-490	Black on white	1.250 (31.8)	0.375 (9.5)	-	1	360	In cartridge
	Fig. 2	M-112-492	B-492	Black on white	1.250 (31.8)	0.375 (9.5)	-	1	360	In cartridge
2 mL	Fig. 33	MC-750-499	B-499	Black on white	0.750 (19.1)	16 ft. (4.8 m)	-	1	1 roll (cont.)	In cartridge
	Fig. 2	M-156-492	B-492	Black on white	1.000 (25.4)	0.900 (22.9)	-	1	180	In cartridge
	Fig. 2	M-49-494-BL	B-494	Black on blue / white	1.000 (25.4)	1.000 (25.4)	-	1	260	In cartridge
	Fig. 2	M-49-494-GN	B-494	Black on green / white	1.000 (25.4)	1.000 (25.4)	-	1	260	In cartridge
	Fig. 2	M-49-494-OR	B-494	Black on orange / white	1.000 (25.4)	1.000 (25.4)	-	1	260	In cartridge
	Fig. 2	M-49-494-PK	B-494	Black on pink / white	1.000 (25.4)	1.000 (25.4)	-	1	260	In cartridge
	Fig. 2	M-49-494-RD	B-494	Black on red / white	1.000 (25.4)	1.000 (25.4)	-	1	260	In cartridge
	Fig. 2	M-49-494-YL	B-494	Black on yellow / white	1.000 (25.4)	1.000 (25.4)	-	1	260	In cartridge
	Fig. 33	MC-1000-499	B-499	Black on white	1.000 (25.4)	16 ft. (4.8 m)	-	1	1 roll (cont.)	In cartridge
	Any top	Fig. 8	M-82-499	B-499	Black on white	0.375 (9.5) dia.	-	-	1	240
Fig. 8		M-82-492	B-492	Black on white	0.375 (9.5) dia.	-	-	1	240	In cartridge
Fig. 8		M-82-499-OR-BK	B-499	Black on orange	0.375 (9.5) dia.	-	-	1	240	In cartridge
Fig. 8		M-82-499-RD-BK	B-499	Black on red	0.375 (9.5) dia.	-	-	1	240	In cartridge
Fig. 8		M-82-499-YL-BK	B-499	Black on yellow	0.375 (9.5) dia.	-	-	1	240	In cartridge
Fig. 8		M-83-499	B-499	Black on white	0.500 (12.7) dia.	-	-	1	240	In cartridge
Fig. 8		M-83-492	B-492	Black on white	0.500 (12.7) dia.	-	-	1	240	In cartridge
Fig. 8		M-83-499-OR-BK	B-499	Black on orange	0.500 (12.7) dia.	-	-	1	240	In cartridge
Fig. 8		M-83-499-RD-BK	B-499	Black on red	0.500 (12.7) dia.	-	-	1	240	In cartridge
Fig. 8		M-83-499-YL-BK	B-499	Black on yellow	0.500 (12.7) dia.	-	-	1	240	In cartridge
BradyPrinter M611 label rolls										
0.5 mL	Fig. 33	M61C-375-499	B-499	White	0.375 (9.5)	30 ft. (9.1 m)	-	1	1 roll (cont.)	M61-R6210
	Fig. 33	M61C-375-461	B-461	White	0.375 (9.5)	50 ft. (15.2 m)	-	1	1 roll (cont.)	M61-R4310
	Fig. 33	M61C-500-499	B-499	White	0.500 (12.7)	30 ft. (9.1 m)	-	1	1 roll (cont.)	M61-R6210
	Fig. 33	M61C-500-492	B-492	White	0.500 (12.7)	50 ft. (15.2 m)	-	1	1 roll (cont.)	M61-R6410
	Fig. 33	M61C-500-461	B-461	White	0.500 (12.7)	50 ft. (15.2 m)	-	1	1 roll (cont.)	M61-R4310
	Fig. 2	M61-130-499	B-499	White	0.825 (21.0)	0.375 (9.5)	-	1	500	M61-R6210

Microcentrifuge and PCR tube labels continue on next page.

*Tube sizes may vary by brand.

Sample collection and processing labels

Microcentrifuge and PCR tube labels (continued)




Figure 1




Figure 2




Figure 5a




Figure 8




Figure 24




Figure 25




Figure 33




Figure 82

Fits tube size*	Figure	Catalog #	Material	Color	Label width A inch (mm)	Label length B inch (mm)	Secondary label C inch (mm)	Labels per row	Labels per roll	Recommended ribbon	
1.5 mL	Fig. 82	M61-47-499	B-499	White	0.500 (12.7)	1.000 (25.4)	-	1	500	M61-R6210	
	Fig. 82	M61-124-492	B-492	White	0.500 (12.7)	1.650 (41.9)	-	1	500	M61-R6410	
	Fig. 24	M61-117-492	B-492	White	0.900 (22.9)	0.500 (12.7)	0.500 (12.7) dia.	1	500	M61-R6410	
	Fig. 24	PTL-117-499	B-499	White	0.900 (22.9)	0.500 (12.7)	0.500 (12.7) dia.	1	500	M61-R6210	
	Fig. 24	M61-98-492	B-492	White	1.000 (25.4)	0.375 (9.5)	0.375 (9.5) dia.	1	500	M61-R6410	
	Fig. 24	M61-98-494-BL	B-494	Blue / white	1.000 (25.4)	0.375 (9.5)	0.375 (9.5) dia.	1	500	M61-R6410	
	Fig. 24	M61-98-494-GN	B-494	Green / white	1.000 (25.4)	0.375 (9.5)	0.375 (9.5) dia.	1	500	M61-R6410	
	Fig. 24	M61-98-494-OR	B-494	Orange / white	1.000 (25.4)	0.375 (9.5)	0.375 (9.5) dia.	1	500	M61-R6410	
	Fig. 24	M61-98-494-PK	B-494	Pink / white	1.000 (25.4)	0.375 (9.5)	0.375 (9.5) dia.	1	500	M61-R6410	
	Fig. 24	M61-98-494-RD	B-494	Red / white	1.000 (25.4)	0.375 (9.5)	0.375 (9.5) dia.	1	500	M61-R6410	
	Fig. 24	M61-98-494-YL	B-494	Yellow / white	1.000 (25.4)	0.375 (9.5)	0.375 (9.5) dia.	1	500	M61-R6410	
	Fig. 24	PTL-98-499	B-499	White	1.000 (25.4)	0.375 (9.5)	0.375 (9.5) dia.	1	500	M61-R6210	
	Fig. 24	M61-98-7425	B-7425	White	1.000 (25.4)	0.375 (9.5)	0.375 (9.5) dia.	1	500	M61-R4310	
	Fig. 2	M61-17-492	B-492	White	1.000 (25.4)	0.500 (12.7)	-	1	250	M61-R6410	
	Fig. 2	M61-17-494-BL	B-494	Blue / white	1.000 (25.4)	0.500 (12.7)	-	1	500	M61-R6410	
	Fig. 2	M61-17-494-GN	B-494	Green / white	1.000 (25.4)	0.500 (12.7)	-	1	500	M61-R6410	
	Fig. 2	M61-17-494-OR	B-494	Orange / white	1.000 (25.4)	0.500 (12.7)	-	1	500	M61-R6410	
	Fig. 2	M61-17-494-PK	B-494	Pink / white	1.000 (25.4)	0.500 (12.7)	-	1	500	M61-R6410	
	Fig. 2	M61-17-494-RD	B-494	Red / white	1.000 (25.4)	0.500 (12.7)	-	1	500	M61-R6410	
	Fig. 2	M61-17-494-YL	B-494	Yellow / white	1.000 (25.4)	0.500 (12.7)	-	1	500	M61-R6410	
	Fig. 2	M61-17-7425	B-7425	White	1.000 (25.4)	0.500 (12.7)	-	1	500	M61-R4310	
	Fig. 24	M61-120-492	B-492	White	1.000 (25.4)	0.500 (12.7)	0.375 (9.5) dia.	1	250	M61-R6410	
	Fig. 24	M61-120-499	B-499	White	1.000 (25.4)	0.500 (12.7)	0.375 (9.5) dia.	1	250	M61-R6210	
Fig. 2	M61-199-7425	B-7425	White	1.300 (33.0)	0.600 (15.2)	-	1	250	M61-R4310		
Fig. 2	PTL-112-490	B-490	White	1.250 (31.8)	0.375 (9.5)	-	1	500	M61-R4310		
2 mL	Fig. 33	M61C-750-499	B-499	White	0.750 (19.1)	30 ft. (9.1 m)	-	1	1 roll (cont.)	M61-R6210	
	Fig. 82	M61-156-492	B-492	White	0.900 (22.9)	1.000 (25.4)	-	1	250	M61-R6410	
	Fig. 2	M61-18-7425	B-7425	White	1.000 (25.4)	0.750 (19.1)	-	1	250	M61-R4310	
	Fig. 2	M61-19-492	B-492	White	1.000 (25.4)	1.000 (25.4)	-	1	250	M61-R6410	
	Fig. 2	M61-19-494-BL	B-494	Blue / white	1.000 (25.4)	1.000 (25.4)	-	1	250	M61-R6410	
	Fig. 2	M61-19-494-GN	B-494	Green / white	1.000 (25.4)	1.000 (25.4)	-	1	250	M61-R6410	
	Fig. 2	M61-19-494-OR	B-494	Orange / white	1.000 (25.4)	1.000 (25.4)	-	1	250	M61-R6410	
	Fig. 2	M61-19-494-PK	B-494	Pink / white	1.000 (25.4)	1.000 (25.4)	-	1	250	M61-R6410	
	Fig. 2	M61-19-494-RD	B-494	Red / white	1.000 (25.4)	1.000 (25.4)	-	1	250	M61-R6410	
	Fig. 2	M61-19-494-YL	B-494	Yellow / white	1.000 (25.4)	1.000 (25.4)	-	1	250	M61-R6410	
	Fig. 33	M61C-1000-499	B-499	White	1.000 (25.4)	30 ft. (9.1 m)	-	1	1 roll (cont.)	M61-R6210	
	Any top	Fig. 8	M61-82-492	B-492	White	0.375 (9.5) dia.	-	-	1	500	M61-R6410
		Fig. 8	PTL-82-499	B-499	White	0.375 (9.5) dia.	-	-	1	500	M61-R6210
Fig. 8		PTL-82-499-YL	B-499	Yellow	0.375 (9.5) dia.	-	-	1	500	M61-R6210	
Fig. 8		M61-82-499-OR	B-499	Orange	0.375 (9.5) dia.	-	-	1	500	M61-R6210	
Fig. 8		M61-82-499-RD	B-499	Red	0.375 (9.5) dia.	-	-	1	500	M61-R6210	
Fig. 8		M61-83-492	B-492	White	0.500 (12.7) dia.	-	-	1	500	M61-R6410	
Fig. 8		PTL-83-499	B-499	White	0.500 (12.7) dia.	-	-	1	500	M61-R6210	
Fig. 8		PTL-83-499-YL	B-499	Yellow	0.500 (12.7) dia.	-	-	1	500	M61-R6210	
Fig. 8		M61-83-499-OR	B-499	Orange	0.500 (12.7) dia.	-	-	1	500	M61-R6210	
Fig. 8	M61-83-499-RD	B-499	Red	0.500 (12.7) dia.	-	-	1	500	M61-R6210		
BradyPrinter i3300 and BBP#33 Printer label rolls											
0.5 mL	Fig. 1	B33-68-490	B-490	White	0.500 (12.7)	0.750 (19.1)	-	4	2,500	B30-R4300	
	Fig. 1	B33-68-499	B-499	White	0.500 (12.7)	0.750 (19.1)	-	4	1,500	B30-R6000	
	Fig. 2	B33-257-7425	B-7425	White	0.825 (21.0)	0.250 (6.4)	-	1	1,500	B30-R4300	
	Fig. 2	B33-125-490	B-490	White	0.900 (22.9)	0.500 (12.7)	-	1	1,650	B30-R4300	
	Fig. 2	B33-125-499	B-499	White	0.900 (22.9)	0.500 (12.7)	-	1	1,500	B30-R6000	
	Fig. 24	B33-126-490	B-490	White	0.900 (22.9)	0.500 (12.7)	0.500 (12.7) dia.	1	2,500	B30-R4300	
Fig. 24	B33-126-499	B-499	White	0.900 (22.9)	0.500 (12.7)	0.500 (12.7) dia.	2	1,500	B30-R4300		

Microcentrifuge and PCR tube labels continue on next page.

*Tube sizes may vary by brand.

Sample collection and processing labels

Microcentrifuge and PCR tube labels (continued)

Fits tube size*	Figure	Catalog #	Material	Color	Label width A inch (mm)	Label length B inch (mm)	Secondary label C inch (mm)	Labels per row	Labels per roll	Recommended ribbon
1.5 mL	Fig. 25	B33-152-492	B-492	White	1.000 (25.4)	0.375 (9.5)	-	2	2,500	B30-R6400
	Fig. 25	B33-152-494-BL	B-494	Blue / white	1.000 (25.4)	0.375 (9.5)	-	2	1,500	B30-R6400
	Fig. 25	B33-152-494-GN	B-494	Green / white	1.000 (25.4)	0.375 (9.5)	-	2	1,500	B30-R6400
	Fig. 25	B33-152-494-OR	B-494	Orange / white	1.000 (25.4)	0.375 (9.5)	-	2	1,500	B30-R6400
	Fig. 25	B33-152-494-PK	B-494	Pink / white	1.000 (25.4)	0.375 (9.5)	-	2	1,500	B30-R6400
	Fig. 25	B33-152-494-RD	B-494	Red / white	1.000 (25.4)	0.375 (9.5)	-	2	1,500	B30-R6400
	Fig. 25	B33-152-494-YL	B-494	Yellow / white	1.000 (25.4)	0.375 (9.5)	-	2	1,500	B30-R6400
	Fig. 25	B33-152-499	B-499	White	1.000 (25.4)	0.375 (9.5)	-	2	1,500	B30-R6000
	Fig. 25	B33-152-7425	B-7425	White	1.000 (25.4)	0.375 (9.5)	-	2	2,000	B30-R4300
	Fig. 24	B33-163-492	B-492	White	1.000 (25.4)	0.375 (9.5)	0.375 (9.5) dia.	1	1,500	B30-R6400
	Fig. 2	B33-163-499	B-499	White	1.000 (25.4)	0.375 (9.5)	0.375 (9.5) dia.	1	1,500	B30-R6000
	Fig. 24	B33-249-7425	B-7425	White	1.000 (25.4)	0.375 (9.5)	0.375 (9.5) dia.	1	1,500	B30-R4300
	Fig. 2	B33-59-7425	B-7425	White	1.000 (25.4)	0.500 (12.7)	-	1	2,000	B30-R4300
	Fig. 24	B33-181-492	B-492	White	1.000 (25.4)	0.500 (12.7)	0.375 (9.5) dia.	1	1,500	B30-R6400
	Fig. 24	B33-250-494-BL	B-494	Blue / white	1.000 (25.4)	0.500 (12.7)	0.375 (9.5) dia.	1	1,500	B30-R6400
	Fig. 24	B33-250-494-GN	B-494	Green / white	1.000 (25.4)	0.500 (12.7)	0.375 (9.5) dia.	1	1,500	B30-R6400
	Fig. 24	B33-250-494-OR	B-494	Orange / white	1.000 (25.4)	0.500 (12.7)	0.375 (9.5) dia.	1	1,500	B30-R6400
	Fig. 24	B33-250-494-PK	B-494	Pink / white	1.000 (25.4)	0.500 (12.7)	0.375 (9.5) dia.	1	1,500	B30-R6400
	Fig. 24	B33-250-494-RD	B-494	Red / white	1.000 (25.4)	0.500 (12.7)	0.375 (9.5) dia.	1	1,500	B30-R6400
	Fig. 24	B33-250-494-YL	B-494	Yellow / white	1.000 (25.4)	0.500 (12.7)	0.375 (9.5) dia.	1	1,500	B30-R6400
Fig. 24	B33-182-492	B-492	White	1.000 (25.4)	0.500 (12.7)	0.440 (11.2) dia.	1	1,500	B30-R6400	
Fig. 24	B33-251-7425	B-7425	White	1.000 (25.4)	0.500 (12.7)	0.440 (11.2) dia.	1	1,500	B30-R4300	
Fig. 2	B33-154-490	B-490	White	1.250 (31.8)	0.375 (9.5)	-	1	1,500	B30-R4300	
Fig. 2	B33-256-7425	B-7425	White	1.300 (33.0)	0.600 (15.2)	-	1	1,500	B30-R4300	
2 mL	Fig. 25	B33-88-7425	B-7425	White	1.000 (25.4)	0.750 (19.1)	-	2	2,000	B30-R4300
	Fig. 2	B33-179-492	B-492	White	1.000 (25.4)	1.000 (25.4)	-	1	1,900	B30-R6400
	Fig. 2	B33-179-494-BL	B-494	Blue / white	1.000 (25.4)	1.000 (25.4)	-	1	1,500	B30-R6400
	Fig. 2	B33-179-494-GN	B-494	Green / white	1.000 (25.4)	1.000 (25.4)	-	1	1,500	B30-R6400
	Fig. 2	B33-179-494-OR	B-494	Orange / white	1.000 (25.4)	1.000 (25.4)	-	1	1,500	B30-R6400
	Fig. 2	B33-179-494-PK	B-494	Pink / white	1.000 (25.4)	1.000 (25.4)	-	1	1,500	B30-R6400
	Fig. 2	B33-179-494-RD	B-494	Red / white	1.000 (25.4)	1.000 (25.4)	-	1	1,500	B30-R6400
	Fig. 2	B33-179-494-YL	B-494	Yellow / white	1.000 (25.4)	1.000 (25.4)	-	1	1,500	B30-R6400
Fig. 2	B33-179-7425	B-7425	White	1.000 (25.4)	1.000 (25.4)	-	1	1,500	B30-R4300	
Any top	Fig. 5a	B33-143-492	B-492	White	0.375 (9.5)	-	-	4	1,500	B30-R6400
	Fig. 5a	B33-143-499	B-499	White	0.375 (9.5)	-	-	4	1,500	B30-R6000
	Fig. 5a	B33-144-490	B-490	White	0.500 (12.7) dia.	-	-	3	2,500	B30-R4300
	Fig. 5a	B33-144-492	B-492	White	0.500 (12.7) dia.	-	-	3	1,500	B30-R6400
	Fig. 5a	B33-144-499	B-499	White	0.500 (12.7) dia.	-	-	3	1,500	B30-R6000
BradyJet J2000 Printer label rolls										
0.5 mL	Fig. 2	J20-257-7425J	B-7425J	White	0.825 (21.0)	0.250 (6.4)	-	1	1,500	J20-CMY
1.5 mL	Fig. 25	J20-152-7425J	B-7425J	White	1.000 (25.4)	0.375 (9.5)	-	2	2,000	J20-CMY
	Fig. 24	J20-249-7425J	B-7425J	White	1.000 (25.4)	0.375 (9.5)	0.375 (9.5) dia.	2	1,500	J20-CMY
	Fig. 2	J20-59-7425J	B-7425J	White	1.000 (25.4)	0.500 (12.7)	-	1	2,000	J20-CMY
	Fig. 24	J20-251-7425J	B-7425J	White	1.000 (25.4)	0.500 (12.7)	0.440 (11.2) dia.	2	1,500	J20-CMY
	Fig. 2	J20-256-7425J	B-7425J	White	1.300 (33.0)	0.600 (15.2)	-	1	1,500	J20-CMY
2 mL	Fig. 25	J20-88-7425J	B-7425J	White	1.000 (25.4)	0.750 (19.1)	-	2	2,000	J20-CMY
Fig. 2	J20-179-7425J	B-7425J	White	1.000 (25.4)	1.000 (25.4)	-	1	1,500	J20-CMY	
BBP®12 Printer and other industry printers using 1-inch small core label rolls										
0.5 mL	Fig. 2	THT-257-7425-2-SC	B-7425	White	0.825 (21.0)	0.250 (6.4)	-	1	2,000	R6413
	Fig. 2	THT-125-499-1.5-SC	B-499	White	0.900 (22.9)	0.500 (12.7)	-	1	1,500	R4913
1.5 mL	Fig. 25	THT-152-492-2.5-SC	B-492	White	1.000 (25.4)	0.375 (9.5)	-	2	2,500	R6413
	Fig. 25	THT-152-499-1.5-SC	B-499	White	1.000 (25.4)	0.375 (9.5)	-	2	1,500	R4913
	Fig. 25	THT-152-7425-2-SC	B-7425	White	1.000 (25.4)	0.375 (9.5)	-	2	2,000	R6413
	Fig. 24	THT-249-492-1.5-SC	B-492	White	1.000 (25.4)	0.375 (9.5)	0.375 (9.5) dia.	1	1,500	R6413
	Fig. 24	THT-249-499-1.5-SC	B-499	White	1.000 (25.4)	0.375 (9.5)	0.375 (9.5) dia.	1	1,500	R4913
	Fig. 24	THT-249-7425-2-SC	B-7425	White	1.000 (25.4)	0.375 (9.5)	0.375 (9.5) dia.	1	2,000	R6413
	Fig. 2	THT-59-492-2.5-SC	B-492	White	1.000 (25.4)	0.500 (12.7)	-	1	2,500	R6413
	Fig. 2	THT-59-7425-2-SC	B-7425	White	1.000 (25.4)	0.500 (12.7)	-	1	2,000	R6413
	Fig. 24	THT-250-492-1.5-SC	B-492	White	1.000 (25.4)	0.500 (12.7)	0.375 (9.5) dia.	1	1,500	R6413
	Fig. 24	THT-251-492-1.5-SC	B-492	White	1.000 (25.4)	0.500 (12.7)	0.440 (11.2) dia.	1	1,500	R6413
Fig. 24	THT-251-7425-2-SC	B-7425	White	1.000 (25.4)	0.500 (12.7)	0.440 (11.2) dia.	1	2,000	R6413	
Fig. 2	THT-256-7425-2-SC	B-7425	White	1.300 (33.0)	0.600 (15.2)	-	1	2,000	R6413	

*Tube sizes may vary by brand.

Microcentrifuge and PCR tube labels continue on next page.

Sample collection and processing labels

Microcentrifuge and PCR tube labels (continued)




Figure 2




Figure 5a




Figure 24




Figure 25

Fits tube size*	Figure	Catalog #	Material	Color	Label width A inch (mm)	Label length B inch (mm)	Secondary label C inch (mm)	Labels per row	Labels per roll	Recommended ribbon
2 mL	Fig. 25	THT-88-7425-2-SC	B-7425	White	1.000 (25.4)	0.750 (19.1)	-	2	2,000	R6413
	Fig. 2	THT-179-492-2.5-SC	B-492	White	1.000 (25.4)	1.000 (25.4)	-	1	2,500	R6413
	Fig. 2	THT-179-7425-2-SC	B-7425	White	1.000 (25.4)	1.000 (25.4)	-	1	2,000	R6413
Any top	Fig. 5a	THT-143-492-1.5-SC	B-492	White	0.375 (9.5) dia.	-	-	4	1,500	R6413
	Fig. 5a	THT-144-492-1.5-SC	B-492	White	0.500 (12.7) dia.	-	-	3	1,500	R6413
	Fig. 5a	THT-144-499-1.5-SC	B-499	White	0.500 (12.7) dia.	-	-	3	1,500	R4913
BradyPrinter i5100, i7100, Brady IP Printers and other industry printers using 3-inch core label rolls										
0.5 mL	Fig. 2	THT-257-7425-2	B-7425	White	0.825 (21.0)	0.250 (6.4)	-	1	2,000	R4302 / IP-R4302
	Fig. 2	THT-125-499-3	B-499	White	0.900 (22.9)	0.500 (12.7)	-	1	3,000	R6000 / IP-R6000
	Fig. 24	THT-126-499-3	B-499	White	0.900 (22.9)	0.500 (12.7)	0.500 (12.7) dia.	1	3,000	R6000 / IP-R6000
1.5 mL	Fig. 25	THT-152-492-3	B-492	White	1.000 (25.4)	0.375 (9.5)	-	2	3,000	R6402 / IP-R6402
	Fig. 25	THT-152-494-BL	B-494	Blue / white	1.000 (25.4)	0.375 (9.5)	-	2	3,000	R6402 / IP-R6402
	Fig. 25	THT-152-494-GN	B-494	Green / white	1.000 (25.4)	0.375 (9.5)	-	2	3,000	R6402 / IP-R6402
	Fig. 25	THT-152-494-OR	B-494	Orange / white	1.000 (25.4)	0.375 (9.5)	-	2	3,000	R6402 / IP-R6402
	Fig. 25	THT-152-494-PK	B-494	Pink / white	1.000 (25.4)	0.375 (9.5)	-	2	3,000	R6402 / IP-R6402
	Fig. 25	THT-152-494-RD	B-494	Red / white	1.000 (25.4)	0.375 (9.5)	-	2	3,000	R6402 / IP-R6402
	Fig. 25	THT-152-494-YL	B-494	Yellow / white	1.000 (25.4)	0.375 (9.5)	-	2	3,000	R6402 / IP-R6402
	Fig. 25	THT-152-499-3	B-499	White	1.000 (25.4)	0.375 (9.5)	-	2	3,000	R6000 / IP-R6000
	Fig. 25	THT-152-7425-2	B-7425	White	1.000 (25.4)	0.375 (9.5)	-	2	2,000	R4302 / IP-R4302
	Fig. 24	THT-163-492-3	B-492	White	1.000 (25.4)	0.375 (9.5)	0.375 (9.5) dia.	1	3,000	R6402 / IP-R6402
	Fig. 24	THT-163-499-3	B-499	White	1.000 (25.4)	0.375 (9.5)	0.375 (9.5) dia.	1	3,000	R6000 / IP-R6000
	Fig. 24	THT-249-7425-2	B-7425	White	1.000 (25.4)	0.375 (9.5)	0.375 (9.5) dia.	1	2,000	R4302 / IP-R4302
	Fig. 2	THT-59-492-10	B-492	White	1.000 (25.4)	0.500 (12.7)	-	1	2 rolls of 5,000	R6400 / IP-R6400
	Fig. 2	THT-59-7425-2	B-7425	White	1.000 (25.4)	0.500 (12.7)	-	1	2,000	R4302 / IP-R4302
	Fig. 24	THT-181-492-3	B-492	White	1.000 (25.4)	0.500 (12.7)	0.375 (9.5) dia.	1	3,000	R6400 / IP-R6400
	Fig. 24	THT-250-494-BL	B-494	Blue / white	1.000 (25.4)	0.500 (12.7)	0.375 (9.5) dia.	1	3,000	R6402 / IP-R6402
	Fig. 24	THT-250-494-GN	B-494	Green / white	1.000 (25.4)	0.500 (12.7)	0.375 (9.5) dia.	1	3,000	R6402 / IP-R6402
	Fig. 24	THT-250-494-OR	B-494	Orange / white	1.000 (25.4)	0.500 (12.7)	0.375 (9.5) dia.	1	3,000	R6402 / IP-R6402
	Fig. 24	THT-250-494-PK	B-494	Pink / white	1.000 (25.4)	0.500 (12.7)	0.375 (9.5) dia.	1	3,000	R6402 / IP-R6402
	Fig. 24	THT-250-494-RD	B-494	Red / white	1.000 (25.4)	0.500 (12.7)	0.375 (9.5) dia.	1	3,000	R6402 / IP-R6402
Fig. 24	THT-250-494-YL	B-494	Yellow / white	1.000 (25.4)	0.500 (12.7)	0.375 (9.5) dia.	1	3,000	R6402 / IP-R6402	
Fig. 24	THT-182-492-3	B-492	White	1.000 (25.4)	0.500 (12.7)	0.440 (11.2) dia.	1	3,000	R6400 / IP-R6400	
Fig. 24	THT-251-7425-2	B-7425	White	1.000 (25.4)	0.500 (12.7)	0.440 (11.2) dia.	1	2,000	R4302 / IP-R4302	
Fig. 2	THT-256-7425-2	B-7425	White	1.300 (33.0)	0.600 (15.2)	-	1	2,000	R4302 / IP-R4302	
2 mL	Fig. 2	THT-180-492-3	B-492	White	0.900 (22.9)	1.000 (25.4)	-	1	3,000	R6400 / IP-R6400
	Fig. 25	THT-88-7425-2	B-7425	White	1.000 (25.4)	0.750 (19.1)	-	2	2,000	R4302 / IP-R4302
	Fig. 2	THT-179-492-3	B-492	White	1.000 (25.4)	1.000 (25.4)	-	1	3,000	R6400 / IP-R6400
	Fig. 2	THT-179-494-BL	B-494	Blue / white	1.000 (25.4)	1.000 (25.4)	-	1	3,000	R6402 / IP-R6402
	Fig. 2	THT-179-494-GN	B-494	Green / white	1.000 (25.4)	1.000 (25.4)	-	1	3,000	R6402 / IP-R6402
	Fig. 2	THT-179-494-OR	B-494	Orange / white	1.000 (25.4)	1.000 (25.4)	-	1	3,000	R6402 / IP-R6402
	Fig. 2	THT-179-494-PK	B-494	Pink / white	1.000 (25.4)	1.000 (25.4)	-	1	3,000	R6402 / IP-R6402
	Fig. 2	THT-179-494-RD	B-494	Red / white	1.000 (25.4)	1.000 (25.4)	-	1	3,000	R6402 / IP-R6402
	Fig. 2	THT-179-494-YL	B-494	Yellow / white	1.000 (25.4)	1.000 (25.4)	-	1	3,000	R6402 / IP-R6402
	Fig. 2	THT-179-7425-2	B-7425	White	1.000 (25.4)	1.000 (25.4)	-	1	2,000	R4302 / IP-R4302
Any top	Fig. 5a	THT-143-492-3	B-492	White	0.375 (9.5) dia.	-	-	3	3,000	R6402 / IP-R6402
	Fig. 5a	THT-143-499-3	B-499	White	0.375 (9.5) dia.	-	-	4	3,000	R6000 / IP-R6000
	Fig. 5a	THT-144-492-3	B-492	White	0.500 (12.7) dia.	-	-	3	3,000	R6402 / IP-R6402
	Fig. 5a	THT-144-499-3	B-499	White	0.500 (12.7) dia.	-	-	3	3,000	R6000 / IP-R6000

Sample collection and processing labels

Self-laminating microcentrifuge and PCR tube labels

- B-461 – Ultra-thin polyester with clear overlamine tail that protects print and allows liquid levels to be visible




Figure 14




Figure 14b




Figure 15



Fits tube size*	Figure	Catalog #	Material	Color	Label width A inch (mm)	Label length B inch (mm)	Print-on area C inch (mm)	Labels per row	Labels per roll	Recommended ribbon
BMP®51 Label Printer cartridges										
0.5 mL	Fig. 15	M-123-461	B-461	Black on white	1.125 (28.6)	0.375 (9.5)	0.750 (19.0)	1	360	In cartridge
	Fig. 15	M-122-461	B-461	Black on white	1.125 (28.6)	0.500 (12.7)	0.750 (19.0)	1	360	In cartridge
	Fig. 15	M-124-461	B-461	Black on white	1.650 (41.9)	0.500 (12.7)	0.750 (19.0)	1	180	In cartridge
	Fig. 15	M-72-461	B-461	Black on white	1.750 (44.5)	0.500 (12.7)	0.750 (19.0)	1	180	In cartridge
1.5 mL	Fig. 15	M-126-461	B-461	Black on white	1.800 (45.7)	0.600 (15.2)	1.000 (25.4)	1	180	In cartridge
	Fig. 15	M-86-461	B-461	Black on white	2.625 (66.7)	0.600 (15.2)	1.000 (25.4)	1	120	In cartridge
2 mL	Fig. 15	M-74-461	B-461	Black on white	2.200 (55.9)	0.500 (12.7)	0.750 (19.0)	1	140	In cartridge
BradyPrinter M611 label rolls										
0.5 mL	Fig. 14	M61-202-461	B-461	White	0.500 (12.7)	1.650 (41.9)	0.750 (19.1)	1	500	M61-R4310
	Fig. 15	M61-123-461	B-461	White	1.125 (28.6)	0.375 (9.5)	0.750 (19.1)	1	500	M61-R4310
	Fig. 15	M61-122-461	B-461	White	1.125 (28.6)	0.500 (12.7)	0.750 (19.1)	1	500	M61-R4310
	Fig. 15	PTL-76-461	B-461	White	1.625 (41.3)	0.600 (15.2)	1.000 (25.4)	1	250	M61-R4310
	Fig. 15	PTL-73-461	B-461	White	1.800 (45.7)	0.375 (9.5)	0.750 (19.1)	1	500	M61-R4310
1.5 mL	Fig. 15	PTL-72-461	B-461	White	1.800 (45.7)	0.500 (12.7)	0.750 (19.1)	1	500	M61-R4310
	Fig. 14	M61-201-461	B-461	White	0.600 (15.2)	1.800 (45.7)	1.000 (25.4)	1	100	M61-R4310
	Fig. 15	PTL-108-461	B-461	White	1.500 (38.1)	0.187 (4.75)	0.700 (17.8)	1	750	M61-R4310
2 mL	Fig. 14	M61-203-461	B-461	White	1.000 (25.4)	1.750 (44.5)	1.000 (25.4)	1	100	M61-R4310
	Fig. 14	PTL-74-461	B-461	White	2.200 (55.9)	0.500 (12.7)	0.750 (19.1)	1	100	M61-R4310
	Fig. 14	PTL-79-461	B-461	White	3.750 (95.3)	0.750 (19.1)	1.000 (25.4)	1	100	M61-R4310
BradyPrinter i3300 and BBP®33 Printer label rolls										
0.5 mL	Fig. 15	B33-123-461	B-461	White	1.800 (45.7)	0.375 (9.5)	0.750 (19.1)	1	1,500	B30-R4300
	Fig. 15	B33-122-461	B-461	White	1.800 (45.7)	0.500 (12.7)	0.750 (19.1)	1	1,500	B30-R4300
	Fig. 14b	B33-127-461	B-461	White	2.200 (55.9)	0.500 (12.7)	0.750 (19.1)	3	1,500	B30-R4300
1.5 mL	Fig. 14b	B33-133-461	B-461	White	1.625 (41.3)	0.600 (15.2)	1.000 (25.4)	3	1,500	B30-R4300
BBP®12 Printer and other industry printers using 1-inch small core label rolls										
0.5 mL	Fig. 15	THT-122-461-1.5-SC	B-461	White	0.500 (12.7)	1.800 (45.7)	0.750 (19.1)	1	1,500	R4313
	Fig. 15	THT-127-461-1.5-SC	B-461	White	0.500 (12.7)	2.200 (55.9)	0.750 (19.1)	1	1,500	R4313
	Fig. 15	THT-123-461-1.5-SC	B-461	White	0.750 (19.1)	1.800 (45.7)	0.750 (19.1)	1	1,500	R4313
1.5 mL	Fig. 15	THT-133-461-1.5-SC	B-461	White	0.600 (15.2)	1.625 (41.3)	1.000 (25.4)	1	1,500	R4313
	Fig. 15	THT-153-461-1.5-SC	B-461	White	0.600 (15.2)	2.625 (66.7)	1.000 (25.4)	1	1,500	R4313
BradyPrinter i5100, i7100, Brady IP Printers and other industry printers using 3-inch core label rolls										
0.5 mL	Fig. 15	THT-123-461-3	B-461	White	0.375 (9.5)	1.800 (45.7)	0.750 (19.1)	1	3,000	R4302 / IP-R4302
	Fig. 15	THT-122-461-3	B-461	White	0.500 (12.7)	1.800 (45.7)	0.750 (19.1)	1	3,000	R4302 / IP-R4302
	Fig. 15	THT-127-461-3	B-461	White	0.500 (12.7)	2.200 (55.9)	0.750 (19.1)	1	3,000	R4300 / IP-R4300
1.5 mL	Fig. 15	THT-133-461-3	B-461	White	0.600 (15.2)	1.625 (41.3)	1.000 (25.4)	1	3,000	R4302 / IP-R4302
	Fig. 15	THT-153-461-3	B-461	White	0.600 (15.2)	2.625 (66.7)	1.000 (25.4)	1	3,000	R4300 / IP-R4300

*Tube sizes may vary by brand.

Bottle, flask and dish labels

Label options range from cost-effective paper to more durable materials with varying adhesive strengths, as well as nylon cloth that survives in harsh environments.

- B-424 – Paper label with good contrast and smear resistance
- B-430 – Clear high-performance polyester label
- B-432 – Clear polyester withstands exposure to solvents and exhibits good adhesion to metal and other surfaces
- B-449 – Removable white polypropylene for temporary labeling applications requiring clean removal. Also ideal for write-on applications
- B-488 – White polyester with moderate chemical resistance and short-term exposure to xylene, DMSO and ethanol
- B-490 – White polyester ideal for flat frozen surfaces. Adheres to frozen tubes and vials when wrapped around itself by at least 0.125 inch
- B-494 – Glossy polyester with a printed color bar that allows for quick visual differentiation between samples. Withstands liquid nitrogen, a variety of chemicals and more. Label dimensions listed below include the color bar area. Vial tops are solid color
- B-498 – Removable white vinyl cloth label for re-identification in general lab applications – can be applied directly on top of previous label
- B-499 – White nylon cloth with permanent adhesive for easy handling and conformability
- B-533 – Removable glossy white polyester for temporary labeling applications requiring clean removal
- B-7425 / B-7425J – Economical white polypropylene label with moderate chemical resistance that is ideal for labeling dry, room-temperature specimen containers. Must be applied to containers prior to being frozen





Figure 1




Figure 2




Figure 3




Figure 33b




Figure 25




Figure 82

Figure	Catalog #	Material	Color	Label width A Inch (mm)	Label length B Inch (mm)	Labels per row	Labels per roll	Recommended ribbon
BMP®21-LAB Label Printer cartridges								
Fig. 33b	M21-250-430	B-430	Black on clear	0.250 (6.4)	21 ft. (6.4 m)	1	1 roll (cont.)	In cartridge
Fig. 33b	M21-250-430-WT-CL	B-430	White on clear	0.250 (6.4)	21 ft. (6.4 m)	1	1 roll (cont.)	In cartridge
Fig. 33b	M21-375-430	B-430	Black on clear	0.375 (9.5)	21 ft. (6.4 m)	1	1 roll (cont.)	In cartridge
Fig. 33b	M21-375-430-WT-CL	B-430	White on clear	0.375 (9.5)	21 ft. (6.4 m)	1	1 roll (cont.)	In cartridge
Fig. 33b	M21-500-430	B-430	Black on clear	0.500 (12.7)	21 ft. (6.4 m)	1	1 roll (cont.)	In cartridge
Fig. 33b	M21-500-430-WT-CL	B-430	White on clear	0.500 (12.7)	21 ft. (6.4 m)	1	1 roll (cont.)	In cartridge
Fig. 33b	M21-750-499	B-499	Black on clear	0.750 (19.1)	16 ft. (4.8 m)	1	1 roll (cont.)	In cartridge
Fig. 33b	M21-750-430	B-430	Black on clear	0.750 (19.1)	21 ft. (6.4 m)	1	1 roll (cont.)	In cartridge
Fig. 33b	M21-750-430-WT-CL	B-430	White on clear	0.750 (19.1)	21 ft. (6.4 m)	1	1 roll (cont.)	In cartridge

Bottle, flask and dish labels continue on next page.

Bottle, flask and dish labels (continued)

Figure	Catalog #	Material	Color	Label width A inch (mm)	Label length B inch (mm)	Labels per row	Labels per roll	Recommended ribbon
BMP®51 Label Printer cartridges								
Fig. 82	M-130-499	B-499	Black on white	0.375 (9.5)	0.825 (21.0)	1	220	In cartridge
Fig. 82	M-131-498	B-498	Black on white	0.500 (12.7)	1.000 (25.4)	1	180	In cartridge
Fig. 82	M-131-499	B-499	Black on white	0.500 (12.7)	1.000 (25.4)	1	180	In cartridge
Fig. 2	M-17-432-CL-BK	B-432	Black on clear	1.000 (25.4)	0.500 (12.7)	1	280	In cartridge
Fig. 2	M-163-498	B-498	Black on white	1.000 (25.4)	0.500 (12.7)	1	180	In cartridge
Fig. 2	M-49-494-BL	B-494	Black on blue / white	1.000 (25.4)	1.000 (25.4)	1	260	In cartridge
Fig. 2	M-49-494-GN	B-494	Black on green / white	1.000 (25.4)	1.000 (25.4)	1	260	In cartridge
Fig. 2	M-49-494-OR	B-494	Black on orange / white	1.000 (25.4)	1.000 (25.4)	1	260	In cartridge
Fig. 2	M-49-494-PK	B-494	Black on pink / white	1.000 (25.4)	1.000 (25.4)	1	260	In cartridge
Fig. 2	M-49-494-RD	B-494	Black on red / white	1.000 (25.4)	1.000 (25.4)	1	260	In cartridge
Fig. 2	M-49-494-YL	B-494	Black on yellow / white	1.000 (25.4)	1.000 (25.4)	1	260	In cartridge
Fig. 2	M-103-498	B-498	Black on white	1.000 (25.4)	1.250 (31.8)	1	160	In cartridge
Fig. 82	M-128-498	B-498	Black on white	1.000 (25.4)	1.900 (48.3)	1	100	In cartridge
Fig. 2	M-132-499	B-499	Black on white	1.500 (38.1)	1.250 (31.8)	1	160	In cartridge
BradyPrinter M611 label rolls								
Fig. 2	M61-130-499	B-499	White	0.825 (21.0)	0.375 (9.5)	1	500	M61-R6210
Fig. 2	M61-18-7425	B-7425	White	1.000 (25.4)	0.750 (19.1)	1	250	M61-R4310
Fig. 2	M61-19-494-BL	B-494	Blue / white	1.000 (25.4)	1.000 (25.4)	1	250	M61-R6410
Fig. 2	M61-19-494-GN	B-494	Green / white	1.000 (25.4)	1.000 (25.4)	1	250	M61-R6410
Fig. 2	M61-19-494-OR	B-494	Orange / white	1.000 (25.4)	1.000 (25.4)	1	250	M61-R6410
Fig. 2	M61-19-494-PK	B-494	Pink / white	1.000 (25.4)	1.000 (25.4)	1	250	M61-R6410
Fig. 2	M61-19-494-RD	B-494	Red / white	1.000 (25.4)	1.000 (25.4)	1	250	M61-R6410
Fig. 2	M61-19-494-YL	B-494	Yellow / white	1.000 (25.4)	1.000 (25.4)	1	250	M61-R6410
Fig. 82	M61-103-498	B-498	White	1.000 (25.4)	1.250 (31.8)	1	250	M61-R6210
Fig. 2	M61-132-499	B-499	White	1.500 (38.1)	1.250 (31.8)	1	250	M61-R6210
Fig. 2	M61-121-7425	B-7425	White	1.900 (48.3)	1.500 (38.1)	1	100	M61-R4310
BradyPrinter i3300 and BBP®33 Printer label rolls								
Fig. 1	B33-68-490	B-490	White	0.500 (12.7)	0.750 (19.1)	4	2,500	B30-R4300
Fig. 1	B33-68-499	B-499	White	0.500 (12.7)	0.750 (19.1)	4	1,500	B30-R6000
Fig. 1	B33-1-498	B-498	White	0.750 (19.1)	0.250 (6.4)	4	5,000	B30-R6200
Fig. 25	B33-152-499	B-499	White	1.000 (25.4)	0.375 (9.5)	2	1,500	B30-R6000
Fig. 1	B33-5-498	B-498	White	1.000 (25.4)	0.500 (12.7)	3	5,000	B30-R6200
Fig. 25	B33-88-498	B-498	White	1.000 (25.4)	0.750 (19.1)	2	2,500	B30-R6200
Fig. 2	B33-179-494-BL	B-494	Blue / white	1.000 (25.4)	1.000 (25.4)	1	1,500	B30-R6400
Fig. 2	B33-179-494-GN	B-494	Green / white	1.000 (25.4)	1.000 (25.4)	1	1,500	B30-R6400
Fig. 2	B33-179-494-OR	B-494	Orange / white	1.000 (25.4)	1.000 (25.4)	1	1,500	B30-R6400
Fig. 2	B33-179-494-PK	B-494	Pink / white	1.000 (25.4)	1.000 (25.4)	1	1,500	B30-R6400
Fig. 2	B33-179-494-RD	B-494	Red / white	1.000 (25.4)	1.000 (25.4)	1	1,500	B30-R6400
Fig. 2	B33-179-494-YL	B-494	Yellow / white	1.000 (25.4)	1.000 (25.4)	1	1,500	B30-R6400
Fig. 2	B33-179-7425	B-7425	White	1.000 (25.4)	1.000 (25.4)	1	1,500	B30-R4300
Fig. 1	B33-89-498	B-498	White	1.000 (25.4)	1.437 (36.5)	3	2,300	B30-R6200
Fig. 2	B33-136-499	B-499	White	1.500 (38.1)	0.750 (19.1)	1	1,500	B30-R6000
Fig. 2	B33-258-7425	B-7425	White	1.500 (38.1)	1.750 (44.5)	1	900	B30-R4300
Fig. 2	B33-17-424	B-424	White	2.000 (50.8)	1.000 (25.4)	1	1,500	B30-R4300
Fig. 2	B33-17-430	B-430	Clear	2.000 (50.8)	1.000 (25.4)	1	1,500	B30-R6000
Fig. 2	B33-137-499	B-499	White	2.000 (50.8)	1.000 (25.4)	1	1,100	B30-R6000
Fig. 2	B33-137-7425	B-7425	White	2.000 (50.8)	1.000 (25.4)	1	1,500	B30-R4300
Fig. 2	B33-18-424	B-424	White	3.000 (76.2)	1.000 (25.4)	1	1,500	B30-R4300
Fig. 2	B33-18-430	B-430	Clear	3.000 (76.2)	1.000 (25.4)	1	1,500	B30-R6000
Fig. 2	B33-19-424	B-424	White	3.000 (76.2)	2.000 (50.8)	1	500	B30-R4300
BradyJet J2000 Printer label rolls								
Fig. 2	J20-179-7425J	B-7425J	White	1.000 (25.4)	1.000 (25.4)	1	1,500	J20-CMY
Fig. 3	J20C-1125-7425J	B-7425J	White	1.125 (28.6)	100 ft. (30.5 m)	1	1 roll (cont.)	J20-CMY
Fig. 2	J20-258-7425J	B-7425J	White	1.500 (38.1)	1.750 (44.5)	1	900	J20-CMY
Fig. 2	J20-137-7425J	B-7425J	White	2.000 (50.8)	1.000 (25.4)	1	1,500	J20-CMY
Fig. 2	J20-260-7425J	B-7425J	White	2.170 (55.1)	3.350 (85.1)	1	340	J20-CMY
Fig. 3	J20C-2250-7425J	B-7425J	White	2.250 (57.2)	100 ft. (30.5 m)	1	1 roll (cont.)	J20-CMY
Fig. 2	J20-261-7425J	B-7425J	White	2.950 (74.9)	4.330 (110.0)	1	270	J20-CMY
Fig. 2	J20-18-7425J	B-7425J	White	3.000 (76.2)	1.000 (25.4)	1	1,000	J20-CMY
Fig. 2	J20-262-7425J	B-7425J	White	4.130 (104.9)	6.100 (154.9)	1	200	J20-CMY

Bottle, flask and dish labels continue on next page.

General labels

Bottle, flask and dish labels (continued)




Figure	Catalog #	Material	Color	Label width A inch (mm)	Label length B inch (mm)	Labels per row	Labels per roll	Recommended ribbon
BBP®12 Printer and other industry printers using 1-inch small core label rolls								
Fig. 25	THT-140-488-1.5-SC	B-488	White	0.900 (22.9)	0.750 (19.1)	2	1,500	R4313
Fig. 2	THT-59-498-2.5-SC	B-498	White	1.000 (25.4)	0.500 (12.7)	1	2,500	R6213
Fig. 2	THT-59-533-2.5-SC	B-533	White	1.000 (25.4)	0.500 (12.7)	1	2,500	R6013
Fig. 2	THT-179-7425-2-SC	B-7425	White	1.000 (25.4)	1.000 (25.4)	1	2,000	R6413
Fig. 2	THT-199-498-3-SC	B-498	White	1.050 (26.7)	0.250 (6.4)	1	3,000	R6213
Fig. 25	THT-6-449-1.5-SC	B-449	White	1.500 (38.1)	0.750 (19.1)	2	1,500	R4313
Fig. 2	THT-136-499-1.5-SC	B-499	White	1.500 (38.1)	0.750 (19.1)	1	1,500	R4913
Fig. 2	THT-258-7425-2-SC	B-7425	White	1.750 (44.5)	1.500 (38.1)	1	2,000	R6413
Fig. 2	THT-17-424-1.5-SC	B-424	White	2.000 (50.8)	1.000 (25.4)	1	1,500	R4313
Fig. 2	THT-17-498-1.5-SC	B-498	White	2.000 (50.8)	1.000 (25.4)	1	1,500	R6213
Fig. 2	THT-137-499-1.5-SC	B-499	White	2.000 (50.8)	1.000 (25.4)	1	1,500	R4913
Fig. 2	THT-17-533-1.5-SC	B-533	White	2.000 (50.8)	1.000 (25.4)	1	1,500	R6013
Fig. 2	THT-137-7425-2-SC	B-7425	White	2.000 (50.8)	1.000 (25.4)	1	2,000	R6413
Fig. 2	THT-18-424-1.5-SC	B-424	White	3.000 (76.2)	1.000 (25.4)	1	1,500	R4313
BradyPrinter i5100, i7100, Brady IP Printers and other industry printers using 3-inch core label rolls								
Fig. 25	THT-140-488-3	B-488	White	0.900 (22.9)	0.750 (19.1)	2	3,000	R4302 / IP-R4302
Fig. 1	THT-5-488-10	B-488	White	1.000 (25.4)	0.500 (12.7)	3	10,000	R4307 / IP-R4307
Fig. 2	THT-59-498-10	B-498	White	1.000 (25.4)	0.500 (12.7)	1	2 rolls of 5,000	R6200 / IP-R6200
Fig. 2	THT-59-533-10	B-533	White	1.000 (25.4)	0.500 (12.7)	1	2 rolls of 5,000	R6000 / IP-R6000
Fig. 2	THT-179-494-BL	B-494	Blue / white	1.000 (25.4)	1.000 (25.4)	1	3,000	R6402 / IP-R6402
Fig. 2	THT-179-494-GN	B-494	Green / white	1.000 (25.4)	1.000 (25.4)	1	3,000	R6402 / IP-R6402
Fig. 2	THT-179-494-OR	B-494	Orange / white	1.000 (25.4)	1.000 (25.4)	1	3,000	R6402 / IP-R6402
Fig. 2	THT-179-494-PK	B-494	Pink / white	1.000 (25.4)	1.000 (25.4)	1	3,000	R6402 / IP-R6402
Fig. 2	THT-179-494-RD	B-494	Red / white	1.000 (25.4)	1.000 (25.4)	1	3,000	R6402 / IP-R6402
Fig. 2	THT-179-494-YL	B-494	Yellow / white	1.000 (25.4)	1.000 (25.4)	1	3,000	R6402 / IP-R6402
Fig. 2	THT-179-7425-2	B-7425	White	1.000 (25.4)	1.000 (25.4)	1	2,000	R4302 / IP-R4302
Fig. 2	THT-199-498-3	B-498	White	1.050 (26.7)	0.250 (6.4)	1	3,000	R6200 / IP-R6200
Fig. 25	THT-6-449-3	B-449	White	1.500 (38.1)	0.750 (19.1)	2	3,000	R4307 / IP-R4307
Fig. 2	THT-136-499-3	B-499	White	1.500 (38.1)	0.750 (19.1)	1	3,000	R6000 / IP-R6000
Fig. 2	THT-258-7425-2	B-7425	White	1.750 (44.5)	1.500 (38.1)	1	2,000	R4302 / IP-R4302
Fig. 2	THT-15-488-2.5	B-488	White	2.000 (50.8)	0.250 (6.4)	1	2,500	R4302 / IP-R4302
Fig. 2	THT-17-424-3	B-424	White	2.000 (50.8)	1.000 (25.4)	1	3,000	R4300 / IP-R4300
Fig. 2	THT-17-430-3	B-430	Clear	2.000 (50.8)	1.000 (25.4)	1	3,000	R6000 / IP-R6000
Fig. 2	THT-17-498-3	B-498	White	2.000 (50.8)	1.000 (25.4)	1	3,000	R6200 / IP-R6200
Fig. 2	THT-137-499-3	B-499	White	2.000 (50.8)	1.000 (25.4)	1	3,000	R6000 / IP-R6000
Fig. 2	THT-17-533-3	B-533	White	2.000 (50.8)	1.000 (25.4)	1	3,000	R6000 / IP-R6000
Fig. 2	THT-137-7425-2	B-7425	White	2.000 (50.8)	1.000 (25.4)	1	2,000	R4302 / IP-R4302
Fig. 2	THT-18-424-3	B-424	White	3.000 (76.2)	1.000 (25.4)	1	3,000	R4307 / IP-R4307
Fig. 2	THT-18-430-3	B-430	Clear	3.000 (76.2)	1.000 (25.4)	1	3,000	R6002 / IP-R6002
Fig. 2	THT-19-430-1	B-430	Clear	3.000 (76.2)	2.000 (50.8)	1	1,000	R6002 / IP-R6002
Fig. 3	THT-21-430	B-430	Clear	3.000 (76.2)	300 ft. (91.4 m)	1	1 roll (cont.)	R6002 / IP-R6002

Water-dissolvable glassware labels

Ideal for temporary applications where easy removal, re-identification or waste reduction are valued. Use in a variety of applications, including vial, tube or glassware labeling, temporary asset tracking and general labeling.

- B-403 – Water-dissolvable paper completely dissolves in warm water within 30 seconds, leaving no adhesive residue




Figure 3




Figure 33




Figure 33b



Figure	Catalog #	Material	Color	Label width A inch (mm)	Label length B inch (mm)	Labels per row	Labels per roll	Recommended ribbon
BMP®21-LAB Label Printer cartridges								
Fig. 33b	M21-750-403	B-403	Black on white	0.750 (19.1)	21 ft. (6.4 m)	1	1 roll (cont.)	In cartridge
BMP®51 Label Printer cartridges								
Fig. 33	MC-750-403	B-403	Black on white	0.750 (19.0)	25 ft. (7.6 m)	1	1 roll (cont.)	In cartridge
Fig. 33	MC-1500-403	B-403	Black on white	1.500 (38.1)	25 ft. (7.6 m)	1	1 roll (cont.)	In cartridge
BradyPrinter M611 label rolls								
Fig. 33	M61C-750-403	B-403	White	0.750 (19.0)	50 ft. (15.2 m)	1	1 roll (cont.)	M61-R4310
Fig. 33	M61C-1500-403	B-403	White	1.500 (38.1)	50 ft. (15.2 m)	1	1 roll (cont.)	M61-R4310
BradyPrinter i3300 and BBP®33 Printer label rolls								
Fig. 3	B30C-1125-403-WT	B-403	White	1.125 (28.6)	100 ft. (30.5 m)	1	1 roll (cont.)	B30-R4300
Fig. 3	B30C-2250-403-WT	B-403	White	2.250 (57.2)	100 ft. (30.5 m)	1	1 roll (cont.)	B30-R4300
Fig. 3	B30C-3000-403-WT	B-403	White	3.000 (76.2)	100 ft. (30.5 m)	1	1 roll (cont.)	B30-R4300
Fig. 3	B30C-4000-403-WT	B-403	White	4.000 (101.6)	100 ft. (30.5 m)	1	1 roll (cont.)	B30-R4300
BradyPrinter i5100, i7100, Brady IP Printers and other industry printers using 3-inch core label rolls								
Fig. 3	THT-500-403	B-403	White	0.500 (12.7)	100 ft. (30.5 m)	1	1 roll (cont.)	R4306 / IP-R4306
Fig. 3	THT-750-403	B-403	White	0.750 (19.1)	100 ft. (30.5 m)	1	1 roll (cont.)	R4306 / IP-R4306
Fig. 3	THT-1000-403	B-403	White	1.000 (25.4)	100 ft. (30.5 m)	1	1 roll (cont.)	R4306 / IP-R4306
Fig. 3	THT-1125-403-WT	B-403	White	1.125 (28.6)	300 ft. (91.4 m)	1	1 roll (cont.)	R4306 / IP-R4306
Fig. 3	THT-1500-403	B-403	White	1.500 (38.1)	100 ft. (30.5 m)	1	1 roll (cont.)	R4302 / IP-R4302
Fig. 3	THT-2000-403	B-403	White	2.000 (50.8)	100 ft. (30.5 m)	1	1 roll (cont.)	R4302 / IP-R4302
Fig. 3	THT-2250-403-WT	B-403	White	2.250 (57.2)	300 ft. (91.4 m)	1	1 roll (cont.)	R4302 / IP-R4302
Fig. 3	THT-3000-403-WT	B-403	White	3.000 (76.2)	300 ft. (91.4 m)	1	1 roll (cont.)	R4300 / IP-R4300
Fig. 3	THT-4000-403-WT	B-403	White	4.000 (101.6)	300 ft. (91.4 m)	1	1 roll (cont.)	R4307 / IP-R4307

Sterilization-indicating labels

Material indicates steam sterilization through a pre-printed header that appears when autoclaved at 250.7° F (121.5° C) for 10 minutes.

- B-7425-AC – Economical white polypropylene label with moderate chemical resistance. These labels must be applied to containers prior to being frozen – freezing may inhibit sterilization indication. Label dimensions listed below include the sterilized header area




Figure 81



Figure	Catalog #	Material	Color	Label width A inch (mm)	Label length B inch (mm)	Labels per row	Labels per roll	Recommended ribbon
BradyPrinter i3300 and BBP®33 Printer label rolls								
Fig. 81	* B33-235-7425-AC	B-7425-AC	White	1.500 (38.1)	1.000 (25.4)	1	1,250	B30-R4300
Fig. 81	* B33-137-7425-AC	B-7425-AC	White	2.000 (50.8)	1.000 (25.4)	1	1,500	B30-R4300
BradyPrinter i5100, i7100, Brady IP Printers and other industry printers using 3-inch core label rolls								
Fig. 81	* THT-235-7425-2-AC	B-7425-AC	White	1.000 (25.4)	1.500 (38.1)	1	2,000	R4306 / IP-R4306
Fig. 81	* THT-137-7425-2-AC	B-7425-AC	White	2.000 (50.8)	1.000 (25.4)	1	2,000	R4302 / IP-R4302

⚠ **WARNING:** Reproductive Harm www.P65Warnings.ca.gov

Lean lab / visual workflow labels

Lean lab and visual workflow labels identify equipment, drawers and cabinets with workflow, instructions and reminders. Bright colors make messages highly visual and clear.

- B-595 – Vinyl labels with aggressive adhesive for general lab labeling
- B-439 – Permanent vinyl designed for use in ambient conditions and limited solvent exposure
- B-2595 – Inkjet-receptive vinyl that is highly conformable for use on rough or curved surfaces
- B-2569 – Inkjet-receptive polyester ideal for smooth surfaces and general lab labeling
- B-7425J – Economical, inkjet-receptive polypropylene with moderate chemical resistance




Figure	Catalog #	Material	Color	Label width A inch (mm)	Label length B inch (mm)	Labels per roll	Recommended ribbon
BMP®21-LAB Label Printer cartridges							
Fig. 33b	M21-250-595-WT	B-595	Black on white	0.250 (6.4)	21 ft. (6.4 m)	1 roll (cont.)	In cartridge
Fig. 33b	M21-375-595-WT	B-595	Black on white	0.375 (9.5)	21 ft. (6.4 m)	1 roll (cont.)	In cartridge
Fig. 33b	M21-500-595-WT	B-595	Black on white	0.500 (12.7)	21 ft. (6.4 m)	1 roll (cont.)	In cartridge
Fig. 33b	M21-750-595-WT	B-595	Black on white	0.750 (19.1)	21 ft. (6.4 m)	1 roll (cont.)	In cartridge
Fig. 33b	M21-250-595-YL	B-595	Black on yellow	0.250 (6.4)	21 ft. (6.4 m)	1 roll (cont.)	In cartridge
Fig. 33b	M21-375-595-YL	B-595	Black on yellow	0.375 (9.5)	21 ft. (6.4 m)	1 roll (cont.)	In cartridge
Fig. 33b	M21-500-595-YL	B-595	Black on yellow	0.500 (12.7)	21 ft. (6.4 m)	1 roll (cont.)	In cartridge
Fig. 33b	M21-750-595-YL	B-595	Black on yellow	0.750 (19.1)	21 ft. (6.4 m)	1 roll (cont.)	In cartridge
Fig. 33b	M21-375-595-RD	B-595	White on red	0.375 (9.5)	21 ft. (6.4 m)	1 roll (cont.)	In cartridge
Fig. 33b	M21-500-595-RD	B-595	White on red	0.500 (12.7)	21 ft. (6.4 m)	1 roll (cont.)	In cartridge
Fig. 33b	M21-750-595-RD	B-595	White on red	0.750 (19.1)	21 ft. (6.4 m)	1 roll (cont.)	In cartridge
Fig. 33b	M21-375-595-BL	B-595	White on blue	0.375 (9.5)	21 ft. (6.4 m)	1 roll (cont.)	In cartridge
Fig. 33b	M21-500-595-BL	B-595	White on blue	0.500 (12.7)	21 ft. (6.4 m)	1 roll (cont.)	In cartridge
Fig. 33b	M21-750-595-BL	B-595	White on blue	0.750 (19.1)	21 ft. (6.4 m)	1 roll (cont.)	In cartridge
Fig. 33b	M21-375-595-GN	B-595	White on green	0.375 (9.5)	21 ft. (6.4 m)	1 roll (cont.)	In cartridge
Fig. 33b	M21-500-595-GN	B-595	White on green	0.500 (12.7)	21 ft. (6.4 m)	1 roll (cont.)	In cartridge
Fig. 33b	M21-750-595-GN	B-595	White on green	0.750 (19.1)	21 ft. (6.4 m)	1 roll (cont.)	In cartridge
BMP®51 Label Printer cartridges							
Fig. 33	MC-250-595-WT-BK	B-595	Black on white	0.250 (6.35)	25 ft. (7.6 m)	1 roll (cont.)	In cartridge
Fig. 33	MC-375-595-WT-BK	B-595	Black on white	0.375 (9.5)	25 ft. (7.6 m)	1 roll (cont.)	In cartridge
Fig. 33	MC-500-595-WT-BK	B-595	Black on white	0.500 (12.70)	25 ft. (7.6 m)	1 roll (cont.)	In cartridge
Fig. 33	MC-750-595-WT-BK	B-595	Black on white	0.750 (19.10)	25 ft. (7.6 m)	1 roll (cont.)	In cartridge
Fig. 33	MC1-1000-595-WT-BK	B-595	Black on white	1.000 (25.40)	25 ft. (7.6 m)	1 roll (cont.)	In cartridge
Fig. 33	MC-1500-595-WT-BK	B-595	Black on white	1.500 (38.10)	25 ft. (7.6 m)	1 roll (cont.)	In cartridge
Fig. 33	MC-250-595-YL-BK	B-595	Black on yellow	0.250 (6.35)	25 ft. (7.6 m)	1 roll (cont.)	In cartridge
Fig. 33	MC-375-595-YL-BK	B-595	Black on yellow	0.375 (9.5)	25 ft. (7.6 m)	1 roll (cont.)	In cartridge
Fig. 33	MC-500-595-YL-BK	B-595	Black on yellow	0.500 (12.70)	25 ft. (7.6 m)	1 roll (cont.)	In cartridge
Fig. 33	MC-750-595-YL-BK	B-595	Black on yellow	0.750 (19.10)	25 ft. (7.6 m)	1 roll (cont.)	In cartridge
Fig. 33	MC1-1000-595-YL-BK	B-595	Black on yellow	1.000 (25.40)	25 ft. (7.6 m)	1 roll (cont.)	In cartridge
Fig. 33	MC-1500-595-YL-BK	B-595	Black on yellow	1.500 (38.10)	25 ft. (7.6 m)	1 roll (cont.)	In cartridge
Fig. 33	MC-375-595-RD-WT	B-595	White on red	0.375 (9.5)	25 ft. (7.6 m)	1 roll (cont.)	In cartridge
Fig. 33	MC-500-595-RD-WT	B-595	White on red	0.500 (12.70)	25 ft. (7.6 m)	1 roll (cont.)	In cartridge
Fig. 33	MC-750-595-RD-WT	B-595	White on red	0.750 (19.10)	25 ft. (7.6 m)	1 roll (cont.)	In cartridge
Fig. 33	MC1-1000-595-RD-WT	B-595	White on red	1.000 (25.40)	25 ft. (7.6 m)	1 roll (cont.)	In cartridge
Fig. 33	MC-1500-595-RD-WT	B-595	White on red	1.500 (38.10)	25 ft. (7.6 m)	1 roll (cont.)	In cartridge

Lean lab / visual workflow labels (continued)

Figure	Catalog #	Material	Color	Label width A inch (mm)	Label length B inch (mm)	Labels per roll	Recommended ribbon
Fig. 33	MC-500-595-BL-WT	B-595	White on blue	0.500 (12.70)	25 ft. (7.6 m)	1 roll (cont.)	In cartridge
Fig. 33	MC-750-595-BL-WT	B-595	White on blue	0.750 (19.10)	25 ft. (7.6 m)	1 roll (cont.)	In cartridge
Fig. 33	MC1-1000-595-BL-WT	B-595	White on blue	1.000 (25.40)	25 ft. (7.6 m)	1 roll (cont.)	In cartridge
Fig. 33	MC-1500-595-BL-WT	B-595	White on blue	1.500 (38.10)	25 ft. (7.6 m)	1 roll (cont.)	In cartridge
Fig. 33	MC-500-595-GN-WT	B-595	White on green	0.500 (12.70)	25 ft. (7.6 m)	1 roll (cont.)	In cartridge
Fig. 33	MC-750-595-GN-WT	B-595	White on green	0.750 (19.10)	25 ft. (7.6 m)	1 roll (cont.)	In cartridge
Fig. 33	MC1-1000-595-GN-WT	B-595	White on green	1.000 (25.40)	25 ft. (7.6 m)	1 roll (cont.)	In cartridge
Fig. 33	MC-1500-595-GN-WT	B-595	White on green	1.500 (38.10)	25 ft. (7.6 m)	1 roll (cont.)	In cartridge
BradyPrinter M611 label rolls							
Fig. 33	M61C-375-439	B-439	White	0.375 (9.5)	50 ft. (15.2 m)	1 roll (cont.)	M61-R6010
Fig. 33	M61C-500-595-WT	B-595	White	0.500 (12.7)	50 ft. (15.2 m)	1 roll (cont.)	M61-R10010
Fig. 33	M61C-1000-595-WT	B-595	White	1.000 (25.4)	50 ft. (15.2 m)	1 roll (cont.)	M61-R10010
Fig. 33	M61C-2000-595-WT	B-595	White	2.000 (50.8)	50 ft. (15.2 m)	1 roll (cont.)	M61-R10010
Fig. 33	M61C-375-439-YL	B-439	Yellow	0.375 (9.5)	50 ft. (15.2 m)	1 roll (cont.)	M61-R6010
Fig. 33	M61C-500-595-YL	B-595	Yellow	0.500 (12.7)	50 ft. (15.2 m)	1 roll (cont.)	M61-R10010
Fig. 33	M61C-1000-595-YL	B-595	Yellow	1.000 (25.4)	50 ft. (15.2 m)	1 roll (cont.)	M61-R10010
Fig. 33	M61C-2000-595-YL	B-595	Yellow	2.000 (50.8)	50 ft. (15.2 m)	1 roll (cont.)	M61-R10010
Fig. 33	M61C-500-595-RD	B-595	Red	0.500 (12.7)	50 ft. (15.2 m)	1 roll (cont.)	M61-R10010-WT
Fig. 33	M61C-1000-595-RD	B-595	Red	1.000 (25.4)	50 ft. (15.2 m)	1 roll (cont.)	M61-R10010-WT
Fig. 33	M61C-2000-595-RD	B-595	Red	2.000 (50.8)	50 ft. (15.2 m)	1 roll (cont.)	M61-R10010-WT
Fig. 33	M61C-500-595-BL	B-595	Blue	0.500 (12.7)	50 ft. (15.2 m)	1 roll (cont.)	M61-R10010-WT
Fig. 33	M61C-1000-595-BL	B-595	Blue	1.000 (25.4)	50 ft. (15.2 m)	1 roll (cont.)	M61-R10010-WT
Fig. 33	M61C-2000-595-BL	B-595	Blue	2.000 (50.8)	50 ft. (15.2 m)	1 roll (cont.)	M61-R10010-WT
Fig. 33	M61C-500-595-GN	B-595	Green	0.500 (12.7)	50 ft. (15.2 m)	1 roll (cont.)	M61-R10010-WT
Fig. 33	M61C-1000-595-GN	B-595	Green	1.000 (25.4)	50 ft. (15.2 m)	1 roll (cont.)	M61-R10010-WT
Fig. 33	M61C-2000-595-GN	B-595	Green	2.000 (50.8)	50 ft. (15.2 m)	1 roll (cont.)	M61-R10010-WT
BradyPrinter i3300 and BBP®33 Printer label rolls							
Fig. 3	B30C-500-595-WT	B-595	White	0.500 (12.7)	100 ft. (30.5 m)	1 roll (cont.)	B30-R10000
Fig. 3	B30C-1125-595-WT	B-595	White	1.125 (28.6)	100 ft. (30.5 m)	1 roll (cont.)	B30-R10000
Fig. 3	B30C-2250-595-WT	B-595	White	2.250 (57.2)	100 ft. (30.5 m)	1 roll (cont.)	B30-R10000
Fig. 3	B30C-3000-595-WT	B-595	White	3.000 (76.2)	100 ft. (30.5 m)	1 roll (cont.)	B30-R10000
Fig. 3	B30C-4000-595-WT	B-595	White	4.000 (101.6)	100 ft. (30.5 m)	1 roll (cont.)	B30-R10000
Fig. 3	B30C-500-595-YL	B-595	Yellow	0.500 (12.7)	100 ft. (30.5 m)	1 roll (cont.)	B30-R10000
Fig. 3	B30C-1125-595-YL	B-595	Yellow	1.125 (28.6)	100 ft. (30.5 m)	1 roll (cont.)	B30-R10000
Fig. 3	B30C-2250-595-YL	B-595	Yellow	2.250 (57.2)	100 ft. (30.5 m)	1 roll (cont.)	B30-R10000
Fig. 3	B30C-3000-595-YL	B-595	Yellow	3.000 (76.2)	100 ft. (30.5 m)	1 roll (cont.)	B30-R10000
Fig. 3	B30C-4000-595-YL	B-595	Yellow	4.000 (101.6)	100 ft. (30.5 m)	1 roll (cont.)	B30-R10000
Fig. 3	B30C-500-595-RD	B-595	Red	0.500 (12.7)	100 ft. (30.5 m)	1 roll (cont.)	B30-R10000-WT
Fig. 3	B30C-1125-595-RD	B-595	Red	1.125 (28.6)	100 ft. (30.5 m)	1 roll (cont.)	B30-R10000-WT
Fig. 3	B30C-2250-595-RD	B-595	Red	2.250 (57.2)	100 ft. (30.5 m)	1 roll (cont.)	B30-R10000-WT
Fig. 3	B30C-4000-595-RD	B-595	Red	4.000 (101.6)	100 ft. (30.5 m)	1 roll (cont.)	B30-R10000-WT
Fig. 3	B30C-500-595-GN	B-595	Green	0.500 (12.7)	100 ft. (30.5 m)	1 roll (cont.)	B30-R10000-WT
Fig. 3	B30C-1125-595-GN	B-595	Green	1.125 (28.6)	100 ft. (30.5 m)	1 roll (cont.)	B30-R10000-WT
Fig. 3	B30C-2250-595-GN	B-595	Green	2.250 (57.2)	100 ft. (30.5 m)	1 roll (cont.)	B30-R10000-WT
Fig. 3	B30C-4000-595-GN	B-595	Green	4.000 (101.6)	100 ft. (30.5 m)	1 roll (cont.)	B30-R10000-WT
BradyJet J2000 Printer label rolls							
Fig. 3	J20C-1125-2569	B-2569	White	1.125 (28.57)	100 ft. (30.5 m)	1 roll (cont.)	J20-CMY
Fig. 3	J20C-1125-2595	B-2595	White	1.125 (28.57)	100 ft. (30.5 m)	1 roll (cont.)	J20-CMY
Fig. 3	J20C-1125-7425J	B-7425J	White	1.125 (28.6)	100 ft. (30.5 m)	1 roll (cont.)	J20-CMY
Fig. 2	J20-260-2569	B-2569	White	2.170 (55.11)	3.350 (85.09)	340	J20-CMY
Fig. 2	J20-260-2595	B-2595	White	2.170 (55.11)	3.350 (85.09)	340	J20-CMY
Fig. 2	J20-260-7425J	B-7425J	White	2.170 (55.1)	3.350 (85.1)	340	J20-CMY
Fig. 3	J20C-2250-2569	B-2569	White	2.250 (57.15)	100 ft. (30.5 m)	1 roll (cont.)	J20-CMY
Fig. 3	J20C-2250-2595	B-2595	White	2.250 (57.15)	100 ft. (30.5 m)	1 roll (cont.)	J20-CMY
Fig. 3	J20C-2250-7425J	B-7425J	White	2.250 (57.15)	100 ft. (30.5 m)	1 roll (cont.)	J20-CMY
Fig. 2	J20-261-2569	B-2569	White	2.950 (74.93)	4.330 (109.98)	270	J20-CMY
Fig. 2	J20-261-2595	B-2595	White	2.950 (74.93)	4.330 (109.98)	270	J20-CMY
Fig. 2	J20-261-7425J	B-7425J	White	2.950 (74.93)	4.330 (109.98)	270	J20-CMY
Fig. 3	J20C-4000-2569	B-2569	White	4.000 (101.6)	100 ft. (30.5 m)	1 roll (cont.)	J20-CMY
Fig. 3	J20C-4000-2595	B-2595	White	4.000 (101.6)	100 ft. (30.5 m)	1 roll (cont.)	J20-CMY
Fig. 2	J20-262-2569	B-2569	White	4.130 (104.90)	6.100 (154.94)	200	J20-CMY
Fig. 2	J20-262-2595	B-2595	White	4.130 (104.90)	6.100 (154.94)	200	J20-CMY
Fig. 2	J20-262-7425J	B-7425J	White	4.130 (104.90)	6.100 (154.94)	200	J20-CMY

Globally Harmonized System (GHS) Labels

Brady carries blank and pre-printed labels to clearly communicate hazards and help keep you compliant with OSHA's Hazard Communication Standards.

- B-569 – Ultra-smooth polyester with permanent adhesive withstands mild chemicals, oil and water
- B-7569 – Vinyl labels with aggressive adhesive withstand exposure to moisture, chemicals and handling
- B-2595 – Inkjet-receptive vinyl that is highly conformable for use on rough or curved surfaces
- B-7425J – Inkjet-receptive polypropylene with moderate chemical resistance




Figure 2




Figure 3



Figure	Catalog #	Material	Color	Label width A inch (mm)	Label length B inch (mm)	Labels per roll	Recommended ribbon
BradyPrinter i3300 and BBP®33 Printer label rolls							
	B30-260-7569-CLP1	B-7569	Red / white	2.170 (55.0)	3.350 (85.0)	340	B30-R10000
	B30-261-7569-CLP1	B-7569	Red / white	2.950 (75.0)	4.330 (110.0)	270	B30-R10000
	B30-262-7569-CLP1	B-7569	Red / white	4.130 (105.0)	6.100 (155.0)	200	B30-R10000
	B30-260-7569-CLP2	B-7569	Red / white	2.170 (55.0)	3.350 (85.0)	340	B30-R10000
	B30-261-7569-CLP2	B-7569	Red / white	2.950 (75.0)	4.330 (110.0)	270	B30-R10000
	B30-262-7569-CLP2	B-7569	Red / white	4.130 (105.0)	6.100 (155.0)	200	B30-R10000
	B30-260-7569-CLP3A	B-7569	Red / white	2.170 (55.0)	3.350 (85.0)	340	B30-R10000
	B30-261-7569-CLP3A	B-7569	Red / white	2.950 (75.0)	4.330 (110.0)	270	B30-R10000
	B30-262-7569-CLP3A	B-7569	Red / white	4.130 (105.0)	6.100 (155.0)	200	B30-R10000
	B30-260-7569-CLP3B	B-7569	Red / white	2.170 (55.0)	3.350 (85.0)	340	B30-R10000
	B30-261-7569-CLP3B	B-7569	Red / white	2.950 (75.0)	4.330 (110.0)	270	B30-R10000
	B30-262-7569-CLP3B	B-7569	Red / white	4.130 (105.0)	6.100 (155.0)	200	B30-R10000
	B30-260-7569-CLP4A	B-7569	Red / white	2.170 (55.0)	3.350 (85.0)	340	B30-R10000
	B30-261-7569-CLP4A	B-7569	Red / white	2.950 (75.0)	4.330 (110.0)	270	B30-R10000
	B30-262-7569-CLP4A	B-7569	Red / white	4.130 (105.0)	6.100 (155.0)	200	B30-R10000
	B30-260-569-CLP4B	B-569	Red / white	2.170 (55.0)	3.350 (85.0)	340	B30-R10000
	B30-260-7569-CLP4B	B-7569	Red / white	2.170 (55.0)	3.350 (85.0)	340	B30-R10000
	B30-261-569-CLP4B	B-569	Red / white	2.950 (75.0)	4.330 (110.0)	270	B30-R10000
	B30-261-7569-CLP4B	B-7569	Red / white	2.950 (75.0)	4.330 (110.0)	270	B30-R10000
	B30-262-569-CLP4B	B-569	Red / white	4.130 (105.0)	6.100 (155.0)	200	B30-R10000
	B30-262-7569-CLP4B	B-7569	Red / white	4.130 (105.0)	6.100 (155.0)	200	B30-R10000
	B30-18-7569-CLP4C	B-7569	Red / white	1.000 (25.4)	3.000 (76.2)	1,000	B30-R10000
BradyJet J2000 Printer label rolls							
Fig. 3	J20C-1125-7425J	B-7425J	White	1.125 (28.57)	100 ft. (30.5 m)	1 roll (cont.)	J20-CMY
Fig. 3	J20C-1125-2595	B-2595	White	1.125 (28.57)	100 ft. (30.5 m)	1 roll (cont.)	J20-CMY
Fig. 2	J20-260-7425J	B-7425J	White	2.170 (55.11)	3.350 (85.09)	340	J20-CMY
Fig. 2	J20-260-2595	B-2595	White	2.170 (55.11)	3.350 (85.09)	340	J20-CMY
Fig. 3	J20C-2250-7425J	B-7425J	White	2.250 (57.15)	100 ft. (30.5 m)	1 roll (cont.)	J20-CMY
Fig. 3	J20C-2250-2595	B-2595	White	2.250 (57.15)	100 ft. (30.5 m)	1 roll (cont.)	J20-CMY
Fig. 2	J20-261-7425J	B-7425J	White	2.950 (74.93)	4.330 (109.98)	270	J20-CMY
Fig. 2	J20-261-2595	B-2595	White	2.950 (74.93)	4.330 (109.98)	270	J20-CMY
Fig. 2	J20-18-7425J	B-7425J	White	3.000 (76.20)	1.000 (25.40)	1,000	J20-CMY
Fig. 2	J20-18-2595	B-2595	White	3.000 (76.20)	1.000 (25.40)	1,000	J20-CMY
Fig. 3	J20C-4000-7425J	B-7425J	White	4.000 (101.6)	100 ft. (30.5 m)	1 roll (cont.)	J20-CMY
Fig. 3	J20C-4000-2595	B-2595	White	4.000 (101.6)	100 ft. (30.5 m)	1 roll (cont.)	J20-CMY
Fig. 2	J20-262-7425J	B-7425J	White	4.130 (104.90)	6.100 (154.94)	200	J20-CMY
Fig. 2	J20-262-2595	B-2595	White	4.130 (104.90)	6.100 (154.94)	200	J20-CMY

Die-cut GHS Pictogram Labels

Use these Die-cut GHS Pictogram Labels to alert employees to the potential dangers of the chemicals in your facility.



Explosion hazard

Explosives, self-reactives, organic peroxides

Catalog #	Size	Material	Qty.
118815	0.625 in.	B-7569 Vinyl	39 / card
118869	0.625 in.	B-235 Paper	39 / card
118816	1 in.	B-7569 Vinyl	16 / card
118870	1 in.	B-235 Paper	16 / card
118817	2 in.	B-7569 Vinyl	4 / card
118871	2 in.	B-235 Paper	4 / card



Flame

Flammable gases, liquids and solids, pyrophorics, self-reactives, self-heating

Catalog #	Size	Material	Qty.
118821	0.625 in.	B-7569 Vinyl	39 / card
118875	0.625 in.	B-235 Paper	39 / card
118822	1 in.	B-7569 Vinyl	16 / card
118876	1 in.	B-235 Paper	16 / card
118823	2 in.	B-7569 Vinyl	4 / card
118877	2 in.	B-235 Paper	4 / card



Flame over circle

Oxidizers, organic peroxides

Catalog #	Size	Material	Qty.
118827	0.625 in.	B-7569 Vinyl	39 / card
118881	0.625 in.	B-235 Paper	39 / card
118828	1 in.	B-7569 Vinyl	16 / card
118882	1 in.	B-235 Paper	16 / card
118829	2 in.	B-7569 Vinyl	4 / card
118883	2 in.	B-235 Paper	4 / card



Gas cylinder

Compressed gases, liquefied gases, dissolved gases

Catalog #	Size	Material	Qty.
118833	0.625 in.	B-7569 Vinyl	39 / card
118887	0.625 in.	B-235 Paper	39 / card
118834	1 in.	B-7569 Vinyl	16 / card
118888	1 in.	B-235 Paper	16 / card
118835	2 in.	B-7569 Vinyl	4 / card
118889	2 in.	B-235 Paper	4 / card



Corrosion

Corrosives to metal, skin corrosion, eye damage

Catalog #	Size	Material	Qty.
118839	0.625 in.	B-7569 Vinyl	39 / card
118893	0.625 in.	B-235 Paper	39 / card
118840	1 in.	B-7569 Vinyl	16 / card
118894	1 in.	B-235 Paper	16 / card
118841	2 in.	B-7569 Vinyl	4 / card
118895	2 in.	B-235 Paper	4 / card



Skull and crossbones

Acute toxicity (severe, fatal)

Catalog #	Size	Material	Qty.
118845	0.625 in.	B-7569 Vinyl	39 / card
118899	0.625 in.	B-235 Paper	39 / card
118846	1 in.	B-7569 Vinyl	16 / card
118900	1 in.	B-235 Paper	16 / card
118847	2 in.	B-7569 Vinyl	4 / card
118901	2 in.	B-235 Paper	4 / card



Health hazard

Carcinogens, respiratory sensitizers, reproductive toxicity, target organ toxicity, germ cell mutagens

Catalog #	Size	Material	Qty.
118851	0.625 in.	B-7569 Vinyl	39 / card
118905	0.625 in.	B-235 Paper	39 / card
118852	1 in.	B-7569 Vinyl	16 / card
118906	1 in.	B-235 Paper	16 / card
118853	2 in.	B-7569 Vinyl	4 / card
118907	2 in.	B-235 Paper	4 / card



Exclamation mark

Irritant, dermal sensitizer, acute toxicity (harmful)

Catalog #	Size	Material	Qty.
118857	0.625 in.	B-7569 Vinyl	39 / card
118911	0.625 in.	B-235 Paper	39 / card
118858	1 in.	B-7569 Vinyl	16 / card
118912	1 in.	B-235 Paper	16 / card
118859	2 in.	B-7569 Vinyl	4 / card
118913	2 in.	B-235 Paper	4 / card



Environment

Aquatic toxicity

Catalog #	Size	Material	Qty.
118863	0.625 in.	B-7569 Vinyl	39 / card
118917	0.625 in.	B-235 Paper	39 / card
118864	1 in.	B-7569 Vinyl	16 / card
118918	1 in.	B-235 Paper	16 / card
118865	2 in.	B-7569 Vinyl	4 / card
118919	2 in.	B-235 Paper	4 / card

Biohazard identification solutions

With the increasing need for emergency preparedness, your facility must be ready to effectively identify and warn of biohazard materials and processes in order to keep patients, employees and the public safe. The most essential step toward safety is the clear and instantly readable identification of all biohazard areas and waste products.

Critical regulations for biohazard identification

OSHA 29 CFR 1910.1030 Bloodborne Pathogens	Labels are required for containers and shipments. Employers to post biohazard signs at the entrance to work areas.
DOT 49 CFR 172.400(b) Labeling Requirements	Labeling is required for infectious substances that meet one or more hazard classifications.
UN2814 Infectious Substances	Ebola virus defined as Category A infectious substance.



Non-Adhesive Biohazard Barricade Tape

- Non-adhesive barricade tape to warn of biohazard areas

Catalog #	Size	Material	Qty.
91447	3 in. x 200 ft.	Polyethylene	1



Toughstripe® Floor Marking Tape

- Highly durable self-sticking floor marking tape to warn of biohazard areas

Catalog #	Size	Material	Qty.
132434	3 in. x 100 ft.	Polyester	1



Toughstripe Floor Marking Sign

- Black areas of floor sign are coated with anti-skid treads that help to prevent slips, trips and falls

Catalog #	Size	Material	Qty.
104506	14 in. x 18 in.	Polyester	1



Biohazard Tape and Labels

- Self-sticking vinyl is ideal for clearly identifying biohazard areas or containers

Catalog #	Size	Material	Qty.
88051	2 in. x 108 ft.	Vinyl	1 roll
LH625E	4 in. x 2.875 in.	Vinyl	25/pk



Waste Disposal Bags

- 3-mil polyethylene bag for disposal of biohazard waste

Catalog #	Size (H x W x D)	Qty.
145586	16 in. x 14 in. x 36 in.	100
145587	30 in. x 14 in. x 39 in.	100



Exposure Control Binder

- Highly visible binder for storing up to 300 pages of your exposure control plan

Catalog #	Size	Qty.
BH2023	1.5-in. ring binder	1



Create custom laboratory signs at BradyID.com

We understand that sometimes you need a sign that is unique to your own lab. With Brady's wide variety of custom sign options, you can create a sign with the exact message you need in the material and size you want.

Laboratory Safety Signs

Brady Safety Signs are used to identify equipment, rooms and materials that contain, or are contaminated with, hazardous agents.



62550	7 in. x 10 in.	FB
62812	10 in. x 14 in.	FB
21800	7 in. x 10 in.	PL
26575	10 in. x 14 in.	PL
83928	3.5 in. x 5 in.	SS
83753	7 in. x 10 in.	SS
46753	10 in. x 14 in.	AL



69768	10 in. x 14 in.	FB
25742	7 in. x 10 in.	PL
25743	10 in. x 14 in.	PL
89170	3.5 in. x 5 in.	SS
85553	7 in. x 10 in.	SS
85554	10 in. x 14 in.	SS



126649	7 in. x 10 in.	PL
126652	10 in. x 14 in.	PL
126650	7 in. x 10 in.	SS
126653	10 in. x 14 in.	SS
126648	7 in. x 10 in.	AL
126651	10 in. x 14 in.	AL



126625	10 in. x 7 in.	PL
126628	14 in. x 10 in.	PL
126626	10 in. x 7 in.	SS
126629	14 in. x 10 in.	SS
126624	10 in. x 7 in.	AL
126627	14 in. x 10 in.	AL




126643	10 in. x 7 in.	PL
126646	14 in. x 10 in.	PL
126644	10 in. x 7 in.	SS
126647	14 in. x 10 in.	SS
126642	10 in. x 7 in.	AL
126645	14 in. x 10 in.	AL



126637	10 in. x 7 in.	PL
126640	14 in. x 10 in.	PL
126638	10 in. x 7 in.	SS
126641	14 in. x 10 in.	SS
126636	10 in. x 7 in.	AL
126639	14 in. x 10 in.	AL



126631	10 in. x 7 in.	PL
126634	14 in. x 10 in.	PL
126632	10 in. x 7 in.	SS
126635	14 in. x 10 in.	SS
126630	10 in. x 7 in.	AL
126633	14 in. x 10 in.	AL



69813	14 in. x 10 in.	FB
89158	10 in. x 7 in.	SS
89159	14 in. x 10 in.	SS



126680	7 in. x 10 in.	SS
126683	10 in. x 14 in.	SS
126679	7 in. x 10 in.	PL
126682	10 in. x 14 in.	PL
126678	7 in. x 10 in.	AL
126681	10 in. x 14 in.	AL



126668	7 in. x 10 in.	SS
126671	10 in. x 14 in.	SS
126667	7 in. x 10 in.	PL
126670	10 in. x 14 in.	PL
126666	7 in. x 10 in.	AL
126669	10 in. x 14 in.	AL



126674	7 in. x 10 in.	SS
126677	10 in. x 14 in.	SS
126673	7 in. x 10 in.	PL
126676	10 in. x 14 in.	PL
126672	7 in. x 10 in.	AL
126675	10 in. x 14 in.	AL



122820	7 in. x 10 in.	FB
70240	10 in. x 14 in.	FB
84351	7 in. x 10 in.	SS
84352	10 in. x 14 in.	SS
22311	7 in. x 10 in.	PL
22312	10 in. x 14 in.	PL
122406	7 in. x 10 in.	AL
40876	10 in. x 14 in.	AL



126158	7 in. x 10 in.	SS
126161	10 in. x 14 in.	SS
126157	7 in. x 10 in.	PL
126160	10 in. x 14 in.	PL
126159	10 in. x 14 in.	AL



70258	10 in. x 14 in.	FB
122650	7 in. x 10 in.	SS
84368	10 in. x 14 in.	SS
22323	7 in. x 10 in.	PL
22324	10 in. x 14 in.	PL
40887	7 in. x 10 in.	AL
40888	10 in. x 14 in.	AL



69097	10 in. x 14 in.	FB
84054	7 in. x 10 in.	SS
84055	10 in. x 14 in.	SS
22082	10 in. x 14 in.	PL
40645	7 in. x 10 in.	AL
40646	10 in. x 14 in.	AL



75639	10 in. x 14 in.	FB
85167	7 in. x 10 in.	SS
85168	10 in. x 14 in.	SS
25654	7 in. x 10 in.	PL
25655	10 in. x 14 in.	PL
43231	7 in. x 10 in.	AL
43232	10 in. x 14 in.	AL



73436	14 in. x 20 in.	FB
85165	7 in. x 10 in.	SS
25652	7 in. x 10 in.	PL
25653	10 in. x 14 in.	PL
43230	10 in. x 14 in.	AL



70350	7 in. x 17 in.	FB
25898	10 in. x 14 in.	PL
43475	10 in. x 14 in.	AL



122826	7 in. x 10 in.	FB
84397	7 in. x 10 in.	SS
84398	10 in. x 14 in.	SS
122519	7 in. x 10 in.	PL
25441	10 in. x 14 in.	PL
43017	7 in. x 10 in.	AL
43018	10 in. x 14 in.	AL



70582	10 in. x 14 in.	FB
88627	7 in. x 10 in.	SS
88628	10 in. x 14 in.	SS
25222	7 in. x 10 in.	PL
25223	10 in. x 14 in.	PL
42799	7 in. x 10 in.	AL
42800	10 in. x 14 in.	AL



84463	7 in. x 10 in.	SS
22355	7 in. x 10 in.	PL
22356	10 in. x 14 in.	PL



84469	7 in. x 10 in.	SS
22362	10 in. x 14 in.	PL
40925	7 in. x 10 in.	AL



129158	7 in. x 10 in.	SS
129161	10 in. x 14 in.	SS
129157	7 in. x 10 in.	PL
129160	10 in. x 14 in.	PL
129156	7 in. x 10 in.	AL
129159	10 in. x 14 in.	AL



88543	7 in. x 10 in.	SS
88544	10 in. x 14 in.	SS
25194	7 in. x 10 in.	PL
25195	10 in. x 14 in.	PL
42771	7 in. x 10 in.	AL
42772	10 in. x 14 in.	AL



47056	7 in. x 10 in.	FB
47094	10 in. x 14 in.	FB
84973	7 in. x 10 in.	SS
22579	7 in. x 10 in.	PL
22580	10 in. x 14 in.	PL
41143	7 in. x 10 in.	AL



126052	10 in. x 14 in.	PL
--------	-----------------	----



74455	10 in. x 14 in.	FB
87782	3.5 in. x 5 in.	SS
85097	7 in. x 10 in.	SS
85098	10 in. x 14 in.	SS
25610	7 in. x 10 in.	PL
25611	10 in. x 14 in.	PL
43187	7 in. x 10 in.	AL
43188	10 in. x 14 in.	AL



30673	7 in. x 10 in.	SS
30670	10 in. x 14 in.	AL



30674	7 in. x 10 in.	SS
30671	10 in. x 14 in.	AL



90933	5 in. x 14 in.	BGPL
-------	----------------	------

Absorbents

Brady SPC absorbents

Be prepared for spills and accidental leaks with absorbent products designed for laboratory work areas. Brady SPC spill control products offer immediate absorption of biohazards and chemicals, therefore limiting your exposure to hazardous conditions while helping you meet OSHA and EPA CFR 1910.40, 112.261 regulations.



High-Visibility Safety Absorbents

- Highly absorbent for all-purpose laboratory applications
- Signage and an absorbent in one unique product

Catalog #	Type	Size	Qty.
With safety print			
CH100	Heavy-wt. pad*	15 in. x 19 in.	100
CH200	Light-wt. pad*	15 in. x 19 in.	200
MXP-C20	Light-wt. pad*	15 in. x 19 in.	5
CH303	Medium-wt. roll*	30 in. x 300 ft.	1
CH30DP	Heavy-wt. roll*	30 in. x 150 ft.	1
CH15DP	Heavy-wt. roll*	15 in. x 150 ft.	1
Without safety print			
CH1212	Medium-wt. pad	12 in. x 12 in.	200
CH1212-50	Medium-wt. pad	12 in. x 12 in.	50

*Perforated



MRO Plus® Absorbents

- All-purpose absorbent for water-, petroleum- and chemical-based fluids
- Perforated to help you reduce waste by using only what you need

Catalog #	Type	Size	Qty.
MRO100	Heavy-wt. pad*	15 in. x 19 in.	100
MRO300	Medium-wt. pad*	15 in. x 19 in.	100
MRO15P	Heavy-wt. roll*	15 in. x 150 ft.	1
MRO350DP	Medium-wt. roll*	30 in. x 150 ft.	1

*Perforated



Universal Plus Absorbents

- Dimpled and perforated green universal absorbent is chemically inert and ideal for aggressive fluids
- 3-ply construction makes it stronger, more durable and reduces lint

Catalog #	Type	Size	Qty.
SPCUN1719	Medium-wt. pad*	15 in. x 19 in.	100
SPCUN1113	Medium-wt. pad	12 in. x 12 in.	200
SPCUN1113-50	Medium-wt. pad	12 in. x 12 in.	50
SRPH200	Heavy-wt. pad	7.5 in. x 10 in.	200
SRPH100	Heavy-wt. pad	7.5 in. x 10 in.	100
SPCUN3800	Medium-wt. roll*	30 in. x 150 ft.	1
SPCUN1019	Medium-wt. roll*	15 in. x 150 ft.	2

*Perforated



Adhesive Absorbent Mats

- Adheres to any surface, reducing the likelihood of slips, trips and falls
- Removes without leaving any residue

Catalog #	Type	Size	Qty.
AD15100	Heavy-wt. roll	15 in. x 100 ft.	1
AD30100	Heavy-wt. roll	30 in. x 100 ft.	1



Spill Response Plus Universal Plus Absorbents

- Small size makes it ideal for use anywhere in your lab
- Chemically inert making it ideal for aggressive fluids

Catalog #	Type	Size	Qty.
SRPH100	Heavy-wt. pad	7.5 in. x 10 in.	100
SRPH200	Heavy-wt. pad	7.5 in. x 10 in.	200
SRPH75P	Heavy-wt. roll*	7.5 in. x 50 in.	1

*Perforated



HAZWIK® SOCs and Pillows

- SOCs easily mold around spills, while pillows are ideal for extra protection during packaging and transport

Catalog #	Type	Size	Qty.
HAZ124	SOC	3 in. x 12 ft.	4
HAZ412	SOC	3 in. x 4 ft.	12
HAZ1818	Pillow	18 in. x 18 in.	16
HAZ1818-2	Pillow	18 in. x 18 in.	8
HAZ99	Pillow	9 in. x 9 in.	32



Economy Spill Kit

- Absorbs up to 5 gallons per kit
- Dimensions: 20 in. L x 16 in. W x 4 in. H
- Contents:
 - 10 pads (15 in. x 19 in.)
 - 2 SOC's (3 in. x 4 ft.)
 - 1 disposal bag
 - 1 pair nitrile gloves
 - 1 instruction sheet

Catalog #	Type	Qty.
SKH-PP†	Chemical	1



2-Gallon Spill Response Plus Spill Kit

- Absorbs up to 4.5 gallons per kit
- Dimensions: 9.75 in. dia. x 10.5 in. H
- Contents:
 - 20 pads (7.5 in. x 10 in.)
 - 1 disposal bag
 - 1 pair nitrile gloves
 - 1 instruction sheet

Catalog #	Type	Qty.
SKH-SRP	Chemical	1



Attack Pac Spill Kit

- Absorbs up to 7 gallons per kit
- Dimensions: 20 in. L x 16 in. W x 4 in. H
- Contents:
 - 15 pads (15 in. x 19 in.)
 - 3 SOC's (3 in. x 4 ft.)
 - 1 disposal bag
 - 1 pair goggles
 - 1 pair nitrile gloves
 - 1 instruction sheet

Catalog #	Type	Qty.
SKH-ATK	Chemical	4



6.5-Gallon Spill Bucket

- Absorbs up to 10 gallons per kit
- Dimensions: 14 in. dia. x 18 in. H
- Contents:
 - 10 pads (15 in. x 19 in.)
 - 3 SOC's (3 in. x 4 ft.)
 - 2 disposal bags
 - 1 pair nitrile gloves
 - 1 instruction sheet

Catalog #	Type	Qty.
SKH-BKT	Chemical	1



Emergency Response Kit

- Absorbs up to 15 gallons per kit
- Dimensions: 25.5 in. L x 17 in. W x 21.5 in. H
- Contents:
 - 32 pads (15 in. x 19 in.)
 - 5 SOC's (3 in. x 12 ft.)
 - 1 pillow (17 in. x 19 in.)
 - 2 disposal bags
 - 1 pair goggles
 - 1 pair nitrile gloves
 - 1 instruction sheet

Catalog #	Type	Qty.
SKH-CFB	Chemical	1



20-Gallon Lab Pack

- Absorbs up to 17 gallons per kit
- Dimensions: 21 in. dia. x 18 in. H
- Contents:
 - 12 pads (15 in. x 19 in.)
 - 3 SOC's (3 in. x 12 ft.)
 - 2 pillows (18 in. x 18 in.)
 - 3 disposal bags
 - 1 pair goggles
 - 1 pair nitrile gloves
 - Emergency Response Guidebook

Catalog #	Type	Qty.
SKH-20	Chemical	1



Slikstopper Drain Seal

- Flexible and chemical-resistant PVC tightly seals drains
- Clean up with soap and water for reuse

Catalog #	Description	Size	Qty.
PVC24	Slikstopper	24 in. x 24 in.	1
PVC18	Slikstopper	18 in. x 18 in.	1



Spill Control Cart

- Secure storage with locking handle (2 keys included), one fixed interior shelf and 3-inch backstop on top surface

Catalog #	Size	Qty.
SC-Cart	36 in. x 18 in. x 38.375 in.	1

Absorbents sold separately.



Spill Control Cabinet

- 2 fully adjustable shelves included
- Durable construction and wall mountable
- Clear window allows for easy inventory

Catalog #	Size	Qty.
SC-Cabinet	19.875 in. x 14.25 in. x 32.75 in.	1

Absorbents sold separately.

† ⚠️ **WARNING:** Cancer and Reproductive Harm www.P65Warnings.ca.gov

Labels by size index

BMP®21-LAB Printer labels by size



Label width (A)	Label height (B)	Material	Catalog #	Pg
0.250 (6.4)	21 ft. (6.4 m)	B-430	M21-250-430	28
0.250 (6.4)	21 ft. (6.4 m)	B-430	M21-250-430-WT-CL	52
0.250 (6.4)	21 ft. (6.4 m)	B-595	M21-250-595-WT	56
0.250 (6.4)	21 ft. (6.4 m)	B-595	M21-250-595-YL	56
0.375 (9.5)	16 ft. (4.8 m)	B-499	M21-375-499	46
0.375 (9.5)	21 ft. (6.4 m)	B-430	M21-375-430	28
0.375 (9.5)	21 ft. (6.4 m)	B-430	M21-375-430-WT-CL	52
0.375 (9.5)	21 ft. (6.4 m)	B-488	M21-375-488	30
0.375 (9.5)	21 ft. (6.4 m)	B-595	M21-375-595-WT	56
0.375 (9.5)	21 ft. (6.4 m)	B-595	M21-375-595-YL	56
0.375 (9.5)	21 ft. (6.4 m)	B-595	M21-375-595-BL	56
0.375 (9.5)	21 ft. (6.4 m)	B-595	M21-375-595-GN	56
0.375 (9.5)	21 ft. (6.4 m)	B-595	M21-375-595-RD	56
0.375 (9.5)	21 ft. (6.4 m)	B-7425	M21-375-7425	28
0.375 (9.5)	21 ft. (6.4 m)	B-461	M21-375-461	34
0.500 (12.7)	16 ft. (4.8 m)	B-499	M21-500-499	32
0.500 (12.7)	21 ft. (6.4 m)	B-430	M21-500-430	26
0.500 (12.7)	21 ft. (6.4 m)	B-430	M21-500-430-WT-CL	52
0.500 (12.7)	21 ft. (6.4 m)	B-488	M21-500-488	30
0.500 (12.7)	21 ft. (6.4 m)	B-595	M21-500-595-WT	56

Label width (A)	Label height (B)	Material	Catalog #	Pg
0.500 (12.7)	21 ft. (6.4 m)	B-595	M21-500-595-YL	56
0.500 (12.7)	21 ft. (6.4 m)	B-595	M21-500-595-BL	56
0.500 (12.7)	21 ft. (6.4 m)	B-595	M21-500-595-GN	56
0.500 (12.7)	21 ft. (6.4 m)	B-595	M21-500-595-RD	56
0.500 (12.7)	21 ft. (6.4 m)	B-7425	M21-500-7425	26
0.500 (12.7)	21 ft. (6.4 m)	B-461	M21-500-461	34
0.750 (19.1)	16 ft. (4.8 m)	B-499	M21-750-499	32
0.750 (19.1)	21 ft. (6.4 m)	B-403	M21-750-403	55
0.750 (19.1)	21 ft. (6.4 m)	B-430	M21-750-430	26
0.750 (19.1)	21 ft. (6.4 m)	B-430	M21-750-430-WT-CL	52
0.750 (19.1)	21 ft. (6.4 m)	B-488	M21-750-488	30
0.750 (19.1)	21 ft. (6.4 m)	B-595	M21-750-595-WT	56
0.750 (19.1)	21 ft. (6.4 m)	B-595	M21-750-595-YL	56
0.750 (19.1)	21 ft. (6.4 m)	B-595	M21-750-595-BL	56
0.750 (19.1)	21 ft. (6.4 m)	B-595	M21-750-595-GN	56
0.750 (19.1)	21 ft. (6.4 m)	B-595	M21-750-595-RD	56
0.750 (19.1)	21 ft. (6.4 m)	B-7425	M21-750-7425	26
0.750 (19.1)	21 ft. (6.4 m)	B-461	M21-750-461	34

Labels by size index

BMP®51 Printer labels by size



Label width (A)	Label height (B)	Secondary label (C)	Material	Catalog #	Pg
0.250 (6.4)	1.100 (27.9)	-	B-427	M-135-427-AW	33
0.250 (6.35)	25 ft. (7.6 m)	-	B-595	MC-250-595-WT-BK	56
0.250 (6.35)	25 ft. (7.6 m)	-	B-595	MC-250-595-YL-BK	56
0.375 (9.5)	0.825 (21.0)	-	B-492	M-130-492	34
0.375 (9.5)	0.825 (21.0)	-	B-499	M-130-499	34
0.375 (9.5) dia.	-	-	B-492	M-82-492	35
0.375 (9.5) dia.	-	-	B-499	M-82-499-OR-BK	35
0.375 (9.5) dia.	-	-	B-499	M-82-499-RD-BK	35
0.375 (9.5) dia.	-	-	B-499	M-82-499	35
0.375 (9.5) dia.	-	-	B-499	M-82-499-YL-BK	35
0.375 (9.5)	16 ft. (4.8 m)	-	B-499	MC-375-499	34
0.375 (9.5)	25 ft. (7.6 m)	-	B-595	MC-375-595-WT-BK	56
0.375 (9.5)	25 ft. (7.6 m)	-	B-595	MC-375-595-YL-BK	56
0.375 (9.5)	25 ft. (7.6 m)	-	B-595	MC-375-595-RD-WT	56
0.375 (9.5)	25 ft. (7.6 m)	-	B-461	MC-375-461	34
0.500 (12.7)	1.000 (25.4)	-	B-492	M-131-492	34
0.500 (12.7)	1.000 (25.4)	-	B-494	M-131-494-BL	34
0.500 (12.7)	1.000 (25.4)	-	B-494	M-131-494-GN	34
0.500 (12.7)	1.000 (25.4)	-	B-494	M-131-494-OR	34
0.500 (12.7)	1.000 (25.4)	-	B-494	M-131-494-PK	34
0.500 (12.7)	1.000 (25.4)	-	B-494	M-131-494-RD	34
0.500 (12.7)	1.000 (25.4)	-	B-494	M-131-494-YL	34
0.500 (12.7)	1.000 (25.4)	-	B-498	M-131-498	53
0.500 (12.7)	1.000 (25.4)	-	B-499	M-131-499	34
0.500 (12.7)	1.650 (41.9)	-	B-490	M-124-490	32
0.500 (12.7)	1.650 (41.9)	-	B-492	M-124-492	34
0.500 (12.7) dia.	-	-	B-492	M-83-492	35
0.500 (12.7) dia.	-	-	B-499	M-83-499-OR-BK	47
0.500 (12.7) dia.	-	-	B-499	M-83-499-RD-BK	47
0.500 (12.7) dia.	-	-	B-499	M-83-499	47
0.500 (12.7) dia.	-	-	B-499	M-83-499-YL-BK	47
0.500 (12.7)	16 ft. (4.8 m)	-	B-499	MC-500-499	34
0.500 (12.7)	25 ft. (7.6 m)	-	B-492	MC-500-492	46
0.500 (12.70)	25 ft. (7.6 m)	-	B-595	MC-500-595-WT-BK	56
0.500 (12.70)	25 ft. (7.6 m)	-	B-595	MC-500-595-YL-BK	56
0.500 (12.70)	25 ft. (7.6 m)	-	B-595	MC-500-595-WT-RD	56
0.500 (12.70)	25 ft. (7.6 m)	-	B-595	MC-500-595-BL-WT	57
0.500 (12.70)	25 ft. (7.6 m)	-	B-595	MC-500-595-GN-WT	57
0.500 (12.70)	25 ft. (7.6 m)	-	B-595	MC-500-595-RD-WT	57
0.500 (12.7)	25 ft. (7.6 m)	-	B-461	MC-500-461	34
0.750 (19.1)	0.500 (12.7)	-	B-498	M-11-498	31
0.750 (19.1)	16 ft. (4.8 m)	-	B-499	MC-750-499	34
0.750 (19.0)	25 ft. (7.6 m)	-	B-403	MC-750-403	55
0.750 (19.10)	25 ft. (7.6 m)	-	B-595	MC-750-595-WT-BK	57
0.750 (19.10)	25 ft. (7.6 m)	-	B-595	MC-750-595-YL-BK	57
0.750 (19.10)	25 ft. (7.6 m)	-	B-595	MC-750-595-WT-RD	57
0.750 (19.10)	25 ft. (7.6 m)	-	B-595	MC-750-595-BL-WT	57
0.750 (19.10)	25 ft. (7.6 m)	-	B-595	MC-750-595-GN-WT	57
0.750 (19.10)	25 ft. (7.6 m)	-	B-595	MC-750-595-RD-WT	57
0.600 (15.2)	1.625 (41.3)	0.375 (9.5) dia.	B-492	M-95-492	32
0.900 (22.9)	0.375 (9.5)	-	B-481	M-195-481	30
0.900 (22.9)	0.750 (19.1)	-	B-481	M-194-481	30
0.900 (22.9)	0.900 (22.9)	-	B-481	M-97-481	30
0.900 (22.9)	0.900 (22.9)	-	B-488	M-97-488	30
1.000 (25.4)	0.375 (9.5)	0.375 (9.5) dia.	B-492	M-118-492	35
1.000 (25.4)	0.375 (9.5)	0.375 (9.5) dia.	B-494	M-118-494-BL	35
1.000 (25.4)	0.375 (9.5)	0.375 (9.5) dia.	B-494	M-118-494-GN	35

BMP51 Printer self-laminating labels by size

Label width (A)	Label height (B)	Print-on area (C)	Material	Catalog #	Pg
1.125 (28.6)	0.375 (9.5)	0.750 (19.0)	B-461	M-123-461	51
1.125 (28.6)	0.500 (12.7)	0.750 (19.0)	B-461	M-122-461	51
1.650 (41.9)	0.500 (12.7)	0.750 (19.0)	B-461	M-124-461	51
1.750 (44.5)	0.500 (12.7)	0.750 (19.0)	B-461	M-72-461	39
1.750 (44.5)	1.000 (25.4)	0.375 (9.5)	B-427	M-133-427	33
1.750 (44.5)	1.000 (25.4)	1.000 (25.4)	B-461	M-125-461	39
1.750 (44.5)	1.500 (38.1)	0.500 (12.7)	B-427	M-134-427	33

Label width (A)	Label height (B)	Secondary label (C)	Material	Catalog #	Pg
1.000 (25.4)	0.375 (9.5)	0.375 (9.5) dia.	B-494	M-118-494-OR	35
1.000 (25.4)	0.375 (9.5)	0.375 (9.5) dia.	B-494	M-118-494-PK	35
1.000 (25.4)	0.375 (9.5)	0.375 (9.5) dia.	B-494	M-118-494-RD	35
1.000 (25.4)	0.375 (9.5)	0.375 (9.5) dia.	B-494	M-118-494-YL	35
1.000 (25.4)	0.375 (9.5)	0.375 (9.5) dia.	B-499	M-118-499	35
1.000 (25.4)	0.375 (9.5)	0.375 (9.5) dia.	B-7425	M-118-7425	35
1.000 (25.4)	0.500 (12.7)	-	B-432	M-17-432-CL-BK	26
1.000 (25.4)	0.500 (12.7)	0.438 (11.1) dia.	B-492	M-117-492	35
1.000 (25.4)	0.500 (12.7)	-	B-498	M-163-498	53
1.000 (25.4)	0.500 (12.7)	0.438 (11.1) dia.	B-499	M-117-499	35
1.000 (25.4)	0.500 (12.7)	0.375 (9.5) dia.	B-492	M-120-492	35
1.000 (25.4)	0.500 (12.7)	0.375 (9.5) dia.	B-494	M-120-494-BL	35
1.000 (25.4)	0.500 (12.7)	0.375 (9.5) dia.	B-494	M-120-494-GN	35
1.000 (25.4)	0.500 (12.7)	0.375 (9.5) dia.	B-494	M-120-494-OR	35
1.000 (25.4)	0.500 (12.7)	0.375 (9.5) dia.	B-494	M-120-494-PK	35
1.000 (25.4)	0.500 (12.7)	0.375 (9.5) dia.	B-494	M-120-494-RD	35
1.000 (25.4)	0.500 (12.7)	0.375 (9.5) dia.	B-494	M-120-494-YL	35
1.000 (25.4)	0.500 (12.7)	0.375 (9.5) dia.	B-499	M-120-499	35
1.000 (25.4)	0.900 (22.9)	-	B-492	M-156-492	35
1.000 (25.4)	1.000 (25.4)	-	B-494	M-49-494-BL	26
1.000 (25.4)	1.000 (25.4)	-	B-494	M-49-494-GN	26
1.000 (25.4)	1.000 (25.4)	-	B-494	M-49-494-OR	26
1.000 (25.4)	1.000 (25.4)	-	B-494	M-49-494-PK	26
1.000 (25.4)	1.000 (25.4)	-	B-494	M-49-494-RD	26
1.000 (25.4)	1.000 (25.4)	-	B-494	M-49-494-YL	26
1.000 (25.4)	1.250 (31.8)	-	B-498	M-103-498	53
1.000 (25.4)	1.900 (48.3)	-	B-432	M-78-432-CL-BK	43
1.000 (25.4)	1.900 (48.3)	-	B-492	M-78-492	43
1.000 (25.40)	1.900 (48.3)	-	B-498	M-128-498	53
1.000 (25.4)	1.900 (48.3)	-	B-499	M-128-499	43
1.000 (25.4)	1.900 (48.3)	-	B-499	MC-1000-499	43
1.000 (25.40)	25 ft. (7.6 m)	-	B-595	MC1-1000-595-WT-BK	57
1.000 (25.40)	25 ft. (7.6 m)	-	B-595	MC1-1000-595-YL-BK	57
1.000 (25.40)	25 ft. (7.6 m)	-	B-595	MC1-1000-595-WT-RD	57
1.000 (25.40)	25 ft. (7.6 m)	-	B-595	MC1-1000-595-BL-WT	57
1.000 (25.40)	25 ft. (7.6 m)	-	B-595	MC1-1000-595-GN-WT	57
1.000 (25.40)	25 ft. (7.6 m)	-	B-595	MC1-1000-595-RD-WT	57
1.000 (25.4)	25 ft. (7.6 m)	-	B-461	MC-1000-461	35
1.250 (31.8)	0.250 (6.4)	-	B-492	M-155-492	33
1.250 (31.8)	0.375 (9.5)	-	B-490	M-112-490	47
1.250 (31.8)	0.375 (9.5)	-	B-492	M-112-492	35
1.250 (31.8)	25 ft. (7.6 m)	-	B-461	MC-1250-461	35
1.500 (38.1)	1.000 (25.4)	-	B-7425	M-31-7425	35
1.500 (38.1)	1.250 (31.8)	-	B-499	M-132-499	32
1.500 (38.1)	2.200 (55.9)	-	B-490	M-113-490	35
1.500 (38.1)	3.750 (95.3)	-	B-490	M-114-490	43
1.500 (38.1)	25 ft. (7.6 m)	-	B-403	MC-1500-403	55
1.500 (38.10)	25 ft. (7.6 m)	-	B-595	MC-1500-595-WT-BK	57
1.500 (38.10)	25 ft. (7.6 m)	-	B-595	MC-1500-595-YL-BK	57
1.500 (38.10)	25 ft. (7.6 m)	-	B-595	MC-1500-595-WT-RD	57
1.500 (38.10)	25 ft. (7.6 m)	-	B-595	MC-1500-595-BL-WT	57
1.500 (38.10)	25 ft. (7.6 m)	-	B-595	MC-1500-595-GN-WT	57
1.500 (38.10)	25 ft. (7.6 m)	-	B-595	MC-1500-595-RD-WT	57
1.750 (44.5)	1.000 (25.4)	-	B-490	M-125-490	35
1.750 (44.5)	1.000 (25.4)	-	B-492	M-125-492	35
1.800 (45.7)	0.600 (15.2)	-	B-490	M-126-490	35
1.900 (48.3)	0.250 (6.4)	-	B-488	M-81-488	26

Labels by size index

Labels by size index

BradyPrinter M611 labels by size



Label width (A)	Label height (B)	Secondary label (C)	Material	Catalog #	Pg
0.250 (6.4)	2.000 (50.8)	-	B-7425	M61-2-7425	28
0.375 (9.5)	0.375 (9.5)	-	B-492	M61-82-492	36
0.375 (9.5)	0.375 (9.5)	-	B-499	M61-82-499-OR	36
0.375 (9.5)	0.375 (9.5)	-	B-499	M61-82-499-RD	36
0.375 (9.5) dia.	-	-	B-499	PTL-82-499	36
0.375 (9.5) dia.	-	-	B-499	PTL-82-499-YL	36
0.375 (9.5)	30 ft. (9.1 m)	-	B-499	M61C-375-499	35
0.375 (9.5)	50 ft. (15.2 m)	-	B-439	M61C-375-439	57
0.375 (9.5)	50 ft. (15.2 m)	-	B-439	M61C-375-439-YL	57
0.375 (9.5)	50 ft. (15.2 m)	-	B-461	M61C-375-461	35
0.500 (12.7)	0.500 (12.7)	-	B-492	M61-83-492	36
0.500 (12.7)	0.500 (12.7)	-	B-499	M61-83-499-OR	36
0.500 (12.7)	0.500 (12.7)	-	B-499	M61-83-499-RD	36
0.500 (12.7)	1.000 (25.4)	-	B-499	M61-47-499	35
0.500 (12.7)	1.650 (41.9)	-	B-490	M61-124-490	32
0.500 (12.7)	1.650 (41.9)	-	B-492	M61-124-492	35
0.500 (12.7) dia.	-	-	B-499	PTL-83-499	36
0.500 (12.7) dia.	-	-	B-499	PTL-83-499-YL	36
0.500 (12.7)	30 ft. (9.1 m)	-	B-499	M61C-500-499	35
0.500 (12.7)	50 ft. (15.2 m)	-	B-492	M61C-500-492	47
0.500 (12.7)	50 ft. (15.2 m)	-	B-595	M61C-500-595-BL	57
0.500 (12.7)	50 ft. (15.2 m)	-	B-595	M61C-500-595-GN	57
0.500 (12.7)	50 ft. (15.2 m)	-	B-595	M61C-500-595-RD	57
0.500 (12.7)	50 ft. (15.2 m)	-	B-595	M61C-500-595-WT	57
0.500 (12.7)	50 ft. (15.2 m)	-	B-595	M61C-500-595-YL	57
0.500 (12.7)	50 ft. (15.2 m)	-	B-461	M61C-500-461	35
0.600 (15.2)	1.625 (41.3)	0.375 (9.5) dia.	B-492	M61-205-492	36
0.750 (19.1)	0.750 (19.1)	-	B-481	M61-194-481	30
0.750 (19.1)	0.900 (22.9)	-	B-488	PTL-80-488	30
0.750 (19.1)	30 ft. (9.1 m)	-	B-499	M61C-750-499	35
0.750 (19.0)	50 ft. (15.2 m)	-	B-403	M61C-750-403	55
0.825 (21.0)	0.375 (9.5)	-	B-492	M61-130-492	35
0.825 (21.0)	0.375 (9.5)	-	B-499	M61-130-499	35
0.900 (22.9)	0.375 (9.5)	-	B-481	M61-195-481	30
0.900 (22.9)	0.500 (12.7)	0.500 (12.7) dia.	B-492	M61-117-492	35
0.900 (22.9)	0.500 (12.7)	0.500 (12.7) dia.	B-499	PTL-117-499	35
0.900 (22.9)	0.900 (22.9)	-	B-481	M61-97-481	30
0.900 (22.9)	0.900 (22.9)	-	B-488	PTL-97-488	30
0.900 (22.9)	1.000 (25.4)	-	B-492	M61-156-492	36
1.000 (25.4)	0.375 (9.5)	0.375 (9.5) dia.	B-492	M61-98-492	35
1.000 (25.4)	0.375 (9.5)	0.375 (9.5) dia.	B-494	M61-98-494-BL	35
1.000 (25.4)	0.375 (9.5)	0.375 (9.5) dia.	B-494	M61-98-494-GN	35
1.000 (25.4)	0.375 (9.5)	0.375 (9.5) dia.	B-494	M61-98-494-OR	35
1.000 (25.4)	0.375 (9.5)	0.375 (9.5) dia.	B-494	M61-98-494-PK	35
1.000 (25.4)	0.375 (9.5)	0.375 (9.5) dia.	B-494	M61-98-494-RD	35
1.000 (25.4)	0.375 (9.5)	0.375 (9.5) dia.	B-494	M61-98-494-YL	35
1.000 (25.4)	0.375 (9.5)	0.375 (9.5) dia.	B-499	PTL-98-499	36

Label width (A)	Label height (B)	Secondary label (C)	Material	Catalog #	Pg
1.000 (25.4)	0.375 (9.5)	0.375 (9.5) dia.	B-7425	M61-98-7425	35
1.000 (25.4)	0.500 (12.7)	-	B-492	M61-17-492	36
1.000 (25.4)	0.500 (12.7)	-	B-494	M61-17-494-BL	36
1.000 (25.4)	0.500 (12.7)	-	B-494	M61-17-494-GN	36
1.000 (25.4)	0.500 (12.7)	-	B-494	M61-17-494-OR	36
1.000 (25.4)	0.500 (12.7)	-	B-494	M61-17-494-PK	36
1.000 (25.4)	0.500 (12.7)	-	B-494	M61-17-494-RD	36
1.000 (25.4)	0.500 (12.7)	-	B-494	M61-17-494-YL	36
1.000 (25.4)	0.500 (12.7)	-	B-7425	M61-17-7425	36
1.000 (25.4)	0.500 (12.7)	0.375 (9.5) dia.	B-492	M61-120-492	36
1.000 (25.4)	0.500 (12.7)	0.375 (9.5) dia.	B-499	M61-120-499	36
1.000 (25.4)	0.750 (19.1)	-	B-7425	M61-18-7425	26
1.000 (25.4)	1.000 (25.4)	-	B-492	M61-19-492	36
1.000 (25.4)	1.000 (25.4)	-	B-494	M61-19-494-BL	26
1.000 (25.4)	1.000 (25.4)	-	B-494	M61-19-494-GN	26
1.000 (25.4)	1.000 (25.4)	-	B-494	M61-19-494-OR	26
1.000 (25.4)	1.000 (25.4)	-	B-494	M61-19-494-PK	26
1.000 (25.4)	1.000 (25.4)	-	B-494	M61-19-494-RD	26
1.000 (25.4)	1.000 (25.4)	-	B-494	M61-19-494-YL	26
1.000 (25.4)	1.250 (31.8)	-	B-498	M61-103-498	53
1.000 (25.4)	1.750 (44.5)	-	B-490	M61-203-490	36
1.000 (25.4)	1.750 (44.5)	-	B-492	M61-203-492	36
1.000 (25.4)	1.900 (48.3)	-	B-492	M61-204-492	44
1.000 (25.4)	30 ft. (9.1 m)	-	B-499	M61C-1000-499	44
1.000 (25.4)	50 ft. (15.2 m)	-	B-595	M61C-1000-595-BL	57
1.000 (25.4)	50 ft. (15.2 m)	-	B-595	M61C-1000-595-GN	57
1.000 (25.4)	50 ft. (15.2 m)	-	B-595	M61C-1000-595-RD	57
1.000 (25.4)	50 ft. (15.2 m)	-	B-595	M61C-1000-595-WT	57
1.000 (25.4)	50 ft. (15.2 m)	-	B-595	M61C-1000-595-YL	57
1.000 (25.4)	50 ft. (15.2 m)	-	B-461	M61C-1000-461	36
1.250 (31.8)	0.375 (9.5)	-	B-490	PTL-112-490	48
1.250 (31.7)	50 ft. (15.2 m)	-	B-461	M61C-1250-461	44
1.000 (30.0)	0.600 (15.2)	-	B-7425	M61-199-7425	36
1.500 (38.1)	1.000 (25.4)	-	B-7425	M61-31-7425	36
1.500 (38.1)	1.250 (31.8)	-	B-499	M61-132-499	32
1.500 (38.1)	2.200 (55.9)	-	B-490	PTL-113-490	36
1.500 (38.1)	3.875 (98.4)	-	B-490	PTL-114-490	44
1.500 (38.1)	50 ft. (15.2 m)	-	B-403	M61C-1500-403	55
1.625 (41.3)	0.600 (15.2)	-	B-490	PTL-76-490	32
1.800 (45.7)	0.600 (15.2)	-	B-490	M61-126-490	36
1.900 (48.3)	1.000 (25.4)	-	B-490	PTL-78-490	44
1.900 (48.3)	1.500 (38.1)	-	B-7425	M61-121-7425	32
2.000 (50.8)	50 ft. (15.2 m)	-	B-595	M61C-2000-595-BL	57
2.000 (50.8)	50 ft. (15.2 m)	-	B-595	M61C-2000-595-GN	57
2.000 (50.8)	50 ft. (15.2 m)	-	B-595	M61C-2000-595-RD	57
2.000 (50.8)	50 ft. (15.2 m)	-	B-595	M61C-2000-595-WT	57
2.000 (50.8)	50 ft. (15.2 m)	-	B-595	M61C-2000-595-YL	57

BradyPrinter M611 self-laminating labels by size

Label width (A)	Label height (B)	Print-on area (C)	Material	Catalog #	Pg
0.500 (12.7)	1.650 (41.9)	0.750 (19.1)	B-461	M61-202-461	51
0.600 (15.2)	1.800 (45.7)	1.000 (25.4)	B-461	M61-201-461	39
0.600 (15.2)	2.625 (66.7)	1.000 (25.4)	B-461	PTL-86-461	39
1.000 (25.4)	1.750 (44.5)	1.000 (25.4)	B-461	M61-203-461	39
1.000 (25.4)	2.625 (66.7)	1.000 (25.4)	B-461	PTL-75-461-YL	39
1.125 (28.6)	0.375 (9.5)	0.750 (19.1)	B-461	M61-123-461	51
1.125 (28.6)	0.500 (12.7)	0.750 (19.1)	B-461	M61-122-461	51
1.500 (38.1)	0.187 (4.75)	0.700 (17.8)	B-461	PTL-108-461	39

Label width (A)	Label height (B)	Print-on area (C)	Material	Catalog #	Pg
1.500 (38.1)	2.500 (63.5)	1.500 (38.1)	B-461	M61-157-461	39
1.625 (41.3)	0.600 (15.2)	1.000 (25.4)	B-461	PTL-76-461	51
1.800 (45.7)	0.375 (9.5)	0.750 (19.1)	B-461	PTL-73-461	39
1.800 (45.7)	0.500 (12.7)	0.750 (19.1)	B-461	PTL-72-461	39
2.200 (55.9)	0.500 (12.7)	0.750 (19.1)	B-461	PTL-74-461	39
2.625 (66.7)	1.000 (25.4)	1.000 (25.4)	B-461	PTL-75-461	39
3.500 (88.9)	1.500 (38.1)	1.250 (31.8)	B-461	PTL-77-461	39
3.750 (95.3)	0.750 (19.1)	1.000 (25.4)	B-461	PTL-79-461	39

BradyPrinter i3300 and BBP®33 Printer labels by size



Label width (A)	Label height (B)	Secondary label (C)	Material	Catalog #	Pg
0.375 (9.5)	-	-	B-492	B33-143-492	49
0.375 (9.5)	-	-	B-499	B33-143-499	49
0.500 (12.7)	0.750 (19.1)	-	B-490	B33-68-490	48
0.500 (12.7)	0.750 (19.1)	-	B-499	B33-68-499	36
0.500 (12.7) dia.	-	-	B-490	B33-144-490	37
0.500 (12.7) dia.	-	-	B-492	B33-144-492	37
0.500 (12.7) dia.	-	-	B-499	B33-144-499	37
0.500 (12.7)	100 ft. (30.5 m)	-	B-595	B30C-500-595-GN	57
0.500 (12.7)	100 ft. (30.5 m)	-	B-595	B30C-500-595-RD	57
0.500 (12.7)	100 ft. (30.5 m)	-	B-595	B30C-500-595-WT	57
0.500 (12.7)	100 ft. (30.5 m)	-	B-595	B30C-500-595-YL	57
0.750 (19.1)	0.250 (6.4)	-	B-498	B33-1-498	53
0.825 (21.0)	0.250 (6.4)	-	B-7425	B33-257-7425	36
0.900 (22.9)	0.375 (9.5)	-	B-481	B33-195-481	30
0.900 (22.9)	0.500 (12.7)	-	B-490	B33-125-490	48
0.900 (22.9)	0.500 (12.7)	-	B-499	B33-125-499	36
0.900 (22.9)	0.500 (12.7)	0.500 (12.7) dia.	B-490	B33-126-490	48
0.900 (22.9)	0.500 (12.7)	0.500 (12.7) dia.	B-499	B33-126-499	36
0.900 (22.9)	0.750 (19.1)	-	B-481	B33-194-481	30
0.900 (22.9)	0.900 (22.9)	-	B-481	B33-141-481	30
0.900 (22.9)	1.000 (25.4)	-	B-492	B33-180-492	37
1.000 (25.4)	0.375 (9.5)	-	B-492	B33-152-492	36
1.000 (25.4)	0.375 (9.5)	-	B-494	B33-152-494-BL	36
1.000 (25.4)	0.375 (9.5)	-	B-494	B33-152-494-GN	36
1.000 (25.4)	0.375 (9.5)	-	B-494	B33-152-494-OR	36
1.000 (25.4)	0.375 (9.5)	-	B-494	B33-152-494-PK	36
1.000 (25.4)	0.375 (9.5)	-	B-494	B33-152-494-RD	36
1.000 (25.4)	0.375 (9.5)	-	B-494	B33-152-494-YL	36
1.000 (25.4)	0.375 (9.5)	-	B-499	B33-152-499	28
1.000 (25.4)	0.375 (9.5)	-	B-7425	B33-152-7425	36
1.000 (25.4)	0.375 (9.5)	0.375 (9.5) dia.	B-492	B33-163-492	49
1.000 (25.4)	0.375 (9.5)	0.375 (9.5) dia.	B-499	B33-163-499	49
1.000 (25.4)	0.375 (9.5)	0.375 (9.5) dia.	B-7425	B33-249-7425	36
1.000 (25.4)	0.500 (12.7)	-	B-492	B33-59-492	36
1.000 (25.4)	0.500 (12.7)	-	B-498	B33-5-498	53
1.000 (25.4)	0.500 (12.7)	-	B-7425	B33-59-7425	28
1.000 (25.4)	0.500 (12.7)	0.375 (9.5) dia.	B-492	B33-181-492	37
1.000 (25.4)	0.500 (12.7)	0.375 (9.5) dia.	B-494	B33-250-494-BL	37
1.000 (25.4)	0.500 (12.7)	0.375 (9.5) dia.	B-494	B33-250-494-GN	37
1.000 (25.4)	0.500 (12.7)	0.375 (9.5) dia.	B-494	B33-250-494-OR	37
1.000 (25.4)	0.500 (12.7)	0.375 (9.5) dia.	B-494	B33-250-494-PK	37
1.000 (25.4)	0.500 (12.7)	0.375 (9.5) dia.	B-494	B33-250-494-RD	37
1.000 (25.4)	0.500 (12.7)	0.375 (9.5) dia.	B-494	B33-250-494-YL	37
1.000 (25.4)	0.500 (12.7)	0.440 (11.2) dia.	B-492	B33-182-492	37
1.000 (25.4)	0.500 (12.7)	0.440 (11.2) dia.	B-7425	B33-251-7425	27
1.000 (25.4)	0.750 (19.1)	-	B-498	B33-88-498	53
1.000 (25.4)	0.750 (19.1)	-	B-7425	B33-88-7425	27
1.000 (25.4)	1.000 (25.4)	-	B-492	B33-179-492	37
1.000 (25.4)	1.000 (25.4)	-	B-494	B33-179-494-BL	27
1.000 (25.4)	1.000 (25.4)	-	B-494	B33-179-494-GN	27
1.000 (25.4)	1.000 (25.4)	-	B-494	B33-179-494-OR	27

Label width (A)	Label height (B)	Secondary label (C)	Material	Catalog #	Pg
1.000 (25.4)	1.000 (25.4)	-	B-494	B33-179-494-PK	27
1.000 (25.4)	1.000 (25.4)	-	B-494	B33-179-494-RD	27
1.000 (25.4)	1.000 (25.4)	-	B-494	B33-179-494-YL	27
1.000 (25.4)	1.000 (25.4)	-	B-7425	B33-179-7425	27
1.000 (25.4)	1.437 (36.5)	-	B-498	B33-89-498	53
1.125 (28.6)	100 ft. (30.5 m)	-	B-403	B30C-1125-403-WT	55
1.125 (28.6)	100 ft. (30.5 m)	-	B-595	B30C-1125-595-GN	57
1.125 (28.6)	100 ft. (30.5 m)	-	B-595	B30C-1125-595-RD	57
1.125 (28.6)	100 ft. (30.5 m)	-	B-595	B30C-1125-595-WT	57
1.125 (28.6)	100 ft. (30.5 m)	-	B-595	B30C-1125-595-YL	57
1.250 (31.8)	0.375 (9.5)	-	B-490	B33-154-490	49
1.300 (33.0)	0.600 (15.2)	-	B-7425	B33-256-7425	27
1.500 (38.1)	0.500 (12.7)	-	B-7425	B33-257-7425	28
1.500 (38.1)	0.750 (19.1)	-	B-499	B33-136-499	32
1.500 (38.1)	0.750 (19.1)	-	B-7425	B33-136-7425	44
1.500 (38.1)	1.000 (25.4)	-	B-7425	B33-235-7425	27
1.500 (38.1)	1.000 (25.4)	-	B-7425-AC	B33-235-7425-AC	55
1.500 (38.1)	1.750 (44.5)	-	B-7425	B33-258-7425	41
1.625 (41.3)	0.600 (15.2)	-	B-490	B33-155-490	32
1.625 (41.3)	0.600 (15.2)	-	B-492	B33-155-492	37
1.625 (41.3)	0.600 (15.2)	0.375 (9.5) dia.	B-492	B33-190-492	32
1.900 (48.3)	1.000 (25.4)	-	B-490	B33-156-490	37
1.900 (48.3)	1.000 (25.4)	-	B-492	B33-156-492	44
2.000 (50.8)	0.250 (6.4)	-	B-490	B33-151-490	28
2.000 (50.8)	0.250 (6.4)	-	B-499	B33-151-499	28
2.000 (50.8)	0.250 (6.4)	-	B-7425	B33-151-7425	28
2.000 (50.8)	0.500 (12.7)	-	B-7425	B33-53-7425	28
2.000 (50.8)	1.000 (25.4)	-	B-424	B33-17-424	40
2.000 (50.8)	1.000 (25.4)	-	B-430	B33-17-430	27
2.000 (50.8)	1.000 (25.4)	-	B-499	B33-137-499	32
2.000 (50.8)	1.000 (25.4)	-	B-7425	B33-137-7425	27
2.000 (50.8)	1.000 (25.4)	-	B-7425-AC	B33-137-7425-AC	55
2.200 (55.9)	1.500 (38.1)	-	B-490	B33-157-490	37
2.250 (57.2)	100 ft. (30.5 m)	-	B-403	B30C-2250-403-WT	55
2.250 (57.2)	100 ft. (30.5 m)	-	B-595	B30C-2250-595-GN	57
2.250 (57.2)	100 ft. (30.5 m)	-	B-595	B30C-2250-595-RD	57
2.250 (57.2)	100 ft. (30.5 m)	-	B-595	B30C-2250-595-WT	57
2.250 (57.2)	100 ft. (30.5 m)	-	B-595	B30C-2250-595-YL	57
3.000 (76.2)	1.000 (25.4)	-	B-424	B33-18-424	40
3.000 (76.2)	1.000 (25.4)	-	B-430	B33-18-430	53
3.000 (76.2)	2.000 (50.8)	-	B-424	B33-19-424	53
3.000 (76.2)	100 ft. (30.5 m)	-	B-403	B30C-3000-403-WT	55
3.000 (76.2)	100 ft. (30.5 m)	-	B-595	B30C-3000-595-WT	57
3.000 (76.2)	100 ft. (30.5 m)	-	B-595	B30C-3000-595-YL	57
3.750 (95.3)	1.500 (38.1)	-	B-490	B33-158-490	44
4.000 (101.6)	100 ft. (30.5 m)	-	B-403	B30C-4000-403-WT	55
4.000 (101.6)	100 ft. (30.5 m)	-	B-595	B30C-4000-595-GN	57
4.000 (101.6)	100 ft. (30.5 m)	-	B-595	B30C-4000-595-RD	57
4.000 (101.6)	100 ft. (30.5 m)	-	B-595	B30C-4000-595-WT	57
4.000 (101.6)	100 ft. (30.5 m)	-	B-595	B30C-4000-595-YL	57

BradyPrinter i3300 and BBP®33 Printer self-laminating labels by size

Label width (A)	Label height (B)	Print-on area (C)	Material	Catalog #	Pg
1.000 (25.4)	5.250 (133.4)	1.500 (38.1)	B-461	B33-139-461	39
1.625 (41.3)	0.600 (15.2)	1.000 (25.4)	B-461	B33-133-461	51
1.750 (44.5)	1.000 (25.4)	0.375 (9.5)	B-427	B33-120-427	33
1.750 (44.5)	1.500 (38.1)	0.500 (12.7)	B-427	B33-121-427	33
1.800 (45.7)	0.375 (9.5)	0.750 (19.1)	B-461	B33-123-461	39

Label width (A)	Label height (B)	Print-on area (C)	Material	Catalog #	Pg
1.800 (45.7)	0.500 (12.7)	0.750 (19.1)	B-461	B33-122-461	39
2.200 (55.9)	0.500 (12.7)	0.750 (19.1)	B-461	B33-127-461	39
2.625 (66.7)	1.000 (25.4)	1.000 (25.4)	B-461	B33-276-461	39
3.500 (88.9)	1.500 (38.1)	1.250 (31.8)	B-461	B33-134-461	45

Labels by size index

BradyJet J2000 Printer labels by size



Label width (A)	Label height (B)	Secondary label (C)	Material	Catalog #	Pg
1.000 (25.4)	0.500 (12.7)	-	B-7425J	J20-59-7425J	29
1.000 (25.4)	0.500 (12.7)	0.440 (11.2) dia.	B-7425J	J20-251-7425J	27
1.000 (25.4)	0.750 (19.1)	-	B-7425J	J20-88-7425J	27
1.000 (25.4)	1.000 (25.4)	-	B-7425J	J20-179-7425J	27
1.125 (28.57)	100 ft. (30.5 m)	-	B-2569	J20C-1125-2569	57
1.125 (28.57)	100 ft. (30.5 m)	-	B-2595	J20C-1125-2595	57
1.125 (28.6)	100 ft. (30.5 m)	-	B-7425J	J20C-1125-7425J	32
1.300 (33.0)	0.600 (15.2)	-	B-7425J	J20-256-7425J	27
1.500 (38.1)	0.500 (12.7)	-	B-7425J	J20-37-7425J	29
1.500 (38.1)	0.750 (19.1)	-	B-7425J	J20-136-7425J	44
1.500 (38.1)	1.000 (25.4)	-	B-7425J	J20-235-7425J	27
1.500 (38.1)	1.750 (44.5)	-	B-7425J	J20-258-7425J	41
2.000 (50.8)	0.250 (6.4)	-	B-7425J	J20-151-7425J	29
2.000 (50.8)	0.500 (12.7)	-	B-7425J	J20-53-7425J	29
2.000 (50.8)	1.000 (25.4)	-	B-7425J	J20-137-7425J	27
2.170 (55.11)	3.350 (85.09)	-	B-2569	J20-260-2569	57

Label width (A)	Label height (B)	Secondary label (C)	Material	Catalog #	Pg
2.170 (55.11)	3.350 (85.09)	-	B-2595	J20-260-2595	57
2.170 (55.1)	3.350 (85.1)	-	B-7425J	J20-260-7425J	53
2.250 (57.15)	100 ft. (30.5 m)	-	B-2569	J20C-2250-2569	57
2.250 (57.15)	100 ft. (30.5 m)	-	B-2595	J20C-2250-2595	57
2.250 (57.2)	100 ft. (30.5 m)	-	B-7425J	J20C-2250-7425J	32
2.950 (74.93)	4.330 (109.98)	-	B-2569	J20-261-2569	57
2.950 (74.93)	4.330 (109.98)	-	B-2595	J20-261-2595	57
2.950 (74.9)	4.330 (110.0)	-	B-7425J	J20-261-7425J	53
3.000 (76.20)	1.000 (25.4)	-	B-2595	J20-18-2595	58
3.000 (76.2)	1.000 (25.4)	-	B-7425J	J20-18-7425J	44
4.000 (101.6)	100 ft. (30.5 m)	-	B-2569	J20C-4000-2569	57
4.000 (101.6)	100 ft. (30.5 m)	-	B-2595	J20C-4000-2595	57
4.000 (101.6)	100 ft. (30.5 m)	-	B-7425J	J20C-4000-7425J	58
4.130 (104.90)	6.100 (154.94)	-	B-2569	J20-262-2569	57
4.130 (104.90)	6.100 (154.94)	-	B-2595	J20-262-2595	57
4.130 (104.9)	6.100 (154.9)	-	B-7425J	J20-262-7425J	53

BBP®12 Printer 1.5-in. core labels by size



Label width (A)	Label height (B)	Secondary label (C)	Material	Catalog #	Pg
0.375 (9.5) dia.	-	-	B-492	THT-143-492-1.5-SC	45
0.500 (12.7) dia.	-	-	B-492	THT-144-492-1.5-SC	37
0.500 (12.7) dia.	-	-	B-499	THT-144-499-1.5-SC	37
0.825 (21.0)	0.250 (6.4)	-	B-7425	THT-257-7425-2-SC	37
0.900 (22.9)	0.500 (12.7)	-	B-499	THT-125-499-1.5-SC	37
0.900 (22.9)	0.750 (19.1)	-	B-481	THT-194-481-2.5-SC	30
0.900 (22.9)	0.750 (19.1)	-	B-488	THT-140-488-1.5-SC	27
0.900 (22.9)	0.900 (22.9)	-	B-481	THT-141-481-1.5-SC	30
0.900 (22.9)	0.900 (22.9)	-	B-481	THT-164-481-2.5-SC	30
0.900 (22.9)	0.900 (22.9)	-	B-481	THT-207-481-2.5-SC	30
0.900 (22.9)	0.900 (22.9)	-	B-488	THT-164-488-1-SC	30
1.000 (25.4)	0.375 (9.5)	-	B-492	THT-152-492-2.5-SC	37
1.000 (25.4)	0.375 (9.5)	-	B-499	THT-152-499-1.5-SC	37
1.000 (25.4)	0.375 (9.5)	-	B-7425	THT-152-7425-2-SC	37
1.000 (25.4)	0.375 (9.5)	0.375 (9.5) dia.	B-492	THT-249-492-1.5-SC	37
1.000 (25.4)	0.375 (9.5)	0.375 (9.5) dia.	B-7425	THT-249-7425-2-SC	37
1.000 (25.4)	0.375 (9.5)	0.375 (9.5) dia.	B-499	THT-249-499-1.5-SC	37
1.000 (25.4)	0.500 (12.7)	-	B-492	THT-59-492-2.5-SC	37
1.000 (25.4)	0.500 (12.7)	-	B-498	THT-59-498-2.5-SC	54
1.000 (25.4)	0.500 (12.7)	-	B-533	THT-59-533-2.5-SC	54
1.000 (25.4)	0.500 (12.7)	-	B-7425	THT-59-7425-2-SC	29
1.000 (25.4)	0.500 (12.7)	0.375 (9.5) dia.	B-492	THT-250-492-1.5-SC	37
1.000 (25.4)	0.500 (12.7)	0.440 (11.2) dia.	B-492	THT-251-492-1.5-SC	37
1.000 (25.4)	0.500 (12.7)	0.440 (11.2) dia.	B-7425	THT-251-7425-2-SC	27
1.000 (25.4)	0.750 (19.1)	-	B-492	THT-88-492-2.5-SC	37
1.000 (25.4)	0.750 (19.1)	-	B-7425	THT-88-7425-2-SC	27

Label width (A)	Label height (B)	Secondary label (C)	Material	Catalog #	Pg
1.000 (25.4)	1.000 (25.4)	-	B-492	THT-179-492-2.5-SC	37
1.000 (25.4)	1.000 (25.4)	-	B-7425	THT-179-7425-2-SC	27
1.050 (26.7)	0.250 (6.4)	-	B-498	THT-199-498-3-SC	54
1.250 (31.8)	0.375 (9.5)	-	B-490	THT-154-490-1.5-SC	37
1.300 (33.0)	0.600 (15.2)	-	B-7425	THT-256-7425-2-SC	27
1.500 (38.1)	0.500 (12.7)	-	B-7425	THT-37-7425-2-SC	29
1.500 (38.1)	0.750 (19.1)	-	B-449	THT-6-449-1.5-SC	54
1.500 (38.1)	0.750 (19.1)	-	B-499	THT-136-499-1.5-SC	41
1.500 (38.1)	0.750 (19.1)	-	B-7425	THT-136-7425-2-SC	45
1.500 (38.1)	1.000 (25.4)	-	B-7425	THT-235-7425-2-SC	27
1.625 (41.3)	0.600 (15.2)	-	B-490	THT-155-490-1.5-SC	33
1.750 (44.5)	1.500 (38.1)	-	B-7425	THT-258-7425-2-SC	41
1.900 (48.3)	1.000 (25.4)	-	B-490	THT-156-490-1.5-SC	33
1.900 (48.3)	1.000 (25.4)	-	B-492	THT-156-492-1.5-SC	33
2.000 (50.8)	0.250 (6.4)	-	B-499	THT-151-499-1.5-SC	42
2.000 (50.8)	0.250 (6.4)	-	B-7425	THT-151-7425-2-SC	29
2.000 (50.8)	0.500 (12.7)	-	B-7425	THT-53-7425-2-SC	29
2.000 (50.8)	1.000 (25.4)	-	B-424	THT-17-424-1.5-SC	40
2.000 (50.8)	1.000 (25.4)	-	B-498	THT-17-498-1.5-SC	54
2.000 (50.8)	1.000 (25.4)	-	B-499	THT-137-499-1.5-SC	33
2.000 (50.8)	1.000 (25.4)	-	B-533	THT-17-533-1.5-SC	54
2.000 (50.8)	1.000 (25.4)	-	B-7425	THT-137-7425-2-SC	27
2.200 (55.9)	1.500 (38.1)	-	B-490	THT-157-490-1.5-SC	37
2.500 (63.5)	0.200 (5.1)	-	B-481	THT-142-481-2.5-SC	30
3.000 (76.2)	1.000 (25.4)	-	B-424	THT-18-424-1.5-SC	40
3.750 (95.3)	1.500 (38.1)	-	B-490	THT-158-490-1.5-SC	45

BBP®12 Printer self-laminating labels by size

Label width (A)	Label height (B)	Print-on area (C)	Material	Catalog #	Pg
0.500 (12.7)	1.800 (45.7)	0.750 (19.1)	B-461	THT-122-461-1.5-SC	39
0.500 (12.7)	2.200 (55.9)	0.750 (19.1)	B-461	THT-127-461-1.5-SC	39
0.600 (15.2)	1.625 (41.3)	1.000 (25.4)	B-461	THT-133-461-1.5-SC	51

Label width (A)	Label height (B)	Print-on area (C)	Material	Catalog #	Pg
0.600 (15.2)	2.625 (66.7)	1.000 (25.4)	B-461	THT-153-461-1.5-SC	39
0.750 (19.1)	1.800 (45.7)	0.750 (19.1)	B-461	THT-123-461-1.5-SC	39

BradyPrinter i5100, i7100 and Brady IP Printer 3-in. core labels by size



Label width (A)	Label height (B)	Secondary label (C)	Material	Catalog #	Pg
0.375 (9.5) dia.	-	-	B-492	THT-143-492-3	50
0.375 (9.5) dia.	-	-	B-499	THT-143-499-3	50
0.500 (12.7) dia.	-	-	B-492	THT-144-492-3	38
0.500 (12.7) dia.	-	-	B-499	THT-144-499-3	38
0.500 (12.7)	100 ft. (30.5 m)	-	B-403	THT-500-403	55
0.750 (19.1)	0.250 (6.4)	-	B-498	THT-252-498-3	31
0.750 (19.1)	0.750 (19.1)	-	B-498	THT-254-498-3	31
0.750 (19.1)	1.500 (38.1)	-	B-427	THT-135-427-3-AW	45
0.750 (19.1)	100 ft. (30.5 m)	-	B-403	THT-750-403	55
0.825 (21.0)	0.250 (6.4)	-	B-7425	THT-257-7425-2	38
0.900 (22.9)	0.375 (9.5)	-	B-481	THT-195-481-3	30
0.900 (22.9)	0.500 (12.7)	-	B-499	THT-125-499-3	33
0.900 (22.9)	0.500 (12.7)	0.500 (12.7) dia.	B-499	THT-126-499-3	50
0.900 (22.9)	0.750 (19.1)	-	B-481	THT-194-481-3	30
0.900 (22.9)	0.750 (19.1)	-	B-488	THT-140-488-3	27
0.900 (22.9)	0.900 (22.9)	-	B-481	THT-141-481-3	30
0.900 (22.9)	0.900 (22.9)	-	B-481	THT-164-481-3	30
0.900 (22.9)	0.900 (22.9)	-	B-488	THT-141-488-3	30
0.900 (22.9)	0.900 (22.9)	-	B-488	THT-164-488-1	30
0.900 (22.9)	1.000 (25.4)	-	B-492	THT-180-492-3	38
1.000 (25.4)	0.375 (9.5)	-	B-492	THT-152-492-3	38
1.000 (25.4)	0.375 (9.5)	-	B-494	THT-152-494-BL	38
1.000 (25.4)	0.375 (9.5)	-	B-494	THT-152-494-GN	38
1.000 (25.4)	0.375 (9.5)	-	B-494	THT-152-494-OR	38
1.000 (25.4)	0.375 (9.5)	-	B-494	THT-152-494-PK	38
1.000 (25.4)	0.375 (9.5)	-	B-494	THT-152-494-RD	38
1.000 (25.4)	0.375 (9.5)	-	B-494	THT-152-494-YL	38
1.000 (25.4)	0.375 (9.5)	-	B-499	THT-152-499-3	38
1.000 (25.4)	0.375 (9.5)	-	B-7425	THT-152-7425-2	38
1.000 (25.4)	0.375 (9.5)	0.375 (9.5) dia.	B-492	THT-163-492-3	50
1.000 (25.4)	0.375 (9.5)	0.375 (9.5) dia.	B-499	THT-163-499-3	50
1.000 (25.4)	0.375 (9.5)	0.375 (9.5) dia.	B-7425	THT-249-7425-2	38
1.000 (25.4)	0.500 (12.7)	-	B-351	THT-5-351-10	41
1.000 (25.4)	0.500 (12.7)	-	B-488	THT-5-488-10	27
1.000 (25.4)	0.500 (12.7)	-	B-492	THT-59-492-10	38
1.000 (25.4)	0.500 (12.7)	-	B-498	THT-59-498-10	54
1.000 (25.4)	0.500 (12.7)	-	B-533	THT-59-533-10	54
1.000 (25.4)	0.500 (12.7)	-	B-7425	THT-59-7425-2	29
1.000 (25.4)	0.500 (12.7)	0.375 (9.5) dia.	B-492	THT-181-492-3	50
1.000 (25.4)	0.500 (12.7)	0.375 (9.5) dia.	B-494	THT-250-494-BL	38
1.000 (25.4)	0.500 (12.7)	0.375 (9.5) dia.	B-494	THT-250-494-GN	38
1.000 (25.4)	0.500 (12.7)	0.375 (9.5) dia.	B-494	THT-250-494-OR	38
1.000 (25.4)	0.500 (12.7)	0.375 (9.5) dia.	B-494	THT-250-494-PK	38
1.000 (25.4)	0.500 (12.7)	0.375 (9.5) dia.	B-494	THT-250-494-RD	38
1.000 (25.4)	0.500 (12.7)	0.375 (9.5) dia.	B-494	THT-250-494-YL	38
1.000 (25.4)	0.500 (12.7)	0.440 (11.2) dia.	B-492	THT-182-492-3	50
1.000 (25.4)	0.500 (12.7)	0.440 (11.2) dia.	B-7425	THT-251-7425-2	27
1.000 (25.4)	0.750 (19.1)	-	B-492	THT-88-492-5	38
1.000 (25.4)	0.750 (19.1)	-	B-7425	THT-88-7425-2	27
1.000 (25.4)	1.000 (25.4)	-	B-492	THT-179-492-3	38
1.000 (25.4)	1.000 (25.4)	-	B-494	THT-179-494-BL	27
1.000 (25.4)	1.000 (25.4)	-	B-494	THT-179-494-GN	27
1.000 (25.4)	1.000 (25.4)	-	B-494	THT-179-494-OR	27
1.000 (25.4)	1.000 (25.4)	-	B-494	THT-179-494-PK	27

Label width (A)	Label height (B)	Secondary label (C)	Material	Catalog #	Pg
1.000 (25.4)	1.000 (25.4)	-	B-494	THT-179-494-RD	27
1.000 (25.4)	1.000 (25.4)	-	B-494	THT-179-494-YL	27
1.000 (25.4)	1.000 (25.4)	-	B-7425	THT-179-7425-2	27
1.000 (25.4)	1.500 (38.1)	-	B-7425-AC	THT-235-7425-2-AC	55
1.000 (25.4)	100 ft. (30.5 m)	-	B-403	THT-1000-403	55
1.050 (26.7)	0.250 (6.4)	-	B-498	THT-199-498-3	54
1.125 (28.6)	300 ft. (91.4 m)	-	B-403	THT-1125-403-WT	55
1.250 (31.8)	0.375 (9.5)	-	B-490	THT-154-490-3	38
1.300 (33.0)	0.600 (15.2)	-	B-7425	THT-256-7425-2	27
1.500 (38.1)	0.500 (12.7)	-	B-351	THT-37-351-10	41
1.500 (38.1)	0.500 (12.7)	-	B-7425	THT-37-7425-2	29
1.500 (38.1)	0.750 (19.1)	-	B-449	THT-6-449-3	54
1.500 (38.1)	0.750 (19.1)	-	B-499	THT-156-499-3	41
1.500 (38.1)	0.750 (19.1)	-	B-7425	THT-136-7425-2	45
1.500 (38.1)	1.000 (25.4)	-	B-7425	THT-235-7425-2	27
1.500 (38.1)	100 ft. (30.5 m)	-	B-403	THT-1500-403	55
1.625 (41.3)	0.600 (15.2)	-	B-490	THT-155-490-3	33
1.625 (41.3)	0.600 (15.2)	-	B-492	THT-155-492-3	38
1.625 (41.3)	0.600 (15.2)	0.375 (9.5) dia.	B-492	THT-190-492-3	33
1.750 (44.5)	1.000 (25.4)	-	B-499	THT-185-499-3	33
1.750 (44.5)	1.500 (38.1)	-	B-7425	THT-258-7425-2	41
1.900 (48.3)	1.000 (25.4)	-	B-490	THT-156-490-3	33
1.900 (48.3)	1.000 (25.4)	-	B-492	THT-156-492-3	33
2.000 (50.8)	0.250 (6.4)	-	B-488	THT-15-488-2.5	27
2.000 (50.8)	0.250 (6.4)	-	B-499	THT-151-499-3	42
2.000 (50.8)	0.250 (6.4)	-	B-7425	THT-151-7425-2	29
2.000 (50.8)	0.500 (12.7)	-	B-351	THT-53-351-3	41
2.000 (50.8)	0.500 (12.7)	-	B-7425	THT-53-7425-2	29
2.000 (50.8)	1.000 (25.4)	-	B-351	THT-17-351-3	41
2.000 (50.8)	1.000 (25.4)	-	B-424	THT-17-424-3	40
2.000 (50.8)	1.000 (25.4)	-	B-430	THT-17-430-3	27
2.000 (50.8)	1.000 (25.4)	-	B-498	THT-17-498-3	54
2.000 (50.8)	1.000 (25.4)	-	B-499	THT-171-499-3	33
2.000 (50.8)	1.000 (25.4)	-	B-533	THT-17-533-3	54
2.000 (50.8)	1.000 (25.4)	-	B-7425	THT-137-7425-2	27
2.000 (50.8)	1.000 (25.4)	-	B-7425-AC	THT-137-7425-2-AC	55
2.000 (50.8)	1.000 (25.4)	-	B-813	DT-17-813-3	40
2.000 (50.8)	1.000 (25.4)	-	B-427	THT-137-427-3-AW	45
2.000 (50.8)	100 ft. (30.5 m)	-	B-403	THT-2000-403	55
2.200 (55.9)	1.500 (38.1)	-	B-490	THT-157-490-3	38
2.250 (57.2)	300 ft. (91.4 m)	-	B-403	THT-2250-403-WT	55
2.500 (63.5)	0.200 (5.1)	-	B-481	THT-142-481-2.5	30
3.000 (76.2)	1.000 (25.4)	-	B-351	THT-18-351-3	41
3.000 (76.2)	1.000 (25.4)	-	B-424	THT-18-424-3	40
3.000 (76.2)	1.000 (25.4)	-	B-430	THT-18-430-3	54
3.000 (76.2)	1.000 (25.4)	-	B-813	DT-18-813-3	40
3.000 (76.2)	2.000 (50.8)	-	B-351	THT-19-351-1	41
3.000 (76.2)	2.000 (50.8)	-	B-430	THT-19-430-1	54
3.000 (76.2)	2.000 (50.8)	-	B-499	THT-19-499-1	33
3.000 (76.2)	300 ft. (91.4 m)	-	B-403	THT-3000-403-WT	55
3.000 (76.2)	300 ft. (91.4 m)	-	B-430	THT-21-430	54
3.750 (95.3)	1.500 (38.1)	-	B-490	THT-158-490-3	45
4.000 (101.6)	300 ft. (91.4 m)	-	B-403	THT-4000-403-WT	55

BradyPrinter i5100, i7100 and Brady IP Printers self-laminating labels by size

Label width (A)	Label height (B)	Print-on area (C)	Material	Catalog #	Pg
0.375 (9.5)	1.800 (45.7)	0.750 (19.1)	B-461	THT-123-461-3	39
0.500 (12.7)	1.800 (45.7)	0.750 (19.1)	B-461	THT-122-461-3	39
0.500 (12.7)	2.200 (55.9)	0.750 (19.1)	B-461	THT-127-461-3	39
0.600 (15.2)	1.625 (41.3)	1.000 (25.4)	B-461	THT-133-461-3	51
0.600 (15.2)	2.625 (66.7)	1.000 (25.4)	B-461	THT-153-461-3	39

Label width (A)	Label height (B)	Print-on area (C)	Material	Catalog #	Pg
1.000 (25.4)	1.000 (25.4)	0.375 (9.5)	B-427	THT-116-427-5	33
1.000 (25.4)	1.750 (44.5)	1.000 (25.4)	B-461	THT-183-461-1	39
1.000 (25.4)	2.625 (66.7)	1.000 (25.4)	B-461	THT-131-461-3	39
1.000 (25.4)	5.250 (133.4)	1.500 (38.1)	B-461	THT-139-461-2	45
1.500 (38.1)	3.500 (88.9)	1.250 (31.8)	B-461	THT-134-461-3	45

Index

Catalog #	Pg	Catalog #	Pg	Catalog #	Pg	Catalog #	Pg	Catalog #	Pg
21800	61	85098	61	122826	61	149044	17	B33-17-430	28
22082	61	85165	61	126052	61	149045	17	B33-17-430	53
22311	61	85167	61	126157	61	149048	17	B33-18-424	40
22312	61	85168	61	126158	61	149050	17	B33-18-424	53
22323	61	85553	61	126159	61	149051	17	B33-18-430	53
22324	61	85554	61	126160	61	149053	17	B33-19-424	53
22355	61	87782	61	126161	61	149054	17	B33-37-7425	28
22356	61	88051	60	126624	61	149055	17	B33-37-7425	32
22362	61	88543	61	126625	61	149056	17	B33-37-7425	37
22579	61	88544	61	126626	61	149069	17	B33-53-7425	28
22580	61	88627	61	126627	61	149208	17	B33-53-7425	42
25194	61	88628	61	126628	61	149452	16	B33-59-492	36
25195	61	89158	61	126629	61	149453	16	B33-59-7425	28
25222	61	89159	61	126630	61	149455	16	B33-59-7425	32
25223	61	89170	61	126631	61	149456	16	B33-59-7425	36
25441	61	90933	61	126632	61	149552	14	B33-59-7425	49
25610	61	91447	60	126633	61	150507	20	B33-68-490	48
25611	61	104506	60	126634	61	150508	20	B33-68-490	53
25652	61	118815	59	126635	61	150639	14	B33-68-499	36
25653	61	118816	59	126636	61	150642	14	B33-68-499	48
25654	61	118817	59	126637	61	150645	14	B33-68-499	53
25655	61	118821	59	126638	61	AD15100	62	B33-88-498	53
25742	61	118822	59	126639	61	AD30100	62	B33-88-7425	27
25743	61	118823	59	126640	61	B30-18-7569-CLP4C	58	B33-88-7425	37
25898	61	118827	59	126641	61	B30-260-7569-CLP1	58	B33-88-7425	49
26575	61	118828	59	126642	61	B30-260-7569-CLP2	58	B33-89-498	53
30670	61	118829	59	126643	61	B30-260-7569-CLP3A	58	B33-120-427	33
30671	61	118833	59	126644	61	B30-260-7569-CLP3B	58	B33-121-427	33
30673	61	118834	59	126645	61	B30-260-7569-CLP4A	58	B33-122-461	39
30674	61	118835	59	126646	61	B30-260-7569-CLP4B	58	B33-122-461	51
40645	61	118839	59	126647	61	B30-261-7569-CLP1	58	B33-123-461	39
40646	61	118840	59	126648	61	B30-261-7569-CLP2	58	B33-123-461	51
40876	61	118841	59	126649	61	B30-261-7569-CLP3A	58	B33-125-490	48
40887	61	118845	59	126650	61	B30-261-7569-CLP3B	58	B33-125-499	36
40888	61	118846	59	126651	61	B30-261-7569-CLP4A	58	B33-125-499	48
40925	61	118847	59	126652	61	B30-261-7569-CLP4B	58	B33-126-490	48
41143	61	118851	59	126653	61	B30-262-7569-CLP1	58	B33-126-499	36
42771	61	118852	59	126666	61	B30-262-7569-CLP2	58	B33-126-499	48
42772	61	118853	59	126667	61	B30-262-7569-CLP3A	58	B33-127-461	39
42799	61	118857	59	126668	61	B30-262-7569-CLP3B	58	B33-127-461	51
42800	61	118858	59	126669	61	B30-262-7569-CLP4A	58	B33-133-461	51
43017	61	118859	59	126670	61	B30-262-7569-CLP4B	58	B33-134-461	45
43018	61	118863	59	126671	61	B30-R4300	14	B33-136-499	32
43187	61	118864	59	126672	61	B30-R6000	14	B33-136-499	44
43188	61	118865	59	126673	61	B30-R6200	14	B33-136-499	53
43230	61	118869	59	126674	61	B30-R6400	14	B33-136-7425	44
43231	61	118870	59	126675	61	B30-R10000	14	B33-137-499	32
43232	61	118871	59	126676	61	B30C-500-595-GN	57	B33-137-499	44
43475	61	118875	59	126677	61	B30C-500-595-RD	57	B33-137-499	53
46753	61	118876	59	126678	61	B30C-500-595-WT	57	B33-137-7425	27
47056	61	118877	59	126679	61	B30C-500-595-YL	57	B33-137-7425	40
47094	61	118881	59	126680	61	B30C-1125-403-WT	55	B33-137-7425	40
62550	61	118882	59	126681	61	B30C-1125-595-GN	57	B33-137-7425	41
62812	61	118883	59	126682	61	B30C-1125-595-RD	57	B33-137-7425	44
69097	61	118887	59	126683	61	B30C-1125-595-WT	57	B33-137-7425	53
69768	61	118888	59	129156	61	B30C-1125-595-YL	57	B33-137-7425-AC	55
69813	61	118889	59	129157	61	B30C-2250-403-WT	55	B33-139-461	39
70240	61	118893	59	129158	61	B30C-2250-595-GN	57	B33-141-481	30
70258	61	118894	59	129159	61	B30C-2250-595-RD	57	B33-143-492	49
70350	61	118895	59	129160	61	B30C-2250-595-WT	57	B33-143-499	49
70582	61	118899	59	129161	61	B30C-2250-595-YL	57	B33-144-490	37
73436	61	118900	59	132434	60	B30C-3000-403-WT	55	B33-144-490	49
74455	61	118901	59	133080	20	B30C-3000-595-WT	57	B33-144-499	37
75639	61	118905	59	133081	20	B30C-3000-595-YL	57	B33-144-492	37
83753	61	118906	59	133082	20	B30C-4000-403-WT	55	B33-144-492	49
83928	61	118907	59	133083	20	B30C-4000-595-GN	57	B33-144-499	49
84054	61	118911	59	133084	20	B30C-4000-595-RD	57	B33-151-490	28
84055	61	118912	59	133085	20	B30C-4000-595-WT	57	B33-151-490	42
84351	61	118913	59	133086	20	B30C-4000-595-YL	57	B33-151-499	28
84352	61	118917	59	133087	20	B31-CCT	14	B33-151-499	42
84368	61	118918	59	133088	20	B31-WIPER	14	B33-151-7425	28
84397	61	118919	59	142267	11	B33-1-498	53	B33-151-7425	42
84398	61	120994	14	142268	11	B33-5-498	53	B33-152-492	36
84463	61	122406	61	143561	20	B33-17-424	40	B33-152-492	49
84469	61	122519	61	145586	60	B33-17-424	40	B33-152-494-BL	36
84973	61	122650	61	145587	60	B33-17-424	53	B33-152-494-BL	49
85097	61	122820	61	149043	17	B33-17-430	27	B33-152-494-GN	36

Catalog #	Pg	Catalog #	Pg	Catalog #	Pg	Catalog #	Pg	Catalog #	Pg
B33-152-494-GN	49	B33-235-7425-AC	55	IP-R4300	16	J20-CMY	15	M-82-499-OR-BK	47
B33-152-494-OR	36	B33-249-7425	36	IP-R4302	16	J20-CUTTER	15	M-82-499-RD-BK	35
B33-152-494-OR	49	B33-249-7425	49	IP-R4306	16	J20-PAD	15	M-82-499-RD-BK	47
B33-152-494-PK	36	B33-250-494-BL	37	IP-R4307	16	J20C-1125-2569	57	M-82-499-YL-BK	35
B33-152-494-PK	49	B33-250-494-BL	49	IP-R6000	16	J20C-1125-2595	57	M-82-499-YL-BK	47
B33-152-494-RD	36	B33-250-494-GN	37	IP-R6002	16	J20C-1125-2595	58	M-83-492	35
B33-152-494-RD	49	B33-250-494-GN	49	IP-R6200	16	J20C-1125-7425J	32	M-83-492	47
B33-152-494-YL	36	B33-250-494-OR	37	IP-R6207	16	J20C-1125-7425J	44	M-83-499	47
B33-152-494-YL	49	B33-250-494-OR	49	IP-R6400	16	J20C-1125-7425J	53	M-83-499-OR-BK	47
B33-152-499	28	B33-250-494-PK	37	IP-R6402	16	J20C-1125-7425J	57	M-83-499-RD-BK	47
B33-152-499	42	B33-250-494-PK	49	J20-18-2595	58	J20C-1125-7425J	58	M-83-499-YL-BK	47
B33-152-499	49	B33-250-494-RD	37	J20-18-7425J	44	J20C-2250-2569	57	M-86-461	39
B33-152-499	53	B33-250-494-RD	49	J20-18-7425J	53	J20C-2250-2595	57	M-86-461	51
B33-152-7425	36	B33-250-494-YL	37	J20-18-7425J	58	J20C-2250-2595	58	M-95-492	32
B33-152-7425	49	B33-250-494-YL	49	J20-37-7425J	29	J20C-2250-7425J	32	M-95-492	43
B33-154-490	49	B33-251-7425	27	J20-37-7425J	32	J20C-2250-7425J	44	M-97-481	30
B33-155-490	32	B33-251-7425	37	J20-37-7425J	37	J20C-2250-7425J	53	M-97-488	30
B33-155-490	44	B33-251-7425	49	J20-53-7425J	29	J20C-2250-7425J	57	M-103-498	53
B33-155-492	37	B33-256-7425	27	J20-53-7425J	42	J20C-2250-7425J	58	M-112-490	47
B33-156-490	37	B33-256-7425	37	J20-59-7425J	29	J20C-4000-2569	57	M-112-492	35
B33-156-490	44	B33-256-7425	49	J20-59-7425J	32	J20C-4000-2595	57	M-112-492	47
B33-156-492	44	B33-257-7425	36	J20-59-7425J	37	J20C-4000-2595	58	M-113-490	35
B33-157-490	37	B33-257-7425	48	J20-59-7425J	49	J20C-4000-7425J	58	M-113-490	43
B33-157-490	44	B33-258-7425	41	J20-88-7425J	27	J2000	15	M-114-490	43
B33-158-490	44	B33-258-7425	44	J20-88-7425J	37	J2000-BWSLAB	15	M-117-492	35
B33-163-492	49	B33-258-7425	53	J20-88-7425J	49	LH625E	61	M-117-492	47
B33-163-499	49	B33-276-461	39	J20-136-7425J	44	M-11-498	31	M-117-499	35
B33-179-492	37	B33-276-461	45	J20-137-7425J	27	M-17-432-CL-BK	26	M-117-499	47
B33-179-492	44	BBP11-U	13	J20-137-7425J	40	M-17-432-CL-BK	28	M-118-492	35
B33-179-492	49	BBP12-NA-PWID	13	J20-137-7425J	49	M-17-432-CL-BK	53	M-118-492	47
B33-179-494-BL	27	BBP12-US	13	J20-137-7425J	41	M-31-7425	35	M-118-494-BL	35
B33-179-494-BL	37	BH2023	60	J20-137-7425J	44	M-49-494-BL	26	M-118-494-BL	47
B33-179-494-BL	44	BMP-UBP	11	J20-137-7425J	53	M-49-494-BL	35	M-118-494-GN	35
B33-179-494-BL	49	BMP-UBP-AC	11	J20-151-7425J	29	M-49-494-BL	43	M-118-494-GN	47
B33-179-494-BL	53	BMP21-LAB	10	J20-151-7425J	42	M-49-494-BL	47	M-118-494-OR	35
B33-179-494-GN	27	BMP21-LAB-BATTCOV	10	J20-152-7425J	37	M-49-494-BL	53	M-118-494-OR	47
B33-179-494-GN	37	BMP21-LAB-KIT1	10	J20-152-7425J	49	M-49-494-GN	26	M-118-494-PK	35
B33-179-494-GN	44	BMP21-MAGNET	10	J20-179-7425J	27	M-49-494-GN	35	M-118-494-PK	47
B33-179-494-GN	49	BMP21-PLUS-BATT	10	J20-179-7425J	37	M-49-494-GN	43	M-118-494-RD	35
B33-179-494-GN	53	BMP21-PLUS-HC	10	J20-179-7425J	41	M-49-494-GN	47	M-118-494-RD	47
B33-179-494-OR	27	BMP51	11	J20-179-7425J	44	M-49-494-GN	53	M-118-494-YL	35
B33-179-494-OR	37	BMP51-KIT-AC	11	J20-179-7425J	49	M-49-494-OR	26	M-118-494-YL	47
B33-179-494-OR	44	BMP51-KIT-LS	11	J20-179-7425J	53	M-49-494-OR	35	M-118-499	35
B33-179-494-OR	49	BMP-HC-2	12	J20-235-7425J	27	M-49-494-OR	43	M-118-499	47
B33-179-494-OR	53	BMP-MAGNET-2	12	J20-235-7425J	32	M-49-494-OR	47	M-118-7425	35
B33-179-494-PK	37	BMP-SC-1	11	J20-249-7425J	37	M-49-494-OR	53	M-118-7425	47
B33-179-494-PK	44	BMP-SC-2	12	J20-249-7425J	49	M-49-494-PK	26	M-119-461	45
B33-179-494-PK	49	BMP-STRAP-1	12	J20-251-7425J	27	M-49-494-PK	35	M-120-492	35
B33-179-494-PK	53	BSM100	62	J20-251-7425J	37	M-49-494-PK	43	M-120-492	47
B33-179-494-RD	27	BSM300	62	J20-251-7425J	49	M-49-494-PK	47	M-120-494-BL	35
B33-179-494-RD	37	BSM3650	62	J20-256-7425J	27	M-49-494-PK	53	M-120-494-BL	47
B33-179-494-RD	44	BWRK-GHS-DWN	19	J20-256-7425J	37	M-49-494-RD	26	M-120-494-GN	35
B33-179-494-RD	49	BWS-LABS-EM	18	J20-256-7425J	49	M-49-494-RD	35	M-120-494-GN	47
B33-179-494-RD	53	BWS-PPS-CD	19	J20-257-7425J	37	M-49-494-RD	43	M-120-494-OR	35
B33-179-494-YL	27	BWS-PPS-CD-VOL	19	J20-257-7425J	49	M-49-494-RD	47	M-120-494-OR	47
B33-179-494-YL	37	BWS-PPS-EM	19	J20-258-7425J	41	M-49-494-RD	53	M-120-494-PK	35
B33-179-494-YL	44	BWS-PPS-EM-VOL	19	J20-258-7425J	44	M-49-494-YL	26	M-120-494-PK	47
B33-179-494-YL	49	BWS-SPS-CD	19	J20-258-7425J	53	M-49-494-YL	35	M-120-494-RD	35
B33-179-494-YL	53	BWS-SPS-EM	19	J20-260-2569	57	M-49-494-YL	43	M-120-494-RD	47
B33-179-7425	27	CH15DP	62	J20-260-2595	57	M-49-494-YL	47	M-120-494-YL	35
B33-179-7425	37	CH100	62	J20-260-2595	58	M-49-494-YL	53	M-120-494-YL	47
B33-179-7425	41	CH1212	62	J20-260-7425J	53	M-72-461	39	M-120-499	35
B33-179-7425	44	CH1212-50	62	J20-260-7425J	57	M-72-461	51	M-120-499	47
B33-179-7425	49	CH30DP	62	J20-260-7425J	58	M-74-461	39	M-122-461	51
B33-179-7425	53	CH200	62	J20-261-2569	57	M-74-461	51	M-123-461	51
B33-180-492	37	CH303	62	J20-261-2595	57	M-75-461	39	M-124-461	51
B33-181-492	37	DISPENSER	21	J20-261-2595	58	M-78-432-CL-BK	43	M-124-490	32
B33-181-492	49	DSP-403KIT	21	J20-261-7425J	53	M-78-492	43	M-124-492	34
B33-182-492	37	DT-17-813-3	40	J20-261-7425J	57	M-81-488	26	M-124-492	47
B33-182-492	49	DT-17-813-3	40	J20-261-7425J	58	M-81-488	28	M-125-461	39
B33-190-492	32	DT-18-813-3	40	J20-262-2569	57	M-81-488	42	M-125-490	35
B33-195-481	30	HAZ99	62	J20-262-2595	57	M-82-492	35	M-125-492	35
B33-194-481	30	HAZ124	62	J20-262-2595	58	M-82-492	47	M-126-461	39
B33-235-7425	27	HAZ412	62	J20-262-7425J	53	M-82-499	35	M-126-461	51
B33-235-7425	32	HAZ1818	62	J20-262-7425J	57	M-82-499	47	M-126-490	35
		HAZ1818-2	62	J20-262-7425J	58	M-82-499-OR-BK	35	M-128-498	53

Index

Catalog #	Pg	Catalog #	Pg	Catalog #	Pg	Catalog #	Pg	Catalog #	Pg
M-128-499	43	M21-500-7425	28	M61-19-494-RD	53	M61-205-492	36	MC1-1000-595-RD-WT	57
M-129-461	45	M21-500-7425	32	M61-19-494-YL	26	M61-R4310	12	MC1-1000-595-WT-BK	57
M-130-492	34	M21-500-7425	34	M61-19-494-YL	36	M61-R6010	12	MC1-1000-595-WT-RD	57
M-130-499	34	M21-500-7425	46	M61-19-494-YL	44	M61-R6210	12	MC1-1000-595-YL-BK	57
M-130-499	46	M21-750-403	55	M61-19-494-YL	48	M61-R6410	12	MC-1250-461	35
M-130-499	53	M21-750-430	26	M61-19-494-YL	53	M61C-375-439	57	MC-1250-461	43
M-131-492	34	M21-750-430	46	M61-31-7425	36	M61C-375-439-YL	57	MC-1500-403	55
M-131-492	47	M21-750-430	52	M61-47-499	35	M61C-375-461	35	MC-1500-595-BL-WT	57
M-131-494-BL	34	M21-750-430-WT-CL	52	M61-47-499	48	M61C-375-461	47	MC-1500-595-GN-WT	57
M-131-494-BL	47	M21-750-461	34	M61-82-492	36	M61C-375-499	35	MC-1500-595-RD-WT	57
M-131-494-GN	34	M21-750-461	46	M61-82-492	48	M61C-375-499	47	MC-1500-595-WT-BK	57
M-131-494-GN	47	M21-750-488	30	M61-82-499-OR	36	M61C-500-461	35	MC-1500-595-WT-RD	57
M-131-494-OR	34	M21-750-499	32	M61-82-499-OR	48	M61C-500-461	47	MC-1500-595-YL-BK	57
M-131-494-OR	47	M21-750-499	43	M61-82-499-RD	36	M61C-500-492	47	MRO15P	62
M-131-494-PK	34	M21-750-499	46	M61-82-499-RD	48	M61C-500-499	35	MRO100	62
M-131-494-PK	47	M21-750-499	52	M61-83-492	36	M61C-500-499	47	MRO300	62
M-131-494-RD	34	M21-750-595-BL	56	M61-83-492	48	M61C-500-595-BL	57	MRO350DP	62
M-131-494-RD	47	M21-750-595-GN	56	M61-83-499-OR	36	M61C-500-595-GN	57	MXP-C20	62
M-131-494-YL	34	M21-750-595-RD	56	M61-83-499-OR	48	M61C-500-595-RD	57	PCK-5	12
M-131-494-YL	47	M21-750-595-WT	56	M61-83-499-RD	36	M61C-500-595-WT	57	PCK-6	14
M-131-498	53	M21-750-595-YL	56	M61-83-499-RD	48	M61C-500-595-YL	57	PCK-6	16
M-131-499	34	M21-750-7425	26	M61-97-481	30	M61C-750-403	55	PCK-6	17
M-131-499	47	M21-750-7425	32	M61-98-492	35	M61C-750-499	35	PTL-72-461	39
M-131-499	53	M21-750-7425	34	M61-98-492	48	M61C-750-499	48	PTL-72-461	51
M-132-499	32	M21-750-7425	46	M61-98-494-BL	35	M61C-1000-461	36	PTL-73-461	39
M-132-499	53	M50-BATT-TRAY	11	M61-98-494-BL	48	M61C-1000-461	44	PTL-73-461	51
M-133-427	33	M50-CABLE	11	M61-98-494-GN	35	M61C-1000-499	44	PTL-74-461	39
M-134-427	33	M50-HC	11	M61-98-494-GN	48	M61C-1000-499	48	PTL-74-461	51
M-135-427-AW	33	M61-2-7425	28	M61-98-494-OR	35	M61C-1000-595-BL	57	PTL-75-461	39
M-156-492	35	M61-2-7425	42	M61-98-494-OR	48	M61C-1000-595-BL	57	PTL-75-461-YL	39
M-156-492	43	M61-17-492	36	M61-98-494-PK	35	M61C-1000-595-RD	57	PTL-76-461	51
M-156-492	47	M61-17-492	48	M61-98-494-PK	48	M61C-1000-595-WT	57	PTL-76-490	32
M-156-492	47	M61-17-494-BL	36	M61-98-494-RD	35	M61C-1000-595-YL	57	PTL-77-461	39
M-163-498	53	M61-17-494-BL	48	M61-98-494-RD	48	M61C-1250-461	44	PTL-77-461	45
M-194-481	30	M61-17-494-GN	36	M61-98-494-YL	35	M61C-1500-403	55	PTL-78-490	44
M-195-481	30	M61-17-494-GN	48	M61-98-494-YL	48	M61C-2000-595-BL	57	PTL-79-461	39
M-AC-110937	10	M61-17-494-OR	36	M61-98-7425	35	M61C-2000-595-GN	57	PTL-79-461	51
M-AC-143110	11	M61-17-494-OR	48	M61-98-7425	48	M61C-2000-595-RD	57	PTL-80-488	30
M-AC-143110	12	M61-17-494-PK	36	M61-103-498	53	M61C-2000-595-WT	57	PTL-82-499	36
M21-250-430	28	M61-17-494-PK	48	M61-117-492	35	M61C-2000-595-YL	57	PTL-82-499	48
M21-250-430	46	M61-17-494-RD	36	M61-117-492	48	M611-AM-BT-W	12	PTL-82-499-YL	36
M21-250-430	52	M61-17-494-RD	48	M61-120-492	36	M611-BC-KIT1	12	PTL-82-499-YL	48
M21-250-430-WT-CL	52	M61-17-494-YL	36	M61-120-492	48	MC-250-595-WT-BK	56	PTL-83-499	36
M21-250-595-WT	56	M61-17-494-YL	48	M61-120-499	36	MC-250-595-YL-BK	56	PTL-83-499	48
M21-250-595-YL	56	M61-17-7425	36	M61-120-499	48	MC-375-461	34	PTL-83-499-YL	36
M21-375-430	28	M61-17-7425	48	M61-121-7425	32	MC-375-499	34	PTL-83-499-YL	48
M21-375-430	46	M61-18-7425	26	M61-121-7425	53	MC-375-499	46	PTL-86-461	39
M21-375-430	52	M61-18-7425	36	M61-122-461	51	MC-375-499	46	PTL-97-488	30
M21-375-430-WT-CL	52	M61-18-7425	48	M61-123-461	51	MC-375-595-RD-WT	56	PTL-98-499	36
M21-375-461	34	M61-18-7425	53	M61-124-490	32	MC-375-595-WT-BK	56	PTL-98-499	48
M21-375-461	46	M61-19-492	36	M61-124-492	35	MC-375-595-YL-BK	56	PTL-108-461	39
M21-375-488	30	M61-19-492	44	M61-124-492	48	MC-500-461	34	PTL-108-461	51
M21-375-499	46	M61-19-492	48	M61-126-490	36	MC-500-461	46	PTL-112-490	48
M21-375-595-BL	56	M61-19-494-BL	26	M61-130-492	35	MC-500-499	34	PTL-113-490	44
M21-375-595-GN	56	M61-19-494-BL	36	M61-130-499	35	MC-500-499	46	PTL-114-490	44
M21-375-595-RD	56	M61-19-494-BL	44	M61-130-499	47	MC-500-595-BL-WT	57	PTL-117-499	35
M21-375-595-WT	56	M61-19-494-BL	48	M61-130-499	53	MC-500-595-GN-WT	57	PTL-117-499	48
M21-375-595-YL	56	M61-19-494-BL	53	M61-132-499	32	MC-500-595-RD-WT	57	PVC18	62
M21-375-7425	28	M61-19-494-GN	26	M61-132-499	53	MC-500-595-WT-BK	56	PVC24	62
M21-375-7425	34	M61-19-494-GN	36	M61-156-492	36	MC-500-595-WT-RD	56	R4300	17
M21-375-7425	46	M61-19-494-GN	44	M61-156-492	44	MC-500-595-YL-BK	56	R4302	17
M21-500-430	26	M61-19-494-GN	48	M61-156-492	48	MC-750-403	55	R4306	17
M21-500-430	28	M61-19-494-GN	53	M61-157-461	39	MC-750-499	34	R4307	17
M21-500-430	46	M61-19-494-OR	26	M61-194-481	30	MC-750-499	47	R4313	13
M21-500-430	52	M61-19-494-OR	36	M61-195-481	30	MC-750-595-BL-WT	57	R4913	13
M21-500-430-WT-CL	52	M61-19-494-OR	44	M61-199-7425	36	MC-750-595-GN-WT	57	R6000	17
M21-500-461	34	M61-19-494-OR	48	M61-199-7425	48	MC-750-595-RD-WT	57	R6002	17
M21-500-461	46	M61-19-494-OR	53	M61-201-461	39	MC-750-595-WT-BK	57	R6013	13
M21-500-488	30	M61-19-494-PK	26	M61-201-461	51	MC-750-595-WT-RD	57	R6200	17
M21-500-499	32	M61-19-494-PK	36	M61-202-461	51	MC-750-595-YL-BK	57	R6207	17
M21-500-499	46	M61-19-494-PK	44	M61-203-461	39	MC-1000-461	35	R6213	13
M21-500-595-BL	56	M61-19-494-PK	48	M61-203-461	51	MC-1000-461	43	R6400	17
M21-500-595-GN	56	M61-19-494-PK	53	M61-203-490	36	MC-1000-499	43	R6402	17
M21-500-595-RD	56	M61-19-494-RD	26	M61-203-490	44	MC-1000-499	47	R6413	13
M21-500-595-WT	56	M61-19-494-RD	36	M61-203-492	36	MC1-1000-595-BL-WT	57	SC-Cabinet	63
M21-500-595-YL	56	M61-19-494-RD	44	M61-204-492	44	MC1-1000-595-GN-WT	57	SC-Cart	63
M21-500-7425	26	M61-19-494-RD	48						

Catalog #	Pg	Catalog #	Pg	Catalog #	Pg	Catalog #	Pg	Catalog #	Pg
SKH-20	63	THT-88-7425-2	38	THT-151-499-1.5-SC	42	THT-179-494-GN	50	THT-251-7425-2	38
SKH-ATK	63	THT-88-7425-2	50	THT-151-499-3	42	THT-179-494-GN	54	THT-251-7425-2	50
SKH-BKT	63	THT-88-7425-2-SC	27	THT-151-7425-2	29	THT-179-494-OR	27	THT-251-7425-2-SC	27
SKH-CFB	63	THT-88-7425-2-SC	37	THT-151-7425-2	42	THT-179-494-OR	38	THT-251-7425-2-SC	37
SKH-PP	63	THT-88-7425-2-SC	50	THT-151-7425-2-SC	29	THT-179-494-OR	45	THT-251-7425-2-SC	49
SKH-SRP	63	THT-116-427-5	33	THT-151-7425-2-SC	42	THT-179-494-OR	50	THT-252-498-3	31
SPCUN1019	62	THT-122-461-1.5-SC	39	THT-152-492-2.5-SC	37	THT-179-494-OR	54	THT-254-498-3	31
SPCUN1113	62	THT-122-461-1.5-SC	51	THT-152-492-2.5-SC	49	THT-179-494-PK	27	THT-256-7425-2	27
SPCUN1113-50	62	THT-122-461-3	39	THT-152-492-3	38	THT-179-494-PK	38	THT-256-7425-2	38
SPCUN1719	62	THT-122-461-3	51	THT-152-492-3	50	THT-179-494-PK	45	THT-256-7425-2	50
SPCUN3800	62	THT-123-461-1.5-SC	39	THT-152-494-BL	38	THT-179-494-PK	50	THT-256-7425-2-SC	27
SRPH75P	62	THT-123-461-1.5-SC	51	THT-152-494-BL	50	THT-179-494-PK	54	THT-256-7425-2-SC	37
SRPH100	62	THT-123-461-3	39	THT-152-494-GN	38	THT-179-494-RD	27	THT-256-7425-2-SC	49
SRPH200	62	THT-123-461-3	51	THT-152-494-GN	50	THT-179-494-RD	38	THT-257-7425-2	38
THT-5-351-10	41	THT-125-499-1.5-SC	37	THT-152-494-OR	38	THT-179-494-RD	45	THT-257-7425-2	50
THT-5-488-10	27	THT-125-499-1.5-SC	49	THT-152-494-OR	50	THT-179-494-RD	50	THT-257-7425-2-SC	37
THT-5-488-10	29	THT-125-499-3	33	THT-152-494-PK	38	THT-179-494-RD	54	THT-257-7425-2-SC	49
THT-5-488-10	54	THT-125-499-3	38	THT-152-494-PK	50	THT-179-494-YL	27	THT-257-7425-2	41
THT-6-449-1.5-SC	54	THT-125-499-3	50	THT-152-494-RD	38	THT-179-494-YL	38	THT-258-7425-2	45
THT-6-449-3	54	THT-126-499-3	50	THT-152-494-RD	50	THT-179-494-YL	45	THT-258-7425-2	54
THT-15-488-2.5	27	THT-127-461-1.5-SC	39	THT-152-494-YL	38	THT-179-494-YL	50	THT-258-7425-2-SC	41
THT-15-488-2.5	29	THT-127-461-1.5-SC	51	THT-152-494-YL	50	THT-179-494-YL	54	THT-258-7425-2-SC	45
THT-15-488-2.5	42	THT-127-461-3	39	THT-152-499-1.5-SC	37	THT-179-7425-2	27	THT-258-7425-2-SC	54
THT-15-488-2.5	54	THT-127-461-3	51	THT-152-499-1.5-SC	49	THT-179-7425-2	38	THT-500-403	55
THT-17-351-3	41	THT-131-461-3	39	THT-152-499-3	38	THT-179-7425-2	41	THT-750-403	55
THT-17-424-1.5-SC	40	THT-133-461-1.5-SC	51	THT-152-499-3	50	THT-179-7425-2	45	THT-1000-403	55
THT-17-424-1.5-SC	54	THT-133-461-3	51	THT-152-7425-2	38	THT-179-7425-2	50	THT-1125-403-WT	55
THT-17-424-3	40	THT-134-461-3	45	THT-152-7425-2	50	THT-179-7425-2	54	THT-1500-403	55
THT-17-424-3	54	THT-135-427-3-AW	45	THT-152-7425-2-SC	37	THT-179-7425-2-SC	27	THT-2000-403	55
THT-17-430-3	27	THT-136-499-1.5-SC	41	THT-152-7425-2-SC	49	THT-179-7425-2-SC	37	THT-2250-403-WT	55
THT-17-430-3	29	THT-136-499-1.5-SC	54	THT-153-461-1.5-SC	39	THT-179-7425-2-SC	41	THT-3000-403-WT	55
THT-17-430-3	54	THT-136-499-3	41	THT-153-461-1.5-SC	51	THT-179-7425-2-SC	45	THT-4000-403-WT	55
THT-17-498-1.5-SC	54	THT-136-499-3	54	THT-153-461-3	39	THT-179-7425-2-SC	50		
THT-17-498-3	54	THT-136-7425-2	45	THT-153-461-3	51	THT-179-7425-2-SC	54		
THT-17-533-1.5-SC	54	THT-136-7425-2-SC	45	THT-154-490-1.5-SC	37	THT-180-492-3	38		
THT-17-533-3	54	THT-137-427-3-AW	45	THT-154-490-3	38	THT-180-492-3	50		
THT-18-351-3	41	THT-137-499-1.5-SC	33	THT-155-490-1.5-SC	33	THT-181-492-3	50		
THT-18-424-1.5-SC	40	THT-137-499-1.5-SC	41	THT-155-490-1.5-SC	45	THT-182-492-3	50		
THT-18-424-1.5-SC	54	THT-137-499-1.5-SC	54	THT-155-490-3	33	THT-183-461-1	39		
THT-18-424-3	40	THT-137-499-3	33	THT-155-490-3	45	THT-185-499-3	33		
THT-18-424-3	54	THT-137-499-3	41	THT-155-492-3	38	THT-190-492-3	33		
THT-18-430-3	54	THT-137-499-3	54	THT-156-490-1.5-SC	33	THT-194-481-2.5-SC	30		
THT-19-351-1	41	THT-137-7425-2	27	THT-156-490-1.5-SC	37	THT-194-481-3	30		
THT-19-430-1	54	THT-137-7425-2	40	THT-156-490-1.5-SC	45	THT-195-481-3	30		
THT-19-499-1	33	THT-137-7425-2	40	THT-156-490-3	33	THT-199-498-3	54		
THT-21-430	54	THT-137-7425-2	41	THT-156-490-3	38	THT-199-498-3-SC	54		
THT-37-351-10	41	THT-137-7425-2	45	THT-156-490-3	45	THT-207-481-2.5-SC	30		
THT-37-7425-2	29	THT-137-7425-2	54	THT-156-492-1.5-SC	33	THT-235-7425-2	27		
THT-37-7425-2	33	THT-137-7425-2-AC	55	THT-156-492-1.5-SC	45	THT-235-7425-2	33		
THT-37-7425-2	38	THT-137-7425-2-SC	27	THT-156-492-3	33	THT-235-7425-2-AC	55		
THT-37-7425-2-SC	29	THT-137-7425-2-SC	40	THT-156-492-3	45	THT-235-7425-2-SC	27		
THT-37-7425-2-SC	32	THT-137-7425-2-SC	40	THT-157-490-1.5-SC	37	THT-235-7425-2-SC	32		
THT-37-7425-2-SC	37	THT-137-7425-2-SC	41	THT-157-490-1.5-SC	45	THT-249-492-1.5-SC	37		
THT-53-351-3	41	THT-137-7425-2-SC	45	THT-157-490-3	38	THT-249-492-1.5-SC	49		
THT-53-7425-2	29	THT-137-7425-2-SC	54	THT-157-490-3	45	THT-249-499-1.5-SC	37		
THT-53-7425-2	42	THT-139-461-2	45	THT-158-490-1.5-SC	45	THT-249-499-1.5-SC	49		
THT-53-7425-2-SC	29	THT-140-488-1.5-SC	27	THT-158-490-3	45	THT-249-7425-2	38		
THT-53-7425-2-SC	42	THT-140-488-1.5-SC	54	THT-163-492-3	50	THT-249-7425-2	50		
THT-59-492-10	38	THT-140-488-3	27	THT-163-499-3	38	THT-249-7425-2-SC	37		
THT-59-492-10	50	THT-140-488-3	54	THT-164-481-2.5-SC	30	THT-249-7425-2-SC	49		
THT-59-492-2.5-SC	37	THT-141-481-1.5-SC	30	THT-164-481-3	30	THT-250-492-1.5-SC	37		
THT-59-492-2.5-SC	49	THT-141-481-3	30	THT-164-488-1	30	THT-250-492-1.5-SC	49		
THT-59-498-10	54	THT-141-488-3	30	THT-164-488-1-SC	30	THT-250-494-BL	38		
THT-59-498-2.5-SC	54	THT-142-481-2.5	30	THT-179-492-2.5-SC	37	THT-250-494-BL	50		
THT-59-533-10	54	THT-142-481-2.5-SC	30	THT-179-492-2.5-SC	45	THT-250-494-GN	38		
THT-59-533-2.5-SC	54	THT-143-492-1.5-SC	45	THT-179-492-2.5-SC	50	THT-250-494-GN	50		
THT-59-7425-2	29	THT-143-492-1.5-SC	50	THT-179-492-3	38	THT-250-494-OR	38		
THT-59-7425-2	33	THT-143-492-3	50	THT-179-492-3	45	THT-250-494-OR	50		
THT-59-7425-2	38	THT-144-492-1.5-SC	37	THT-179-492-3	50	THT-250-494-PK	38		
THT-59-7425-2	50	THT-144-492-1.5-SC	50	THT-179-494-BL	27	THT-250-494-PK	50		
THT-59-7425-2-SC	29	THT-144-492-3	38	THT-179-494-BL	38	THT-250-494-RD	38		
THT-59-7425-2-SC	32	THT-144-492-3	50	THT-179-494-BL	45	THT-250-494-RD	50		
THT-59-7425-2-SC	37	THT-143-499-3	50	THT-179-494-BL	50	THT-250-494-YL	38		
THT-59-7425-2-SC	49	THT-144-499-1.5-SC	37	THT-179-494-BL	54	THT-250-494-YL	50		
THT-88-492-2.5-SC	37	THT-144-499-1.5-SC	50	THT-179-494-GN	27	THT-251-492-1.5-SC	37		
THT-88-492-5	38	THT-144-499-3	38	THT-179-494-GN	38	THT-251-492-1.5-SC	49		
THT-88-7425-2	27	THT-144-499-3	50	THT-179-494-GN	45	THT-251-7425-2	27		

A

Absorbents	62
Autoclave-indicating labels	55

B

Barcode scanners	20
BBP12 Printer	13
BBP12 Ribbons	13
Benchtop printers	13
Biohazard tape	60
Blood tube labels	40
BMP21-LAB Printer	10
BMP51 Printer	11
Bottle, flask and dish labels	52
Brady Workstation Software	18
BradyJet J2000 Printer	15
BradyPrinter i3300 Printer	14
BradyPrinter i3300 Ribbons	14
BradyPrinter i5100 Printer	16
BradyPrinter i5100 Ribbons	16
BradyPrinter i7100 Printer	17
BradyPrinter i7100 Ribbons	17
BradyPrinter M611 Printer	12
BradyPrinter M611 Ribbons	12

C

Cell culture flask labels	26
Centrifuge and conical tube labels	43
Centrifuge and conical tube self-lam labels	45
Cryo box and freezer rack labels	32
Cryo cane labels	33
Cryo tube labels	34
Cryo tube self-laminating labels	39
Cryopreservation labels	32

D

Dispenser	21
-----------------	----

G

General lab labels	52
GHS labels	58
Glassware labels	52

H

Histology labels	30
------------------------	----

L

Lab diagram guide.....	6
Label dispenser.....	21
Label material comparison guide.....	25
Lean lab / visual workflow labels.....	56

M

Material comparison guide.....	25
Microbiology labels	26
Microcentrifuge labels.....	46
Microcentrifuge self-laminating labels	51
Microplate labels	42
Microscope slide labels	30

P

PCR tube labels.....	46
PCR tube self-laminating labels	51
Petri dish labels	28
Portable printers.....	10
Printer comparison guide.....	8

S

Sample bag labels	41
Sample collection and processing labels	40
Scanners	20
Signs.....	61
Slide labels	30
Software.....	18
Specimen cup labels.....	40
Spill kits	63
Straw and rod labels	33

T

Tamper-resistant seals.....	41
Tissue cassette labels.....	31

V

Visual workflow labels	56
------------------------------	----

W

Water-dissolvable labels	55
Well labels	42

We identify and protect premises, products and people.

For more information or to learn more about our innovative products, visit BradyID.com or call 1-888-272-3946.



Y4178816 MS4178816

USA
Customer Service: 1-888-272-3946
Inside Sales: 1-888-311-0775
BradyID.com

Canada
Customer Service: 1-800-263-6179
BradyCanada.ca

Mexico
Customer Service: 1-800-262-7777
Inside Sales: 1-800-262-7777 ext 177
BradyLatinAmerica.com

