

# Leonardo Mechanography Kiosk

Evaluation results  
in less than  
2 minutes

# Father & Son




# Complete Self Guided System




- In addition to the two Leonardo Mechanograph GRFP PRO force plates, the Kiosk Bundle contains a 23" touch panel terminal and an RFID reader, which allow easy use by the end user without the need for additional staff.

# The Modular System for Sports and Sports Research

- Flexible and mobile system for use in sports, studies or for advanced users. Force plate for spatially resolved measurement of ground reaction forces. Two separated platforms for analysis of asymmetries in movement. Versatile use due to two separate, portable force plates that can be arranged in different orders



# Choice of 6 Different Test Parameters



**NOVO TEC MEDICAL**  
Get fit. Stay fit.

## Assessment Selection

Assessment List

- Fitness (Push-Ups) (BT, Jmp, Tap, PuT)**
- Fitness (Hoppings) (BT, Jmp, Tap, 1LH)
- Athletic (Hoppings) (Jmp, Tap, PuT, 1LH)
- Athletic (GF & Hoppings) (Jmp, Tap, GF, 1LH)
- Best Ager (Grip Force) (BT, Jmp, Tap, GF)
- Geriatric (GF & CRT) (BT, Jmp, GF, CRT)

**Start ▶**

BT  
Jmp  
Tap  
PuT

TOP 20

# Leonardo Mechanography Athletic (Imp, Tap, GF, 1LH)




2 Leg jump  
max power test



Tapping  
20 seconds




Grip Force




One Legged Hopping  
6 seconds straight knee

# Results After 2 Minutes of Testing



# Time Line Trending Results!





For the Athlete and the Researcher!

StimDesigns LLC  
32683 Coast Ridge Dr  
Carmel, CA 93923  
831-224-0982

