


# Accurate Detection of Bacterial Vaginosis and Trichomonas in Minutes



Diagnosing Bacterial Vaginosis and *Trichomonas* can be difficult, time consuming and labor intensive. OSOM® BVBLUE® and OSOM® Trichomonas Tests transform the diagnostic process by minimizing the need to search for clue cells or use a microscope to find a motile organism.

## OSOM® BVBLUE® Test

- Detects bacterial vaginosis without any equipment
- Sensitivity: 92.8%; Specificity: 98.0% compared to Gram staining & Nugent's score\*
- 1-minute hands-on time, 10 minute read time

## OSOM® Trichomonas Test

- Detects *Trichomonas vaginalis* without any equipment
- Sensitivity: 96%; Specificity: 95.0% compared to Wet Mount\*
- 1-minute hands-on time, 10 minute read time



Visit [osomtests.com](http://osomtests.com) for more information or call 800-332-1042

\*See package insert for details

**SEKISUI**  
DIAGNOSTICS  
Because every result matters™